

A METHOD AND APPARATUS FOR AUTOMATIC FOCUSING  
OF A CONFOCAL LASER MICROSCOPE

5

CROSS-REFERENCE TO ATTACHED APPENDICES

Appendices A-E (2 sheets of 158 frames) that are attached herewith, are parts of the present disclosure and are incorporated herein by reference in their entirety. Appendices A and B are a listing of routines that are used in microprocessors of coarse and fine Z controllers in two alternative embodiments of a microscope system. Appendix C is a listing of routines that are used in a programmer to program a Programmable Logic Device (PLD). Appendix D is a listing of components used in one embodiment of a microscope system. Appendix E is a listing of routines that implement computations for various auto-focus operations.

Appendix F is a listing of host workstation software that sets up various parameters for use in the area scan method. Appendix G, entitled "Ultrapointe Model 1000 Laser Imaging System Users' Manual" attached herewith is part of present disclosure and is incorporated herein by reference in its entirety.

NOTICE OF COPYRIGHT RIGHTS

A portion of the disclosure of this patent document contains material which is subject to copyright protection. The copyright owner has no objection to the facsimile reproduction by anyone of the patent document or the patent disclosure, as it appears in the Patent and Trademark Office patent files or records, but otherwise reserves all copyright rights whatsoever.

FIELD OF THE INVENTION

The present invention relates generally to a method and apparatus for automatic focusing. More particularly, the present invention relates to automatic focusing of a confocal microscope.

BACKGROUND OF THE INVENTION

Certain prior art auto-focus microscopes generate an electronic focus signal, that is related to the distance of the target from an objective lens of the microscope, to determine the position at which the target is in focus (henceforth "focus position").

In other auto-focus microscopes, a white light camera is used to detect a maximum contrast in the camera picture to indicate when the target is in focus. Measuring the picture contrast requires a large amount of computer processing. This method is therefore slow and complicated. Additionally, this method fails if the target has no contrast (i.e., if the target is flat and has a uniform color, or is mirror-like).

The above method can fail depending on the amount of noise in the electronic focus signal, such as the noise introduced by vibration of the target during relative movement between the target and the microscope. Such noise can result in detection of a false focus position when a microprocessor uses a low power objective lens. A low power objective lens is a lens that has a low power, such as 1.5X or 5X, as opposed to a high power objective lens that has a high power, such as 100X or 200X.

Use of a low power objective lens can result in a large depth of focus, for example, greater than half the range of motion of the target with respect to the microscope. A large depth of focus combined with noise can result in an electronic focus signal that does not have a single distinct well defined peak that is

essential for proper operation of certain conventional microscopes. In such a case, a microscope can use a false peak to move the target to a position at a distance from the focus position, thereby requiring manual focusing.

In certain conventional microscopes, an electronic focus signal that is used for automatic focusing is generated from an illumination spot held stationary on the target. The illumination spot could be broad band white light or monochromatic laser light. Local height variations in the surface of the target can cause the focused condition for the entire field of view to be less than optimal if the spot at which focusing is performed is significantly higher or lower than the elevation of the remainder of the target's surface in the field of view. Therefore in targets, such as semiconductor substrates, that consist of highly regular orthogonal lines, for example, row and column bit lines of a DRAM, the microscope can produce different results depending on whether a portion of a DRAM line or a portion of the substrate between two DRAM lines lies within the spot. If the spot is on a portion of the substrate in a trench, the microscope focuses on the bottom of the trench. If the spot is on the top of a levee, the microscope focuses on the top of the levee.

However, it is desirable that the microscope focus on the same layer all the time, for example, on the top layer, regardless of what layer happened to be within the spot, and regardless of the materials on the top layer. In the multi-layer structure of a substrate, the different layers can be of different materials with different reflectivities. In a typical wafer, the top layer is a transparent passivation layer, and a predetermined layer that is typically desired is under the top layer. There is a need to focus on such a

predetermined layer even if the predetermined layer does not have the highest reflectivity.

Speed, accuracy and repeatability are additional desirable characteristics for focusing a microscope.

5 The time required to automatically focus a microscope determines the number of targets that can be viewed in a given time period, or, equivalently, determines the amount of time required to view a given number of targets, and thereby determines the cost associated  
10 with viewing each sample.

Certain prior art microscopes utilize auto-focus optics that are separate from the imaging optics of the microscope. In such microscopes, any "drift" between the auto-focus optics and the imaging optics results in  
15 loss of auto-focusing accuracy and repeatability.

It is therefore desirable to have a microscope that can automatically focus on a predetermined layer of a target quickly, accurately and repeatably.

## 20 SUMMARY OF THE INVENTION

In accordance with the present invention a confocal microscope system (henceforth "microscope system") uses a median point method for a coarse autofocus operation and an area scan method for a fine  
25 autofocus operation to provide accurate, repeatable and high-speed automatic focusing on any predetermined layer of a target. During an auto-focus operation, such as a coarse auto-focus operation or a fine auto-focus operation, the microscope system moves a target  
30 with respect to an objective lens, while an electronic focus signal is measured. In a coarse auto-focus operation, the value of the electronic focus signal is recorded periodically at large distances between a target's elevation, as compared to smaller  
35 distances used in a fine auto-focus operation.

The electronic focus signal generated by the



microscope system is a novel signal that has a magnitude proportional to the amount of light reflected by the target. The electronic focus signal reaches a maximum (sometimes referred to as "peak") when the target is in focus. The electronic focus signal is used to control an auto-focus operation in which the target is moved to its focus position (an elevation at which the target is in focus) with respect to the objective lens of a confocal microscope. Such an electronic focus signal provides a very precise indication of the focus position of the target. In one embodiment, the electronic focus signal is also used to generate an image of the target. Since the same signal (electronic focus signal) is used for both focusing and imaging, a good image signal is obtained after the auto-focus operation. The narrow peak also allows the microscope system to discriminate between semi-transparent layers of the target.

In a coarse auto-focus operation, when using an objective lens having low power (e.g. 1.5x or 5x), the microscope system moves the target in a start-up move to a predetermined starting position so that the direction of the focus position is in a first predetermined direction (e.g. positive Z-axis direction) with respect to the starting position. Then the microscope system moves the target in a first pass (a movement during which electronic focus signal values are recorded), in the first direction through a predetermined first distance (e.g. one-third of the complete range of motion) and after completion of the target's movement calculates a first estimate of the focus position by the median point method. In the median point method, a host workstation calculates the sum of the recorded values and determines the target's position at which half of this sum was exceeded, to be a first estimate of the focus position. After

calculating an estimate of the focus position, the microscope system moves the target to this focus position estimate.

5 Depending on the amount of jitter, or lack thereof in the movement of the target, additional or fewer estimates of the focus position can be calculated prior to movement of the target to the focus position. The number of estimates are kept to a minimum, to get the highest speed possible.

10 In one embodiment, after a first pass in a coarse auto-focus operation, the microscope system makes an optional second pass (similar to the first pass), to move the target through a predetermined second distance that includes positions on both sides of the focus  
15 position's first estimate. During the second pass, the microscope system again records the values of the electronic focus signal, and to calculates a second estimate of the focus position, again using the median point method. Although the second pass is optional, a  
20 coarse auto-focus operation for a low power objective lens that includes two passes has the advantage of accurate and repeatable focusing even in the presence of noise that typically causes conventional microscopes to malfunction. Also, such an auto-focus operation  
25 does not require a threshold and allows focus to be found quickly over a large range, but with high accuracy from the second pass. Finally, the median point method is less susceptible to stage backlash and stage motion irregularity than other conventional  
30 methods.

An auto-focus operation can be implemented with any one of three methods for generating the electronic focus signal: (a) spot method, (b) line scan method or (c) area scan method, in which a microscope system  
35 measures the value of an electronic focus signal for each elevation of the target (a) as the illumination

spot is held stationary, (b) as the illumination spot is scanned along a line, or (c) as the illumination spot is scanned in an area in the field of view respectively. In a line scan method and an area scan method, the electronic focus signal can be measured at discrete spots or alternatively measured continuously over an infinite number of spots. In each of these three methods, the microscope system moves the target through several elevations, and determines the elevation that generates the peak (e.g. largest value) of the electronic focus signal.

Automatic focusing can be based on a line scan, so that the electronic focus signal is generated for several spots along a line segment on the target and the microscope automatically focuses on the brightest feature of the target that is in that line segment.

Automatic focusing can be based on an area scan, using an area peak detector (sometimes referred to as "peak detector") that generates and a microprocessor that records the largest value of the electronic focus signal as the illumination spot scans an area.

The microprocessor computes an estimate of the focus position by using the recorded values of the electronic focus signal for each target elevation. In one embodiment, during a pass in a fine auto-focus operation, the microprocessor estimates the focus position to be the target elevation at which the largest recorded value (sometimes referred to as "peak value") of the electronic focus signal was generated, and the microscope system moves the target to focus on the layer ("brightest layer") that generated this largest recorded value.

A predetermined layer that the user is interested in focusing on can be reached using an offset from the brightest layer. Use of the brightest layer and the offset results in an accurate and repeatable auto-focus

operation that is dependent only on the brightest layer's reflectivity and is not influenced by the shape, location, or extent of other features on a target. Such an auto-focus operation is also not  
5 dependent on prior knowledge of the reflectivity of the predetermined layer.

The invention will be more fully understood in light of the following drawings taken together with the detailed description.

10

#### BRIEF DESCRIPTION OF THE DRAWINGS

FIG. 1 is a simplified block diagram that illustrates a confocal microscope system according to an embodiment of the invention.

15

FIGS. 2A-2C show a target of a confocal microscope below the focus position, at the focus position and above the focus position, respectively, illustrating, at each position, the pattern of light reflected from the target.

20

FIG. 3A illustrates the relation between the magnitude of an idealized electronic focus signal of a confocal microscope and the position of a target with respect to an objective lens as the target is moved along the Z-axis.

25

FIG. 3B illustrates an idealized electronic focus signal after the auto-focus system has been initialized.

FIG. 3C illustrates an actual electronic focus signal for a microscope system that uses a low power  
30 objective lens.

FIG. 4 is a block diagram of one embodiment of a Z-axis controller used to control a fine Z-stage and to provide feedback to a coarse Z-stage.

FIGS. 5A-5E are schematic diagrams of the  
35 embodiment of the Z-axis controller of FIG. 4.

FIG. 6 is a graphic representation of three passes

performed during a coarse auto-focus operation according to the invention.

FIG. 7 is a graphic representation illustrating one embodiment of a fine auto-focus operation.

5 FIG. 8A is a graphic representation illustrating another embodiment of a fine auto-focus operation.

FIGS. 8B and 8C illustrate the movement of a target along the Z axis of FIG. 8A with respect to time in two alternative embodiments.

10 FIG. 9 is an isometric view of a fine Z-stage, FIG. 10 is a cross sectional view of the fine Z-stage of FIG. 9.

FIG. 11 is an exploded isometric view of a piezoelectric element, sensor and bottom plate of the  
15 fine Z-stage of FIG. 9.

FIG. 12 is a block diagram of another embodiment of a Z-axis controller used to control a fine Z-stage and to provide feedback to a coarse Z-stage.

FIGS. 13A, 13B and 13C are timing diagrams  
20 illustrating various signals in one embodiment of a confocal microscope system.

FIGS. 14-21 are schematic diagrams of the embodiment of the Z-axis controller of FIG. 12.

## 25 DETAILED DESCRIPTION

FIG. 1 is a simplified block diagram of a confocal microscope system 100 described briefly below and described in more detail in commonly owned, copending U.S. Patent Application, Serial No. 08/080,014,  
30 entitled "Laser Imaging System for Inspection and Analysis of Sub-Micron Particles", filed by Bruce W. Worster et al, on June 17, 1993, that is incorporated herein by reference in its entirety.

Laser 102 (FIG. 1) generates a laser beam 123I  
35 that is transmitted through a polarizing beam splitter 104, a converging lens 130, a spatial filter 131 (such

as a pinhole), a collimator lens 132, a quarter waveplate 133, reflected from an X-mirror 106 and a Y-mirror 108, and transmitted through an objective lens 110 to the surface of a target 112. In one embodiment, laser 102 is a conventional argon-ion laser, however, other types of lasers may be utilized in alternate embodiments. Target 112 is an object, such as a semiconductor wafer, that is to be viewed using microscope system 100. X-mirror 106 and Y-mirror 108 are each rotatable about an axis such that the illumination spot created by incident laser beam 123I can be moved along an X-axis and a Y-axis, respectively, of target 112. Laser 102, beam splitter 104, converging lens 130, spatial filter 131, collimator lens 132, quarter waveplate 133, X-mirror 106, Y-mirror 108 and objective lens 110 are each conventional structures that are well known by those skilled in the art of confocal microscopes.

Laser beam 123I reflects from the surface of target 112 in a pattern (illustrated in FIGs. 2A-2C) that is dependent upon the distance of objective lens 110 from target 112. FIGs. 2A-2C show target 112 below focus position 203, at focus position 203 and above focus position 203, respectively. When target 112 is positioned below or above focus position 203, respectively, a small percentage of light 123R from incident laser beam 123I that is originally transmitted through objective lens 110 is reflected back through objective lens 110 in a coherent manner. The amount of reflected light 123R is illustrated by the graph of the electronic focus signal as a function of position in FIG. 3A. There is a gradual transition in the amount of reflected light (depending on the target's position) between FIGs. 2A-2C. As is well known in microscope engineering art, the amount of reflected light also depends on, for example, reflectivity of target, angle

of target and angle of the chuck that supports the target.

Microscope system 100 includes a spatial filter 131 with a pin hole that permits passage of incident laser beam 123I and some portion of reflected light 123R, the portion depending on the position of target 112. Spatial filter 131 blocks off the rest of reflected light 123R as illustrated in FIGs. 2A and 2C. However, when target 112 is positioned at focus position 203 (FIG. 2B), a maximum amount of light from incident laser beam 123I is reflected (light 123R) and transmitted back through objective lens 110 and the pinhole along a conjugate path of microscope system 100.

Reflected light 123R passes through objective lens 110, is reflected by Y-mirror 108 (FIG. 1), X-mirror 106 and reflects off beam splitter 104 to a photodetector 114. Photodetector 114 is a device, such as a photo-multiplier tube (PMT) or photo-diode, that generates an electronic focus signal 115 that is an analog signal that has a magnitude (voltage in one embodiment) proportional to the intensity of reflected light 123R as measured by photodetector 114. The photodetector's gain is adjustable to accommodate different laser power, laser wavelength and target reflectivity. As is well known in microscope engineering, when the target is in focus, the value of the detected electronic focus signal is proportional to:  $\text{PMT's GAIN} * \text{TARGET'S REFLECTIVITY} * \text{LASER'S POWER}$ . The PMT's gain is adjusted empirically, as discussed below.

Electronic focus signal 115 (FIG. 1) is provided to host workstation 116 and to Z-axis controller 118. Z-axis controller 118 is directly coupled to fine Z-stage 120 and is indirectly coupled to coarse Z-stage 122 through host workstation 116 and coarse Z-axis

controller 117. Coarse Z-stage 122 uses a motor, such as a stepper motor, to move target 112 relative to objective lens 110 through a relatively large range of motion (e.g. a large percentage of total possible movement) along the Z-axis, in an operation referred to as "coarse auto-focus operation". In one embodiment, coarse Z-stage 122 moves target 112 through 1000 microns in a start-up move of a coarse auto-focus operation to a starting position 516 (FIG. 6) from an initial position of target prior to the auto-focus operation (the total possible movement being approximately 6000  $\mu\text{m}$  microns).

In one embodiment of the present invention, coarse Z-axis controller 117 is a conventional stepper motor controller available as part number 310MX3 from New England Affiliated Technology (NEAT) of 620 Essex St., Lawrence, MA 01841. Coarse Z-stage 122 can be driven by a conventional stepper motor (not shown) such as the Vexta C5858-9012 available from Oriental Motor of 16-17 Veno 6-Chome Taito-Ku, Tokyo, Japan. In this embodiment, coarse Z-stage 122 is Part # 1930237 available from NEAT (above).

As explained in more detail below, fine Z-stage 120 uses a piezoelectrically driven element to move target 112 through a relatively small range of motion (e.g. a small percentage of total possible movement) along the Z-axis, in an operation referred to as "fine auto-focus operation." In one embodiment, fine Z-stage 120 moves target 112 through 50 microns in a fine auto-focus operation (the total possible movement being approximately 6000  $\mu\text{m}$  microns). Although the invention is described as having a movable target 112 and a stationary objective lens 110, target 112 can be held stationary while objective lens 110 is moved.

FIG. 3A is an idealized graph of the magnitude (sometimes referred to as "strength" or "intensity") of



an electronic focus signal 115A as a function of the position, e.g. elevation of target 112 along the Z-axis. Electronic focus signal 115A has a magnitude that is theoretically proportional to a sinc-squared function  $((\sin(x)/x)^2)$ , having a full width half max measurement 305A that varies based on the numerical aperture of objective lens 110 and the wavelength of laser beam 123. The full width half max measurement 305A is the width, along the Z-axis, between two points 308A and 308B at which the magnitude of electronic focus signal 115A is at half of its maximum magnitude 301A at focus position 203. As illustrated in FIG. 3A, electronic focus signal 115A has a value of  $(A1+A2)/2$  for the full width half max measurement 305A, wherein A1 is the background value (described below) and A2 is the maximum value for focus position 203. The relationship of the numerical aperture and the wavelength to the magnitude of the electronic focus signal is well known to a person skilled in the art of confocal microscopes.

In FIG. 3A, focus position 203 at the center of main lobe 307A is a distinct position, as shown by a sharp, well defined peak 301A in the magnitude of electronic focus signal 115A. Electronic focus signal 115A also exhibits side lobes such as lobes 306a-306b. Depth of focus (sometimes referred to as "width of focus") 302A spans a range along the Z-axis in which the magnitude of electronic focus signal 115 is greater than a background value 303. Background value 303 is a small value of electronic focus signal 115A that is non-zero and results from leakage currents and a small amount of background light that reaches photodetector 114.

Depth of focus 302A becomes smaller as the numerical aperture of objective lens 110 increases or as the wavelength of laser beam 123 decreases, as is

well known. In the following discussion, objective lens 110 has a power of 100X and numerical aperture of 0.95, and laser beam 123 has a wavelength of 488 nm in one embodiment. In this embodiment, electronic focus  
5 signal 115A has a full width half max measurement 305A of approximately 0.5 microns.

Microscope system 100 automatically moves target 112 to an estimate of focus position 203 that results from one or more auto-focus operations as described  
10 below. To initialize microscope system 100 for a coarse auto-focus operation, the gain of photodetector 114 is increased to an empirically predetermined autofocus gain. Also, an empirically predetermined zero position offset 310 is applied to electronic focus  
15 signal 115 such that background value 303 and side lobes 306a-306b are effectively eliminated from electronic focus signal 115A. FIG. 3B illustrates electronic focus signal 115A after microscope system 100 has been initialized.

20 The auto-focus gain of photodetector 114 (FIG. 1) and zero position offset 310 (FIG. 3B) are empirically determined by performing a series of auto focus operations using different targets 112 to cover as wide a range of targets as typically used by the user (e.g.  
25 an aluminum target and a semiconductor target). Zero position offset 310 is selected to ensure that background value 303 and side lobes 306a-306b are below zero position 320. A small positive value greater than zero position 320 is then selected for use as threshold  
30 value 304. Once threshold value 304 has been selected, the threshold value 304 can be used for subsequent auto-focus operations on a wide variety of targets.

When the magnitude of electronic focus signal 115A exceeds a predetermined threshold value 304, a near-  
35 focused condition is said to exist. During a near-focused condition, target 112 is relatively close to

focus position 203 (illustrated by positions 308a and 308b of FIG. 3A). The gain of photodetector 114 and the zero position offset 310 are then maintained during the coarse auto-focus operation.

5        Certain photodetectors 114, particularly photo multiplier tubes (PMTs), can be damaged when exposed to high signals (laser signals having a magnitude greater than magnitude that PMTs are designed for) that result from increasing the gain of photodetector 114. Such  
10       damage is not a problem in the current invention because PMT control circuitry has an overload sensing circuit that automatically reduces PMT's gain to zero by a conventional method, before damage occurs.

FIG. 3C illustrates variations in magnitude of an  
15       actual electronic focus signal 115C. Magnitude of electronic focus signal 115C is irregular due to noise caused by, for example, vibration between target 112 and objective lens 110 during movement of target 112. In FIG. 3C, objective lens 110 has low power (for  
20       example magnification of 1.5X or 5X). The low power of objective lens 110 results in a large depth of focus, such as distance 305C (e.g. 30-100 microns) between elevations 335A and 335B for the full width half max measurement. The large depth of focus and the  
25       irregular waveform of electronic focus signal 115C can result in a false focus position, such as one of position 333A or 333B, if peak detection is used to estimate focus position 203.

Use of a threshold, as described above in  
30       reference to FIGs. 3A and 3B, when electronic focus signal 115C is generated by a low power objective lens has the disadvantage that a small change in threshold, (e.g. from 0.5 to 0.6) results in a large change in the target's elevation along the Z axis (e.g. from points  
35       335A and 335B to points 336A and 336B) at which the threshold is exceeded. In such a case, focus position

203 of idealized main lobe 307 is estimated in one embodiment without using a threshold, by summing up the magnitude of electronic focus signal 115C at several elevations of target 112 in the range of motion along the Z-axis, in a median point method that is described below.

FIG. 4 is a block diagram of Z-axis controller 118 that controls fine Z-stage 120 and also provides feedback used to control coarse Z-stage 122. Operation of some structures in FIG. 4 is described below, while operation of the rest is described in reference to FIGs. 5A-5E and 6-8.

Within Z-axis controller 118 (FIG. 4), electronic focus signal 115 is transmitted to a first input terminal of comparator 401 and to an input terminal of an analog to digital converter (ADC) 407. ADC 407 measures the magnitude of electronic focus signal 115 when commanded by microprocessor 403, for example eighty (80) times a second. A digital output signal (on 8 lines in one embodiment) of ADC 407 drives microprocessor 403. Comparator 401 sets a flip-flop 402 when the magnitude of electronic focus signal 115 exceeds a threshold signal at the comparator's second input terminal that is driven by a digital to analog converter (DAC) 404. DAC 404 in turn receives a threshold value from microprocessor 403.

The output terminal of comparator 401 (FIG. 4) is coupled to the set terminal of latching flip flop 402. The Q output terminal of flip-flop 402 is coupled to an input terminal of status register 405.

An output terminal of control register 406 is coupled to the reset terminal of flip flop 402. Control register 406 resets flip-flop 402 on command from microprocessor 403.

Microprocessor 403 has terminals coupled to status register 405, control register 406, ADC 407 and host

work station 116. Microprocessor 403 is also coupled to position control register 408. The output of position control register 408 is transmitted through DAC 409, summing node 410, integrator 420, and amplifier 411 to provide a control voltage signal to a piezoelectric element 1130 of fine Z-stage 120. Summing node 410 also receives a feedback signal (on line 435) that is linearly proportional to the position of target 112 along the Z-axis i.e. distance from a proximity sensor 1135 of fine Z-stage 120. The feedback signal is used to linearize the behavior of a piezoelectric element 1130 (FIG. 11).

FIGs. 5A-5E are schematic diagrams of one embodiment of Z-axis controller 118 of FIG. 4. Similar elements in FIGs. 4 and 5A-5E are labelled with the same number. Ratings of various components illustrated in FIGs. 5A-5E are listed in microfiche appendix D (at pages 121-126) for one embodiment.

Central processing unit (CPU) 5000 (FIG. 5D) of microprocessor 403 transmits and receives information through bus transceiver 5034 to 8-bit data bus 5006. CPU 5000 is, for example, a TP-RS485 twisted pair control module, model number 55050-00, available from Echelon of 4015 Miranda Avenue, Palo Alto, CA 94304. Bus transceiver 5034 provides additional drive capability to CPU 5000. Bus transceiver 5034 is a well known device, available for example as part number 74ALS245, from Texas Instruments (TI) of 7839 Churchill Way, Dallas, TX 75251.

Address register 5036 receives addressing information from CPU 5000 through data bus 5006. This addressing information determines the register or device within Z-axis controller 118 that the microprocessor 403 accesses. Address register 5036 is available, for example, from TI (above) as part number 74ALS573.

The output of address register 5036 is provided to address decoders 5038 and 5040. Address decoders 5038 and 5040 decode the addressing information and generate signals that enable a register or device within Z-axis controller 118 that microprocessor 403 accesses. Address decoders 5038 and 5040 are available for example from TI (above) as part numbers 74ALS138 and 74LS139, respectively. Microprocessor 403 communicates with position control register 408, status register 405, control register 406, DAC 404 and ADC 407 using 8-bit data bus 5006.

Registers 5001 and 5003 (FIG. 5A) within position control register 408 receive positioning information from microprocessor 403 on data bus 5006. Registers 5001-5004 are well known in the art of designing microprocessor based systems. Registers 5001 and 5003 are available for example from TI (above) as part number 74ALS574. Registers 5002 and 5004 are available for example from TI (above) as part number 74ALS273. Eight bit words on data bus 5006 are transmitted through registers 5001-5004 of position control register 408 to provide a 12-bit input to 12-bit DAC unit 5008 of DAC 409.

DAC unit 5008 is a conventional digital to analog converter (DAC), known in the electronics art, and available for example from Analog Devices (AD) of 1 Technology Way, Norwood, MA 02062, as part number AD7541. The remaining ancillary elements of DAC 409 including operational amplifier 5009 and the illustrated resistors, capacitors and diodes are conventional elements commonly seen in conventional circuits using DAC 409, as is well known to a person skilled in the art of electronic engineering. Operational amplifier 5009 is available for example from AD (above) as part number OP-177E. DAC 409 provides an analog output signal on lead 5010.

As shown in FIG. 5B, lead 5010 is connected to one input terminal of operational amplifier 5012 of summing node 410. Operational amplifier 5012 is available for example as part number AD712 from AD (above). The  
5 signal on the other input terminal operational amplifier 5012 is derived from the position feedback signal (on line 435 in FIG. 4) of sensor 1135 (FIG. 11) in fine Z-stage 120. Sensor 1135 provides an input signal to operational amplifier 5018. Operational  
10 amplifier 5018 is available for example from AD (above) as part number AD712.

Various illustrated resistors and capacitors, such as R7, C8, R8 that are coupled to operational amplifier 5018 create a conventional buffer. Components rated as  
15 DNL (Do Not Load) in FIGs. 5A-5E and FIGs. 14-21 are not used (i.e. not part of the circuit). The output signal of operational amplifier 5018 is provided to the other input terminal of operational amplifier 5012. The output terminal of summing node 410 is coupled to  
20 the input terminal of integrator 420. Integrator 420 includes an operational amplifier such as part number AD712, available from AD (above). The associated resistors, capacitors and diodes such as R19, C20, R20, R21, C25, D5 of integrator 420 are known to one skilled  
25 in the art of designing electronics.

The output signal of integrator 420 is provided to notch filter 5013, which includes two resistors (e.g. R24, R23) and two capacitors (e.g. C22, C27). The  
output signal of notch filter 5013 is provided to  
30 operational amplifier 5014, that is available, for example, from AD (above) as part number AD712. The output signal of operational amplifier 5014 is provided to the input terminal of amplifier 411.

Amplifier 411 is a conventional amplifier that  
35 includes an operational amplifier available from for example Apex Microtechnology of 5980 N. Shannon Rd.,

Tucson, AZ 85741 as part number PA85. The illustrated diodes, resistors and capacitors such as R22, D3, D4, D8, C34, R14, C14, R5, R6, C7, R2, C1, R4 of amplifier 411 are all known in the electronics art. The output  
5 signal of amplifier 411 is provided to piezoelectric element 1130 within fine Z-stage 120.

As shown in FIG. 5C, electronic focus signal 115 is provided to ADC 407. Electronic focus signal 115 is routed through multiplexer 5022 to operational

10 amplifier 5020. Multiplexer 5022 is a conventional part available, for example from Siliconix of 2201 Laurelwood Road, Santa Clara, CA 95058 as part number DG 211. Operational amplifier 5020 is available from for example AD (above) as part number AD843.

15 Operational amplifier 5020 buffers electronic focus signal 115. The output of operational amplifier 5020 is provided to an input of ADC unit 5021. ADC unit 5021 is a conventional part available for example as part number AD7575 from AD (above). The other devices  
20 C17, C18, D2, R26, C23, C51, C52, U27 coupled to ADC unit 5021, as illustrated in FIG. 5C, are known in the art of electronics. In response to electronic focus signal 115, ADC unit 5021 supplies an 8-bit digital signal representative of electronic focus signal 115.

25 The 8-bit digital signal of ADC unit 5021 is provided to microprocessor 403 on data bus 5006.

FIG. 5C also illustrates flip flop 402. Flip flop 402 is programmed as one of the devices present within programmable logic device (PLD) 5023. PLD 5023 is  
30 available for example from Lattice Semiconductor of 5555 NE Moore Ct., Hillsboro, OR 97124, as part number GAL20RA10. The input signals to PLD 5023 include: a set input signal from comparator 401 and a reset input signal from control register 406. PLD 5023 processes  
35 these input signals and generates a Q output signal representing the output signal at Q output terminal of



flip flop 402. This Q output signal is provided to status register 405. (PLD 5023 also has input signals and output signals unrelated to auto-focus operations.) One embodiment of instructions to a programmer for programming PLD 5023 are listed microfiche appendix C that is incorporated herein by reference in its entirety. Instructions in microfiche appendix C can instruct programmer model no. BP-1200 available from BP Microsystems, Inc. of 1000 N. Post Oak Road, Houston, TX 77055.

FIG. 5D illustrates status register 405. Status register 405 is a conventional register available for example as part number 74ALS541 from TI (above). As previously discussed, status register 405 receives the Q output signal of flip flop 402 from PLD 5023 (status register 405 also receives other information unrelated to auto-focus operations). The 8-bit output signal of status register 405 is provided to data bus 5006 such that microprocessor 403 can detect when flip flop 402 is set.

Control register 406 (FIG. 5D) receives an 8-bit input signal from microprocessor 403 on data bus 5006. Control register 406 is available for example from TI (above) as part number 74ALS273. An output terminal of control register 406 is coupled to PLD 5023, such that a signal from control register 406 can reset flip flop 402.

DAC 404 also receives an 8-bit input signal from microprocessor 403 on data bus 5006. This 8-bit input signal is transmitted through register 5007 (available for example from TI as part number 74ALS574) to conventional DAC unit 5011 (available for example from National Semiconductor as part number DAC0808). DAC unit 5011 converts the incoming 8-bit signal into an analog output signal. This analog output signal is provided to an input terminal of operational amplifier

5017 (available for example from AD as part number AD712). The output signal of operational amplifier 5017 is provided to an input terminal of comparator 401. Electronic focus signal 115 is provided to the other input terminal of comparator 401. Comparator 401 includes comparator unit 5019, available for example from National Semiconductor of 2900 Semiconductor Drive, Santa Clara, CA 95062, as part number LM311. The output signal of comparator 401 is provided to flip flop 402.

FIG. 5E illustrates conventional structures used in the power supply connections 5041, 5042 that supply power to fine Z-stage controller 118. Analog/digital grounding structure 5043 is used to connect analog and digital grounds on fine Z-stage controller 118.

In one embodiment of the present invention, a coarse auto-focus operation is performed as follows. Host work station 116 (FIG. 1), through coarse Z-stage controller 117, instructs coarse Z-stage 122 to move target 112 in a start-up move (not shown) downward through a startup predetermined distance (e.g. 1000 microns) along the Z-axis from the current position of target 112. If there is less distance from the start-up predetermined distance between the current position and a lowermost position (not shown) along the Z-axis, then target 112 is moved to the lowermost position along the Z-axis. The lowermost position is that position along the Z-axis below which target 112 cannot be moved by any mechanism in microscope system 100. This startup pass ensures that target 112 is positioned at a starting position 516 below focus position 203. For example, as all focus positions are within a limited distance ( $\approx 200 \mu\text{m}$  in one embodiment) of an upper travel limit, a position e.g.  $1000 \mu\text{m}$  below any initial position (not shown) guarantees that starting position 516 is below focus position 203. Travel

limits are two positions of a target that are farthest from each other, such as a lower most position (not shown) and first safe operating position 520 of FIG. 6. If the first time a target is loaded, coarse Z stage 122 is at its lowermost position, there is no startup pass. A startup pass also ensures that target 112 is positioned at least a minimum distance below focus position 203, thereby allowing coarse Z-stage 122 to achieve consistent start-up characteristics (i.e., velocity, acceleration, etc.) before target 112 encounters a focused condition.

Host work station 116 instructs microprocessor 403 to send a digital signal to DAC 404 such that the output signal of DAC 404 has a voltage level corresponding to threshold value 304 that was determined during initialization of microscope system 100. In one embodiment, the input signal to DAC 404 is a fixed digital input signal having a value of 24 (out of a range of values between 0 and 255). In response to another signal generated by host workstation 116, control register 406 transmits a reset signal to reset flip-flop 402 to its initial state (e.g., a logic "0").

FIG. 6 is a graphic representation of three coarse passes 513, 514 and 515 that are performed during a coarse auto-focus operation according to one embodiment of the invention. A number of coarse passes other than three can be used in other embodiments. The vertical axis in FIG. 6 illustrates the position (e.g. elevation) of target 112 (FIG. 1) along the Z-axis. The horizontal axis in FIG. 6 illustrates the magnitude of the electronic focus signal 115.

To perform a first coarse pass, host work station 116 (FIG. 1) instructs coarse Z-stage controller 117 to move (illustrated by arrow 512) coarse Z-stage 122, and thereby target 112, from starting position 516 in positive Z direction 533 through a first safe operating

distance 512 to a first safe operating position 520. If no focus condition is detected during a first coarse pass, target 112 is moved to first safe operating position 520. If a focus position is detected, the movement of target 112 is illustrated by the three coarse passes, 513, 514 and 515.

First safe operating distance 512 is selected so that there is no chance that target 112 will contact objective lens 110 that is located at position 517. Any amount of clearance between first safe operating position 520 and position 517 can be selected. In one embodiment, the focus position of each objective lens is measured and travel limits are programmed (e.g. soft coded into a software configuration table) so that no contact occurs between target 112 and objective lens 110. All entries in a configuration table are defined by a user during installation of microscope system 100. In particular, travel limits for movement of the target are set by manual focusing on a target. In one embodiment, first safe operating position 520 is empirically set at 50 microns above focus position 203 (i.e., 50 microns above the focal point of objective lens 110) based on various parameters such as thickness of target and the working distance (e.g. 270 microns) as described below. In one embodiment, first safe operating distance 512, is selected from the range of 1000-6500 microns, depending on the initial position of target 112, i.e., before target 112 is moved to starting position 516.

The working distance of objective lens 110 is the distance from the objective lens's position 517 (FIG. 6) to focus position 203. In various embodiments, the working distance varies from 270 microns to several millimeters, depending on the numerical aperture of objective lens 110.

Safe operating position 520 depends on many

factors, such as flatness of the target, flatness of the XY stage on which the target is supported, position of fine Z stage 120, thickness of the target and repeatability of a "home" position of coarse Z stage 122. Safe operating position 520 is empirically selected to be above a typical focus position 203 for a typical target, but below the position at which the target touches objective lens 110. Safe operating position 520 is selected to be sufficiently away from focus position 203 to ensure enough travel to obtain a well defined peak in the magnitude of electronic focus signal 115.

In a threshold method, during first coarse pass 513, comparator 401 continuously compares incoming electronic focus signal 115 with threshold value 304 received from DAC 404. Because target 112 is initially out of focus, electronic focus signal 115 is initially less than threshold value 304. Under these conditions, the output signal of comparator 401 has a positive voltage. As target 112 approaches focus position 203, the magnitude of electronic focus signal 115 increases. When the magnitude of electronic focus signal 115 exceeds threshold value 304, the output signal of comparator 401 transitions to a negative voltage, thereby setting flip flop 402.

Once flip flop 402 (FIG. 4) is set, the voltage at the Q output terminal transitions to a logic high state in this embodiment. Such a logic high state at the Q output is transmitted to status register 405, causing a bit within status register 405 to change value. Microprocessor 403, which continuously monitors status register 405, thereby detects that flip flop 402 has latched. Upon detecting this latched condition, microprocessor 403 signals host work station 116. In response, host work station 116 instructs coarse Z-axis controller 117 to stop coarse Z-stage 122, and thereby

stop movement of target 112.

Because electronic focus signal 115 is an analog signal and flip flop 402 latches when the magnitude of electronic focus signal 115 exceeds threshold value 304, the possibility of missing a focused condition is eliminated if the peak magnitude is greater than the threshold value. There are no discrete sampling periods during which the focused condition may be missed. Thus, flip flop 402 latches even for a focus signal 115 having a narrow depth of focus 302, as illustrated in FIG. 3A. Consequently, target 112 can be moved upward at a much faster velocity than was possible with auto-focus microscopes of the prior art.

In one embodiment, the average velocity of target 112 during first coarse pass 513 is dependent on first safe operating distance 512. Target 112 is moved at an average velocity that enables target 112 to move through first safe operating distance 512 in approximately one second. However, the maximum average velocity is approximately 3000 microns per second in this embodiment. This allows first coarse pass 513 to be completed in approximately 1 second.

If a focus condition is not detected before target 112 completes movement through first safe operating distance 512, host work station 116 instructs coarse 2-stage controller 117 to stop target 112 when first safe operating position 520 is reached. A focus position can remain undetected depending on the parameters used to detect focus position 203. For example, when a threshold value is larger than the peak magnitude of electronic focus signal 115, the output of comparator 401 does not transition to a negative voltage even as target 112 passes through focus position 203.

Because of the high velocity (e.g. 3000 microns per second) at which target 112 is moved during first coarse pass 513 and the significant amount of time

(e.g. 20-60 milliseconds) required to stop the movement of target 112 after a focused condition is detected, target 112 moves ("target overshoot") to first stopping position 505 (FIG. 6) above focus position 203 at the time that target 112 comes to rest. This target overshoot is illustrated in FIG. 6, which shows that a focused condition is detected at first trip position 504, and that target 112 comes to rest at first stopping position 505. Thus, it is necessary to perform one or more additional passes to position target 112 closer to focus position 203.

The distance required to stop target 112 during first coarse pass 513 is dependent upon threshold value 304 (because a lower threshold value 304 causes flip flop 402 to latch sooner), the system gain (because a higher system gain allows for earlier detection of a focused condition), the time required for microprocessor 403 to detect that flip flop 402 has latched, the time required for microprocessor 403 to communicate this information to host work station 116, the time required for host work station 116 to issue a command to stop coarse Z-stage 122, and the velocity of coarse Z-stage 122 at first trip position 504. Such dependencies are well known to a person skilled in engineering.

Host work station 116, through microprocessor 403, then instructs control register 406 to transmit another reset signal to clear flip flop 402 (FIG. 4) before second coarse pass 514 is begun. This reset signal allows flip flop 402 to detect the next time electronic focus signal 115 exceeds threshold value 304.

The steps performed during second coarse pass 514 are similar or identical to those of first coarse pass 513. Host work station 116 instructs coarse Z-stage 122 (through coarse Z-stage controller 117) to move target 112 in negative Z direction 532 at a velocity

that is slower than the velocity of target 112 during first coarse pass 513. This motion continues until a near-focused condition is detected at second trip position 506, causing flip flop 402 to latch, or until  
5 target 112 moves through a predetermined second safe operating distance 507 to second safe operating position 522. This second safe operating distance 507 is selected to assure that target 112 travels through focus position 203, based on the width of electronic  
10 focus signal 115 and the range of anticipated overshoot of target 112 during first coarse pass 513. Second operating distance 507 is selected greater than the overshoot from first coarse pass 513, plus the width of electronic focus signal 115.

15 In one embodiment, second safe operating distance 507 is 1200 microns, unless target 112 was stopped less than 50 microns from first safe operating position 520. If target 112 was stopped less than 50 microns from first safe operating position 520, the second safe  
20 operating distance 507 is smaller, e.g. 200 microns in this embodiment. The latter second safe operating distance is smaller because target 112 is decelerating near the end of first safe operating distance 512 so that it is easier to stop target 112, thereby  
25 decreasing the amount of target overshoot so that target 112 is stopped closer to focus position 203 than would otherwise be the case. Because target 112 is closer to focus position 203, a smaller second safe operating distance 507 can be used.

30 In one embodiment, it takes approximately 1 second for target 112 to pass through the entire second safe operating distance 507. The average velocity of target 112 during second coarse pass 514 therefore is approximately 200 microns per second or 1200 microns  
35 per second, depending upon second safe operating distance 507.



While traversing second safe operating distance 507, target 112 encounters a focused condition at second trip position 506. When electronic focus signal 115 exceeds threshold value 304, flip flop 402 latches and target 112 is stopped in the same manner as in first coarse pass 513. Target overshoot results in target 112 coming to rest below first trip position 504. In FIG. 6, target 112 stops at second stopping position 508 and a microscope system 100 estimates the location of focus position 203 based on overshoot and width of focus.

One or more additional passes can be used to make a better estimate (than the previous estimate) of focus position 203, depending on the resolution of ADC 407.

The steps performed during third coarse pass 515 are also similar or identical to the steps described above for second coarse pass 514 and first coarse pass 513. Host work station 116, through microprocessor 403, instructs coarse Z-stage controller 117 to move target 112 in positive Z direction 533 at a velocity slower than the velocity of target 112 during second coarse pass 514. This motion continues until a near-focused condition 509 causes flip flop 402 to latch or until target 112 moves through a predetermined third safe operating distance 510 to a third safe operating position 524. Again, third safe operating distance 510 is selected to assure that target 112 travels through focus position 203 and is dependent on the amount of target overshoot associated with second coarse pass 514.

In one embodiment, the third safe operating distance 510 is 140 microns. For target 112 to take approximately one second to pass through the entire third safe operating distance 510, the average velocity of target 112 during the third coarse pass 515 is approximately 140 microns per second in one embodiment.

Prior to traversing the entire third safe operating distance 510, target 112 again encounters a near-focused condition at first trip position 504. When the electronic focus signal 115 exceeds threshold value 304, flip flop 402, having been reset, latches and target 112 is stopped in the same manner as in coarse passes 513 and 514.

Because of the relatively low velocity of coarse Z-stage 122 during third coarse pass 515, target 112 comes to rest at a third stopping position 511 that can either be above or below focus position 203. After third coarse pass 515, the third stopping position 511 of target 112 is within approximately  $\pm 10 \mu\text{m}$  microns of focus position 203. To achieve this degree of accuracy, all of the factors contributing to target overshoot are considered to empirically determine the amount of time required to stop target 112 upon encountering a focused condition during the third coarse pass 515. Given the time required to stop target 112, the velocity of target 112 during third coarse pass 515 is selected to assure that third stopping position 511 of target 112 is within approximately  $\pm 10 \mu\text{m}$  microns of focus position 203.

If full width half max measurement 305 of electronic focus signal 115 is sufficiently wide, the coarse auto-focus operation described above is modified. Although there is no clear boundary that determines when an electronic focus signal 115 is "sufficiently wide", an electronic focus signal 115 created with a laser beam 123 having a wavelength of 488 nm and an objective lens 110 having a numerical aperture of 0.13 or lower is considered "wide," and has full width half max measurement 305 of approximately 29 microns. Electronic focus signal 115 exhibits a wider full width half max measurement 305 as the magnification and numerical aperture of objective lens

110 decreases or as the wavelength of laser beam 123 increases. For example, use of a low power (e.g. 1.5x or 5x) objective lens results in a large depth of focus, as illustrated by electronic focus signal 115C (FIG. 3C).

When attempting a first coarse pass 513 on a "wide" electronic focus signal 115, first stopping position 505 of target 112 is relatively close to the second trip position 506. In certain cases, first stopping position 505 of target 112 is at a position where the value of electronic focus signal 115 exceeds the threshold value 304. That is, there is not enough target overshoot to guarantee that target 112 "escapes" the electronic focus signal 115. Target 112 may "escape" if as a result of first coarse pass 513, target 112 is close to second trip position 506 and is unable to obtain a high velocity before flip-flop 402 latches during second coarse pass 514. As a result, second stopping position 508 of target 112 can be well short of peak 301. This condition causes the coarse auto-focus operation described above to miss focus position 203.

In one embodiment, to perform a coarse auto-focus operation on a "wide" electronic focus signal 115 (FIG. 3C), target 112 (FIG. 1) is moved in a start-up move by 2000 microns (rather than 1000 microns described above) below the target's initial position (not shown) to starting position 516 (FIG. 6). If target 112 is within 2000 microns of the lowermost position (not shown) along the Z-axis, target 112 is moved this lowermost position. This starting position can provide target 112 with an additional distance in which to accelerate before reaching first trip position 504. This acceleration and starting position ensure that target 112 overshoots second trip position 506. The rest of the characteristics of first coarse pass 513

for wide electronic focus signal 115C are identical to those described above.

During second coarse pass 514 for wide electronic focus signal 115C, target 112 is moved downward along the Z-axis through second safe operating distance 507. The second safe operating distance 507 and the velocity of target 112 during the second coarse pass 514 is determined in the manner previously described (i.e. based on a first estimate determined by the threshold method described above). However, rather than monitoring the status of flip-flop 402 during second coarse pass 514, host work station 116 instructs microprocessor 403 to record the values of the electronic focus signal output by ADC 407 repeatedly at a predetermined interval (e.g. 80 times a second) while coarse Z-stage 122 is moving target 112 at an approximately constant velocity.

The values recorded by microprocessor 403 therefore roughly correspond to magnitude of electronic focus signal 115 at regular distances along the Z-axis for the entire predetermined second safe operating distance 507. Based on these recorded values, focus position 203 is calculated by a median point method described below. Because movement of coarse Z-stage 122 generally induces noise in the electronic focus signal 115, the values recorded by microprocessor 403 can exhibit peaks at positions other than focus position 203 (as illustrated by FIG. 3C). Microprocessor 403 therefore filters the noise to obtain a better estimate of focus position 203 in the "median point" method.

In the median point method host workstation 116 calculates the sum of the values previously recorded by microprocessor 403 and determines the elevation along the Z-axis at which half of this sum was reached during second coarse pass 514. Host workstation 116 then

issues a command to move target 112 in the positive Z-direction to this position. The median point method assumes that noise in electronic focus signal 115 is due to irregularity of movement of coarse Z-stage 122 and that this irregularity occurs randomly and is  
5 equally probable on either side of the focus position 203 along the Z-axis.

In another embodiment of a microscope system 100 that uses a low power objective lens, during first  
10 coarse pass 513, instead of monitoring the status of flip-flop 402 according to the threshold method, microprocessor 403 records the values of electronic focus signal 115C (FIG. 3C) in the manner described above for second coarse pass 514 for use in the median  
15 point method. Therefore, during first coarse pass 513, coarse Z-stage 122 (FIG. 1) moves target 112 through a predetermined first safe operating distance 512 that is a safe travel limit of the movement of target 112. In this embodiment, starting position 516 is a  
20 predetermined distance of 1000 microns below the upper travel limit set by first safe operating position 520. Such a starting position 516 that is independent of the target's initial position avoids the inherent uncertainty in the location of the focus position  
25 relative to starting position if the target is moved by a predetermined distance from the target's initial position. At the end of first coarse pass 513, host work station 116 calculates a first estimate of focus position 203 by using the median point method described  
30 above.

Then host workstation 116 determines starting and stopping positions of second coarse pass 514. Starting position 531 is at a distance D (of 120 microns in one embodiment), above first estimate of focus position 203  
35 while stopping position 522 is at the distance D after first estimated focus position 203. Distance D is

chosen based on the desired precision and the accuracy of first estimate of focus position 203. If 80 values are collected during second coarse pass 514, a 240 micron range of movement yields a resolution of 4  
5 microns.

In this embodiment, host workstation 116 uses values recorded during first coarse pass 513, to calculate (as described below) and set an auto-focus gain (e.g. gain of a PMT), so that, the intensity of  
10 the electronic focus signal during second coarse pass 514 is optimal. After the auto-focus gain is set, host workstation 116 instructs coarse Z stage 122 to move to second coarse pass starting position 531. After coarse stage 122 has moved to starting position 531, host  
15 workstation 116 instructs microprocessor 403 to start recording the values of electronic focus signal 115, and instructs coarse Z stage 122 to move to stopping position 522.

The calculation for optimal auto-focus gain  
20 depends on the specific sensor, and in one embodiment is as follows:

if "peak intensity" of the electronic focus signal is below optimal: then new auto-focus gain=old auto-focus  
25 gain +  $\frac{\log(195)}{\log(\text{peak intensity})} * 9$ ;

if the peak intensity is above optimal, then new auto-focus gain=old auto-focus gain - (peak intensity - 195)/10,

where sensor gain is normalized to a range of 1 to 100, 1 being minimum gain, electronic focus signal intensity  
30 being normalized to the range of 0 to 225, 0 being dark. The peak intensity (e.g. 195) of the electronic focus signal is at maximum fraction (e.g. 3/4) of the maximum permitted intensity (e.g. 255), in one  
35 embodiment.

A maximum fraction (e.g. 3/4) is empirically

chosen so that electronic focus signal 115 is not too high (saturation) at focus position 203 and not too low (indistinguishable from noise) at positions other than focus position 203. Other maximum fractions, such as 2/3 or 4/5 can also be chosen, and the exact fraction that is chosen is not a critical aspect of this invention. If for some reason, the peak intensity of electronic focus signal is less than a given minimum, the minimum is used in the calculation. In one embodiment, 15 is the minimum intensity.

In one embodiment, first safe operating distance 512 is 1000 microns, yielding a 12 micron resolution for 80 samples of electronic focus signal 115. In this embodiment, second safe operating distance 537 is 240 microns, that yields a 4 micron resolution for 80 samples of electronic focus signal 115.

Appendix E illustrates one embodiment of software code that implements in host workstation 116 various computations for the auto-focus operations described above for use with a low power objective lens. Routine lonui\_StageAF in host workstation 116 initializes variables and then invokes routine lon\_APIQuery that determines the characteristics of objective lens 110. Then host workstation 116 determines a travel limit fZLimUm and then finds out magnification iObjPwr of objective lens 110. Then host workstation 116 instructs fine Z-stage 120 to move to a predetermined middle position HDWR\_FASTZ\_MID\_POS.

Host workstation 116 then saves values of several parameters currently in use, such as the intensity, system gain, zero value and Y axis scan amplitude. Then host workstation 116 sets a new zero value and a new system gain and invokes routine lonuiLoPwrStgAF for a low power objective lens or alternatively, performs various steps described below for a high power objective lens.

Routine lonuiLoPwrStgAF initializes variables, such as variable dHalfRange2ndPass and variable fZPos and then calls routine lonuiMoveStageZUm that moves target 112 to 1,000 microns below the uppermost travel limit set by first safe operating position 520 in one embodiment to starting position 516 (FIG. 6). Then host workstation 116 calls routine mot\_SetCurrAxisSpeed that sets the target speed for travelling for one second. Next, host workstation 116 calls routine lonuiMoveStgAndReadInten to move coarse Z-stage 122 to a stopping position 522 and during this movement reads the intensity (i.e. magnitude) of electronic focus signal 115. Routine lonuiMoveStgAndReadInten returns with an estimate of focus position 203 and the peak intensity. Based on the peak focus signal intensity, host workstation 116 adjusts auto-focus gain for optimal focus signal intensity in the next pass. If the peak intensity is lower than optimal, auto-focus gain is increased, if peak intensity is higher than optimal, auto-focus gain is decreased.

Routine zr\_SetTargetLaserIFrImgInten sets the auto-focus gain based on the intensity. Then host workstation 116 reads the current position by calling routine stg\_ReadStageZUm.

After optimal auto-focus gain is set, host workstation 116 then determines target 112's position in second coarse pass 514 as variable fTarget. Variable fTarget is normally the Z position of target 1128 for the first estimate of focus position 203. If, during first coarse pass 513 electronic focus signal 115's intensity is all zero, or if photodetector 114 is overloaded, then variable fTarget is set to 400 micron away from coarse Z stage 122's position. Since actual focus position 203 can be at some distance away from the first estimate of focus position 203 the second coarse pass 514 starts at some distance before the



first estimate and ends at some distance after the first estimate.

The distance before and after the first estimate is chosen to be several times the depth-of-focus of objective lens 110 (in this embodiment 120 microns) and is contained in variable dHalfRange2ndPass. If after first coarse pass 513, coarse Z stage 122 is at a position greater than distance dHalfRange2ndPass away from estimate of focus position 203, then host workstation 116 instructs coarse Z stage 122 to move to a starting position at a distance dHalfRange2ndPass before first estimate of focus position, and second coarse pass 514 starts from this position. Otherwise, second coarse 514 starts from wherever first coarse pass 513 ends. The stopping position of second coarse pass 514 is calculated to be two times dHalfRange2ndPass away from starting position 514.

Then host workstation 116 changes the speed so as to move through second safe operating distance 537 within one second, by using routine mot\_SetCurrAxisSpeed. Host workstation 116 then moves target 112 through second safe operating distance 537 by calling routine lonuiMoveStgAndReadInten, and computes the second estimate of focus position 203.

After computing host workstation 116 moves target 112 to the second estimate by calling routine lonuiMoveStageZUm.

After one of the previously described coarse auto-focus operations is completed, an optional fine auto-focus operation can be performed, if electronic focus signal 115 has a sharp peak, as illustrated in FIG. 3A. The zero position offset 310 and photodetector gain illustrated in FIG. 3B are not utilized during a fine auto-focus operation. During a fine auto-focus operation, target 112 is moved by a fine Z-stage 120 that includes a piezoelectric element 1130. Fine Z-

stage 120 is described in detail later.

Before a fine auto-focus operation is performed, target 112 is positioned by a coarse auto-focus operation (above) such that focus position 203 is within the operating range of fine Z-stage 120 (FIG. 1). In one embodiment, fine Z-stage 120 has an operating range of 50 microns along the Z-axis. Consequently, by performing the coarse auto-focus operation described above, in which target 112 is positioned within +/- 10 microns of the focus position 203, target 112 is positioned such that focus position 203 is within the operating range of fine Z-stage 120.

FIG. 7 is a graphic representation of one embodiment (henceforth "first" embodiment) of a fine auto-focus operation as encoded in routine AFocusServo on page 68 of microfiche appendix A. Although FIG. 7 illustrates two fine passes, any number of passes other than two can be used in other embodiments. The vertical axis in FIG. 7 illustrates the position of target 112 along the Z-axis. The horizontal axis in FIG. 7 illustrates the magnitude of electronic focus signal 115.

Prior to a first fine pass 601, microprocessor 403 sends a zeroing signal to position control register 408 (FIG. 4). This signal is transmitted through DAC 409, summing node 410 and amplifier 411 to piezoelectric element 1130 (FIG. 11) in fine Z-stage 120. Fine Z-stage 120, which was positioned in the middle of its operating range during the coarse auto-focus operation, moves to the lowermost position 619 (FIG. 7) of its operating range (in a startup pass similar to that described above for a coarse auto-focus operation) in response to the zeroing signal. Microprocessor 403 then instructs control register 406 to transmit a reset signal to clear flip flop 402. This instruction allows flip flop 402 to detect when electronic focus signal

115 exceeds threshold value 304.

To begin first fine pass 601, microprocessor 403 transmits a series of signals to position control register 408, thereby causing fine Z-stage 120 to move target 112 in positive Z direction 631 at a relatively high velocity. In one embodiment, this velocity is approximately 75 microns per second.

While it is necessary to precisely control the movement of target 112 during first fine pass 601, piezoelectric element 1130 (FIG. 11) of fine Z-stage 120 has a slightly non-linear position response to the voltage supplied by amplifier 411 (FIG. 4). To correct for this non-linear characteristic, a proximity sensor 1135 (FIG. 11) in fine Z-stage 120 produces an electrical feedback signal that is transmitted to summing node 410 (FIG. 4) and subtracted from the output signal of DAC 409 to create an error signal. When integrator 420 receives any non-zero error signal, integrator 420 generates an output signal that forces the error signal to zero. In this manner, integrator 420 compensates for the non-linear response of the piezoelectric element 1130, thereby allowing for linear control of the target's position by piezoelectric element 1130.

First fine pass 601 uses the threshold method so that before target 112 reaches the top of the operating range 620 of fine Z-stage 120, electronic focus signal 115 exceeds the threshold value 304 at a first trip position 604, because, as described above, the coarse auto-focus operation ensures that focus position 203 lies within the range of fine Z-stage 120 (FIG. 7). At this time, flip flop 402 latches, thereby enabling a bit in status register 405. Microprocessor 403, which continuously monitors status register 405, detects that the change in status of flip flop 402 and signals position control register 408 to stop movement of fine

Z-stage 120, and thereby stop movement of target 112. Target 112 comes to rest at a first stopping position 605 which can either be above or below focus position 203.

5           Consequently, prior to performing a second fine pass 615, microprocessor 403 instructs fine Z-stage 120 to reposition target 112 a predetermined fixed distance 614 in negative Z direction 632 such that a second  
10           stopping position 608 of target 112 is below focus position 203. Factors that must be considered when selecting distance 614 are similar or identical to the factors described above regarding the distance required to stop target 112 during first coarse pass 513. In one embodiment, distance 614 is 3.6 microns.

15           After target 112 has been repositioned at second stopping position 608, a second fine pass 615 is performed by moving target 112 in positive Z direction 631 at a relatively low velocity (e.g. 7.5 microns per second). During second fine pass 615, microprocessor  
20           403 monitors the output of ADC 407, rather than the status of flip flop 402. The output of ADC 407 is a digital representation of electronic focus signal 115. In this manner, microprocessor 403 measures the value of electronic focus signal 115 as target 112 moves  
25           along the Z-axis.

          Second fine pass 615 uses a peak detection method in which software within microprocessor 403 maintains an updated record of the maximum value output by ADC 407 and the position of target 112 at this maximum  
30           output value. Upon detecting an increasing value of electronic focus signal 115, followed by a decreasing value of electronic focus signal 115, microprocessor 403 instructs fine Z-stage 120 to stop the motion of target 112 at third stopping position 611. Third  
35           stopping position 611 is located above focus position 203 because of target overshoot which occurs for

reasons similar to those described above.

Microprocessor 403 then instructs fine Z-stage 120 to move target 112 in negative Z direction 632 to position 612 where the maximum value output by ADC 407 was detected.

In one embodiment, position 612 is typically within one tenth of a micron of focus position 203. This accuracy is determined by the velocity of target 112 and the number of bits used in ADC 407. The slower the velocity of target 112 during second fine pass 615 and the greater the number of bits used in ADC 407, the closer position 612 will be to focus position 203. The entire fine auto-focus operation is completed in approximately one second in one embodiment. The completion time is dependent upon the location of focus position 203 within the range of fine Z-stage 120. Computer code used to perform an auto-focus operation in accordance with one embodiment of fine auto-focus method is set forth in Appendix A. The computer code of Appendix A is written in neuron C language, which requires an ECHELON compiler, available from Echelon Corp. of 4015 Miranda Ave., Palo Alto, CA 94304.

FIG. 8A is a graphic representation of another embodiment (henceforth "second" embodiment) of a fine auto-focus operation. Although FIG. 8A illustrates two fine passes, a number of passes other than two can be used in other embodiments. The vertical axis in FIG. 8A illustrates the position (e.g. elevation) of target 112 along the Z-axis. The horizontal axis in FIG. 8A illustrates the magnitude (i.e. strength) of electronic focus signal 115. In first fine pass 750, target 112 is moved to each of positions 701-732. In second fine pass 850, target 112 is moved to each of positions 801-832.

Prior to first fine pass 750, microprocessor 403 sends a zeroing signal to position control register 408

(FIG. 4). This signal is transmitted through DAC 409, integrator 420, summing node 410 and amplifier 411 to piezoelectric element 1130 of fine Z-stage 120. Fine Z-stage 120, which was positioned in the middle of its operating range before a coarse auto-focus operation (e.g., Step No. 2048 in FIG. 8A), moves to the bottom of its operating range (i.e., Step No. 0 in FIG. 8A) in response to the zeroing signal, as illustrated by FIG. 8B.

First fine pass 750 uses a peak detection method. During first fine pass 750, target 112 is moved in the positive direction, e.g. upward through the full range of motion of fine Z-stage 120 (i.e., 50 microns). As shown in FIG. 8A, this range is divided into 4096 steps. Other numbers of steps can be used in other embodiments. Microprocessor 403 can position target 112 at any one of these steps by sending a digital word to DAC 409 (FIG. 4). During first fine pass 750, microprocessor 403 sequentially provides 32 digital words to DAC 409, causing fine Z-stage 120 to sequentially move target 112 to each of 32 positions 701-732. Each of the 32 positions 701-732 are separated by 128 steps (approximately 1.56 microns in one embodiment). At each of the 32 positions 701-732, the output voltage of ADC 407 (corresponding to the electronic focus signal 115) is repeatedly measured and digitally filtered (low-pass) by microprocessor 403 to obtain a single value for each of the 32 positions 701-732. Microprocessor 403 saves the peak value of the output voltage and the position at which the peak value occurred.

Because the 32 positions 701-732 of first fine pass 750 are spaced 1.56 microns apart, and electronic focus signal 115 has a depth of focus 302 of approximately 2.54 microns, a focus position is not missed. A relative peak value of the electronic focus

signal is found within the 32 positions 701-732 as long as there is enough system gain to distinguish electronic focus signal 115 from background value 303 and enough system gain such that electronic focus signal 115 can be detected by ADC 407. System gain is determined by the laser power, the gain of photodetector 114, and the reflectivity of the sample. The present invention therefore has an advantage over prior art auto-focus microscopes, which require a much larger system gain in order for an auto-focus operation to be performed.

The position 718 at which the peak value is detected during first fine pass 750 (illustrated as Step No. 2304 in FIG. 8A) becomes the center of the second fine pass 850. Microprocessor 403 instructs fine Z-stage 120 to move target 112 to 128 steps below position 718 (i.e., to Step No. 2176). In other embodiments, a different number of steps can be used.

Next, microprocessor 403 sequentially provides 32 digital words to DAC 409, causing fine Z-stage 120 to sequentially move target 112 upward to each of 32 positions 801-832. Each of 32 positions 801-832 are separated by 8 steps (approximately 0.0977 microns). The total distance of second fine pass is 3.125 microns. At each of 32 positions 801-832, output voltage of ADC 407 (corresponding to electronic focus signal 115) is repeatedly measured and digitally filtered (low-pass) by microprocessor 403 to obtain a single value for each of 32 positions 801-832. Microprocessor 403 saves the peak value of the output voltage and the position at which the peak value occurred.

In yet another embodiment (henceforth "third" embodiment), first fine pass 750 is performed as described above for the second embodiment and second fine pass 850 is performed by positioning target 112 at

128 steps above position 718 sample value was detected during first fine pass 750 (e.g., at Step No. 2432 in FIG. 8A and at position 852 in FIG. 8C).

Microprocessor 403 then sequentially provides 32  
5 digital words to DAC 409, causing the fine Z-stage 120 to sequentially move target 112 in the same direction (e.g. downward) through 32 positions 832-801 in FIG. 8A. (See also FIG. 8C.) Moving in the same direction in second fine pass 850 avoids the large acceleration  
10 caused by direction reversal at position 851 (FIG. 8B).

After passing through the entire range (3.125 micron in one embodiment) of second fine pass 850, microprocessor 403 instructs fine Z-stage 120 to position target 112 at position 812 at which a peak  
15 value is detected during second fine pass 850 (illustrated as Step No. 2272 in FIG. 8A). At the end of second fine pass 850, target 112 is positioned adjacent to actual focus position 203, i.e., within half the length between each of the 32 positions of the  
20 second fine pass 850 (within at least 0.0488 microns in one embodiment).

An advantage of a second fine pass (and additional fine passes) that target 112 is positioned quickly and reliably without relying on a threshold value. In one  
25 embodiment, fine Z-stage 120 positions target 112 sixty-seven times (including sixty-four positions at which measurements are taken and three positions at which measurements are not taken) to perform the fine auto-focus operation in approximately 0.9 seconds.  
30 Computer code used to perform one embodiment an auto-focus operation for a fine pass is set forth in microfiche appendix A.

A preferred embodiment of fine Z-stage 120 is shown in FIGs. 9-11. Appendix D at pages 127 and 128  
35 lists various parts used in the embodiment shown in FIGs. 9-11. Referring to FIG. 9, fine Z-stage 120



includes a square bottom plate 1010, four stationary bars 1020 fixedly connected along the edges of the bottom plate 1010, four rotating bars 1030 pivotally connected to the stationary bars 1020 such that each rotating bar 1030 is connected to one stationary bar 1020, and a square top plate 1050 pivotally connected along its edges to the four rotating bars 1030.

A first set of horizontally-disposed flexures 1040 is connected between upper surfaces of the stationary bars 1020 and the rotating bars 1030, and a second set of vertically-disposed flexures 1060 is connected between side surfaces of the rotating bars 1030 and the edges of the top plate 1050. In addition, the fine Z-stage 120 includes a piezoelectric actuator mechanism 1100 disposed in a space formed between the top plate 1050 and the bottom plate 1010. Finally, an optional biasing spring 1140 is connected between the top plate 1150 and the bottom plate 1110 for biasing the top plate 1150 toward the bottom plate 1110.

Referring to FIG. 10, the bottom plate 1010 is preferably a flat aluminum sheet 8 by 8 inches wide and 0.375 inch thick. As shown in FIG. 10, the bottom plate 1010 includes a receiving hole 1011 within which is located a pin 1012 for securing a first end of the optional biasing spring 1140. The bottom plate 1010 also includes an upper surface 1013.

Referring back to FIG. 9, the stationary bars 1020 are preferably aluminum bars 5.5 inches long, 0.65 inch high and 0.5 inch wide. The stationary bars 1020 are connected to the upper surface 1013 of the bottom plate 1010 using screws. The stationary bars 1020 are formed into a square frame and located along the outer edges of the bottom plate 1020. Referring again to FIG. 10, each stationary bar 1020 includes an upper surface 1021. A lip 1022 is formed along an outer edge of the upper surface 1021 of each stationary bar 1020.

The rotating bars 1030 are preferably aluminum bars which are 5.5 inches long, 0.5 inch high and 0.6 inch wide. The rotating bars 1030 are pivotally connected to the stationary bars 1020, each rotating  
5 bar 1030 being connected to one stationary bar 1020. Each rotating bar 1030 includes an upper surface 1031 and an inner side surface 1032. A lip 1033 is formed along a lower edge of the inner side surface 1032 of each rotating bar 1030.

10 As shown in FIG. 10, thin horizontally-disposed flexures 1040 are connected between the stationary bars 1020 and the rotating bars 1030. Each flexure 1040 is a thin sheet of 303 stainless steel which is 1 inch long, 1 inch wide and approximately 0.01 inch thick.  
15 Each flexure 1040 has one portion connected to the upper surface 1021 of a stationary bar 1020 by a first fixture 1041, a second portion connected to the upper surface 1031 of a rotating bar 1030 by a second fixture 1042, and a small pivot portion 1043 located between  
20 the stationary bars 1020 and the rotating bars 1030. Two flexures 1040 are connected between each stationary bar 1020 and its associated rotating bar 1030, thereby restricting each rotating bar 1030 to pivot around the pivot portion 1043 such that the rotating bar 1030  
25 remains in a parallel relationship with its associated stationary bar 1020.

The top plate 1050 is preferably a flat aluminum sheet 5.5 by 5.5 inches wide and 0.5 inch thick. Referring to FIG. 10, the top plate 1050 includes a  
30 receiving hole 1051 within which is located a pin 1052 for securing a second end of the optional biasing spring 1140 (discussed below). In addition, the top plate 1050 defines threaded holes 1053 (only one shown) for receiving preload screws 1054 and 1055 (see FIG.  
35 9). Finally, the top plate 1050 includes side surfaces 1056.

As shown in FIG. 10, thin vertically-disposed flexures 1060 are connected between the top plate 1050 and the rotating bars 1030. Similar to the above-described horizontally-disposed flexures 1040, each

5 vertically-disposed flexure 1060 is a thin sheet of 303 stainless steel which is 1 inch long, 1 inch wide and 0.01 inch thick. Each vertically-disposed flexure 1060 has one portion connected to the inner side surface 1032 of a rotating bar 1030 by a fixture 1061, a second

10 portion connected to a side surface 1056 of top plate 1050 by a fixture 1062, and a pivot portion 1063 located between the rotating bars 1030 and the top plate 1050. Two flexures 1060 are connected between each rotating bar 1030 and one edge of the top plate

15 1050, thereby restricting the top plate 1050 to pivot with respect to the rotating bars 1030 such that the top plate 1050 remains in a parallel relationship with the rotating bars 1030.

Referring to FIGS. 10 and 11, the piezoelectric actuator mechanism 1100 includes a retainer block 1110

20 connected to the bottom plate 1010, a rotating block 1120 integrally and pivotally connected to the retainer block 1110, and an extending portion 1114 integrally and pivotally connected to the retainer block 1110.

25 Further, a piezoelectric element 1130 is received in the retainer block 1110 and has a free end contacting the rotating block 1120. Finally, a sensor 1135 is received in the retainer block 1110 and generates a signal corresponding to an amount of rotation of the

30 extending portion 1114.

Referring to FIG. 11, the retainer block 1110 is formed from 7075-T6 high-strength aluminum alloy and is approximately 1.1 inch long, 0.6 inch wide and 0.1 inch thick. The retainer block 1110 defines a first

35 through-hole 1111 for receiving the piezoelectric element 1130. A stainless steel plate 1112 is

connected by fasteners 1113 to retain the piezoelectric element 1130 within the through-hole 1111. In addition, the retainer block 1110 defines a second through-hole 1117 for receiving the sensor 1135. The  
5 retainer block 1110 also includes a hole 1118 for receiving a sensor lock screw 1119 which contacts and secures the sensor 1135.

As shown in FIG. 11, an extended portion 1114 is formed from 7075-T6 high-strength aluminum alloy and is  
10 integrally connected to the upper surface of the retainer block 1110 by a thin flexure. The extended portion 1114 includes a socket 1115 for receiving a 0.125 inch diameter steel ball 1116 which contacts the top plate 1150. The extended portion 1114 is preloaded  
15 downward by the preload screw pressing against the steel ball 1116 (see FIG. 9). The sensor 1135 transmits a signal representing a distance between the sensor and a side wall of the extended portion 1114 which varies in response to the rotation of the  
20 extended portion 1114. The signal is used to determine the vertical displacement of the top plate 1050, as discussed above.

As shown in FIGs. 10 and 11, the rotating block 1120 is a prism formed from 7075-T6 high-strength  
25 aluminum alloy and is approximately 1 inch long. As shown in FIG. 10, the rotating block 1120 includes a vertical side wall 1121, an upper wall 1122 and a diagonal wall 1123. The rotating block 1120 is integrally connected to the upper surface of the  
30 retainer block 1110 by a flexure 1124. The side wall 1121 defines a socket 1125 for receiving a 0.125 inch diameter steel ball 1126 which contacts an end of the piezoelectric element 1130. In addition, the upper wall 1122 defines a second socket 1127 for receiving  
35 another 0.125 inch diameter steel ball 1128 which contacts the preload screw 1154 mounted in the top

plate 1050. It is noted that the steel ball 1128 is located further from the flexure 1124 than the steel ball 1116, which is mounted on the extended portion 1114.

5       The piezoelectric element 1130 is a cylindrical unit housed in the through-hole 1111 of the retainer block 1110 such that movement of a first end is prevented by the plate 1112, and a second end contacts the vertical side wall 1121 of the rotating block 1120  
10       through the ball 1126. The piezoelectric element 1130 is connected to amplifier 411, as discussed above. A preferred piezoelectric element is sold by Physic Instrument of Waldbronn, Germany under model number P830.20.

15       The sensor 1135 is also a cylindrical unit housed in the through-hole 1117 and spaced a predetermined distance from the vertical side wall of the extended portion 1114. The sensor 1135 is connected to summing node 410, as discussed above. Once the sensor 1135 is  
20       mounted a predetermined distance from the side wall of the extended portion 1114, the sensor lock screw 1118 is tightened against the side of the sensor to prevent movement of the sensor 1135 within the through-hole 1117. A preferred sensor is sold by Kaman  
25       Instrumentation of Colorado Springs, CO under model number SMU 9000-15N.

      Finally, an optional biasing spring 1140 may be connected between the pin 1012 formed in the bottom plate 1010 and the pin 1052 connected to the top plate  
30       1050.

      In operation, when no actuating voltage is applied to the piezoelectric element 1130, the optional biasing spring 1140 pulls the top plate toward the bottom plate until the top plate abuts and rests against the ball  
35       1116. In this position, the upper surface 1021 of the stationary bars 1020 and the upper surface 1031 of the

rotating bars 1030 are aligned such that the horizontal flexures 1040 are substantially planar. In addition, the inner side surfaces 1032 of the rotating bars 1030 and the side surfaces 1055 of the top plate 1050 are aligned such that the vertical flexures 1060 are substantially planar.

Upon application of an actuating voltage across the piezoelectric element 1130, the piezoelectric element 1130 presses against the side wall 1121 of the rotating block 1120 through the ball 1126, thereby causing the rotating block 1120 to rotate about a pivot portion 1062 of the flexure 1060 connecting the rotating block 1120 to the retainer block 1110. As the rotating block 1120 is rotated away from the retainer block 1110, the ball 1128 presses upward on the preload screw 1154, causing the top plate 1050 to move upward.

Upward movement of the top plate 1050 causes a rotation of the rotating bars 1030. As the rotating bars 1030 are rotated, the horizontally disposed flexures 1040 restrain the rotating bars 1030 such that they remain parallel with their associated stationary bars 1020. In addition, the vertically disposed flexures 1060 cause the top plate 1050 to remain parallel with each of the rotating bars 1030. As a result, the rotating bars 1030 act as torsion bars which prevent unwanted rotation or translation of the top plate 1050, thereby resulting in the upper surface of the top plate 1050 remaining parallel with the bottom plate 1010.

Returning to FIG. 1, the mirror control 124 of FIG. 1 is used to rotate X-mirror 106 and Y-mirror 108 such that laser beam 123 can scan more than a single point on target 112 while performing the coarse and/or fine auto-focus operations. Thus, if X-mirror 106 is rotated while Y-mirror 108 is held still, laser beam 123 will trace a line along the X-axis on target 112.

Similarly, X-mirror 106 can be held still while Y-mirror 108 is rotated, thereby tracing a line along the Y-axis on target 112.

In one embodiment of the present invention, the  
5 coarse and fine auto-focus operations are performed by scanning a line, rather than a spot, on target 112. In this embodiment, the line scan is performed at a frequency of approximately 8 Khz while target 112 is held stationary at one of several elevations in the  
10 range of movement along the Z-axis.

For an area scan method, X-mirror 106 and Y-mirror 108 can both be rotated to trace either a small area or selected parts of a larger area in the X-Y plane of target 112. The small area can be an area in the  
15 center of the field of view, while the larger area can be the entire field of view. During the area scan, microscope system 100 covers various features in the field of view in a path similar to a raster path of a television tube. To generate an electronic focus  
20 signal at a current Z-axis elevation, an area peak detector 1210 records and a microprocessor 403 reads the value generated by the highest amount of reflected light 123R that occurs during the area scan.

By using a line scan method or an area scan method  
25 and averaging the result of the scan to create electronic focus signal 115, the estimate of focus position 203 during a coarse auto-focus operation or a fine auto-focus operation becomes less sensitive to local height variations on the surface of target 112.

30 In a preferred embodiment, the area scan method is used for a fine auto-focus operation, although an area scan can also be used for a coarse auto-focus operation. In this embodiment, area peak detector 1210 (FIG. 12) is synchronously reset at the start of the  
35 area scan (anywhere on the area can be used as the start, as long as the reset is synchronous). FIG. 12

is similar or identical to FIG. 4, except for certain components that are described below. Area peak detector 1210 acquires or "loads up" during the rest of the area scan as follows. Area peak detector 1210 has  
5 a storage capacitor C67 (FIG. 14) that is discharged when area peak detector 1210 is reset. The voltage on storage capacitor C67 is continuously compared to the input voltage of op-amp U26. When op-amp U26's input voltage exceeds the stored voltage of capacitor C67,  
10 capacitor C67 is charged up to equal the input voltage. This process happens continuously, so that storage capacitor C67 records the highest input voltage to op-amp U26 no matter how short its duration, (within bandwidth limitations) that occurs after the reset.  
15 The only way that storage capacitor C67 discharges is through the reset mechanism or through leakage currents.

Just before being reset for the next area scan, a sampling ADC 5021 digitizes the output of area peak  
20 detector 1210. The output of ADC 5021 represents electronic focus signal 115 for the current elevation of target 112. ADC 5021 has a track-and-hold buffer in front (in one embodiment buffer is part of the architecture of ADC 5021), so that acquiring the output  
25 of area peak detector 1210 is virtually instantaneous. The ADC can digitize a voltage at the output of area peak detector 1210 at  $0.386 \text{ V}/\mu\text{s}$  without error, that translates to an acquisition time of about 25 ns. Area peak detector 1210 can be reset very quickly, using for  
30 example, a 500 nanosecond pulse. The time during which output of area peak detector 1210 is digitized and area peak detector 1210 reset is an extremely small portion of the overall area scan time so that only a very small portion of the area (e.g. 0.15%) is ignored while area  
35 peak detector 1210 is reset.

Normally an X-axis laser scanner (not shown) (also



referred to as "line scanner") that is used in microscope system 100 is resonant, while a Y-axis scanner (not shown), is a closed-loop galvo scanner that follows a "sawtooth" waveform (relatively slow scan followed by a fast flyback) during normal imaging that is unrelated to an auto-focus operation.

For an area scan used to find focus position 203, the Y-axis scanner (also referred to as "page scanner") follows a triangular wave profile 1301 (FIG. 13A) that is different from the "sawtooth" waveform at a much higher scan rate (e.g. 125 Hz) than normal (e.g. 13 Hz). The high scan rate allows microscope system 100 to trace a raster (e.g. sine wave) path in the area in the field of view quickly (8msec in one embodiment), albeit at less resolution (e.g. 32 lines/frame). At the end of an auto-focus operation, once target 112 has been positioned at focus position 203, the page scanner is returned to its normal slow sawtooth waveform.

Buffer 1230 connected to the output terminals of ADC 5021 allows signals on ADC 5021's output terminals to be turned on and held on without interfering with the operation of data bus 5006. This is required because ADC 5021 used in this embodiment cannot be configured to do a conversion without turning on signals at its output terminals. Conversion is initiated by hardware at a very precise predetermined time, at the end of an area scan. Signals that result from conversion remain active at the output terminals of ADC 5021 until microprocessor 403 detects that a scan is complete, reads ADC 5021 through buffer 1230 by asserting signal R\_ADC\* and resets ADC 5021 by pulsing signal RST-ZSYNC\*.

FIG. 13B illustrates events after the rising edge of signal ZSYNC. Signal ZSYNC rises at the end of each page scan to indicate completion of a scan cycle. The time scale in FIG. 13B is different from that of FIG.

13A. Signal ZSYNCD\*, that is a saved version of signal ZSYNC, goes low in response to a rising edge in signal ZSYNC. At the same time, signal CNVRT\* (FIG. 13B) also goes low, initiating an analog to digital conversion and causing signal ADC\_BSY\* (FIG. 13B) to go low.  
Signal BUSY1\* and signal BUSY2\* are delayed versions of signal ADC\_BSY\*, and are used to generate an active high pulse in signal RST\_PKDET (FIG. 13B) after signal ADC\_BSY\* returns high, to reset peak detector 1220 after ADC's conversion is completed. So ADC 5021 automatically starts a conversion when signal ZSYNC goes high, and area peak detector 1210 is automatically reset when the conversion is complete. Because there is a T/H in ADC 5021, the peak detector 1220 can be reset after the T/H has acquired, and before the conversion is complete.

Signal CNVRT\* stays low until signal RST\_ZSYNC\* falls, which keeps the conversion result present at the output terminals of ADC 5021. Software in microprocessor 5000 reads ADC 5021 (via signal R\_ADC\*) before resetting signal ZSYNCD\*. Signal ZSYNCD\* is an input signal to microprocessor 5000 and indicates that a frame is complete.

In one embodiment, there is a relatively long and unpredictable latency in the response of host workstation 116 to signal ZSYNCD\*. Hence most of the control signals (e.g. signals RST\_PKDET and CNVRT\*) for area peak detector 1210 and ADC 5021 are generated in hardware.

Area peak detector 1210 consists primarily of op-amps U26 and U27 (FIG. 14). The "hold" capacitor for area peak detector 1210 is C67. When the input voltage at pin 3 of op-amp U26 is greater than the hold voltage at capacitor C67, D9 is reverse-biased, so op-amp U26 has no feedback (at least momentarily). The output of op-amp U26 will therefore rise until op-amp U27

responds to the increased output of op-amp U26 and provides feedback via resistor R47. In this manner, when the input voltage to op-amp U26 is greater than the hold voltage, op-amp U26 drives the hold voltage higher to equal the input voltage. When the input voltage to op-amp U26 is less than the hold voltage, a diode D12 is reverse-biased and op-amp U26 receives a feedback signal through diode D9. Because diode D12 is reverse-biased, changes in the input voltage at op-amp U26 do not affect the hold voltage when op-amp U26's input voltage is less than the hold voltage.

Op-amp U27 buffers the hold voltage to minimize the current drawn from hold capacitor C67 (and hence the droop rate), and to provide feedback to op-amp U26. Transistor Q2 resets hold capacitor 67. Transistor Q1 guarantees that the input voltage of op-amp U26 is less than the hold voltage during reset (the hold voltage is zero during reset), so that transistor U26 does not slew positive during the reset operation. When the input voltage at U26 pin 3 is greater than the hold voltage at C67, D9 is reverse biased, so U26 has no feedback (at least momentarily). U26's output rises (slews positive) at the maximum rate possible for the particular op-amp (its "slow side") until U27 responds to the increased output of U26 and provides feedback via R47.

Diode D10 and resistor R46 establish a small negative voltage (approx. -0.4V in one embodiment) for clamping the input signal during reset. Resistors R53 and R47 slow down area peak detector 1210, to reduce overshoot that is common in conventional area peak detectors. Resistor R53 limits the rate at which capacitor C67 can be charged, resulting in an acquisition bandwidth of 2.3 MHz in one embodiment. Resistor R47 works with parasitic capacitance of U26 and D9 to form a feedback filter for composite op-amp

U26/U27. Resistor R47 acts to stabilize (reduce overshoot) of the composite op-amp. Both resistors R53 and R47 are determined empirically. Even with resistors R53 and R46, area peak detector 1210 provides an extremely fast response, and the combination of resistor R52 and capacitor C64 acts as a filter to limit the bandwidth of the input signal, so that high frequency noise does not pass through to area peak detector 1210.

An auto-focus routine AFPeakSync (page 98 of microfiche appendix B) that performs an autofocus operation during using an area scan method by a microscope system 100 is similar to auto-focus routine AFFast (page 97 of microfiche appendix B) described above (FIG. 8B). As seen in FIG. 8B, in a first fine pass, auto-focus routine AFFast steps from bottom to top of the fine Z-axis range of movement in large steps, then returns to a position below the coarse step that had the highest focus signal.

In a second fine pass, auto-focus routine AFFast then steps upward again, but this time in small steps (e.g. 0.098  $\mu\text{m}$ ). Auto-focus routine AFFast causes a large acceleration due to direction reversal just before the second fine pass by stepping to a position just below the estimated focus position (FIG. 8B). [Larger accelerations imply larger position errors in fine Z-stage 120 and are preferably avoided to improve positioning fidelity.]

Auto-focus routine AFPeakSync also uses two fine passes: a first fine pass and a second fine pass. The first fine pass of routine AFPeakSync (henceforth "AFPeakSync first fine pass") differs from the coarse pass of the coarse Z-stage in the manner described below. The AFPeakSync first fine pass operates similar to the coarse pass of auto-focus routine AFFast. There is no synchronization, and auto-focus routine

AFPeakSync steps and measures as quickly as possible, using large steps (e.g. 1.56 micrometer).

Software (Appendix F) in host workstation 116 sets up the page scanner (not shown) to perform the area scan method during the AFPeakSync first fine pass.  
5 Routine lonuiSuperfineAF sets up all parameters for use in the area scan method. In routine lonuiSuperFineAF, function lonui\_Get\_IndexFrPixZoomEnum returns with the current X and Y scanner parameters, such as amplitude  
10 of oscillation. After the scanner parameters are set, host workstation 116 instructs fine Z-axis controller 118 to execute a fine auto-focus operation. After the fine auto-focus operation is done, host workstation 116 returns scanner parameters to their original values.

15 In alternative embodiments, the page scanner can execute a line scan method or a spot method during the AFPeakSync first fine pass. The AFPeakSync first fine pass is merely used to move target 112 close enough to focus position 112 so that a AFPeakSync second fine  
20 pass can correct residual errors.

Auto-focus routine AFPeakSync avoids the large acceleration of auto-focus routine AFFast in a second fine pass by stepping to a position above the estimated focus position at the end of a first fine pass and  
25 continuing in small steps for the second fine pass in the same direction, as illustrated by FIG 8C.

In one embodiment, a second fine pass in auto-focus routine AFPeakSync covers twice the range of a fine pass in auto-focus routine AFFast, which reduces  
30 sensitivity to manufacturing variations and improves performance of the autofocus operation. The increased range of the AFPeakSync second fine pass compensates for any position errors in the AFPeakSync first fine pass that could result in a AFPeakSync step with the  
35 highest focus strength being off by one rough step.

Auto-focus routine AFPeakSync spends two page

scanner cycles at each position of the AFPeakSync second fine pass. The first cycle allows the fine Z-axis (signal ACTUAL Z POSITION in FIG. 13C) to settle out, and the second cycle acquires the peak magnitude of electronic focus signal 115 from the page scan at the current elevation. Another embodiment of an auto-focus routine can spend only one page scanner cycle by allowing the fine Z movement to settle during a trace 1310 and by acquiring the peak during the retrace 1320.

10       Auto-focus routine AFPeakSync executes in about 1.5 seconds, compared to 0.9 seconds for auto-focus routine AFFast because of the time required to perform an area scan at each of several elevations of target 112. In spite of such an increased execution time, a  
15       microscope system that uses auto-focus routine AFPeakSync is faster overall as compared to a microscope system that uses auto-focus routine AFFast, because with increased reliability of the auto-focus operation using routine AFPeakSync, manual focus  
20       adjustments are needed only rarely.

      The first fine pass of auto-focus routine AFPeakSync steps as quickly as possible irrespective of where the scanner is during sampling. In the second fine pass, however, auto-focus routine AfPeakSync runs  
25       the page scanner at a high rate (>100Hz) and allows one page scanner cycle for fine Z-stage 120 to settle, before doing data acquisition on the next cycle, as illustrated in FIG. 13C. Auto-focus routine AFPeakSync therefore spends two page scanner cycles at each Z  
30       position (e.g. elevation). Furthermore, in area peak detector mode, an analog-to-digital converter 5021 is automatically started when a rising edge in signal ZSYNC occurs, and area peak detector 1210 is automatically reset. So, auto-focus routine AFPeakSync  
35       causes host microprocessor 5000 to wait for a rising edge in signal ZSYNC and then read the peak magnitude

detected from the previous page scan from analog-to-digital converter 5021. Resetting signal ZSYNC also resets analog-to-digital converter 5021 and so ADC 5021 is read before resetting signal ZSYNC.

5        In peak detector mode, analog-to-digital converter 5021 is controlled only by a rising edge in signal ZSYNC. An asynchronous conversion is not possible without resetting the peak detector mode bit. Peak  
10       detector 1210 obtains a representation of the highest intensity that occurred in a video frame at a given Z elevation. Peak detector 1210 allows microscope system  
100 to repeatably focus on the layer with the highest reflectivity.

      Microscope system 100 can use programmable offset  
15       FocOff to allow the user to select any predetermined layer to focus on. In one embodiment, microscope system 100 initially focuses on a layer ("brightest layer") of target 112 that generates the largest electronic focus signal. The user can adjust the  
20       position of target 112 to bring a predetermined layer into focus. Microscope system 100 records the target's adjustment as an offset from the brightest layer. In subsequent auto-focus operations, microscope system 100 focuses on the brightest layer and then automatically  
25       moves target 112 through the adjustment offset that was specified by the user, to focus on the predetermined layer.

      As the relative reflectivities of the layers in a wafer are repeatable and independent of the  
30       illumination power, the relative positions of a predetermined layer chosen by a user is at a repeatable offset from the brightest layer. So microscope system 100 can repeatably focus on any predetermined layer after a single adjustment. Microscope system 100  
35       focuses on such a predetermined layer as long as the layer structure stays constant, for example for array

type structures, such as RAM arrays. If, after selecting a predetermined layer, areas off a RAM array are to be imaged, a simple re-adjustment allows the user to continue.

5 In one embodiment, the function to specify the offset is implemented as a system function that can be mapped by a user to any function key of host workstation 116. Once mapped to a function key, the offset function can be accessed by the user from the  
10 keyboard at any time. When so accessed, host workstation 116 compares the current target position with the focus position estimate of last auto-focus operation, calculates the difference and stores the difference as the offset for finding the predetermined  
15 layer in future auto-focus operations.

The offset is limited to +/- 5 microns, to make it easier to recover from an improperly set offset (a typical semiconductor topology is less than 1  $\mu\text{m}$ ) and only fine Z-stage 120 is moved during offset  
20 calculation. The offset is calculated assuming coarse Z-stage 122 has not moved since last auto-focus operation.

While the present invention has been described in connection with specific embodiments, variations on  
25 these embodiments will be obvious to those having ordinary skill in the art. For example, target 112 may be moved by means other than a stepper motor or a piezoelectric element, such as a linear voice coil motor or an electrostrictive actuator. Furthermore,  
30 target 112 can initially be positioned above the focus position, with the first pass in any of the embodiments beginning by moving target 112 in a negative Z direction. In addition, although the present invention was described in connection with a microscope that  
35 reflects a maximum intensity to the photodetector during a focused condition, it is clear that the



invention may be modified to operate with a microscope that reflects a minimum intensity to the photodetector during a focused condition.

Furthermore, although the invention, as described,  
5 utilizes a laser beam 123 to perform both the auto-focusing and imaging operations, it is understood that a confocal laser optical system could be used to perform an auto-focusing operation for a non-confocal microscope system which utilizes only a white light  
10 source to perform the imaging operation. Such an application is advantageous because a non-confocal microscope system utilizing a white light imaging source results in a focus signal which is sinc function, rather than a sinc squared function. Because  
15 the sinc function does not exhibit a peak which is as pronounced as the sinc squared function, it is more difficult to determine the focus position using a focus signal generated by a non-confocal, white light optical system. Therefore, using a confocal laser optical  
20 system to generate the focus signal used to perform the auto-focus operation in a white light microscope results in a superior electronic focus signal, thereby allowing the auto-focus operation to be performed with greater precision. Confocal white light can also be  
25 used, with scanning replaced by detector array, such as a camera.

Although a certain number and types of passes are described above, other number and types of passes can also be used. For example a microscope system 100 can  
30 perform a single pass in a coarse auto-focus operation using a median point method and adjust the gain of photodetector 114 followed by two passes in a fine auto-focus operation to estimate focus position 203.

Various embodiments of the invention described  
35 above are encompassed by the attached claims.

## APPENDIX A

```
5 // Fast Z-Axis Controller (P/N 000491) Test Program

// Source File: Fastz.nc
// Author: Chris F.
// Target Board: 000491, Fast Z-Axis Controller
10 // Target Socket: M1-U2
// Ultrapointe Corp.
// PROPRIETARY! ALL RIGHTS RESERVED
// **** FOR INTERNAL USE ONLY ****
//

15 // The FastZ PCB controls a piezo-driven stage with analog position
// feedback. The position command
// is a digital 12-bit word written by the Neuron, and is double-
// buffered in front of the position command DAC. The position loop
20 // is closed in analog hardware. The position command can be updated
// in one of two modes: in sync mode, the Neuron writes to the DAC
// buffer, and the DAC receives the updated command when a ZSYNC
// pulse is recieved from the page scanner controller (P/N 000134);
// in async mode, the Neuron writes to the DAC buffer and then forces
25 // the command into the DAC.
// For autofocus purposes, there is a latching video comparator on
// the board, with a programmable threshold (similar to that on the
// PMT Amplifier (P/N 000276). There is also a fast 8-bit ADC for use
// by autofocus routines involving the stepper motor Z-axis and the
30 // fast (piezo) Z-axis. This same ADC can be used to read the high
// voltage driving the piezo crystal, and to read back the analog
// position feedback.
// Finally, an emergency retract input is available that forces the
// position command to zero. This feature is not used at present.
35 //
// There is a custom bus structure on the FastZ PCB that allows
```

```
// Neuron to address many devices with an 8-bit bi-directional bus,  
// yet uses only the general purpose I/O lines of the device. This  
// allows a general purpose Neuron module to be used because the main  
// data/address buses of the Neuron are not required. In this scheme,  
5 // I/O bits 0-7 are used as the data bus, bit 8 is a read/not-write  
// line, bit 9 is an address/not-data line, and bit 10 is a low-true  
// strobe line. In general, the Neuron will write an address to the  
// address register, then read or write the intended data value; the  
// strobe line executes the intended operation after address and data  
10 // have been appropriately set up.  
// The following rules are adhered to by the hardware and should  
// be respected by software:  
// 1. No one drives the data bus unless STB* (low-true STroBe) is  
// asserted. The only exception is that the Neuron drives the  
15 // bus with an address value without asserting STB* when  
// ADDR/DATA* is set high.  
// 2. Setting the ADDR/DATA* (ADDReSS/not-DATA) line forces a  
// write operation (the RD/WR* line is don't care when  
// ADDR/DATA* = 1).  
20 // 3. The address latch is transparent, and controlled by ADDR/DATA*,  
// so no strobe is required to write the address.  
// 4. Because of rule 1, no one should read the bus until STB* is  
// low.  
// This bus arrangement means a lot of bit twiddling to do a read or  
25 // write, but if the address does not change between operations, it  
// can be speeded up considerably. FastRead and FastWrite routines are  
// provided that do not set up addresses are provided for this purpose.
```

```
#include <control.h>
```

30

```

//*****
//***** DEFINITIONS *****
//*****

// Constants

5
// Mode Commands:
#define SELF_TEST      0x00      // commands from host
#define AUTO_FOCUS     0x10
#define READ_VIDEO     0x11
10 #define READ_FDBK     0x12
#define READ_HV        0x13
#define SET_FOR_VOL    0x14
#define MONITOR_FOCUS  0x15
#define ABORT_CMD      0x16
15 #define RD_STATUS     0x17
#define AF_FULLRANGE   0x18      // af function from 0 to 4095 tvl
#define AF_PMTVALS     0x19      // new af with pmt reading storage tvl
#define AF_FAST_T      0x1A      // fast af algorithm for testing tvl
#define AF_HALFRANGE   0x1B      // af function from 1024 to 3036 tvl
20 #define AFFR_FAST     0x1C      // fast af function from 0 to 4095 tvl
#define AF_FAST        0x1D      // fast af algorithm tvl
#define NULL_MODE      0xFF      // used to prevent repeated mode cmds

// status codes
25 #define SUCCESS       0
// autofocus responses:
#define NO_FOCUS_FF    1          // focus FF never tripped
#define TOO_FAR        2          // never saw increasing signal
#define NOT_FAR_ENUF    3          // never saw decreasing signal
30 #define END_OF_TRAVEL 4          // success, but at end of travel

// volume mode responses:
#define RDY_FOR_VOL     5
#define VOLUME_COMPL    6
#define UNDERFLOW      7

```

```

// coarse autofocus responses:
#define COMPARATOR_READY      8
#define THRESH_TOO_LOW       9
#define COMP_TRIPPED          0x0A
5  #define ABORTED             0x0B
// misc
#define CMD_UNKNOWN           0xFF

// Register addresses: note that two MSBs of address are mux
10 // select bits.
// Write-Only addresses:
#define CTRL_REG              0           // Control register address
#define THRESH_DAC            1           // Autofocus threshold DAC
#define POS_CMD_LO            2           // Position command lo byte
15 #define POS_CMD_HI          3           // Position command hi byte
// Read-Only addresses:
#define STS_REG               0           // Status Register
#define ADC_ADDR              1           // ADC

20 // mux addresses: logical OR with register addresses
#define SEL_VIDEO              0           // video for autofocus
#define SEL_FDBK               0x40       // analog position feedback
#define SEL_HV                 0x80       // high voltage output

25 // An 'n' suffix on a signal name indicates low true
// control register bits:
#define SYNC_MODE              0x01       // set to one for sync DAC update
#define ASYNC_TRIG             0x02       // pulse hi to update DAC in async mode
#define RST_ZSYNCn             0x04       // pulse low to reset ZSYNC flip flop
30 #define RST_FOCUSn           0x08       // pulse low to reset FOCUS flip flop
#define RST_RETRACTn           0x10       // pulse low to reset EMERG RETRACT
#define RST_INTEGn             0x20       // hold low to turn off position loop
// integrator
#define SPARE_TEST             0x40

```

```

#define TEST_LEDn      0x80  // set low to illuminate TEST LED

// status register bits:
#define ZSYNCDn        0x01  // low means ZSYNC has occurred
5  #define FOCUSDn      0x02  // low means FOCUS has occurred
#define RETRACTDn      0x04  // low means EMERG RETRACT has occurred
#define ADC_BUSYn      0x08  // low means conversion in process
#define TEST_IN        0x20
#define LED_TST_IN     0x40
10  #define TEST_JMPRn   0x80  // low means test jumper is in place

// control line defintions:
#define READ            1      // for RD_WRn line
#define WRITE          0
15  #define ADDR        1      // for ADDR_DATAn line
#define DATA          0

// autofocus constants:
#define MAX_Z           4095  // max position DAC range
20  #define AF_FAST_DELTA 10   // step size for fast autofocus sweep
#define AF_SLOW_DELTA   1     // step size for final autofocus
#define Z_INSURANCE     300   // distant to step down before fine AF
#define MIN_DELTA       15    // min delta from ADC to indicate slopes

25  // intial values
#define ZCTRL_INIT      0x22;
#define INIT_THRESH     0x16;

30  /*******
    /******* PRAGMAS *****/
    /*******
    #pragma net_buf_in_count 7
    #pragma net_buf_out_count 7

```

```

//*****
//***** I/O DECLARATIONS *****
//*****

5  IO_0 output  byte  WRITE_DATA;      // parallel port
   IO_0 input   byte  READ_DATA;

   IO_8 output  bit   RD_WRn;          // ReaD/not-WRite
   IO_9 output  bit   ADDR_DATAn;     // ADDRess/not-DATA
10  IO_10 output bit   STBn;           // low-true STroBe

//*****
//***** NETWORK VARIABLES *****
15 //*****

   network output polled unsigned SW_Rev = SWR;    // software version
   network output unsigned ResetSts;              // for indicating reset
   network input  unsigned OpMode=0;              // for remote proc. calls
   network output unsigned OpStatus;              // response to calls
20  network input  unsigned OpParArray[2];         // Operation parameters
   network output unsigned OpStsArray[2];         // Extra Status info

   network input  unsigned AFocusThresh=16;        // 0-255 for Autofocus threshold
   network input  unsigned long ZPosCmd=0;         // Position Command 0-4095
25  network input  unsigned Z_Control;             // Control byte
   network input  unsigned MuxSel=0;              // ADC mux selector
   network input  unsigned ZStep=4;               // step size for vol acquire
   network input  unsigned NumFrames=64;          // # frames for volume

30  network output polled unsigned Z_Status;        // Status byte
   network output polled unsigned long PosStatus;  // more status

struct PMTARRAY          // define structure to hold pmt values   tvt
{

```

```

    unsigned    pmt[16];
};

// declare 16 network variables as structure of PMTARRAY type
5 // to hold 256 pmt readings (one every 16 steps of the fast z stage tv

network output polled struct PMTARRAY pmtvals1;
network output polled struct PMTARRAY pmtvals2;
network output polled struct PMTARRAY pmtvals3;
10 network output polled struct PMTARRAY pmtvals4;
network output polled struct PMTARRAY pmtvals5;
network output polled struct PMTARRAY pmtvals6;
network output polled struct PMTARRAY pmtvals7;
network output polled struct PMTARRAY pmtvals8;
15 network output polled struct PMTARRAY pmtvals9;
network output polled struct PMTARRAY pmtvals10;
network output polled struct PMTARRAY pmtvals11;
network output polled struct PMTARRAY pmtvals12;
far network output polled struct PMTARRAY pmtvals13;
20 far network output polled struct PMTARRAY pmtvals14;
far network output polled struct PMTARRAY pmtvals15;
far network output polled struct PMTARRAY pmtvals16;

//*****
25 //***** GLOBAL VARIABLES *****
//*****

//*****
30 //***** FUNCTION DEFINITIONS *****
//*****

unsigned short ReadPort (unsigned short TargetAddr)
// read a port, including set up address
{

```



```
unsigned short result;
```

```

    io_set_direction (WRITE_DATA, IO_DIR_OUT);
    io_out (ADDR_DATAn, ADDR);           // set ADDR_DATAn to ADDR
5   io_out (WRITE_DATA, TargetAddr | MuxSel); // write address
    io_out (ADDR_DATAn, DATA);         // set ADDR_DATAn to DATA
    io_set_direction (READ_DATA, IO_DIR_IN);
    io_out (RD_WRn, READ);              // set RD_WRn for READ
    io_out (STBn, 0);                  // assert STBn
10  result = io_in (READ_DATA);         // read data
    io_out (STBn, 1);                  // de-assert STBn
    return (result);                   // return result
}
```

15

```
void WritePort (unsigned short TargetAddr, unsigned short WriteData)
```

```
// write to a port, including set up address
```

```
{
```

```

20  io_set_direction (WRITE_DATA, IO_DIR_OUT);
    io_out (ADDR_DATAn, ADDR);           // set ADDR_DATAn to ADDR
    io_out (WRITE_DATA, TargetAddr | MuxSel); // write address
    io_out (ADDR_DATAn, DATA);         // set ADDR_DATAn to DATA
    // no need to change port direction for write
25  io_out (RD_WRn, WRITE);              // set RD_WRn for WRITE
    io_out (WRITE_DATA, WriteData);     // set up data
    io_out (STBn, 0);                  // assert STBn
    io_out (STBn, 1);                  // de-assert STBn
}
```

30

```
unsigned FastRead (void)
```

```
// read a port, without address set up
```

```
{
```

```

    unsigned result;

    io_out (STBn, 0);                // assert STBn
    result = io_in (READ_DATA);      // read data
5   io_out (STBn, 1);                // de-assert STBn
    return (result);                 // return result
}

10 void FastWrite (unsigned short WriteData)
    // write to a port, without address set up
    {

    io_out (WRITE_DATA, WriteData); // set up data
15   io_out (STBn, 0);                // assert STBn
    io_out (STBn, 1);                // de-assert STBn
    }

20 void SetAddressRead (unsigned short TargetAddr)
    // Setup address for fast read
    {

    io_set_direction (WRITE_DATA, IO_DIR_OUT);
25   io_out (ADDR_DATAn, ADDR);        // set ADDR_DATAn to ADDR
    io_out (WRITE_DATA, TargetAddr | MuxSel); // write address
    io_out (ADDR_DATAn, DATA);        // set ADDR_DATAn to DATA
    io_set_direction (READ_DATA, IO_DIR_IN);
    io_out (RD_WRn, READ);             // set RD_WRn for READ
30   }

void SetAddressWrite (unsigned short TargetAddr)
    // Setup address for fast read

```

```

{

    io_set_direction (WRITE_DATA, IO_DIR_OUT);
    io_out (ADDR_DATAn, ADDR);           // set ADDR_DATAn to ADDR
5   io_out (WRITE_DATA, TargetAddr | MuxSel); // write address
    io_out (ADDR_DATAn, DATA);         // set ADDR_DATAn to DATA
    // no need to change port direction for write
    io_out (RD_WRn, WRITE);             // set RD_WRn for WRITE
}

10

void SetCtrl (unsigned NewCtrl)
// sets the control word
{
15   Z_Control = NewCtrl;
    WritePort (CTRL_REG, NewCtrl);
}

20 void MoveZ (unsigned long target)
// move the Z-axis to the current commanded position (ZPosCmd)
{
    if (target > MAX_Z)
        target = MAX_Z;
25   ZPosCmd = target;
    WritePort (POS_CMD_LO, (unsigned int)(ZPosCmd & 0xFF)); // set position
    WritePort (POS_CMD_HI, (unsigned int)(ZPosCmd >> 8) & 0xFF);
    SetCtrl (Z_Control & (~ASYNC_TRIG) & (~SYNC_MODE)); // clock DAC
    SetCtrl (Z_Control | ASYNC_TRIG);
30 }

```

## ORIGINAL ALGORITHM

```

unsigned AFocusServo (unsigned Multiplier)
// Parameter is a multiplier for the speed of autofocus, and should
// be changed as a function of the numerical aperature of the
5 // objective, or as a function of the depth of focus. The depth of
// focus is inversely proportional to the square of the numerical
// aperature. For an NA of 1, the multiplier is 1, for a NA of
// 0.5, the multiplier should be  $1/((0.5)^2) = 4$ , etc.
//
10 // The autofocus routines use a latching focus comparator with a
// programmable threshold (the 'focus flip-flop'), and a fast 8-bit
// ADC. Both of these devices monitor the video signal from the
// PMT preamplifier, but the focus flip flop uses a high-bandwidth
// video signal, and the ADC uses a low bandwidth version. The
15 // reason for this is that the focus flip flop must detect the
// focus condition, no matter how short in time it persists, but
// the ADC is used for the final autofocus operations and should
// deliver an average value of the video signal, representing the
// video level over a spacial 'patch' of the target image. Since the
20 // ZSYNC signal is also present on the FastZ board, it is possible
// to sample the ADC synchronous with the video frame and thereby
// obtain a repeatable autofocus. There is no attempt at this time
// to do that - we merely rely on the bandwidth of the video input
// to the ADC being matched to the video frame rate. Of course, the
25 // video frame rate is adjustable, but the ADC bandwidth is not, so
// an approximation is the best we can do.
//
// Autofocus for the Fast Z-axis works as follows:
// 1. The stage is moved to its extreme lower limit (ZPosCmd = 0).
30 // 2. The autofocus threshold is appropriately set. 1/16 of full
// scale might be a good starting value; we expect the autozero
// routine to set the 'black' level to approx. 1/64 of full
// scale. Note that the PMT gain must be appropriately set
// for the incident laser power (and wavelength) and type of

```

```

//    substrate being scanned. An autogain cycle is expected after
//    autofocus completes.
//    3. The focus flip flop is cleared.
//    4. The stage is moved up a step at a time (but rapidly) while
5    //    the focus flip flop is monitored.
//    5. When the focus flip flop trips, the stage is stopped. We expect
//    that stage may have to move up further to attain peak focus,
//    but depending on many factors (e.g. servo bandwidth, stage speed,
//    the rate at which the focus FF is sampled), the stage could
10   //    have moved beyond focus (this is almosta certain for the
//    stepper motor autofocus).
//    6. The stage is therefore moved back down a fixed amount to
//    ensure it is below the focus position.
//    7. The stage is then moved up slowly while the video signal is
15   //    monitored via the ADC.
//    8. When a local peak is detected above a second (software)
//    threshold, the movement stops and the stage returns to the peak.
//    9. If a peak is found, and the stage ends up at a point other
//    than the end points, the autofocus operation is considered
20   //    successful. The final commanded position and success/failure
//    are status variables.
{
    unsigned short BaseVal, LastVal, NextVal, BestVal; // ADC values
    unsigned long BestPos;           // position of peak signal
25   unsigned long Move_Delta;
    short Delta;
    boolean PkFound;
    unsigned Stat;

30   if (Multiplier == 0)
        Multiplier = 1;
    SetCtrl (Z_Control & (~ SYNC_MODE)); // async mode
    // move to the extreme lower limit
    MoveZ (0);

```

```

delay(1977);                      // wait 50 ms
// clear the focus FF
SetCtrl (Z_Control & (~RST_FOCUSn));
SetCtrl (Z_Control | RST_FOCUSn);
5  Move_Delta = AF_FAST_DELTA * Multiplier;
// we assume the autofocus flip flop threshold has been set
if ((ReadPort (STS_REG) & FOCUSDn) != 0)
{
    // if focus FF not tripped immediatley
    // note that this loop should look at RETRACTDn, if it were used
10  while (((ReadPort (STS_REG) & FOCUSDn) != 0) &&
        ((ZPosCmd + Move_Delta) < MAX_Z))
    {
        MoveZ (ZPosCmd + Move_Delta);
        watchdog_update();
15  }
    if ((ZPosCmd + Move_Delta) < MAX_Z) // if focus FF tripped
    {
        if (ZPosCmd - Z_INSURANCE >= 0)
            MoveZ (ZPosCmd - Z_INSURANCE);    // back off
20  else
            MoveZ (0);
        delay (4000);                      // wait 100 ms
        // now move up looking for peak signal
        PkFound = FALSE;
25  MuxSel = SEL_VIDEO;
        NextVal = ReadPort (ADC_ADDR);
        LastVal = NextVal;
        BaseVal = NextVal;
        BestVal = NextVal;
30  BestPos = ZPosCmd;
        Move_Delta = AF_SLOW_DELTA * Multiplier;
        // start by moving up until get an increasing signal
        while ((PkFound==FALSE) && ((ZPosCmd + Move_Delta)<MAX_Z))
        {

```

```

        MoveZ (ZPosCmd + Move_Delta);
        LastVal = NextVal;
        NextVal = ReadPort (ADC_ADDR);
        if (NextVal > BestVal)
5           {
                BestPos = ZPosCmd;
                BestVal = NextVal;
            }
        Delta = NextVal - BaseVal;
10        if (Delta > MIN_DELTA)
                PkFound = TRUE;
                watchdog_update();
            }
        if (PkFound == TRUE)          // if did find a rising slope
15        {
                // now keep going until negative slope
                BaseVal = BestVal;
                PkFound = FALSE;
                while ((PkFound == FALSE) && ((ZPosCmd + Move_Delta) <
20        MAX_Z))
                {
                        MoveZ (ZPosCmd + Move_Delta);
                        LastVal = NextVal;
                        NextVal = ReadPort (ADC_ADDR);
25                        if (NextVal > BestVal)
                                {
                                        BestVal = NextVal;
                                        BaseVal = NextVal;
                                        BestPos = ZPosCmd;
30                                }
                        Delta = BaseVal - NextVal;
                        if (Delta > MIN_DELTA)
                                PkFound = TRUE;
                                watchdog_update();

```

```

    }
    if (PkFound == TRUE)
    {
        if ((ZPosCmd + Move_Delta) < MAX_Z)
5         Stat = SUCCESS;
        else
            Stat = END_OF_TRAVEL;
    }
    else
10     Stat = NOT_FAR_ENUF;
    }
    else
        Stat = TOO_FAR;
        MoveZ (BestPos);
15     }
    else
        Stat = NO_FOCUS_FF;
    }
    else
20     Stat = THRESH_TOO_LOW;
    return (Stat);
}

unsigned long AFFullRange(void) // tvf
25 // Function to step through entire fast z stage range at 16 counts per step
// Stage is returned to fast z position which gives the highest PMT value.
{
    unsigned short NextVal, BestVal; // ADC values
    unsigned long BestPos; // position of peak signal
30 unsigned long step; // z stage step loop counter
    static int FocOff; // Focus offset

    // Set FocOff
    FocOff = OpParArray[0];

```



```

SetCtrl (Z_Control & (~SYNC_MODE)); // async mode
// move to the extreme lower limit
MoveZ (0); // Start at the bottom
delay(1977); // wait 50 ms

5 // step through entire Z range checking for peak PMT and storing values

NextVal=ReadPort(ADC_ADDR); // Store current PMT value
BestVal=NextVal; // Clear Largest value
10 BestPos=0; // Clear Best Position

for (step=0; step <= 4096; step += 16) // Entire range in steps of 16
{
    if(step == 0) step = 1;
    15 MoveZ(step-1); // Move 0 to 4095
    watchdog_update();
    NextVal = ReadPort(ADC_ADDR); // Store current PMT value

    if (NextVal > BestVal) // If Current value is larger
    { // than the best (largest) value
        20 BestPos=step; // this is the best position
        BestVal=NextVal; // Update best value
    }
}

25 OpStatus = BestVal; // send the largest value out
MoveZ(BestPos-1-FocOff); // move to the best position
return BestPos-1-FocOff; // return the best positon
}

30

```

#### NEW ALGORITHM

```

unsigned long AFFast(void) // tvf
// Fast Autofocus Function: This funciton does a coarse then a fine autofocus
//

```

```

{
    unsigned short NextVal, BestVal;    // ADC values
    unsigned long BestPos;              // position of peak signal
    unsigned long AvgVal;               // average value storage
5   unsigned long FineCenter;          // center of fine focus range
    unsigned long Move_Delta;           // steps per moove
    unsigned long step;                 // z stage step loop counter
    unsigned m;                         // array member
    unsigned l;                         // loop counter
10  static int FocOff;                 // Focus offset
    unsigned int StepSize;              // step size during coarse focus

    // Set FocOff and StepSize
    FocOff= OpParArray[0];
15  StepSize = 128;

    SetCtrl (Z_Control & (~SYNC_MODE)); // async mode
    // move to the extreme lower limit
    MoveZ (0);                          // Start at the bottom
20  delay(1977);                       // wait 50 ms

    // ***** coarse autofocus *****
    // step through entire Z range with large steps
    // checking for peak PMT and storing values
25  AvgVal=0;
    MuxSel = SEL_VIDEO;
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
30  AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    NextVal=AvgVal/4;
    BestVal=NextVal;                    // Clear Largest value
    BestPos=0;                         // Clear Best Position

```

```

for (step=0; step <= 4096; step += StepSize) // Entire range in steps of 1
{
    if(step == 0) step = 1;
    MoveZ(step-1); // Move 0 to 4095
5    watchdog_update();
    AvgVal=0;
    MuxSel = SEL_VIDEO;
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
10    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    NextVal=AvgVal/4;

    if (NextVal > BestVal) // If Current value is larger
15    { // than the best (largest) value
        BestPos=step; // this is the best position
        BestVal=NextVal; // Update best value
    }
}
20

// Move to start value for fine autofocus and wait
if (BestPos > 129)
    MoveZ (BestPos-128-1); // Start at the bottom
else
25    MoveZ (1);
    delay(395); // wait 10 ms

// Set up comparison values
    AvgVal=0;
30    MuxSel = SEL_VIDEO;
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value

```

```

    NextVal = AvgVal/4;
    BestVal = NextVal;                                // Clear Largest value

    // ***** fine autofocus *****
5    // using position from coarse focus
    FineCenter = BestPos; // Set center of focus to the best coarse focus pos.
    if(FineCenter < 129) // Make sure it's valid (at least 129)
        FineCenter = 129;
    for (step = FineCenter-128; step <= FineCenter+128; step += 8) // Entire r
10    {
        MoveZ(step-1);                                // Move through fine af range
        watchdog_update();
        AvgVal = 0;
        MuxSel = SEL_VIDEO;
15    AvgVal = AvgVal + ReadPort(ADC_ADDR); // Store current PMT value
        AvgVal = AvgVal + ReadPort(ADC_ADDR); // Store current PMT value
        AvgVal = AvgVal + ReadPort(ADC_ADDR); // Store current PMT value
        AvgVal = AvgVal + ReadPort(ADC_ADDR); // Store current PMT value
        NextVal = AvgVal/4;
20
        if (NextVal > BestVal) // If Current value is larger
        { // than the best (largest) value
            BestPos = step; // this is the best position
            BestVal = NextVal; // Update best value
25    }
    }

    OpStatus = BestVal; // send the largest value out
    MoveZ(BestPos-1-FocOff); // move to the best position
30    return BestPos-1-FocOff; // return the best position
}

```

```

unsigned long AFFastTest(void) // tvf

```

```

// Fast Autofocus Function: This function does a coarse then a fine autofocus.

```

```

// It is a test version of AFFast w/ variable averaging and step size and
// recording of PMT Values in network variables.

{
    PROFILE ALGORITHM
    unsigned short NextVal, BestVal;    // ADC values
5    unsigned long AvgVal;              // Average PMT Value
    unsigned long BestPos;              // position of peak signal
    unsigned long FineCenter;           // center of fine focus range
    unsigned long Move_Delta;           // steps per moove
    unsigned long step;                 // z stage step loop counter
10    unsigned m;                       // array member
    unsigned l;                         // loop counter
    unsigned a;                         // averaging setting
    static int FocOff;                  // Focus offset
    unsigned int StepSize;              // step size during coarse focus
15
    // Set FocOff and StepSize
    a = OpParArray[0];
    StepSize = OpParArray[1];

20    // Make sure StepSize is good (> 7)
    if (StepSize < 8)    // must be at least 8
        StepSize = 8;

    // Make sure averaging is good (> 3)
25    if (a < 4)          // must be at least 4
        a = 4;

    SetCtrl (Z_Control & (~SYNC_MODE)); // async mode
    // move to the extreme lower limit

30    MoveZ (0);           // Start at the bottom
    delay(1977);          // wait 50 ms

    // ***** coarse autofocus *****
    // step through entire Z range with large steps

```

```

// checking for peak PMT and storing values
    AvgVal=0;
    for (l=0;l<a;l++)
    {
5      MuxSel = SEL_VIDEO;
      AvgVal=AvgVal+ReadPort(ADC_ADDR);    // Store current PMT value
    }
    NextVal=AvgVal/a;
    BestVal=NextVal;                // Clear Largest value
10   BestPos=0;                    // Clear Best Position

    for (step=0; step <= 4096; step += StepSize) // Entire range in steps of 1
    {
        if(step == 0) step = 1;
15     MoveZ(step-1);              // Move 0 to 4095
        watchdog_update();
        AvgVal=0;
        for(l=0;l<a;l++)
        {
20         MuxSel = SEL_VIDEO;
            AvgVal = AvgVal+ReadPort(ADC_ADDR);    // Get current PMT value
        }
        NextVal=AvgVal/a;

25         if (NextVal > BestVal)                // If Current value is larger
        {                                       // than the best (largest) value
            BestPos=step;                      // this is the best position
            BestVal=NextVal;                   // Update best value
30     }

    /* Write PMT Values to nv's */

    if (StepSize == 128)

```

```

    {
        if (step < 2048) // use structure 1
        {
            m = step/StepSize;
5           pmtvals1.pmt[m] = NextVal;
        }

        else if ((step >= 2048) && (step < 4096)) // use structure 2
        {
10          m = (step-2048)/StepSize;
            pmtvals2.pmt[m]=NextVal;
        }
    }

15 }

// Move to start value for fine autofocus and wait
if (BestPos > 129)
    MoveZ (BestPos-128-1); // Start at the bottom
20 else
    MoveZ (1);
    delay(395); // wait 10 ms

// Set up comparison values
25 AvgVal=0;
    for(l=0;l<a;l++)
    {
        MuxSel = SEL_VIDEO;
        AvgVal = AvgVal+ReadPort(ADC_ADDR); // Get current PMT value
30    }
    NextVal=AvgVal/a;
    BestVal=NextVal; // Clear Largest value

```

```

// ***** fine autofocus *****
// using position from coarse focus
FineCenter=BestPos; // Set center of focus to the best coarse focus pos.
if(FineCenter<129) // Make sure it's valid (at least 129
5   FineCenter=129;
for (step = FineCenter-128; step <= FineCenter+128; step +=8)
// go through fine autofocus range in steps of 8
{
    MoveZ(step-1); // Move through fine af range
10   watchdog_update();
    AvgVal=0;
    for(l=0;l<a;l++)
    {
        MuxSel = SEL_VIDEO;
15   AvgVal = AvgVal+ReadPort(ADC_ADDR); // Get current PMT value
    }
    NextVal=AvgVal/a;

    if (NextVal > BestVal) // If Current value is larger
20   { // than the best (largest) value
        BestPos=step; // this is the best position
        BestVal=NextVal; // Update best value
    }

25   /* Write PMT Values to nv's */

    if (step < FineCenter) // use structure 3
    {
        m = 16-((FineCenter-step)/8);
30   pmtvals3.pmt[m] = NextVal;
    }

    else if (step >= FineCenter) // use structure 4
    {

```



```

        m = (step-FineCenter)/8;
        pmtvals4.pmt[m]=NextVal;
    }

5    }

    OpStatus = BestVal;           // send the largest value out
    MoveZ(BestPos-1-FocOff);      // move to the best position
    return BestPos-1-FocOff;      // return the best position
10 }

```

## OTHER FUNCTIONS

```

unsigned long AFPMTVals(void)    // tvt
// Function to step through entire fast z stage range at 16 counts per step
// storing PMT values at each step in network variable structures (16 structure
15 // of 16 element arrays). Stage is returned to fast z position which gives the
// highest PMT value.
{
    unsigned short NextVal, BestVal;    // ADC values
    unsigned long BestPos;              // position of peak signal
20    unsigned long Move_Delta;          // steps per moove
    unsigned long step;                 // z stage step loop counter
    unsigned m;                         // array member
    static int FocOff;                  // Focus offset
    // short Delta;
25    // unsigned Stat;

    SetCtrl (Z_Control & (~SYNC_MODE)); // async mode
    // move to the extreme lower limit
    MoveZ (0);
30    delay(1977);                       // wait 50 ms

    // step through Z range checking for peak PMT and storing values

```

```

NextVal=ReadPort(ADC_ADDR);
BestVal=NextVal;
BestPos=0;
*****
5   for (step=0; step <= 4096; step += 16)
    {
        if(step == 0) step = 1;
        MoveZ(step-1);
        watchdog_update();
10   NextVal = ReadPort(ADC_ADDR);

        if (NextVal > BestVal)
        {
            BestPos=step;
15   BestVal=NextVal;
        }

        /* Write PMT Values to nv's */

20   if (step < 256) // use structure 1
        {
            m = step/16;
            pmtvals1.pmt[m] = NextVal;
        }

25   else if ((step >= 256) && (step < 512)) // use structure 2
        {
            m = (step-256)/16;
            pmtvals2.pmt[m]=NextVal;
30   }

        else if ((step >= 512) && (step < 768)) // use structure 3
        {
            m= (step-512)/16;

```

```
        pmtvals3.pmt[m]=NextVal;
    }

    else if ((step >= 768) && (step < 1024)) // use structure 4
5    {
        m = (step - 768) / 16;
        pmtvals4.pmt[m]=NextVal;
    }

10    else if ((step >= 1024) && (step < 1280)) // use structure 5
    {
        m = (step - 1024) / 16;
        pmtvals5.pmt[m]=NextVal;
    }

15    else if ((step >= 1280) && (step < 1536)) // use structure 6
    {
        m = (step - 1280) / 16;
        pmtvals6.pmt[m]=NextVal;
20    }

    else if ((step >= 1536) && (step < 1792)) // use structure 7
    {
        m = (step - 1536) / 16;
25    pmtvals7.pmt[m]=NextVal;
    }

    else if ((step >= 1793) && (step < 2048)) // use structure 8
    {
30    m = (step - 1793) / 16;
        pmtvals8.pmt[m]=NextVal;
    }

    else if ((step >= 2048) && (step < 2304)) // use structure 9
```

```

    {
        m = (step - 2048) / 16;
        pmtvals9.pmt[m]=NextVal;
    }
5
    else if ((step >= 2304) && (step < 2560)) // use structure 10
    {
        m = (step - 2304) / 16;
        pmtvals10.pmt[m]=NextVal;
10    }

    else if ((step >= 2560) && (step < 2816)) // use structure 11
    {
        m = (step - 2560) / 16;
15    pmtvals11.pmt[m]=NextVal;
    }

    else if ((step >= 2816) && (step < 3072)) // use structure 12
    {
20    m = (step - 2816) / 16;
        pmtvals12.pmt[m]=NextVal;
    }

    else if ((step >= 3072) && (step < 3328)) // use structure 13
25    {
        m = (step - 3072) / 16;
        pmtvals13.pmt[m]=NextVal;
    }

    else if ((step >= 3328) && (step < 3584)) // use structure 14
30    {
        m = (step - 3328) / 16;
        pmtvals14.pmt[m]=NextVal;
    }
}
```

```

        else if ((step >= 3584) && (step < 3840)) // use structure 15
        {
            m = (step - 3584) / 16;
            pmtvals15.pmt[m]=NextVal;
5          }

        else // if ((step >= 3840) && (step < 4096)) // use structure 16
        {
            m = (step - 3840) / 16;
10         pmtvals16.pmt[m]=NextVal;
        }

    }
15

    OpStatus = BestVal;           // send the largest value out
    MoveZ(BestPos-1-FocOff);      // move to the best position
    return BestPos-1-FocOff;      // return the best position
    }
20

unsigned CheckParms(void)
{
    signed long FinalPos;
25    unsigned short RetVal;

    FinalPos = ZPosCmd - (NumFrames-1)*ZStep;
    if (FinalPos < 0)
        RetVal = UNDERFLOW;
30    else
        RetVal = RDY_FOR_VOL;
    return (RetVal);
}

```

```

unsigned SelfTest (void)
{
    unsigned AllErrors, OldCtrl;

5    AllErrors = 0;
    OldCtrl = Z_Control;
    // turn off two test bits
    SetCtrl ( Z_Control & (~(SPARE_TEST | TEST_LEDn)));
    if ((ReadPort (STS_REG) & (TEST_IN | LED_TST_IN)) != 0)
10    AllErrors = AllErrors | 0x1;
    SetCtrl ( Z_Control | SPARE_TEST);
    if ((ReadPort (STS_REG) & TEST_IN) != TEST_IN)
        AllErrors = AllErrors | 0x2;
    SetCtrl ( (Z_Control | TEST_LEDn) & (~SPARE_TEST));
15    if ((ReadPort (STS_REG) & (TEST_IN | LED_TST_IN)) != LED_TST_IN)
        AllErrors = AllErrors | 0x4;

    SetCtrl (OldCtrl);
    return (AllErrors);
20 }

//*****
25 //***** WHEN STATEMENTS *****
//*****

when (reset)
{
    unsigned long i;

30

    ResetSts = 0;                // indicate reset to host
    Z_Control = ZCTRL_INIT;
    WritePort (CTRL_REG, Z_Control);
    AFocusThresh = INIT_THRESH;

```

```

MoveZ (0);
MuxSel = SEL_VIDEO;

if (((Z_Status = ReadPort (STS_REG)) & TEST_JMPRn) == 0) // if test mode
5   {
    // exercise everything
    SetCtrl (Z_Control);
    while (TRUE)
    // until next reset
    {
    for (i=0; i<4096; i+=512)
10   {
        MoveZ (i);
        ReadPort (STS_REG);
        delay(787);
        watchdog_update();
15   }

    for (i=4095; i>=512; i-=512)
    {
        MoveZ (i);
20   ReadPort (STS_REG);
        delay(787);
        watchdog_update();
    }
25 }
else
{
    SetCtrl (TEST_LEDn | ASYNC_TRIG);
    // set control bits: async mode,
    // reset all FF's, reset integ,
30   // turn off test LED

    WritePort (THRESH_DAC, AFocusThresh);
    SetCtrl (TEST_LEDn | RST_RETRACTn | RST_FOCUSn | RST_INTEGn
|
    RST_ZSYNCn | ASYNC_TRIG); // release resets except integ

```

```

        for (i=ZPosCmd; i > 1; i--)
            MoveZ (ZPosCmd);
    }
    /*
5      unsigned long i;
      ResetSts = 0;           ?? indicate reset to host
      Z_Control = ZCTRL_INIT;
      WritePort (CTRL_REG, Z_Control);
      AFocusThresh = INIT_THRESH;
10     MoveZ (0);
      MuxSel = SEL_VIDEO;

      if (((Z_Status = ReadPort (STS_REG)) & TEST_JMPRn) == 0) ?? if test mode
      {
          ?? exercise everything
15     SetCtrl (0);
      while (TRUE)           ?? until next reset
      {
          watchdog_update();
          SetCtrl (~Z_Control);    ?? toggle ctrl bits
20     if (Z_Control == 0)
      {
          SetCtrl (ASYNC_TRIG);    ?? write to DAC
          SetCtrl (0);
      }
25     WritePort (THRESH_DAC, AFocusThresh);    ?? set threshold
      AFocusThresh++;
      ZPosCmd++;
      WritePort (POS_CMD_LO, (unsigned int)((ZPosCmd % 4) & 0xFF));
      WritePort (POS_CMD_HI, (unsigned int)((ZPosCmd % 4) >> 8) & 0xFF);
30     ReadPort (STS_REG);
      Z_Status = ReadPort (ADC_ADDR);
      post_events();
    }
  }

```



```

else
{
    SetCtrl (TEST_LEDn | ASYNC_TRIG);    ?? set control bits: async mode,
                                         ?? reset all FF's, reset integ,
5                                         ?? turn off test LED
    WritePort (THRESH_DAC, AFocusThresh);
    SetCtrl (TEST_LEDn | RST_RETRACTn | RST_FOCUSn | RST_INTEGn
|
    RST_ZSYNCn | ASYNC_TRIG);    ?? release resets except integ
10 }
    OpStatus = 5;

    */

15 }

when (nv_update_occurs (AFocusThresh))
// Update the auto focus threshold DAC
{
20 WritePort (THRESH_DAC, AFocusThresh);
}

when (nv_update_occurs (ZPosCmd))
25 // update the position command DAC
{
    if (ZPosCmd > MAX_Z)
        ZPosCmd = MAX_Z;
    MoveZ (ZPosCmd);
30 }

when (nv_update_occurs (Z_Control))
// update the control word
{

```

```

WritePort (CTRL_REG, Z_Control);
}

// *****
5 // ***** remote procedure calls *****
// *****
// when (nv_update_occurs) is not used here so that OpMode can
// be monitored within the when clauses, specifically to break
// out of the MONITOR_FOCUS mode. OpMode is changed to NULL_MODE
10 // at the end of each clause to ensure that they are only invoked
// once.

when (OpMode == SELF_TEST)
{
15   OpStatus = SelfTest();
   OpMode = NULL_MODE;
}

when (OpMode == AUTO_FOCUS)
20   // Perform autofocus function)
   {
   OpStatus = AFocusServo (OpParArray[0]);
   PosStatus = ZPosCmd;
   OpMode = NULL_MODE;
25   }

when (OpMode == AF_FULLRANGE)
{
   // Perform full range autofocus function
30   MuxSel = SEL_VIDEO;
   PosStatus = AFFullRange();
   OpMode = NULL_MODE;
}

```

```
when (OpMode == AF_FAST)
{
    // Perform fast autofocus function
    MuxSel = SEL_VIDEO;
5    PosStatus = AFFast();
    OpMode = NULL_MODE;
}

10  when (OpMode == AF_FAST_T)
    {
        // Perform fast autofocus function
        MuxSel = SEL_VIDEO;
        PosStatus = AFFastTest();
15    OpMode = NULL_MODE;
    }

when (OpMode == AF_PMTVALS)
{
20    // Perform autofocus function which writes PMT values
    MuxSel = SEL_VIDEO;
    PosStatus = AFPMTVals();
    OpMode = NULL_MODE;
}

25  when (OpMode == READ_VIDEO)
    // Read ADC on stated mux channel)
    {
        MuxSel = SEL_VIDEO;
30    OpStatus = ReadPort (ADC_ADDR);
        OpMode = NULL_MODE;
    }

when (OpMode == READ_HV)
```

```

    // Read ADC on stated mux channel)
    {
        MuxSel = SEL_HV;
        OpStatus = ReadPort (ADC_ADDR);
5      OpMode = NULL_MODE;
    }

    when (OpMode == READ_FDBK)
        // Read ADC on stated mux channel)
10     {
        MuxSel = SEL_FDBK;
        OpStatus = ReadPort (ADC_ADDR);
        OpMode = NULL_MODE;
    }

15     when (OpMode == SET_FOR_VOL)
        // setup for volume acquisition)
        {
            unsigned FrameCtr;

20         if ((OpStatus = CheckParms()) == RDY_FOR_VOL)
            {
                // clear ZSYNC FF, wait for syncs and step when they come
                // assumes have moved to start position (closest to objective)
                // assumes ZSYNCs have been stopped, so clear ZSYNC FF and wait
25         SetCtrl (Z_Control & (~RST_ZSYNCn));           // clear ZSYNC
                SetCtrl (Z_Control | SYNC_MODE | RST_ZSYNCn); // go to sync
            mode
                // load DAC register with current position so it stays put
                // on first ZSYNC
30         WritePort (POS_CMD_LO, (unsigned int)(ZPosCmd & 0xFF));
                WritePort (POS_CMD_HI, (unsigned int)(ZPosCmd >> 8) & 0xFF);
                // alert host that ready for volume
                OpStatus = RDY_FOR_VOL;

```

```

    post_events();
    // take the volume)
    for (FrameCtr = 0; FrameCtr < NumFrames; FrameCtr++)
    {
        5      // wait for ZSYNC
        while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
            ;
        // reset the ZSYNC flip flop
        SetCtrl (Z_Control & (~RST_ZSYNCDn));
        10      SetCtrl (Z_Control | RST_ZSYNCDn);
        // step down
        ZPosCmd = ZPosCmd - ZStep;
        WritePort (POS_CMD_LO, (unsigned int)(ZPosCmd & 0xFF));
        WritePort (POS_CMD_HI, (unsigned int)(ZPosCmd >> 8) & 0xFF);
        15      watchdog_update();          // prevent reset
    }
    SetCtrl (Z_Control & (~SYNC_MODE)); // async mode
    OpStatus = VOLUME_COMPL;
    }
    20      PosStatus = ZPosCmd;
    OpMode = NULL_MODE;
    }

    when (OpMode == MONITOR_FOCUS)
    25      // Clear focus flipflop and wait for trip
    {
        // we assume the autofocus flip flop threshold has been set
        // note that this loop should look at RETRACTDn, if it were used
        // clear the focus FF
        30      SetCtrl (Z_Control & (~RST_FOCUSn));
        SetCtrl (Z_Control | RST_FOCUSn);
        if ((ReadPort (STS_REG) & FOCUSDn) == 0) // if trip immediately
            OpStatus = THRESH_TOO_LOW;
        else

```

```
OpStatus = COMPARATOR_READY;
post_events();
    // now hang out and wait for tripped
while (((ReadPort (STS_REG) & FOCUSDn) != 0) &&
5      (OpMode == MONITOR_FOCUS))
    {
        post_events();
    }
if (OpMode == MONITOR_FOCUS)
10   OpStatus = COMP_TRIPPED;
else
    OpStatus = ABORTED;
    post_events();
    OpMode = NULL_MODE;
15 }

when (OpMode == RD_STATUS)
    // read status reg
    {
20   OpStatus = ReadPort (STS_REG);
    OpMode = NULL_MODE;
    }

when (OpMode == ABORT_CMD)
25   // this shouldn't happen unless we are in another op mode
    {
        OpStatus = ABORTED;
        OpMode = NULL_MODE;
    }
```

```

unsigned ReadVidLoop(void)
{
    unsigned  l,a;          // loop and average counter
    unsigned  MaxVid;       // Maximum Video Value
5    unsigned  CurVid;      // Current Video Value
    unsigned long AvgVid;   // Average Video
    int       MaxStatus;    // Satus of MaxVal
    unsigned  Avg;          // Averaging value
    unsigned  Count;        // loop count
10
    // InMaximum to starting value

    OpStatus 0;             // Clear OpStatus
    OpStsArray[0] 0;        // Clear OpParArray[0]
15    OpStsArray[1] 0;      // Clear OpParArray[1]
    Count OpParArray[0];    // Set loop count
    Avg OpParArray[1];      // Set averaging

    if ((Avg < 1) || (Avg > 255))
20    Avg 1;                 // Make sure averaging is in range

    if ((Count < 64) || (Count > 255))
        Count 64;          // Make sure count is in range

25

    // PMT Read Loop
    for (l 0;l < Count;l++)
        // for 64 counts
30    {
        AvgVid 0;

        for (a 0; a < Avg; a++)
            AvgVid  AvgVid + ReadPort(A C A R); // Get Current Video Level

```

```

CurVid  AvgVid/a;          // averaged over a

if (CurVid > MaxVid)      // If Current level is greater
{
5   MaxVid CurVid;          // than max, reset max.
   MaxStatus 1;            // status of maximum.
}

// Store Video Level

10  if (l < ARRAYSIZE)
    {
        pmtvals1.pmt[l%ARRAYSIZE] CurVid;    // Store Current Video
        if (MaxStatus 1)                      // if this is max
15     {
            OpStsArray[0] 0;                  // Set Par0 to array
            OpStsArray[1] l%ARRAYSIZE;        // Set Par1 to element
        }
    }
20  else if ((l > ARRAYSIZE) && (l < ARRAYSIZE*2))
    {
        pmtvals2.pmt[l%ARRAYSIZE] CurVid;    // Store Current Video
        if (MaxStatus 1)                      // if this is max
        {
25     OpStsArray[0] 1;                  // Set Par0 to array
            OpStsArray[1] l%ARRAYSIZE;        // Set Par1 to element
        }
    }

30  else if ((l > ARRAYSIZE*2) && (l < ARRAYSIZE*3))
    {
        pmtvals3.pmt[l%ARRAYSIZE] CurVid;    // Store Current Video
        if (MaxStatus 1)                      // if this is max
        {

```



```

        OpStsArray[0] 2;           // Set Par0 to array
        OpStsArray[1] 1%ARRAYSIZE; // Set Par1 to element
    }
}
5   else if ((l > ARRAYSIZE*3) && (l < ARRAYSIZE*4))
    {
        pmtvals4.pmt[1%ARRAYSIZE] CurVid; // Store Current Video
        if (MaxStatus 1) // if this is max
        {
10      OpStsArray[0] 3;           // Set Par0 to array
        OpStsArray[1] 1%ARRAYSIZE; // Set Par1 to element
        }
    }
    else if ((l > ARRAYSIZE*4) && (l < ARRAYSIZE*5))
15  {
        pmtvals5.pmt[1%ARRAYSIZE] CurVid; // Store Current Video
        if (MaxStatus 1) // if this is max
        {
            OpStsArray[0] 4;           // Set Par0 to array
20      OpStsArray[1] 1%ARRAYSIZE; // Set Par1 to element
        }
    }
    else if ((l > ARRAYSIZE*5) && (l < ARRAYSIZE*6))
    {
25      pmtvals6.pmt[1%ARRAYSIZE] CurVid; // Store Current Video
        if (MaxStatus 1) // if this is max
        {
            OpStsArray[0] 5;           // Set Par0 to array
            OpStsArray[1] 1%ARRAYSIZE; // Set Par1 to element
30      }
        }
    }
    else if ((l > ARRAYSIZE*6) && (l < ARRAYSIZE*7))
    {
        pmtvals7.pmt[1%ARRAYSIZE] CurVid; // Store Current Video

```

```

        if (MaxStatus 1)                // if this is max
        {
            OpStsArray[0] 6;              // Set Par0 to array
            OpStsArray[1] 1%ARRAYSIZE;    // Set Par1 to element
5      }
    }
    else if ((l > ARRAYSIZE*7) && (l < ARRAYSIZE*8))
    {
        pmtvals8.pmt[1%ARRAYSIZE] CurVid; // Store Current Video
10    if (MaxStatus 1)                    // if this is max
        {
            OpStsArray[0] 7;              // Set Par0 to array
            OpStsArray[1] 1%ARRAYSIZE;    // Set Par1 to element
15    }
    }
    else if ((l > ARRAYSIZE*8) && (l < ARRAYSIZE*9))
    {
        pmtvals9.pmt[1%ARRAYSIZE] CurVid; // Store Current Video
        if (MaxStatus 1)                  // if this is max
20    {
            OpStsArray[0] 8;              // Set Par0 to array
            OpStsArray[1] 1%ARRAYSIZE;    // Set Par1 to element
        }
    }
25    else if ((l > ARRAYSIZE*9) && (l < ARRAYSIZE*10))
    {
        pmtvals10.pmt[1%ARRAYSIZE] CurVid; // Store Current Video
        if (MaxStatus 1)                  // if this is max
        {
30    OpStsArray[0] 9;                    // Set Par0 to array
            OpStsArray[1] 1%ARRAYSIZE;    // Set Par1 to element
        }
    }
    }
    else if ((l > ARRAYSIZE*10) && (l < ARRAYSIZE*11))

```

```

{
    pmtvals11.pmt[1%ARRAYSIZE] CurVid;    // Store Current Video
    if (MaxStatus 1)                        // if this is max
    {
5        OpStsArray[0] 10;                // Set Par0 to array
        OpStsArray[1] 1%ARRAYSIZE;        // Set Par1 to element
    }
}
else if ((l > ARRAYSIZE*11) && (l < ARRAYSIZE*12))
10 {
    pmtvals12.pmt[1%ARRAYSIZE] CurVid;    // Store Current Video
    if (MaxStatus 1)                        // if this is max
    {
        OpStsArray[0] 11;                // Set Par0 to array
15        OpStsArray[1] 1%ARRAYSIZE;        // Set Par1 to element
    }
}
else if ((l > ARRAYSIZE*12) && (l < ARRAYSIZE*13))
{
20    pmtvals13.pmt[1%ARRAYSIZE] CurVid;    // Store Current Video
    if (MaxStatus 1)                        // if this is max
    {
        OpStsArray[0] 12;                // Set Par0 to array
        OpStsArray[1] 1%ARRAYSIZE;        // Set Par1 to element
25    }
}
else if ((l > ARRAYSIZE*13) && (l < ARRAYSIZE*14))
{
    pmtvals14.pmt[1%ARRAYSIZE] CurVid;    // Store Current Video
30    if (MaxStatus 1)                        // if this is max
    {
        OpStsArray[0] 13;                // Set Par0 to array
        OpStsArray[1] 1%ARRAYSIZE;        // Set Par1 to element
    }
}

```

```
    }
    else if ((l > ARRAYSIZE*14) && (l < ARRAYSIZE*15))
    {
        pmtvals15.pmt[l%ARRAYSIZE] CurVid;    // Store Current Video
5      if (MaxStatus 1)                        // if this is max
        {
            OpStsArray[0] 14;                // Set Par0 to array
            OpStsArray[1] l%ARRAYSIZE;        // Set Par1 to element
        }
10     }
    else // if ((l > ARRAYSIZE*15) && (l < ARRAYSIZE*16))
    {
        pmtvals16.pmt[l%ARRAYSIZE] CurVid;    // Store Current Video
        if (MaxStatus 1)                      // if this is max
15     {
            OpStsArray[0] 15;                // Set Par0 to array
            OpStsArray[1] l%ARRAYSIZE;        // Set Par1 to element
        }
    }
20
    MaxStatus 0;
    watchdog update();
}
return MaxVid;
25 }
□
```

5

**A METHOD AND APPARATUS FOR PERFORMING  
AN AUTOMATIC FOCUS OPERATION**

10

**Christopher R. Fairley  
Timothy V. Thompson  
Ken K. Lee**

15

20

**APPENDIX B**

25

**M-2464-1P US**

```

// Fast Z-Axis Controller

// Source File:   Fastz.nc
// Author:       Chris F.
5 // Target Board: 0001003, Fast Z-Axis Controller
// Target Socket: M1-U2
// Ultrapointe Corp.
// PROPRIETARY! ALL RIGHTS RESERVED
// **** FOR INTERNAL USE ONLY ****
10 //
// Revision History:
//
//

15 // !!!!***** REVISION NUMBER *****!!!!
#define SWR 0x21 // software revision 2.1
//
//
//

20 // The FastZ PCB controls a piezo-driven stage with analog position
// feedback. The position command is a digital 14-bit word
// (12 bits used) written by the Neuron, and is double-
// buffered in front of the position command DAC. The position loop
25 // is closed in analog hardware. The position command can be updated
// in one of two modes: in sync mode, the Neuron writes to the DAC
// buffer, and the DAC receives the updated command when a ZSYNC
// pulse is recieved from the page scanner controller (P/N 000134);
// in async mode, the Neuron writes to the DAC buffer and the DAC is
30 // updated when hi byte is written.
// For autofocus purposes, there is a latching video comparator on
// the board, with a programmable threshold (similar to that on the
// PMT Amplifier (P/N 000276). There is also a fast 8-bit ADC for use
// by autofocus routines involving the stepper motor Z-axis and the

```

```
// fast (piezo) Z-axis. This same ADC can be used to read the analog
// position feedback or a peak detector gated by ZSYNC.
```

```
//
```

```
// There is a custom bus structure on the FastZ PCB that allows
```

```
5 // Neuron to address many devices with an 8-bit bi-directional bus,
// yet uses only the general purpose I/O lines of the device. This
// allows a general purpose Neuron module to be used because the main
// data/address buses of the Neuron are not required. In this scheme,
// I/O bits 0-7 are used as the data bus, bit 8 is a read/not-write
10 // line, bit 9 is an address/not-data line, and bit 10 is a low-true
// strobe line. In general, the Neuron will write an address to the
// address register, then read or write the intended data value; the
// strobe line executes the intended operation after address and data
// have been appropriately set up.
```

```
15 // The following rules are adhered to by the hardware and should
// be respected by software:
```

```
// 1. No one drives the data bus unless STB* (low-true STroBe) is
// asserted. The only exception is that the Neuron drives the
// bus with an address value without asserting STB* when
```

```
20 // ADDR/DATA* is set high.
```

```
// 2. Setting the ADDR/DATA* (ADDRess/not-DATA) line forces a
// write operation (the RD/WR* line is don't care when
// ADDR/DATA* = 1).
```

```
// 3. The address latch is transparent, and controlled by ADDR/DATA*,
25 // so no strobe is required to write the address.
```

```
// 4. Because of rule 1, no one should read the bus until STB* is
// low.
```

```
// This bus arrangement means a lot of bit twiddling to do a read or
// write, but if the address does not change between operations, it
```

```
: 30 // can be speeded up considerably. FastRead and FastWrite routines are
// provided that do not set up addresses are provided for this purpose.
```

```
#include <control.h>
```

```

//*****
//***** DEFINITIONS *****
//*****

// Constants

5
// Mode Commands:
#define SELF_TEST    0x00    // commands from host
#define READ_VIDEO   0x11
#define READ_FDBK    0x12
10 #define READ_HV     0x13
#define SET_FOR_VOL  0x14
#define MONITOR_FOCUS 0x15
#define ABORT_CMD    0x16
#define RD_STATUS    0x17
15 #define AF_PMTVALS  0x19    // new af with pmt reading storage tvt
#define AF_PEAK_SYNC  0x1C    // patch scan autofocus w/ peak detect & sync
#define AF_FAST       0x1D    // fast af algorithm tvt
#define RD_VID_LP     0x1E    // funcion to read PMT 64 times in 1 sec.
#define NULL_MODE     0xFF    // used to prevent repeated mode cmds

20
// status codes
#define SUCCESS       0

// volume mode responses:
25 #define RDY_FOR_VOL  5
#define VOLUME_COMPL  6
#define UNDERFLOW     7

// coarse autofocus responses:
30 #define COMPARATOR_READY 8
#define THRESH_TOO_LOW  9
#define COMP_TRIPPED    0x0A
#define ABORTED         0x0B

```



```

// misc
#define CMD_UNKNOWN 0xFF
#define DAC_SHIFT 2 // shift value for Z Command DAC

5 // Register addresses: note that two MSBs of address are mux
  // select bits.
  // Write-Only addresses:
  #define CTRL_REG 0 // Control register address
  #define THRESH_DAC 1 // Autofocus threshold DAC
10 #define POS_CMD_LO 2 // Position command lo byte
    #define POS_CMD_HI 3 // Position command hi byte

  // Read-Only addresses:
  #define STS_REG 0 // Status Register
15 #define ADC_ADDR 1 // ADC
    #define RST_ZSYNCn 0x02 // pulse low to reset ZSYNC flip flop

  // mux addresses: logical OR with register addresses
  #define SEL_VIDEO 0 // video for autofocus
20 #define SEL_FDBK 0x40 // analog position feedback
    #define SEL_HV 0x80 // high voltage output
    #define SEL_PKDET 0x80 // peak detector replaces HV input

  // An 'n' suffix on a signal name indicates low true
25 // control register bits:
    #define SYNC_MODE 0x01 // set to one for sync DAC update
    #define SPARE_TEST1 0x04 // control reg test line 1
    #define RST_FOCUSn 0x08 // pulse low to reset FOCUS flip flop
    #define SPARE_TEST2 0x10 // control reg test line 2
30 #define RST_INTEGn 0x20 // hold low to turn off position loop
    // integrator

    #define PKDET_MODE 0x40 // spare now allocated to peak detector AF
    #define TEST_LEDn 0x80 // set low to illuminate TEST LED

```

```

// status register bits:
#define ZSYNCDn    0x01    // low means ZSYNC has occurred
#define FOCUSDn    0x02    // low means FOCUS has occurred
#define TEST_IN1    0x04    // first test line
5  #define ADC_BUSYn  0x08    // low means conversion in process
#define TEST_IN2    0x10    // second test line
#define PKDET_IN    0x20    // test peak det line
#define LED_TST_IN  0x40    // test led line
#define TEST_JMPRn  0x80    // low means test jumper is in place
10

// control line definitions:
#define READ        1        // for RD_WRn line
#define WRITE        0
#define ADDR        1        // for ADDR_DATAn line
15 #define DATA        0

// autofocus constants:
/* #define MAX_Z      4095    // max position DAC range
#define AF_FAST_DELTA  10    // step size for fast autofocus sweep
20 #define AF_SLOW_DELTA  1    // step size for final autofocus
#define ARRAYSIZE      16    // size of arrays for PMT values
#define COARSE_STEP     64    // step sizes for coarse and fine portions
#define FINE_STEP        8    // of autofocus
#define FINE_RANGE      128    // 1/2 of fine range
25 #define DAC_RANGE      4095 // range of position cmd DAC shifted 2 bits

// initial values
#define ZCTRL_INIT  0x20;
        //!SYNC_MODE,!TESTn,!PKDET_MODE,!RST_FOCUSn,RST_INTEGN
30 #define INIT_THRESH  0x16;

```

```

//*****
//***** PRAGMAS *****
//*****

#pragma net_buf_in_count 7
5  #pragma net_buf_out_count 7

//*****
//***** I/O DECLARATIONS *****
10 //*****

IO_0 output byte WRITE_DATA; // parallel port
IO_0 input byte READ_DATA;

IO_8 output bit RD_WRn; // ReaD/not-WRite
15 IO_9 output bit ADDR_DATAn; // ADDRess/not-DATA
IO_10 output bit STBn=1; // low-true STroBe

//*****
20 //***** NETWORK VARIABLES *****
//*****

network output polled unsigned SW_Rev = SWR; // software version
network output unsigned ResetSts; // for indicating reset
network input unsigned OpMode=NULL_MODE; // for remote proc. calls
25 network output unsigned OpStatus; // response to calls
network input unsigned OpParArray[2]; // Operation parameters
network output unsigned OpStsArray[2]; // Extra Status info

network input unsigned AFocusThresh; // 0-255 for Autofocus threshold
: 30 network input unsigned long ZPosCmd; // Position Command 0-4095
network input unsigned Z_Control; // Control byte
network input unsigned MuxSel; // ADC mux selector
network input unsigned ZStep; // step size for vol acquire
network input unsigned NumFrames; // # frames for volume

```

```
network output polled unsigned Z_Status;      // Status byte
network output polled unsigned long PosStatus; // more status

struct PMTARRAY          // define structure to hold pmt values   tvl
5  {
    unsigned   pmt[16];
};

// declare 16 network variables as structure of PMTARRAY type
10 // to hold 256 pmt readings (one every 16 steps of the fast z stage tvl

network output polled struct PMTARRAY pmtvals1;
network output polled struct PMTARRAY pmtvals2;
network output polled struct PMTARRAY pmtvals3;
15 network output polled struct PMTARRAY pmtvals4;
network output polled struct PMTARRAY pmtvals5;
network output polled struct PMTARRAY pmtvals6;
network output polled struct PMTARRAY pmtvals7;
network output polled struct PMTARRAY pmtvals8;
20 network output polled struct PMTARRAY pmtvals9;
network output polled struct PMTARRAY pmtvals10;
network output polled struct PMTARRAY pmtvals11;
network output polled struct PMTARRAY pmtvals12;
far network output polled struct PMTARRAY pmtvals13;
25 far network output polled struct PMTARRAY pmtvals14;
far network output polled struct PMTARRAY pmtvals15;
far network output polled struct PMTARRAY pmtvals16;

//*****
30 //***** GLOBAL VARIABLES *****
//*****
```

```

/*****
/***** FUNCTION DEFINITIONS *****/
/*****

```

```

unsigned short ReadPort (unsigned short TargetAddr)

```

```

5 // read a port, including set up address

```

```

{

```

```

    unsigned short result;

```

```

    io_set_direction (WRITE_DATA, IO_DIR_OUT);

```

```

10 io_out (ADDR_DATAn, ADDR);          // set ADDR_DATAn to ADDR

```

```

    io_out (WRITE_DATA, TargetAddr | MuxSel);    // write address

```

```

    io_out (ADDR_DATAn, DATA);          // set ADDR_DATAn to DATA

```

```

    io_set_direction (READ_DATA, IO_DIR_IN);

```

```

    io_out (RD_WRn, READ);              // set RD_WRn for READ

```

```

15 io_out (STBn, 0);                    // assert STBn

```

```

    result = io_in (READ_DATA);          // read data

```

```

    io_out (STBn, 1);                    // de-assert STBn

```

```

    return (result);                     // return result

```

```

}

```

```

20

```

```

void WritePort (unsigned short TargetAddr, unsigned short WriteData)

```

```

// write to a port, including set up address

```

```

{

```

```

25

```

```

    io_set_direction (WRITE_DATA, IO_DIR_OUT);

```

```

    io_out (ADDR_DATAn, ADDR);          // set ADDR_DATAn to ADDR

```

```

    io_out (WRITE_DATA, TargetAddr | MuxSel);    // write address

```

```

    io_out (ADDR_DATAn, DATA);          // set ADDR_DATAn to DATA

```

```

: 30 // no need to change port direction for write

```

```

    io_out (RD_WRn, WRITE);              // set RD_WRn for WRITE

```

```

    io_out (WRITE_DATA, WriteData);        // set up data

```

```

    io_out (STBn, 0);                    // assert STBn

```

```

    io_out (STBn, 1);                    // de-assert STBn

```

```

    }

```

```

unsigned FastRead (void)

```

```

5  // read a port, without address set up
    {
        unsigned result;

        io_out (STBn, 0);           // assert STBn
10  result = io_in (READ_DATA);     // read data
        io_out (STBn, 1);           // de-assert STBn
        return (result);            // return result
    }

```

```

15
void FastWrite (unsigned short WriteData)

```

```

// write to a port, without address set up
    {

20  io_out (WRITE_DATA, WriteData); // set up data
        io_out (STBn, 0);           // assert STBn
        io_out (STBn, 1);           // de-assert STBn
    }

```

```

25
void SetAddressRead (unsigned short TargetAddr)

```

```

// Setup address for fast read
    {

30  io_set_direction (WRITE_DATA, IO_DIR_OUT);
        io_out (ADDR_DATAn, ADDR); // set ADDR_DATAn to ADDR
        io_out (WRITE_DATA, TargetAddr | MuxSel); // write address
        io_out (ADDR_DATAn, DATA); // set ADDR_DATAn to DATA
        io_set_direction (READ_DATA, IO_DIR_IN);

```

```

io_out (RD_WRn, READ);          // set RD_WRn for READ
}

```

```

5 void SetAddressWrite (unsigned short TargetAddr)
  // Setup address for fast read
  {

io_set_direction (WRITE_DATA, IO_DIR_OUT);
10 io_out (ADDR_DATAn, ADDR);      // set ADDR_DATAn to ADDR
io_out (WRITE_DATA, TargetAddr | MuxSel); // write address
io_out (ADDR_DATAn, DATA);      // set ADDR_DATAn to DATA
// no need to change port direction for write
io_out (RD_WRn, WRITE);          // set RD_WRn for WRITE
15 }

```

```

void SetCtrl (unsigned NewCtrl)
// sets the control word
20 {
  Z_Control = NewCtrl;
  WritePort (CTRL_REG, NewCtrl);
}

```

```

25 void MoveZ (unsigned long target)
  // move the Z-axis to the current commanded position (ZPosCmd)
  {
    if (target > DAC_RANGE)
; 30 target = DAC_RANGE;
    ZPosCmd = (target << DAC_SHIFT);
    WritePort (POS_CMD_LO, (unsigned)(ZPosCmd & 0xFF)); // set position
    WritePort (POS_CMD_HI, (unsigned)((ZPosCmd >> 8) & 0xFF));
    ZPosCmd = target;

```

```

}
```

```

void ResetZSync(void)
{

```

```

5   io_set_direction (WRITE_DATA, IO_DIR_OUT);    // set port to out
    io_out (ADDR_DATAn, ADDR);                    // set ADDR_DATAn to ADDR
    io_out (WRITE_DATA, (~RST_ZSYNCn) | MuxSel);  // write address
    io_out(RD_WRn, READ);                          // set read line
    io_out (STBn, 0);                              // assert STBn
10  io_out (STBn, 1);                              // de-assert STBn
}

```

```

unsigned long AFFast(void) // tvf

```

```

// Fast Autofocus Function: This function does a coarse then a fine autofocus

```

```

15 {
    unsigned short NextVal, BestVal; // ADC values
    unsigned long BestPos;           // position of peak signal
    unsigned long AvgVal;            // average value storage
    unsigned long FineCenter;        // center of fine focus range
20  unsigned long Move_Delta;         // steps per move
    unsigned long step;              // z stage step loop counter
    unsigned m;                      // array member
    unsigned l;                      // loop counter
    static int FocOff;               // Focus offset
25  unsigned int StepSize;           // step size during coarse focus

```

```

// Set FocOff and StepSize

```

```

    FocOff= OpParArray[0];

```

```

    StepSize = 128;

```

```

30

```

```

    SetCtrl (Z_Control & (~SYNC_MODE)); // async mode

```

```

// move to the extreme lower limit

```

```

    MoveZ (0); // Start at the bottom

```

```

    delay(1977); // wait 50 ms

```



```
// ***** coarse autofocus *****
```

```
// step through entire Z range with large steps
```

```
// checking for peak PMT and storing values
```

```
AvgVal=0;
```

```
5 MuxSel = SEL_VIDEO;
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
10 NextVal=(short)(AvgVal/4);
```

```
BestVal=NextVal; // Clear Largest value
```

```
BestPos=0; // Clear Best Position
```

```
for (step=0; step <= 4096; step +=StepSize) // Entire range in steps of 128
```

```
15 {
```

```
if(step == 0) step = 1;
```

```
MoveZ(step-1); // Move 0 to 4095
```

```
watchdog_update();
```

```
AvgVal=0;
```

```
20 MuxSel = SEL_VIDEO;
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
```

```
25 NextVal=(short)(AvgVal/4);
```

```
if (NextVal > BestVal) // If Current value is larger
```

```
{ // than the best (largest) value
```

```
BestPos=step; // this is the best position
```

```
: 30 BestVal=NextVal; // Update best value
```

```
}
```

```
}
```

```
// Move to start value for fine autofocus and wait
```

```

    if (BestPos > 129)
        MoveZ (BestPos-128-1); // Start at the bottom
    else
        MoveZ (1);
5    delay(395);                // wait 10 ms

// Set up comparison values
    AvgVal=0;
    MuxSel = SEL_VIDEO;
10    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
    NextVal=(short)(AvgVal/4);
15    BestVal=NextVal;          // Clear Largest value

// ***** fine autofocus *****
// using position from coarse focus
    FineCenter=BestPos; // Set center of focus to the best coarse focus pos.
20    if(FineCenter < 129) // Make sure it's valid (at least 129
        FineCenter=129;
    for (step = FineCenter-128; step <= FineCenter+128; step +=8)
    {
        MoveZ(step-1);          // Move through fine af range
25        watchdog_update();
        AvgVal=0;
        MuxSel = SEL_VIDEO;
        AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
        AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
30        AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
        AvgVal=AvgVal+ReadPort(ADC_ADDR); // Store current PMT value
        NextVal=(short)(AvgVal/4);

        if (NextVal > BestVal)    // If Current value is larger

```

```

        {
            // than the best (largest) value
            BestPos=step;           // this is the best position
            BestVal=NextVal;        // Update best value
        }
5      }

    OpStatus = BestVal;           // send the largest value out
    MoveZ(BestPos-1-FocOff);       // move to the best position
    return BestPos-1-FocOff;       // return the best position
10  }

unsigned long AFPeakSync(void)    // tvf
// Fast Autofocus Function with peak detection and synchronous operation
// Note that in synchronous mode, when a write to the position command
15 // registers is executed, the stage does not receive the command until the
// next ZSYNC, and the stage does not settle until a few msec beyond that.
// Do not use the 'MoveZ' routine while in synchronous mode because it
// reverts to asynchronous mode.
// The coarse pass of this Af routine operates in asynchronous mode,
20 // stepping as quickly as possible, and not paying attention to where
// the scanners are during sampling. In the fine pass, however, the page
// scanner is run at a high
// rate (> 100Hz) and allow one page scanner cycle for the FastZ to settle
// before doing data acquisition on the next. The routine therefore
25 // spends two page scanner cycles at each Z position. Furthermore, in
// peak detector mode, an ADC conversion is automatically started when
// ZSYNC occurs, and the peak detector is automatically reset. So when
// the routine waits for ZSYNC and then reads the ADC, the ADC reading
// applies to the peak detected signal from the previous page scan.
: 30 // Also note that resetting ZSYNC also resets the ADC, so the ADC should
// be read before resetting ZSYNC. In peak detector mode, the ADC is
// controlled ONLY by ZSYNC; there is no way to do an asynchronous
// conversion without resetting the peak detector mode bit. The peak
// detector is used to obtain a representation of the highest intensity

```

```

// that occurred in a video frame at a given Z elevation. The intention
// is to repeatedly focus on the layer with the highest reflectivity, and
// use programmable offsets (FocOff) so that the user can select any layer
// to focus on.
5 // One final note: this routine will time out if the page scanner
// rate is not fast enough, because 'watchdog_update' is used sparingly.
{
    unsigned short NextVal, BestVal; // ADC values
    unsigned long BestPos;           // position of peak signal
10 unsigned long step;               // z stage step loop counter
    unsigned long RunningStart, AcquirePos, ZPosShift;
    boolean MoveDown;
    unsigned long RUN_START;
    unsigned m;                      // array member
15 unsigned l;                      // loop counter
    static int FocOff;               // Focus offset

    FocOff = OpParArray[0];
    RUN_START = 0;
20 // move to the extreme lower limit
    MoveZ (0);                       // Start at the bottom
    delay (1977);                    // 50ms delay

    // ***** coarse autofocus *****
25 // step through entire Z range with large steps
    // checking for peak PMT and storing values
    MuxSel = SEL_VIDEO;
    NextVal = 0;
    BestVal=0;                       // Clear Largest value
30 BestPos=0;                       // Clear Best Position

    // assumes have moved to start position, use async mode for coarse pass
    SetCtrl (Z_Control & (~SYNC_MODE) & (~PKDET_MODE));

```

```

for (step = 0; step <= DAC_RANGE + 1; step += COARSE_STEP)
{
    if (step >= DAC_RANGE + 1) step = DAC_RANGE;
    MoveZ (step);
5    NextVal = ReadPort (ADC_ADDR);
    if (NextVal > BestVal)          // If Current value is larger
    {                               // than the best (largest) value
        BestPos=step;              // this is the best position
        BestVal=NextVal;           // Update best value
10    }
}

```

```

OpStsArray[0] = (short)(BestPos & 0xFF); // return value for diagnostics
OpStsArray[1] = (short)((BestPos & 0xFF00) >> 8);
15 // Move to start value for fine autofocus
    // top two possible positions from coarse pass are special
    // and so is bottom position
    if (BestPos >= DAC_RANGE + 1 - FINE_RANGE) { // top 2 positions
        RunningStart = DAC_RANGE + 1 - (2 * FINE_RANGE) - RUN_START;
20    AcquirePos = DAC_RANGE + 1 - (2 * FINE_RANGE);
        MoveDown = FALSE; }
    else if (BestPos < FINE_RANGE) {           // bottom position
        RunningStart = (2 * FINE_RANGE) + RUN_START;
        AcquirePos = (2 * FINE_RANGE);
25    MoveDown = TRUE; }                      // other
    else {
        RunningStart = BestPos + FINE_RANGE + RUN_START;
        AcquirePos = BestPos + FINE_RANGE;
        MoveDown = TRUE; }
30    MoveZ (RunningStart);

```

```

// fine pass in peak detect/sync mode

```

```

MuxSel = SEL_PKDET;          // move to sync, pkdet modes

```

```

SetCtrl (Z_Control | SYNC_MODE | PKDET_MODE );
ResetZSync();
while ((ReadPort(STS_REG) & ZSYNCDn) != 0) // wait for sync
;
5   ResetZSync();
    // Set up comparison values
    NextVal=0;
    BestVal=0; // clear Largest value

10  // ***** fine autofocus *****
    // using position from coarse focus
    if (MoveDown == TRUE) // downward pass
    {
        for (step = RunningStart; // running start
15         step > AcquirePos;
            step -= FINE_STEP)
        {
            // wait for ZSYNC (dummy sync while stage in transit)
            while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
20         ;
            // note that the following move happens on the next ZSYNC
            ZPosShift = step << DAC_SHIFT;
            WritePort (POS_CMD_LO, (unsigned int)(ZPosShift & 0xFF));
            WritePort (POS_CMD_HI, (unsigned int)(ZPosShift >> 8) & 0xFF);
25         // reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
            ResetZSync();
            // wait for sync (peak detector is acquiring)
            while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
            ;
30         // reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
            ResetZSync();
        }
        for (step = AcquirePos; // fine pass data acq
            (step >= AcquirePos - (2 * FINE_RANGE)) &&

```

```

        (step < DAC_RANGE);      // because unsigned
        step -= FINE_STEP)
    {
        // wait for ZSYNC (dummy sync while stage in transit)
5       while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
            ;

        // note that the following move happens on the next ZSYNC
        ZPosShift = step << DAC_SHIFT;
        WritePort (POS_CMD_LO, (unsigned int)(ZPosShift & 0xFF));
10      WritePort (POS_CMD_HI, (unsigned int)(ZPosShift >> 8) & 0xFF);
        // reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
        ResetZSync();

        // wait for sync (peak detector is acquiring)
        while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
15          ;

        NextVal = ReadPort(ADC_ADDR);
        if (NextVal > BestVal)      // If Current value is larger
        {                          // than the best (largest) value
            BestPos = step + FINE_STEP;  // this is the best position
20          BestVal = NextVal;        // Update best value
        }

        // reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
        ResetZSync();
    }

25  // now process the last step (because of pipeline)
    while ((ReadPort(STS_REG) & ZSYNCDn) != 0) // wait for last step
        ;                                     // to settle
    ResetZSync();
    while ((ReadPort(STS_REG) & ZSYNCDn) != 0) // wait for peak det
30      ;                                     // to load up
    NextVal = ReadPort(ADC_ADDR);
    if (NextVal > BestVal)      // If Current value is larger
    {                          // than the best (largest) value
        BestPos = step + FINE_STEP;  // this is the best position

```

```

        if (BestPos >= DAC_RANGE + 1) BestPos = 0;
        BestVal = NextVal;          // Update best value
    }
}

5   else                                // upward pass
    {
        for (step = RunningStart;      // running start
            step < AcquirePos;
            step += FINE_STEP)
10  {
        // wait for ZSYNC (dummy sync while stage in transit)
        while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
            ;

        // note that the following move happens on the next ZSYNC
15  ZPosShift = step << DAC_SHIFT;
        WritePort (POS_CMD_LO, (unsigned int)(ZPosShift & 0xFF));
        WritePort (POS_CMD_HI, (unsigned int)(ZPosShift >> 8) & 0xFF);
        // reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
        ResetZSync();
20  // wait for sync (peak detector is acquiring)
        while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
            ;

        // reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
        ResetZSync();
25  }

        for (step = AcquirePos;        // fine pass data acq
            step <= AcquirePos + (2 * FINE_RANGE);
            step += FINE_STEP)
        {
30  if (step >= (DAC_RANGE + 1)) step = DAC_RANGE;
        // wait for ZSYNC (dummy sync while stage in transit)
        while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
            ;

        // note that the following move happens on the next ZSYNC

```



```

ZPosShift = step < < DAC_SHIFT;
WritePort (POS_CMD_LO, (unsigned int)(ZPosShift & 0xFF));
WritePort (POS_CMD_HI, (unsigned int)(ZPosShift > > 8) & 0xFF);
// reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
5   ResetZSync();
    // wait for sync (peak detector is acquiring)
    while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
        ;
    NextVal = ReadPort(ADC_ADDR);
10   if (NextVal > BestVal)          // If Current value is larger
        {                          // than the best (largest) value
            BestPos = step - FINE_STEP;    // this is the best position
            BestVal = NextVal;             // Update best value
        }
15   // reset the ZSYNC flip flop for next time (also resets ADC CNVRT)
        ResetZSync();
    }
    // now process the last step (because of pipeline)
    while ((ReadPort(STS_REG) & ZSYNCDn) != 0) // wait for last step
20   ;                                  // to settle
        ResetZSync();
    while ((ReadPort(STS_REG) & ZSYNCDn) != 0) // wait for peak det
        ;                                  // to load up
    NextVal = ReadPort(ADC_ADDR);
25   if (NextVal > BestVal)          // If Current value is larger
        {                          // than the best (largest) value
            BestPos = step - FINE_STEP;    // this is the best position
            if (BestPos >= DAC_RANGE+1) BestPos = DAC_RANGE;
            BestVal = NextVal;             // Update best value
30   }
    }

ResetZSync();
SetCtrl (Z_Control & (~PKDET_MODE));

```

```

    OpStatus = BestVal;          // send the largest value out
    ZPosShift = (BestPos - FocOff) << DAC_SHIFT; // move to the best position
    WritePort (POS_CMD_LO, (unsigned int)(ZPosShift & 0xFF));
    WritePort (POS_CMD_HI, (unsigned int)(ZPosShift >> 8) & 0xFF);
5   return (BestPos - FocOff);    // return the best position
}

unsigned long AFPMTVals(void) // tvr
// Function to step through entire fast z stage range at 16 counts per step
10 // storing PMT values at each step in network variable structures (16 structures
// of 16 element arrays). Stage is returned to fast z position which gives the
// highest PMT value.
{
    unsigned short NextVal, BestVal; // ADC values
15   unsigned long BestPos;          // position of peak signal
    unsigned long Move_Delta;      // steps per moove
    unsigned long step;            // z stage step loop counter
    unsigned m;                   // array member
    static int FocOff;             // Focus offset
20
    SetCtrl (Z_Control & (~SYNC_MODE)); // async mode
    // move to the extreme lower limit
    MoveZ (0);
    delay(1977);                  // wait 50 ms
25
    // step through entire Z range checking for peak PMT and storing values

    NextVal=ReadPort(ADC_ADDR);
    BestVal=NextVal;
30   BestPos=0;

    for (step=0; step <= 4096; step += 16)
    {
        if(step == 0) step = 1;

```

```

        MoveZ(step-1);
        watchdog_update();
        NextVal = ReadPort(ADC_ADDR);

5      if (NextVal > BestVal)
        {
            BestPos=step;
            BestVal=NextVal;
        }

10     /* Write PMT Values to nv's */

        if (step < 256) // use structure 1
        {
15         m = (short)(step/16);
            pmtvals1.pmt[m] = NextVal;
        }

        else if ((step >= 256) && (step < 512)) // use structure 2
20     {
            m = (short)((step-256)/16);
            pmtvals2.pmt[m]=NextVal;
        }

        else if ((step >= 512) && (step < 768)) // use structure 3
25     {
            m = (short)((step-512)/16);
            pmtvals3.pmt[m]=NextVal;
        }

        else if ((step >= 768) && (step < 1024)) // use structure 4
30     {
            m = (short)((step - 768) / 16);
            pmtvals4.pmt[m]=NextVal;
        }

```

```
    }  
  
    else if ((step >= 1024) && (step < 1280)) // use structure 5  
    {  
5      m = (short)((step - 1024) / 16);  
      pmtvals5.pmt[m]=NextVal;  
    }  
  
    else if ((step >= 1280) && (step < 1536)) // use structure 6  
10   {  
      m = (short)((step - 1280) / 16);  
      pmtvals6.pmt[m]=NextVal;  
    }  
  
    else if ((step >= 1536) && (step < 1792)) // use structure 7  
15   {  
      m = (short)((step - 1536) / 16);  
      pmtvals7.pmt[m]=NextVal;  
    }  
  
    else if ((step >= 1793) && (step < 2048)) // use structure 8  
20   {  
      m = (short)((step - 1793) / 16);  
      pmtvals8.pmt[m]=NextVal;  
25   }  
  
    else if ((step >= 2048) && (step < 2304)) // use structure 9  
    {  
      m = (short)((step - 2048) / 16);  
30   pmtvals9.pmt[m]=NextVal;  
    }  
  
    else if ((step >= 2304) && (step < 2560)) // use structure 10  
    {
```

```

        m = (short)((step - 2304) / 16);
        pmtvals10.pmt[m]=NextVal;
    }

5      else if ((step >= 2560) && (step < 2816)) // use structure 11
    {
        m = (short)((step - 2560) / 16);
        pmtvals11.pmt[m]=NextVal;
    }

10     else if ((step >= 2816) && (step < 3072)) // use structure 12
    {
        m = (short)((step - 2816) / 16);
        pmtvals12.pmt[m]=NextVal;
15     }

    else if ((step >= 3072) && (step < 3328)) // use structure 13
    {
        m = (short)((step - 3072) / 16);
20     pmtvals13.pmt[m]=NextVal;
    }

    else if ((step >= 3328) && (step < 3584)) // use structure 14
    {
25     m = (short)((step - 3328) / 16);
        pmtvals14.pmt[m]=NextVal;
    }

    else if ((step >= 3584) && (step < 3840)) // use structure 15
30     {
        m = (short)((step - 3584) / 16);
        pmtvals15.pmt[m]=NextVal;
    }

```

```

        else // if ((step >= 3840) && (step < 4096)) // use structure 16
        {
            m = (short)((step - 3840) / 16);
            pmtvals16.pmt[m]=NextVal;
5          }

        }

10      OpStatus = BestVal;          // send the largest value out
      MoveZ(BestPos-1-FocOff);      // move to the best position
      return BestPos-1-FocOff;      // return the best position
    }

15  unsigned ReadVidLoop(void)
    {
        // This functions reads the pmt video value 'a' times to find an average at
        // points in the Fine Z travel with 'l' step size. a and l values are
20      // updated by network variables to give the desired time for the loop
        // which is about 1 second in the LIS 1000

        unsigned   l,a;          // loop and average counter
        unsigned   MaxVid;       // Maximum Video Value
25      unsigned   CurVid;       // Current Video Value
        unsigned long AvgVid;    // Average Video
        int        MaxStatus;    // Satus of MaxVal
        unsigned   Avg;          // Averaging value
        unsigned   Count;       // loop count
30

        // Initialize
        MuxSel=SEL_VIDEO;        // Select Video;
        SetCtrl(Z_Control & (~SYNC_MODE) & (~PKDET_MODE));
        CurVid=ReadPort(ADC_ADDR); // Get Current Video Value

```

```

    MaxVid=CurVid;          // Set Maximum to starting value
    OpStatus=0;              // Clear OpStatus
    OpStsArray[0]=0;         // Clear OpParArray[0]
    OpStsArray[1]=0;         // Clear OpParArray[1]
5    Count=OpParArray[0];     // Set loop count
    Avg=OpParArray[1];        // Set averaging

    if (Avg == 0)
        Avg=1;               // Make sure averaging is in range
10
    if (Count < 64)
        Count=64;            // Make sure count is in range

15
    // PMT Read Loop
    for (l=0;l <=Count;l++)
        // for 64 counts
    {
20        AvgVid=0;

        for (a=0; a < Avg; a++)
            AvgVid = AvgVid + ReadPort(ADC_ADDR); // Get Current Video Level
        CurVid = (short)(AvgVid/a);              // averaged over a
25
        if (CurVid > MaxVid)                    // If Current level is greater
        {
            MaxVid=CurVid;                      // than max, reset max.
            MaxStatus=1;                          // status of maximum.
30        }

        // Store Video Level

    if (l < ARRAYSIZE)

```

```

{
    pmtvals1.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
    if (MaxStatus==1)                      // if this is max
    {
5       OpStsArray[0]=0;                    // Set Par0 to array
        OpStsArray[1]=l%ARRAYSIZE;        // Set Par1 to element
    }
}
else if ((l >= ARRAYSIZE) && (l < ARRAYSIZE*2))
10 {
    pmtvals2.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
    if (MaxStatus==1)                      // if this is max
    {
        OpStsArray[0]=1;                  // Set Par0 to array
15     OpStsArray[1]=l%ARRAYSIZE;          // Set Par1 to element
    }
}

else if ((l >= ARRAYSIZE*2) && (l < ARRAYSIZE*3))
20 {
    pmtvals3.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
    if (MaxStatus==1)                      // if this is max
    {
        OpStsArray[0]=2;                  // Set Par0 to array
25     OpStsArray[1]=l%ARRAYSIZE;          // Set Par1 to element
    }
}

else if ((l >= ARRAYSIZE*3) && (l < ARRAYSIZE*4))
{
30     pmtvals4.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
    if (MaxStatus==1)                      // if this is max
    {
        OpStsArray[0]=3;                  // Set Par0 to array
        OpStsArray[1]=l%ARRAYSIZE;        // Set Par1 to element
    }
}

```



```

    }
  }
  else if ((l >= ARRAYSIZE*4) && (l < ARRAYSIZE*5))
  {
5      pmtvals5.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
      if (MaxStatus == 1)                      // if this is max
      {
          OpStsArray[0]=4;                    // Set Par0 to array
          OpStsArray[1]=1%ARRAYSIZE;          // Set Par1 to element
10      }
      }
  else if ((l >= ARRAYSIZE*5) && (l < ARRAYSIZE*6))
  {
      pmtvals6.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
15      if (MaxStatus == 1)                      // if this is max
      {
          OpStsArray[0]=5;                    // Set Par0 to array
          OpStsArray[1]=1%ARRAYSIZE;          // Set Par1 to element
20      }
      }
  else if ((l >= ARRAYSIZE*6) && (l < ARRAYSIZE*7))
  {
      pmtvals7.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
      if (MaxStatus == 1)                      // if this is max
25      {
          OpStsArray[0]=6;                    // Set Par0 to array
          OpStsArray[1]=1%ARRAYSIZE;          // Set Par1 to element
      }
      }
  else if ((l >= ARRAYSIZE*7) && (l < ARRAYSIZE*8))
30  {
      pmtvals8.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
      if (MaxStatus == 1)                      // if this is max
      {

```

```

        OpStsArray[0]=7;           // Set Par0 to array
        OpStsArray[1]=1%ARRAYSIZE; // Set Par1 to element
    }
}
5  else if ((l >= ARRAYSIZE*8) && (l < ARRAYSIZE*9))
    {
        pmtvals9.pmt[1%ARRAYSIZE]=CurVid; // Store Current Video
        if (MaxStatus == 1)                // if this is max
        {
10         OpStsArray[0]=8;           // Set Par0 to array
            OpStsArray[1]=1%ARRAYSIZE; // Set Par1 to element
        }
    }
    else if ((l >= ARRAYSIZE*9) && (l < ARRAYSIZE*10))
15  {
        pmtvals10.pmt[1%ARRAYSIZE]=CurVid; // Store Current Video
        if (MaxStatus == 1)                // if this is max
        {
            OpStsArray[0]=9;           // Set Par0 to array
20         OpStsArray[1]=1%ARRAYSIZE; // Set Par1 to element
        }
    }
    else if ((l >= ARRAYSIZE*10) && (l < ARRAYSIZE*11))
    {
25         pmtvals11.pmt[1%ARRAYSIZE]=CurVid; // Store Current Video
        if (MaxStatus == 1)                // if this is max
        {
            OpStsArray[0]=10;          // Set Par0 to array
            OpStsArray[1]=1%ARRAYSIZE; // Set Par1 to element
30         }
        }
    else if ((l >= ARRAYSIZE*11) && (l < ARRAYSIZE*12))
    {
        pmtvals12.pmt[1%ARRAYSIZE]=CurVid; // Store Current Video

```

```

        if (MaxStatus == 1)                // if this is max
        {
            OpStsArray[0]=11;                // Set Par0 to array
            OpStsArray[1]=1%ARRAYSIZE;        // Set Par1 to element
5      }
    }
    else if ((l >= ARRAYSIZE*12) && (l < ARRAYSIZE*13))
    {
        pmtvals13.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
10     if (MaxStatus == 1)                // if this is max
        {
            OpStsArray[0]=12;                // Set Par0 to array
            OpStsArray[1]=1%ARRAYSIZE;        // Set Par1 to element
15     }
    }
    else if ((l >= ARRAYSIZE*13) && (l < ARRAYSIZE*14))
    {
        pmtvals14.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
20     if (MaxStatus == 1)                // if this is max
        {
            OpStsArray[0]=13;                // Set Par0 to array
            OpStsArray[1]=1%ARRAYSIZE;        // Set Par1 to element
        }
    }
25     else if ((l >= ARRAYSIZE*14) && (l < ARRAYSIZE*15))
    {
        pmtvals15.pmt[l%ARRAYSIZE]=CurVid;    // Store Current Video
        if (MaxStatus == 1)                // if this is max
        {
30         OpStsArray[0]=14;                // Set Par0 to array
            OpStsArray[1]=1%ARRAYSIZE;        // Set Par1 to element
        }
    }
    else // if ((l >= ARRAYSIZE*15) && (l < ARRAYSIZE*16))

```

```

    {
        pmtvals16.pmt[1%ARRAYSIZE]=CurVid;    // Store Current Video
        if (MaxStatus == 1)                      // if this is max
        {
5           OpStsArray[0]=15;                    // Set Par0 to array
            OpStsArray[1]=1%ARRAYSIZE;          // Set Par1 to element
        }
    }

10     MaxStatus=0;
        watchdog_update();
    }

    return MaxVid;

15 }

unsigned CheckParms(void)
    {
20     signed long FinalPos;
        unsigned short RetVal;

        FinalPos = ZPosCmd - (NumFrames-1)*ZStep;
        if (FinalPos < 0)
25         RetVal = UNDERFLOW;
        else
            RetVal = RDY_FOR_VOL;
        return (RetVal);
    }

30 unsigned SelfTest (void)
    {
        unsigned AllErrors, OldCtrl;

```

```

AllErrors = 0;
OldCtrl = Z_Control;
// turn all test bits off
SetCtrl ( Z_Control & (~(SPARE_TEST1 | SPARE_TEST2 | TEST_LEDn |
5                               PKDET_MODE)));
// test
if ((ReadPort (STS_REG) & (TEST_IN1 | TEST_IN2 | PKDET_MODE |
                               LED_TST_IN)) != 0)
    AllErrors = AllErrors | 0x1;
10 // turn spare test 1 on
SetCtrl ( Z_Control | SPARE_TEST1);
// test
if ((ReadPort (STS_REG) & TEST_IN1) != TEST_IN1)
    AllErrors = AllErrors | 0x2;
15 // turn spare test 1 off, spare test 2 on
SetCtrl ( (Z_Control | SPARE_TEST2) & (~ SPARE_TEST1));
// test
if ((ReadPort (STS_REG) & (TEST_IN1 | (TEST_IN2))) != TEST_IN2)
    AllErrors = AllErrors | 0x4;
20 // turn spare test 2 off, peak det on
SetCtrl ( (Z_Control | PKDET_MODE) & (~ SPARE_TEST2));
// test
if ((ReadPort (STS_REG) & (TEST_IN1 | TEST_IN2 | PKDET_IN))
    != PKDET_IN)
25     AllErrors = AllErrors | 0x8;
// turn peak det off, led on
SetCtrl ( (Z_Control | TEST_LEDn) & (~ PKDET_MODE));
// test
if ((ReadPort (STS_REG) & (TEST_IN1 | TEST_IN2 | PKDET_IN |
30                               LED_TST_IN)) != LED_TST_IN)
    AllErrors = AllErrors | 0x10;
// set control to original value
SetCtrl (OldCtrl);
return (AllErrors);

```

```

    }

    /*******
    /******* WHEN STATEMENTS *****/
5  /*******
    when (reset)
    {
        unsigned long i,Zmove,Thresh;
        short dummy;

10
        ResetSts = 0;           // indicate reset to host
        SetCtrl(ZCTRL_INIT);
        AFocusThresh = INIT_THRESH;
        MoveZ (0);
15
        MuxSel = SEL_VIDEO;

        if (((Z_Status = ReadPort (STS_REG)) & TEST_JMPRn) == 0) // if test mode
        {
            // exercise everything (SYNC_MODE,RST_FOCUS*,RST_INTEG*,
20
                                   PKDET_MODE
            // W_CNTL*, W_THRESH*, W_DACLO*, W_DACHI*, R_STS*, R_ADC*,
                                   RST_ZSYNC*)

            SetCtrl (0);
25
            while (TRUE)           // until next reset
            {
                watchdog_update();
                SetCtrl(~Z_Control); // toggle control bits
                dummy=ReadPort (STS_REG);
30
                Z_Status = ReadPort(ADC_ADDR);
                for(i=0; i<DAC_RANGE; i+=32)
                {
                    //AFocusThresh++;
                    //ZPosCmd++;

```

```

        Zmove=i<<DAC_SHIFT;
        WritePort(POS_CMD_LO, ((unsigned int)(Zmove & 0xFF))); //
        WritePort(POS_CMD_HI, ((unsigned int)((Zmove >> 8) & 0xFF)));
                                                    //& 0xFF
5        WritePort(THRESH_DAC, (unsigned int)(i>>4)); // set threshold
        post_events();
    }
    for(i=DAC_RANGE; i>=65; i-=32)
    {
10        //AFocusThresh++;
        //ZPosCmd++;
        Zmove=i<<DAC_SHIFT; //
        WritePort(POS_CMD_LO, ((unsigned int)(Zmove)& 0xFF));
        WritePort(POS_CMD_HI, ((unsigned int)((Zmove >> 8) & 0xFF))); //
15        WritePort(THRESH_DAC, (unsigned int)(i>>4)); // set threshold
        post_events();
    }
}
}
}
20 else
{
    SetCtrl (TEST_LEDn); // set control bits: async mode,
                        // reset all FF's, reset integ,
                        // turn off test LED
25    WritePort (THRESH_DAC, AFocusThresh);
    SetCtrl (TEST_LEDn | RST_FOCUSn | RST_INTEGn);
}

}

30
when (nv_update_occurs (AFocusThresh))
// Update the auto focus threshold DAC
{
    WritePort (THRESH_DAC, AFocusThresh);

```

```

    }

    when (nv_update_occurs (ZPosCmd))
5   // update the position command DAC
    {
        if (ZPosCmd > DAC_RANGE)
            ZPosCmd = DAC_RANGE;
        SetCtrl(Z_Control & (~ SYNC_MODE)); //async mode
10   MoveZ (ZPosCmd);
    }

    when (nv_update_occurs (Z_Control))
15   // update the control word
    {
        WritePort (CTRL_REG, Z_Control);
    }

20
    // *****
    // ***** remote procedure calls *****
    // *****
    // when (nv_update_occurs) is not used here so that OpMode can
25 // be monitored within the when clauses, specifically to break
    // out of the MONITOR_FOCUS mode. OpMode is changed to NULL_MODE
    // at the end of each clause to ensure that they are only invoked
    // once.

30   when (OpMode == SELF_TEST)
    {
        OpStatus = SelfTest();
        OpMode = NULL_MODE;
    }

```



```
when (OpMode == AF_FAST)
{
    // Perform fast autofocus function
    MuxSel = SEL_VIDEO;
5    PosStatus = AFFast();
    OpMode = NULL_MODE;
}

when (OpMode == AF_PEAK_SYNC)
10 {
    // Perform fast autofocus function
    MuxSel = SEL_VIDEO;
    PosStatus = AFPeakSync();
    OpMode = NULL_MODE;
15 }

when (OpMode == AF_PMTVALS)
{
    // Perform autofocus function which writes PMT values
20 MuxSel = SEL_VIDEO;
    PosStatus = AFPMTVals();
    OpMode = NULL_MODE;
}

25 when (OpMode == READ_VIDEO)
    // Read ADC on stated mux channel)
    {
        MuxSel = SEL_VIDEO;
        OpStatus = ReadPort (ADC_ADDR);
30 OpMode = NULL_MODE;
    }

when (OpMode == RD_VID_LP)
    // call function to read PMT values
```

```

    {
        MuxSel = SEL_VIDEO;
        OpStatus = ReadVidLoop();
        OpMode = NULL_MODE;
5      }

    when (OpMode == READ_HV)
        // Read ADC on stated mux channel)
        {
            MuxSel = SEL_HV;
10       OpStatus = ReadPort (ADC_ADDR);
            OpMode = NULL_MODE;
        }

    when (OpMode == READ_FDBK)
15     // Read ADC on stated mux channel)
        {
            MuxSel = SEL_FDBK;
            OpStatus = ReadPort (ADC_ADDR);
            OpMode = NULL_MODE;
20     }

    when (OpMode == SET_FOR_VOL)
        // setup for volume acquisition)
        {
25     unsigned FrameCtr;
        unsigned long ZPosShift;

        if ((OpStatus = CheckParms()) == RDY_FOR_VOL)
            {
30         // clear ZSYNC FF, wait for syncs and step when they come
            // assumes have moved to start position (closest to objective)
            // assumes ZSYNCs have been stopped, so clear ZSYNC FF and wait

            //SetCtrl (Z_Control & (~RST_ZSYNCn));           // clear ZSYNC

```

```

ResetZSync();
SetCtrl (Z_Control | SYNC_MODE); // go to sync mode
// load DAC register with current position so it stays put
// on first ZSYNC
5   ZPosShift = ZPosCmd << DAC_SHIFT;
    WritePort (POS_CMD_LO, (unsigned int)(ZPosShift & 0xFF));
    WritePort (POS_CMD_HI, (unsigned int)(ZPosShift >> 8) & 0xFF);
    // alert host that ready for volume
    OpStatus = RDY_FOR_VOL;
10  post_events();
    // take the volume)
    for (FrameCtr = 0; FrameCtr < NumFrames; FrameCtr++)
    {
        // wait for ZSYNC
15  while ((ReadPort(STS_REG) & ZSYNCDn) != 0)
        ;
        // reset the ZSYNC flip flop
        ResetZSync();
        // step down
20  ZPosCmd = (ZPosCmd - ZStep);
        ZPosShift = ZPosCmd << DAC_SHIFT;
        WritePort (POS_CMD_LO, (unsigned int)(ZPosShift & 0xFF));
        WritePort (POS_CMD_HI, (unsigned int)(ZPosShift >> 8) & 0xFF);
        watchdog_update(); // prevent reset
25  }
    SetCtrl (Z_Control & (~SYNC_MODE)); // async mode
    OpStatus = VOLUME_COMPL;
    }
    PosStatus = ZPosCmd;
30  OpMode = NULL_MODE;
    }

when (OpMode == MONITOR_FOCUS)
    // Clear focus flipflop and wait for trip

```

```

    {
        // we assume the autofocus flip flop threshold has been set
        // note that this loop should look at RETRACTDn, if it were used
        // clear the focus FF
5      SetCtrl (Z_Control & (~RST_FOCUSn));
        SetCtrl (Z_Control | RST_FOCUSn);
        if ((ReadPort (STS_REG) & FOCUSDn) == 0) // if trip immediately
            OpStatus = THRESH_TOO_LOW;
        else
10      OpStatus = COMPARATOR_READY;
        post_events();
        // now hang out and wait for tripped
        while (((ReadPort (STS_REG) & FOCUSDn) != 0) &&
            (OpMode == MONITOR_FOCUS))
15      {
            post_events();
        }
        if (OpMode == MONITOR_FOCUS)
            OpStatus = COMP_TRIPPED;
20      else
            OpStatus = ABORTED;
        post_events();
        OpMode = NULL_MODE;
    }
25
    when (OpMode == RD_STATUS)
        // read status reg
        {
            OpStatus = ReadPort (STS_REG);
30      OpMode = NULL_MODE;
        }

    when (OpMode == ABORT_CMD)
        // this shouldn't happen unless we are in another op mode

```

```
{  
  OpStatus = ABORTED;  
  OpMode = NULL_MODE;  
}  
□
```

5

**A METHOD AND APPARATUS FOR PERFORMING  
AN AUTOMATIC FOCUS OPERATION**

10

15

Christopher R. Fairley  
Timothy V. Thompson  
Ken K. Lee

20

**APPENDIX C**

25

**M-2464-1P US**

```

module fastz_02
title 'Fastz ctrl pld
Chris Fairley Ultrapointe Corp. San Jose CA'

```

```

5  " FASTZ_02.abl
   " Fast-Z axis logic
   " Ultrapointe Corp.
   " !!!! PROPRIETARY !!!!
   " FOR INTERNAL USE ONLY
10 "
   "
   " Target Device: GAL20RA10
   " Target Socket: U29 of FAST Z-AXIS CTRLLR
   "
15 " Revised:
   "

```

```

fastz_02 device 'P20RA10';
PLOADn      pin 1;
20 CLK4MHZ    pin 2;
R_ADCn      pin 3;
ZSYNCn      pin 4;
RST_ZSYNCn  pin 5;
TRIPPEDn    pin 6;
25 RST_FOCUSn pin 7;
READ_WRITEn pin 8;
ADDR_DATAn  pin 9;
SYNC_MODE   pin 10;
W_DACHIn    pin 11;
30 OUT_En    pin 13;
ADC_BSYn    pin 14;
PKDET_MODE  pin 15;
BUSY2n      pin 16;
BUSY2n      istype 'reg_D, pos';
35 RST_PKDET  pin 17;
RST_PKDET   istype 'reg_D, pos';

```

	LD_DACn	pin 18;
	LD_DACn	istype 'com, pos';
	BUSY1n	pin 19;
	BUSY1n	istype 'reg_D, pos';
5	XCVR_RWn	pin 20;
	XCVR_RWn	istype 'com, pos';
	CNVRTn	pin 21;
	CNVRTn	istype 'com, pos';
	ZSYNCDn	pin 22;
10	ZSYNCDn	istype 'reg_D, pos';
	FOCUSDn	pin 23;
	FOCUSDn	istype 'reg_D, pos';
	H, L, X, C	= 1, 0, .X., .C.;
15	RWn	= READ_WRITEn;
	ADn	= ADDR_DATAn;
	ZSn	= ZSYNCn;
	RZSn	= RST_ZSYNCn;
	SMOD	= SYNC_MODE;
20	WDACn	= W_DACHIn;
	PKMOD	= PKDET_MODE;
	RADCn	= R_ADCn;
	ZSDn	= ZSYNCDn;
	LDACn	= LD_DACn;
25	CVRTn	= CNVRTn;
	AA,BB,FF	= TRIPPEDn,RST_FOCUSn,PLOADn;
	CC,DD,EE	= H,H,H;

## equations

30	* LD DAC on high byte write in async mode,	
	* when ZSYNC happens in sync mode	
	LD_DACn	= !( (SYNC_MODE & !ZSYNCDn) #
		(!SYNC_MODE & !W_DACHIn) );
	LD_DACn.oe	= H;
35		



```

" always writing in address mode, follow R/W in data mode
XCVR_RWn      = !((!ADDR_DATAn & !READ_WRITE) #
                ADDR_DATAn);
XCVR_RWn.oe    = H;

5
" ADC convert when read ADC in normal mode.
" assert and hold convert after ZSYNC happens
CNVRTn        = ! ( (!PKDET_MODE & !R_ADCn) #
                    (PKDET_MODE & !ZSYNCDn & RST_ZSYNCDn) #
                    (PKDET_MODE & !CNVRTn & RST_ZSYNCDn) );
10
CNVRTn.oe      = H;

" edge not required for FOCUSDn, so if already in focus
" after RST_FOCUSn, will set anyway without edge.
15
FOCUSDn       := H;
FOCUSDn.oe     = H;
FOCUSDn.clk    = H;
FOCUSDn.ap     = !TRIPPEDn;
FOCUSDn.ar     = !RST_FOCUSn;
20

" ZSYNCDn edges are reliable and may want to reset ZSYNC
" while ZSYNCDn is low (wide pulse), so use CLK input
ZSYNCDn       := L;
ZSYNCDn.oe     = H;
25
ZSYNCDn.clk    = !ZSYNCDn;
ZSYNCDn.ap     = L;
ZSYNCDn.ar     = !RST_ZSYNCDn;

" delayed versions of ADC_BSYn for pulse shaping
30
BUSY1n        := ADC_BSYn;
BUSY1n.oe      = H;
BUSY1n.clk     = CLK4MHZ;
BUSY1n.ap      = L;
BUSY1n.ar      = L;
35

```

```

        BUSY2n      := BUSY1n;
        BUSY2n.oe    = H;
        BUSY2n.clk   = CLK4MHZ;
        BUSY2n.ap     = L;
5         BUSY2n.ar    = L;

        " reset peak detector as soon as conversion complete
        RST_PKDET    := ADC_BSYn & !BUSY2n;
        RST_PKDET.oe  = H;
10        RST_PKDET.clk = CLK4MHZ;
        RST_PKDET.ap   = L;
        RST_PKDET.ar   = L;

15        " ZSYNCDn, FOCUSDn:
        test_vectors
        (( FF, ZSYNCDn, RST_ZSYNCDn, TRIPPEDn, RST_FOCUSn ] -> [ ZSYNCDn, FOCUSDn ])
        [ EE, H , L , H , L ] -> [ H , H ];
        [ EE, H , L , H , H ] -> [ H , H ];
20        [ EE, H , H , H , H ] -> [ H , H ];
        [ EE, L , H , H , H ] -> [ L , H ];
        [ EE, L , H , H , H ] -> [ L , H ];
        [ EE, H , H , H , H ] -> [ L , H ];
        [ EE, H , L , H , H ] -> [ H , H ];
25        [ EE, H , H , H , H ] -> [ H , H ];
        [ EE, L , H , H , H ] -> [ L , H ];
        [ EE, H , L , H , H ] -> [ H , H ];
        [ EE, H , H , L , H ] -> [ H , L ];
        [ EE, H , H , L , H ] -> [ H , L ];
30        [ EE, H , H , H , H ] -> [ H , L ];
        [ EE, H , H , H , L ] -> [ H , H ];
        [ EE, H , H , H , L ] -> [ H , H ];
        [ EE, H , H , H , H ] -> [ H , H ];
        [ EE, H , H , L , H ] -> [ H , L ];
35        [ EE, H , H , H , L ] -> [ H , H ];

```

```

" XCVR_RWn:
test_vectors
  ([ FF, AA, BB, ADn, RWn ] -> [XCVR_RWn])
    [ EE, CC, DD, H , H ] -> [ L  ];
5    [ EE, CC, DD, H , L ] -> [ L  ];
    [ EE, CC, DD, L , H ] -> [ H  ];
    [ EE, CC, DD, L , L ] -> [ L  ];

" LD_DACn:
10 test_vectors
  ([FF, AA, BB, ZSn, RZSn, SMOD, WDACn, PKMOD, RADCN] -> [ZSDn, LDACn, CVRTn])
    [ EE, CC, DD, H , L , L , H , L , H ] -> [ H , H , H ];
    [ EE, CC, DD, H , H , L , H , L , H ] -> [ H , H , H ];
    [ EE, CC, DD, H , H , L , L , L , H ] -> [ H , L , H ];
15 [ EE, CC, DD, L , H , L , H , L , H ] -> [ L , H , H ];
    [ EE, CC, DD, H , H , L , L , L , H ] -> [ L , L , H ];
    [ EE, CC, DD, H , L , H , L , L , H ] -> [ H , H , H ];

    [ EE, CC, DD, H , H , H , H , L , H ] -> [ H , H , H ];
20 [ EE, CC, DD, L , H , H , H , L , H ] -> [ L , L , H ];
    [ EE, CC, DD, L , H , H , H , L , H ] -> [ L , L , H ];
    [ EE, CC, DD, H , L , H , H , L , H ] -> [ H , H , H ];
    [ EE, CC, DD, H , H , H , H , L , H ] -> [ H , H , H ];
    [ EE, CC, DD, L , H , H , H , L , H ] -> [ L , L , H ];
25 [ EE, CC, DD, H , H , H , H , L , H ] -> [ L , L , H ];
    [ EE, CC, DD, H , H , L , H , L , H ] -> [ L , H , H ];
    [ EE, CC, DD, H , H , H , H , L , H ] -> [ L , L , H ];
    [ EE, CC, DD, H , L , H , H , L , H ] -> [ H , H , H ];

" CNVRTn:
30 "[FF, BB, ZSn, RZSn, SMOD, WDACn, PKMOD, RADCN] -> [ ZSDn, LDACn, CVRTn ]
    [ EE, CC, DD, H , H , L , H , L , H ] -> [ H , H , H ];
    [ EE, CC, DD, H , H , L , H , L , L ] -> [ H , H , L ];
    [ EE, CC, DD, H , H , L , H , L , L ] -> [ H , H , L ];
    [ EE, CC, DD, H , H , L , H , L , H ] -> [ H , H , H ];
35 [ EE, CC, DD, H , H , L , H , H , H ] -> [ H , H , H ];
    [ EE, CC, DD, L , H , L , H , H , H ] -> [ L , H , L ];

```

```

[ EE,CC,DD, H , H , L , H , H , H ]-> [ L , H , L ];
[ EE,CC,DD, H , H , L , H , H , H ]-> [ L , H , L ];
[ EE,CC,DD, H , H , L , H , L , H ]-> [ L , H , H ];
[ EE,CC,DD, H , H , L , H , H , H ]-> [ L , H , L ];
5 [ EE,CC,DD, H , L , L , H , H , H ]-> [ H , H , H ];
[ EE,CC,DD, H , H , L , H , H , H ]-> [ H , H , H ];
[ EE,CC,DD, H , H , L , H , H , H ]-> [ H , H , H ];

```

```

" RST_PKDET:

```

```

10 test_vectors

```

```

([ FF, AA,BB, CLK4MHZ, ADC_BSYn ]-> [ BUSY1n, BUSY2n, RST_PKDET ])

```

```

[ EE, CC,DD, C , H ]-> [ H , H , L ];
[ EE, CC,DD, C , H ]-> [ H , H , L ];
[ EE, CC,DD, C , H ]-> [ H , H , L ];
15 [ EE, CC,DD, C , H ]-> [ H , H , L ];
[ EE, CC,DD, C , H ]-> [ H , H , L ];
[ EE, CC,DD, C , H ]-> [ H , H , L ];
[ EE, CC,DD, C , H ]-> [ H , H , L ];
[ EE, CC,DD, C , H ]-> [ H , H , L ];
20 [ EE, CC,DD, C , H ]-> [ H , H , L ];
[ EE, CC,DD, C , H ]-> [ H , H , L ];

```

```

end fastz_02

```

5

**A METHOD AND APPARATUS FOR PERFORMING  
AN AUTOMATIC FOCUS OPERATION**

10

15

**Christopher R. Fairley  
Timothy V. Thompson  
Ken K. Lee**

20

**APPENDIX D**

25

**M-2464-1P US**

## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 1 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

Alternate BOM

5	Level-I0-Part Number	Description	QTY/Per	QTY--Effective--UM	ITM	
			BAS	FROM	TO	#ON Process-Op
						DRW
	0 000491	PCB,FAST-Z AXIS CNTRLR ASSY				EA
10	1 000061	CONN, RJ-45 PHONE JACK J4 J5	2.	42393	999999	EA 0
	1 000066	DIODE, SA5.0A, TRANSORB D12	1.	70192	999999	EA 72
	1 000067	DIODE, LED, RED, T1-1/4 D1 D11	2.	50193	999999	EA 0
15	1 000068	CAP, 0.047 UF, 50V, CER, AXIAL36. C11 C12 C13 C16 C18 C19 C2 C21 C24 C26 C29 C3 C30 C31 C32 C33 C35 C36	42493	999999	EA	0
20		C37 C38 C39 C4 C40 C41 C42 C43 C44 C45 C46 C49 C5 C52 C53 C54 C57 C9				
	1 000069	CAP, 470 PF, 50V, CER, AXIAL C6	1.	70192	999999	EA 17
25	1 000073	CONN, SOC, 18 PIN, FEM J6	1.	70192	999999	EA 27
	1 000074	CONN, SOC, 6 PIN, FEM J3	1.	70192	999999	EA 70
	1 000078	CAP, 33 UF, 16V, LYTIC, RAD 3. C50 C55 C56	70192	999999	EA	42
30	1 000079	SOC, 24 PIN, DIP, 0.3" CTRS (U26)	1.	70192	999999	EA 71
	1 000101	IC, DS1233 ECONO RESET U2	1.	70192	999999	EA 33
35	1 000105	IC, LM7805, REG, 5V U27	1.	70192	999999	EA 38

## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 2 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

## Alternate BOM

5	Level-I0-Part Number	Description	QTY/Per	QTY--Effective--UM BAS FROM TO	ITM #ON Process-Op DRW
	0 000491	PCB,FAST-Z AXIS CNTRLR ASSY		EA	
10	1 000148	IC,74ALS245,OCTAL BUS XCVR 1. U17	70192	999999 EA	3
	1 000152	IC,74ALS574,OCTAL FLIP-FLOP 3. U12 U19 U7	70192	999999 EA	7
	1 000165	RES, 75 OHM, 1%, 1/4W,MF 1. R37	70192	999999 EA	60
15	1 000166	RES, 100 OHM, 1%, 1/4W,MF 6. R16 R18 R25 R38 R39 R6	42493	999999 EA	0
	1 000168	RES, 261 OHM, 1%, 1/4W,MF 2. R3 R43	70192	999999 EA	44
20	1 000171	RES, 20K OHM, 1%, 1/4W,MF 1. R29	70192	999999 EA	61
	1 000172	RES, 10K OHM, 1%, 1/4W,MF 8. R1 R10 R11 R13 R15 R28 R36 R9	70192	999999 EA	48
25	1 000173	RES, 5.11K OHM, 1%, 1/4W,MF 1. R30	70192	999999 EA	62
	1 000174	RES, 1K OHM, 1%, 1/4W,MF 4. R19 R35 R41 R42	70192	999999 EA	45
	1 000176	IC, AD712JN, DUAL OP AMP 3. U11 U22 U6	70192	999999 EA	13
30	1 000179	CONN,SMB,JACK,PCB,RT ANGLE,75 2. J7 J8	70192	999999 EA	66
	1 000194	CAP,1000PF,50V,CER,AXIAL 2. C15 C48	70192	999999 EA	18
35	1 000197	CAP,33PF,100V,CER,AXIAL 2. C10 C28	70192	999999 EA	25
	1 000199	TESTPOINT,PCB,STR,YEL 12.	70192	999999 EA	67

## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 3 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

Alternate BOM

5	Level-I0-Part Number	Description	QTY/Per QTY--Effective--UM			ITM
			BAS	FROM	TO	#ON Process-Op
						DRW
	0 000491	PCB,FAST-Z AXIS CNTRLR ASSY				EA
		TP1 TP10 TP11 TP12 TP2 TP3				
10		TP4 TP5 TP6 TP7 TP8 TP9				
	1 000246	SW SPST PCB NO MOM 6 MM 2.	70192	999999	EA	43
		S1 S2				
	1 000279	MODULE,LON XCVR,TP RS-485 1.	70192	999999	EA	68 M1
	1 000281	DIODE, 1N5711, SCHOTTKY 6.	70192	999999	EA	9
15		D10 D3 D4 D6 D7 D9				
	1 000283	CAP,1UF,50V,CER,AXIAL 1.	70192	999999	EA	26
		C25				
	1 000289	IC,AD843,OP-AMP,FET INPUT 1.	70192	999999	EA	16
		U14				
20	1 000292	IC,LM336-2.5,VOLTAGE REF 1.	42393	999999	EA	0
		U1				
	1 000295	CAP,6.8UF,20V,TANT,RADIAL 2.	70192	999999	EA	41
		C17 C51				
	1 000300	RES,1M OHM,1%,1/4W,MF 2.	42393	999999	EA	0
25		R14 R27				
	1 000302	RES,100K OHM,1%,1/4W,MF 2.	70192	999999	EA	46
		R26 R40				
	1 000305	RES,12.1K OHM,1%,1/4W,MF 1.	42393	999999	EA	0
		R17				
30	1 000355	SCR,PAN,PHL,SS,6-32X1/4,SEM 6.	50193	999999	EA	0
		(M1)				
	1 000400	CONN,3 PIN,HDR,PCB,RT ANGLE 1.	70192	999999	EA	28
		J1				
	1 000401	CONN,4 PIN,HDR,PCB,RT ANGLE 1.	70192	999999	EA	29
35		J10				
	1 000402	CONN,5 PIN,HDR,PCB,RT ANGLE 1.	70192	999999	EA	30



## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 4 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

## Alternate BOM

5	Level-I0-Part Number	Description	QTY/Per	QTY--Effective--UM	ITM	BAS	FROM	TO	#ON	Process-Op
	0	000491	PCB,FAST-Z AXIS CNTRLR ASSY			EA				
		J2								
10	1	000491SD	PCB,FAST-Z AXIS CNTRLR SCHEM 0.	50193	999999	EA			0	
	1	000493	PCB,FAST-Z AXIS CNTRLR FAB 1.	70192	999999	EA			74	
	1	000495	IC,PA85,POWER OP-AMP	1.	70192	999999	EA		40	
		U8								
	1	000500	IC,AD7575,8-BIT ADC,5 USEC	1.	70192	999999	EA		15	
15		U9 1	000501	IC,AD7541,12-BIT DAC	1.	70192				
	999999	EA	14							
		U4								
	1	000502	IC,DAC0808,8-BIT DAC	1.	70192	999999	EA		31	U15
	1	000504	IC,DG211,QUAD ANALOG SWITCH	1.	70192	999999	EA		32	
20		U13								
	1	000506	CONN,8 PIN,HDR,PCB,STR,MOLEX	1.	70192	999999	EA		39	
		J9								
	1	000507	CAP,68PF,200V,5%,CER,RAD	2.	22294	999999	EA		0	
		C14	C7							
25	1	000509	CAP,0.1UF,200V,20%,CER,RAD	2.	70192	999999	EA		21	
		C1	C34							
	1	000537	IC,74LS14,HEX SCHMITT TRGR INV	1.	70192	999999	EA		1	
		U18								
	1	000538	IC,74ALS138,3-8 DECODER	1.	70192	999999	EA		2	
30		U25								
	1	000539	IC,74ALS273,OCTAL FLIP-FLOP	3.	70192	999999	EA		4	
		U10	U24	U5						
	1	000540	IC,74ALS541,OCTAL BUFFER	1.	70192	999999	EA		5	
		U20								
35	1	000541	IC,74ALS573,OCTAL XPRNT LATCH	1.	70192	999999	EA		6	
		U16								

-157-

RECTIFIED SHEET (RULE 91)

## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 5 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

## Alternate BOM

5	Level-I/O-Part Number	Description	QTY/Per	QTY--Effective--UM	ITM	#ON Process-Op
			BAS	FROM	TO	DRW
	0 000491	PCB.FAST-Z AXIS CNTRLR ASSY			EA	
	1 000542	IC,74ALS139,DUAL 2-4 DECODER 1.	1.	70192	999999	EA 8
10		U23				
	1 000544	IC,AD589,1.23V REF	1.	70192	999999	EA 12
		D2				
	1 000545	CAP,100PF,100V,CER,AXIAL 1.	1.	70192	999999	EA 24
		C23				
15	1 000546	IC,LM311,COMPARATOR	1.	70192	999999	EA 36
		U21				
	1 000547	RES,0 OHM,JUMPER	2.	120193	999999	EA 0
		R4 R7				
	1 000548	RES,16.2K OHM,1%,1/4W,MF 1.	1.	70192	999999	EA 49
20		R21				
	1 000549	RES,1.62K OHM,1%,1/4W,MF 1.	1.	70192	999999	EA 50
		R12				
	1 000550	RES,68.1K OHM,1%,1/4W,MF 1.	1.	70192	999999	EA 51
		R24				
25	1 000551	RES,162K OHM,1%,1/4W,MF 1.	1.	70192	999999	EA 53
		R20				
	1 000552	RES,3.83K OHM,1%,1/4W,MF 1.	1.	70192	999999	EA 54
		R23				
	1 000554	RES,10.0 OHM,1%,1/4W,MF 1.	1.	70192	999999	EA 56
30		R2				
	1 000555	RES,8.25 OHM,1%,1/4W,MF 1.	1.	70192	999999	EA 57
		R5				
	1 000556	RES,2.49K OHM,1%,1/4W,MF 2.	2.	70192	999999	EA 58
		R31 R32				
35	1 000557	RES,1.21K OHM,1%,1/4W,MF 2.	2.	70192	999999	EA 63

## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 6 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

## Alternate BOM

5	Level-I0-Part Number	Description	QTY/Per	QTY--Effective--UM			ITM
				BAS	FROM	TO	
							#ON Process-Op DRW
	0 000491	PCB,FAST-Z AXIS CNTRLR ASSY					EA
		R33 R34					
10	1 000562	CAP,0.1UF,50V,POLY,RAD	2.	70192	999999	EA	19
		C22 C27					
	1 000653	NUT,HEX,SS,6-32,W/LOCK WSHR	2.	42393	999999	EA	0
	1 000728	IC,OP-177GP,OP-AMP	1.	42393	999999	EA	0
		U3					
15	1 000760	RES,49.9K,OHM,1%,1/4W,MF	1.	42393	999999	EA	0
		R22					
	1 000803	SCR,BUT,SOC,SS,6-32X3/8	2.	42393	999999	EA	0
	1 000868	DIODE,1N5235A,ZENER,6.8V	1.	42393	999999	EA	0
		D5					
20	1 000872	RES,562 OHM,1%,1/4W,MF	1.	42393	999999	EA	0
		R44					
	1 000874	HEATSINK,TO-3 XSTR	1.	42393	999999	EA	0
		(U8)					
	1 006008	RES,POT,10K OHM,1T,3/8 SQ	1.	22294	999999	EA	0
25	1 011012	DIODE,SA150,TRANSIENT SUPPRESS	1.	80193	999999	EA	0
		D8					
	1 012022	STANDOFF HEX AL F/F,6-32X1/2	3.	50193	999999	EA	0
		(M1)					
	1 012048	CONN,2 PIN SHUNT JUMPER	1.	120193	999999	EA	0
30		(JP1)					
	1 012057	CONN,2 PIN HDR, STRAIGHT	1.	120193	999999	EA	0
		JP1					

## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 1 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

## Alternate BOM

5	Level-I0-Part Number	Description	QTY/Per	QTY-Effective-UM			ITM
				BAS	FROM	TO	
							#ON Process-Op
							DRW
	0	000625	ASSY,X-Y-Z-AF STAGE			EA	
10	1	000250	STAGE SYSTEM X Y Z	1.	70192	999999	EA 1
	1	000253	PZT MTG AFSTAGE	1.	70192	999999	EA 5
	1	000254	FLEXURE AFSTAGE	16.	70192	999999	EA 14
	1	000256	PZT RETAINER AFSTAGE	1.	70192	999999	EA 6
	1	000257	ROTATING BAR AFSTAGE	4.	70192	999999	EA 10
15	1	000258	FLEXURE CLAMP AFSTAGE	32.	70192	999999	EA 15
	1	000259	PRELOAD SCREW AFSTAGE	2.	70192	999999	EA 20
	1	000260	STATIONARY BAR AFSTAGE	4.	70192	999999	EA 9
	1	000261	TOP PLATE AFSTAGE	1.	70192	999999	EA 11
	1	000263	PZT 30 MICRON 100VDC	1.	70192	999999	EA 7
20	1	000264	COVER AFSTAGE	1.	70192	999999	EA 12
	1	000266	VACUUM PLATE,AF STAGE	1.	70192	999999	EA 13
	1	000344	SCR,CAP,SOC,SS,10-32X5/8	10.	70192	999999	EA 33
	1	000357	SCR,CAP,SOC,SS,6-32X3/8	2.	31894	999999	EA 0
	1	000418	SENSOR,PROX,50-150UM,KAMAN	1.	70192	999999	EA 3
25	1	000516	PLATE,ZSTAGE ADAPTOR	1.	70192	999999	EA 2
	1	000530	WAFER CHUCK	1.	70192	999999	EA 18
	1	000581	SCR,CAP,SOC,SS,10-32 X 3/8	4.	31894	999999	EA 0
	1	000614	SCR,PAN,SLT,SS,8-32X1/2,W/SEAL	4.	70192	999999	EA 19
	1	000780	HOSE BARB,,10-32 TO 1/16 HOSE	3.	50193	999999	EA 0
30	1	000781	SOLENOID VALVE,5 PORT,BASE MTG	3.	90193	999999	EA 0
	1	000784	HOSE,1/16ID X 1/8OD,8X,URTHN	3.	50193	999999	EA 0
	1	000798	SCR,BUT,SOC,SS,4-40X3/8	13.	70192	999999	EA 34

## ULTRAPOINTE CORPORATION

SINGLE LEVEL EXPLOSION C/N - 220 CHRISF PAGE 2 RKPSTD01

Format Type Product Structure Fields Type of Components Input and Output

## Alternate BOM

5	Level-I0-Part Number	Description	QTY/Per	QTY--Effective--UM	ITM	#ON Process-Op
			BAS	FROM	TO	DRW
	0 000625	ASSY,X-Y-Z-AF STAGE			EA	
10	1 000799	SCR,SET,SOC,SS,4-40X9/32,NYTIP 1.	70192	999999	EA	35
	1 000805	SCR,CAP,SOC,SS,8-32X3/8	4.	70192	999999	EA 28
	1 000813	SCR,CAP,SOC,SS,10-32X3/4	4.	31894	999999	EA 0
	1 000853	BALL,1/8 DIA,SS	3.	70192	999999	EA 8
	1 000854	SCR,FLT,SOC,SS,6-32X1/2	64.	31894	999999	EA 0
15	1 000858	SCR,SET,SOC,SS,8-32X3/4,CUP PT 4.	31894	999999	EA	0
	1 000860	O-RING,1/16WX1/160DX11/640D NO 3.	70192	999999	EA	24
	1 012050	SCR,FLT HD, PHLPS,SS,#4X0.5 LG 4.	90193	999999	EA	0
	1 012058	CONN,AC PWR INLET,250VAC,10 A 1.	11894	999999	EA	0

5

**A METHOD AND APPARATUS FOR PERFORMING  
AN AUTOMATIC FOCUS OPERATION**

10

**Christopher R. Fairley  
Timothy V. Thompson  
Ken K. Lee**

15

20

**APPENDIX E**

25

**M-2464-1P US**

```

/*****
 * function description:
 *
 * Some functions are not listed
 * These functions are described here
 *****/
*/

```

```

10  cfg_GetLockRdConfig
    Get the pointer to system configuration data structure

    da_GlobalLonPointer
    Get the pointer to system status data structure

15  lon_APIQuery
    Query value of a specified hardware sybsystem parameter

    lon_APICmd
    Set a specified hardware sybsystem parameter to specified value
20

    lonui_SetFastZPos
    Set the fine Z stage to specified position

    lonuiLocalFineAF
25    Do fine stage auto-focus

    util_sginap
    Causes the program to sleep( wait) for specified number of 10 milli-seconds.

30  XtMalloc
    Allocated system memory

    XtFree
    Release previously allocated memory.
35

    longmi_UpdateNV
    Update specified entry in hardware status table

```

lonimg\_QueryNetVar

Query and get specified value(8bit) of hardware subsystem parameter.

5 lonimg\_QueryNetVar16

Query and get specified value(16 bit) of hardware subsystem parameter.

LonTimeBlockReadOneAndDecodeMsg

10 Wait for a message from hardware subsystem, when the message is received,  
decode it to get the message's subsystem parameter ID, and its value.

/\*  
\*\*\*\*\*

\* (c) Copyright 1992,1993,1994 Ultrapointe Corp.

15 \* All rights reserved

\*

\* Auto focus

\*

\*\*\*\*\*

20 \*/

#define WORKING\_DIST\_UM 800

#define WORKING\_DIST\_10X\_UM 1600

25 #define BACKUP\_BEFORE\_AF\_UM 2000 /\* if closer than this backdown first \*/

#define AF2\_INTEN\_DEC .85 /\* if PMT reset, lower by this amount \*/

#define AF2\_INTEN\_DEC\_NORMAL .95

#define PMT\_CLAMPED 1 /\* pmt latched overload \*/

30 static UI\_LON\_IMAGE UIImage;

#define LASER\_I\_PERCENT\_CHG\_PWR 70 /\* above this intensity, change power \*/

#define MAX\_LASER\_PWR 125

#define MIN\_PMT\_GAIN 245

35 #define MAX\_PMT\_ZERO 132

#define MIN\_PMT\_ZERO 70

/\* laser pwr: 0-255 = 0mw-50mw \*/



```

#define MAX_LASER_PWR          125

#define NUM_NV_IN_NODE          50
5  #define GENERIC_OPMODE_NV_ADDR      2
#define GENERIC_OPSTATUS_NV_ADDR  3

#define OPMODE_TIMEOUT_MS          8000
#define OPMODE_TIMEOUT_TICK        (OPMODE_TIMEOUT_MS * CLK_TCK /1000)
10 #define AF_READ_MS                20
#define AF_READ_TICK                (AF_READ_MS * CLK_TCK /1000)
#define VOL_CHK_LOOP_TIME          30 /* 10ms */
#define C_AF_COMPARATOR_RDY        8
#define C_AF_THRES_TOO_LO         9
15 #define C_AF_COMP_TRIP            0xa
#define C_AF_ABORT                 0xb

#define PMT_CLAMPED 1           /* pmt latched overload */
#define LASER_PWR_SLIDER_RANGE    100
20

typedef struct CALLBK_STRUCT_tag
{
    void (*CallBkPtr)(void *, int, int, int); /* FctPtr to callback */
    void *ParmPtr; /* ParameterPtr for callback */
25 void *NextCBPtr;
} CALLBK_STRUCT;

typedef struct HARDWR_IMAGE_tag
{
30 char *NVName;
    int iNVValue;
    char chArray[LON_ARRAY_DATA_SZ]; /* for FastZ intensity array */
    int iNumByte; /* 1 or 2 or 16byte NV */
35 CALLBK_STRUCT *CBStructPtr; /* start of list callback struct */
    NETVAR_DIRECT NVDirection; /*
    BOOL UpdatedFlg;

```

```
} HARDWR_IMAGE_ENTRY :
```

```
typedef struct LON_NODE_IMAGE_tag
```

```
{
```

```
5 char *NodeName;
```

```
HARDWR_IMAGE_ENTRY NVImg[ NUM_NV_IN_NODE];
```

```
} LON_NODE_IMAGE;
```

```
10 /******
```

```
* struct for the hardware image table
```

```
* items in the table are sequential, keyed to the
```

```
* nodeID and NV number
```

```
*
```

```
15 /******
```

```
*/
```

```
static LON_NODE_IMAGE HardWrImgTable[] = {
```

```
{ "dummy" ,
```

```
{ /* 0, node 0 is not used in LON-talk */
```

```
20 { "dummy", 0,{0}, -1,NULL, FALSE, FALSE, },
```

```
} },
```

```
{ "LaserPwrControl",
```

```
{ /* 1, Laser power controller */
```

```
25 /* 0*/ { "SoftWr Ver", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Laser SW Rev,
```

```
*/
```

```
/* 1*/ { "ResetStatus",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Laser ResetStat, */
```

```
/* 2*/ { "ModeCmd", 0,{0}, 1, NULL, eNetVarFromSGI,FALSE }, /* Laser Mode,
```

```
cmd */
```

```
30 /* 3*/ { "LaserStatus",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Laser Status, */
```

```
/* 4*/ { "LaserOn/Off",HDWR_LASER_OFF, {0}, 1, NULL, eNetVarFromSGI,FALSE },
```

```
/* Laser On/Off = 1/0, cmd */
```

```
/* 5*/ { "Run/Idle", HDWR_LASER_IDLE, {0}, 1, NULL, eNetVarFromSGI,FALSE
```

```
}, /* Laser Run/Idle = 1/0,cmd */
```

```
35 /* 6*/ { "Light/Curr", HDWR_LASER_LIGHT, {0}, 1, NULL, eNetVarFromSGI,FALSE
```

```
}, /* Laser ModeSel, cmd light/current 1/0 */
```

```
/* 7*/ { "CurrentCmd", HDWR_DEF_LSR_CUR, {0}, 1, NULL, eNetVarFromSGI,FALSE
```

```

    }, /* Laser Current,      cmd 0-255 */
    /* 8*/ { "PowerCmd",      HDWR_DEF_LSR_PWR, {0}, 1, NULL, eNetVarFromSGI, FALSE
    }, /* Laser Power,        cmd 0-255 */
    /* 9*/ { "WhiteLtpwr",    HDWR_DEF_LSR_WL_PWR, {0}, 1, NULL, eNetVarFromSGI,
5  FALSE }, /* WhiteL Pwr,    cmd 0-255 */
    /*10*/ { "SpareCmd",      0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* SpareCmd cmd
    0-255 */
    /*11*/ { "InterLock",     0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Laserinterlock
    open/close= ??????? */
10  /*12*/ { "CurrCheck",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
    LaserCurrCheck DAC      0-255 */
    /*13*/ { "PwrCheck",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O LaserPwrCheck
    DAC      0-255 */
    /*14*/ { "WLPwrCheck",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
15  WhiteLpwrCheck DAC      0-255 */
    /*15*/ { "SpareCheck",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SpareCheck
    DAC      0-255 */
    /*16*/ { "HeadTempMon",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
    LaserHeadTempMonitor    0-255 */
20  /*17*/ { "SpareTemp",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SpareTemp
    0-255 */
    /*18*/ { "LaserPwrMon",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
    LaserPwrMonitor          0-255 */
    /*19*/ { "CurrMonitor",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
25  laserCurrMonitor          0-255 */
    { NULL,                    0,{0}, -1, NULL, FALSE, FALSE },
    } },
    { "CameraFilterWheel",
    {
    /* 2, camera filter wheel */
30  /* 0*/ { "SoftWr Ver",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev,*/
    /* 1*/ { "ResetStatus",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
    /* 2*/ { "ModeCmd",        0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* Mode, cmd */
    /* 3*/ { "Status",         0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
    /* 4*/ { "PositionCmd",HDWR_LSR_CAM_APASS,{0}, 1, NULL, eNetVarFromSGI, FALSE
35  }, /* wheel position      1-6 */
    /* 5*/ { "HomeStatus",     0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O home
    sensor status on/off= ???? */

```

```

/* 6*/ { "PosStatus",      1,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O wheel
pos status      1-6          */
/* 7*/ { "StepCnt", 1,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O Step Count
??? */
5 /* 8*/ { "StepDelay".600,{0}, 2. NULL, eNetVarFromSGI.FALSE },      /*O Step
delay      ???          */
      { NULL,      0,{0}, -1. NULL, FALSE, FALSE },
    } },
    { "LaserLineFilter",
10 {
      /* 3, laser line filter wheel */
/* 0*/ { "SoftWr Ver", 0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O SW Rev, */
/* 1*/ { "ResetStatus",0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O ResetStat, */

/* 2*/ { "ModeCmd", 0,{0}, 1. NULL, eNetVarFromSGI.FALSE },      /* Mode,
15 cmd          */
/* 3*/ { "Status", 0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O Status, */
/* 4*/ { "PositionCmd".HDWR_LSR_LINE_ALL,{0}, 1. NULL, eNetVarFromSGI.FALSE },/*
wheel position      1-6          */
/* 5*/ { "HomeStatus", 0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O home sensor
20 status on/off= ??? */
/* 6*/ { "PosStatus",      1,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O wheel pos status
      1-6 */
/* 7*/ { "StepCnt", 1,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O Step Count
??? */
25 /* 8*/ { "StepDelay".600,{0}, 2. NULL, eNetVarFromSGI.FALSE },      /*O Step
delay      ???          */
      { NULL,      0,{0}, -1. NULL, FALSE, FALSE },
    } },
    { "LaserAttenFilter",
30 {
      /* 4, laser attenuator filter wheel */
/* 0*/ { "SoftWr Ver", 0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O SW Rev, */
/* 1*/ { "ResetStatus",0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O ResetStat, */
/* 2*/ { "ModeCmd", 0,{0}, 1. NULL, eNetVarFromSGI.FALSE },      /* Mode, cmd */
/* 3*/ { "Status", 0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O Status, */
35 /* 4*/ { "PositionCmd".HDWR_LSR_ATTEN_HIG,{0}, 1. NULL, eNetVarFromSGI.FALSE
      }, /* wheel position      1-6          */
/* 5*/ { "HomeStatus", 0,{0}, 1. NULL, eNetVarToSGI, FALSE },      /*O home sensor

```

```

status on/off= ??? */
/* 6*/ { "PosStatus",      1,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O wheel pos status
1-6 */
/* 7*/ { "StepCnt", 1,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Step Count   ??? */
5 /* 8*/ { "StepDelay".600,{0}, 2, NULL, eNetVarFromSGI.FALSE }, /*O Step delay??? */
    { NULL,      0,{0}, -1, NULL, FALSE, FALSE },
    } },
    { "Bright/DarkField",
    {
10 /* 0*/ { "SoftWr Ver",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev,   */
    /* 1*/ { "ResetStatus".0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
    /* 2*/ { "ModeCmd",      0,{0}, 1, NULL, eNetVarFromSGI.FALSE }, /* Mode, cmd */
    /* 3*/ { "Status",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status,  */
    /* 4*/ { "PositionCmd".HDWR_BRDRK_FLD_BRI,{0}, 1, NULL, eNetVarFromSGI.FALSE
15 }, /* pos select bright/dark= ??? */
    /* 5*/ { "PosStatus",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /* pos status
    bright/dark= ??? */
    /* 6*/ { "Sw1Status",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O sw1 status
    on/off = ??? */
20 /* 7*/ { "Sw2Status",      1,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O sw2 statuson/off
    = ??? */
    /* 8*/ { "StepCnt", 1,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Step Count   ??? */
    { NULL,      0,{0}, -1, NULL, FALSE, FALSE },
    } },
25 { "Page Scanner",
    {
    /* 6, page scanner */
    /* 0*/ { "SoftWr Ver",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev,*/
    /* 1*/ { "ResetStatus".0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
    /* 2*/ { "ModeCmd",      0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* Mode, cmd
30 */
    /* 3*/ { "Status",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
    /* 4*/ { "VertSync",      HDWR_DEF_PAGE_VERT, {0}, 2, NULL, eNetVarFromSGI,
    FALSE }, /* vertSync      0-65536 */
    /* 5*/ { "TraceTiming".HDWR_DEF_PAGE_TRACE, {0}, 2, NULL, eNetVarFromSGI,
35 FALSE }, /* TraceTiming      0-65536 */
    /* 6*/ { "RetraceTime".HDWR_DEF_PAGE_RETRACE,{0}, 2, NULL, eNetVarFromSGI,
    FALSE }, /* RetraceTiming    0-65536 */

```

```

/* 7*/ { "StartAmpOff", HDWR_DEF_PAGE_OFFSET, {0}, 2, NULL, eNetVarFromSGI,
FALSE }, /* StartAmplitude-offset 0-65536 */
/* 8*/ { "Increment", HDWR_DEF_PAGE_INCREM, {0}, 2, NULL, eNetVarFromSGI,
FALSE }, /* Incr Amplitude 0-65536 */
5 /* 9*/ { "Decrement", HDWR_DEF_PAGE_DECREM, {0}, 2, NULL, eNetVarFromSGI,
FALSE }, /* Decr Amplitude 0-65536 */
/*10*/ { "CtlRegValue", HDWR_DEF_PAGE_CTRL_REG, {0}, 1, NULL,
eNetVarFromSGI, FALSE }, /* ControlRegValue 0-255 */
/*11*/ { "XAxisAmpl", HDWR_DEF_PAGE_XAMP, {0}, 2, NULL,
10 eNetVarFromSGI, FALSE }, /* X-axis Amplitude 0-4095 */
/*12*/ { "YAxisAmpl", HDWR_DEF_PAGE_YAMP, {0}, 2, NULL,
eNetVarFromSGI, FALSE }, /* Y-axis Amplitude 0-4095 */
/*13*/ { "XAxisPhase", HDWR_DEF_PAGE_XPHASE, {0}, 2, NULL,
eNetVarFromSGI, FALSE }, /* X-axisPhase 0-4095 */
15 { NULL, 0,{0}, -1, NULL, eNetVarFromSGI, FALSE },
}},
{ "FastZ Ctrl",
{
/* 7, fast-z-axis-controller */
/* 0*/ { "SoftWr Ver", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev, */
20 /* 1*/ { "ResetStatus",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
/* 2*/ { "OpMode", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* OpMode,
cmd */
/* 3*/ { "OpStatus", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O OpStatus, */
/* 4*/ { "OpParArr1", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* cmd, */
25 /* 5*/ { "OpParArr2", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* cmd, */
/* 6*/ { "OpStsArr1", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
/* 7*/ { "OpStsArr2", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
/* 8*/ { "AFocusThresh", HDWR_DEF_Z_AFTHRESH, {0}, 1, NULL, eNetVarFromSGI,
FALSE }, /* AFocusThresh, cmd */
30 /* 9*/ { "ZPositionCmd", HDWR_DEF_ZPOSCMD, {0}, 2, NULL, eNetVarFromSGI, FALSE
}, /* ZPosCmd, cmd */
/*10*/ { "ZControl", HDWR_DEF_ZCONTROL, {0}, 1, NULL, eNetVarFromSGI,
FALSE }, /* ZControl, cmd */
/*11*/ { "MuxSelect", HDWR_DEF_ZMUXSEL, {0}, 1, NULL, eNetVarFromSGI,
35 FALSE }, /* MuxSel, cmd */
/*12*/ { "ZStep", HDWR_DEF_ZSTEP, {0}, 1, NULL, eNetVarFromSGI, FALSE }, /*
ZStep, cmd */

```

```

/*13*/ { "NumFrames", HDWR_DEF_ZNUMFRAME,{0}.1, NULL, eNetVarFromSGI,
FALSE }, /* NumFrames, cmd */
/*14*/ { "ZStatus", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ZStatus, */
/*15*/ { "PosStatus", 0,{0}, 2, NULL, eNetVarToSGI, FALSE }, /*O PositionStatus
5 */
/*16*/ { "Intensity0", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*17*/ { "Intensity1", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*18*/ { "Intensity2", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*19*/ { "Intensity3", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
10 /*20*/ { "Intensity4", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*21*/ { "Intensity5", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*22*/ { "Intensity6", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*23*/ { "Intensity7", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*24*/ { "Intensity8", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
15 /*25*/ { "Intensity9", 0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*26*/ { "Intensity10",0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*27*/ { "Intensity11",0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*28*/ { "Intensity12",0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*29*/ { "Intensity13",0,{0}, 16, NULL, eNetVarToSGI, FALSE },
20 /*30*/ { "Intensity14",0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*31*/ { "Intensity15",0,{0}, 16, NULL, eNetVarToSGI, FALSE },
/*32*/ { "SelfTest", FASTZ_OPCMD_SELFTEST, {0}.1, NULL, eOpModeCmd, FALSE },
/* SelfTest, cmd */
/*33*/ { "AutoFocus", FASTZ_OPCMD_AUTOFOCUS,{0}.1, NULL, eOpModeCmd,
25 FALSE }, /* AutoFocus, cmd */
/*34*/ { "ReadVideo", FASTZ_OPCMD_READVIDEO,{0}.1, NULL, eOpModeCmd,
FALSE }, /* ReadVideo, cmd */
/*35*/ { "ReadFeedBk", FASTZ_OPCMD_READFDBK, {0}.1, NULL, eOpModeCmd,
FALSE }, /* ReadFeedBk, cmd */
30 /*36*/ { "ReadHV", FASTZ_OPCMD_READHV, {0}.1, NULL, eOpModeCmd, FALSE
}, /* ReadHV, cmd */
/*37*/ { "SetForVol", FASTZ_OPCMD_SETFORVOL,{0}.1, NULL, eOpModeCmd,
FALSE }, /* SetForVol, cmd */
/*38*/ { "MonitorFoc", FASTZ_OPCMD_MONITORFOC,{0}.1, NULL, eOpModeCmd,
35 FALSE }, /* monitor focus.cmd */
/*39*/ { "Abort", FASTZ_OPCMD_ABORT, {0}.1, NULL, eOpModeCmd, FALSE }, /*
Abort, cmd */

```

```

/*40*/ { "ReadStatus", FASTZ_OPCMD_RDSTATUS, {0},1, NULL, eOpModeCmd,
FALSE }, /* ReadStatus, cmd */
/*41*/ { "AutoFocus0", FASTZ_OPCMD_FAST_AF0, {0},1, NULL, eOpModeCmd, FALSE
}, /* ReadStatus, cmd */
5 /*42*/ { "AutoFocus1", FASTZ_OPCMD_FAST_AF1, {0},1, NULL, eOpModeCmd, FALSE
}, /* ReadStatus, cmd */
/*43*/ { "AutoFocus2", FASTZ_OPCMD_FAST_AF2, {0},1, NULL, eOpModeCmd, FALSE
}, /* ReadStatus, cmd */
/*44*/ { "AutoFocus3", FASTZ_OPCMD_FAST_AF3, {0},1, NULL, eOpModeCmd, FALSE
10 }, /* ReadStatus, cmd */
/*45*/ { "AutoFocus4", FASTZ_OPCMD_FAST_AF4, {0},1, NULL, eOpModeCmd, FALSE
}, /* ReadStatus, cmd */
/*46*/ { "AutoFocus5", FASTZ_OPCMD_FAST_AF5, {0},1, NULL, eOpModeCmd, FALSE
}, /* ReadStatus, cmd */
15 /*47*/ { "Rd Intensity",FASTZ_OPCMD_RD_INTEN, {0},1, NULL, eOpModeCmd, FALSE },
/* StgAF read inten */
{ NULL, 0,{0}, -1, NULL, eOpModeCmd, FALSE },
}},
{ "PMT Amp",
20 { /* 8, PMT preamp controller */
/* 0*/ { "SoftWr Ver", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev, */
/* 1*/ { "ResetStatus", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
/* 2*/ { "ModeCmd", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* Mode, cmd */
/* 3*/ { "Status", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
25 /* 4*/ { "OpParArr1", HDWR_DEF_PMT_AZFCT, {0}, 1, NULL,
eNetVarFromSGI, FALSE }, /* cmd, */
/* 5*/ { "OpParArr2", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* cmd, */
/* 6*/ { "OpStsArr1", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
/* 7*/ { "OpStsArr2", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
30 /* 8*/ { "ResetFocus", HDWR_PMT_RST_FOCUS_OFF, {0}, 1, NULL,
eNetVarFromSGI, FALSE }, /* resetFocus, */
/* 9*/ { "ResetPeakDet", HDWR_PMT_RST_PKDET_OFF, {0}, 1, NULL,
eNetVarFromSGI, FALSE }, /* resetPeakDet, */
/*10*/ { "FocThresStat", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
35 FocusThresholdExceed */
/*11*/ { "PMTOverloaded",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O overload, */
/*12*/ { "PMTClamped", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O clamped, */

```



```

/*13*/ { "AFThresCmd", HDWR_DEF_PMT_AFTHRESH, {0}, 1, NULL,
eNetVarFromSGI.FALSE }, /* set threshold */
/*14*/ { "PMTVoltCmd", HDWR_DEF_PMT_PMTVOLT, {0}, 1, NULL,
eNetVarFromSGI.FALSE }, /* set PMT voltage, */
5 /*15*/ { "AutoZeroCmd", HDWR_DEF_PMT_AZERO, {0}, 1, NULL,
eNetVarFromSGI.FALSE }, /* start autozero, */
/*16*/ { "AZeroGainK1", HDWR_DEF_PMT_AZEROGAINK1,{0}, 2, NULL,
eNetVarFromSGI.FALSE }, /*O AutoZeroGainK1, */
/*17*/ { "AZeroGainK2", HDWR_DEF_PMT_AZEROGAINK2,{0}, 2, NULL,
10 eNetVarFromSGI.FALSE }, /*O AutoZeroGainK2, */
/*18*/ { "AZeroGainK3", HDWR_DEF_PMT_AZEROGAINK3,{0}, 2, NULL,
eNetVarFromSGI.FALSE }, /*O AutoZeroGainK3, */
/*19*/ { "AvgAutoGainK1",HDWR_DEF_PMT_AVGAGAINK1, {0}, 2, NULL,
eNetVarFromSGI.FALSE }, /*O AvgAutoGainK1, */
15 /*20*/ { "AveAutoGainK2",HDWR_DEF_PMT_AVGAGAINK2, {0}, 2, NULL,
eNetVarFromSGI.FALSE }, /*O AvgAutoGainK2, */
/*21*/ { "AveAutoGainK3",HDWR_DEF_PMT_AVGAGAINK3, {0}, 2, NULL,
eNetVarFromSGI.FALSE }, /*O AvgAutoGainK3, */
/*22*/ { "SelfTest", PMT_OPCMD_SELFTEST, {0}, 1, NULL, eOpModeCmd,
20 FALSE }, /* SelfTest, cmd */
/*23*/ { "AutoZero", PMT_OPCMD_AUTOZERO, {0}, 1, NULL, eOpModeCmd,
FALSE }, /* AutoZero, cmd */
/*24*/ { "PkAutoGain", PMT_OPCMD_PK_AGAIN, {0}, 1, NULL, eOpModeCmd,
FALSE }, /* PeakAutoGain, cmd */
25 /*25*/ { "AvgAutoGain", PMT_OPCMD_AVG_AGAIN, {0}, 1, NULL, eOpModeCmd,
FALSE }, /* AvgAutoGain, cmd */
/*26*/ { "ReadAvgVideo", PMT_OPCMD_RD_AVG_VIDEO, {0}, 1, NULL, eOpModeCmd,
FALSE }, /* ReadAvgVideo, cmd */
/*27*/ { "ReadPkVideo", PMT_OPCMD_RD_PK_VIDEO, {0}, 1, NULL, eOpModeCmd,
30 FALSE }, /* ReadPeakVideo, cmd */
/*28*/ { "ReadPMTHV", PMT_OPCMD_RD_PMT_HV, {0}, 1, NULL, eOpModeCmd,
FALSE }, /* ReadPMTHV, cmd */
/*29*/ { "ReadLaserPwr", PMT_OPCMD_RD_LSR_PWR, {0}, 1, NULL, eOpModeCmd,
FALSE }, /* ReadLaserPwr, cmd */
35 /*30*/ { "ReadVideoZero",PMT_OPCMD_RD_VIDEOZERO, {0}, 1, NULL, eOpModeCmd,
FALSE }, /* ReadVideoZero, cmd */
/*31*/ { "ReadFocus", PMT_OPCMD_RD_FOCUS, {0}, 1, NULL, eOpModeCmd,

```

```

FALSE },          /* ReadFocus,  cmd */
/*32*/ { "ReadPmtOvld", PMT_OPCMD_RD_PMTOVLD, {0}, 1, NULL, eOpModeCmd,
FALSE },          /* ReadPmtOvld, cmd */
        { NULL,      0, {0}, -1, NULL, FALSE, FALSE },
5    } },
    { "Turret/OpticsCtrl",
    {
        /* 9, Turret/Optics controller */
        /* 0*/ { "SoftWr Ver", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev, */
        /* 1*/ { "ResetStatus",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
10    /* 2*/ { "ModeCmd", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* Mode, cmd */
        /* 3*/ { "Status", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
        /* 4*/ { "PositionCmd",HDWR_TURRET_HOME_OBJ, {0}, 1, NULL, eNetVarFromSGI,
FALSE },          /* TurretPosCommand 1-5 */
        /* 5*/ { "WL/Laser Sw",HDWR_LSR_WHLT_WHITE, {0}, 1, NULL, eNetVarFromSGI,
15    FALSE },          /* ModeSw, White/Laser 1/0 */
        /* 6*/ { "Shutter O/C",HDWR_LSR_SHUTTER_CLS, {0}, 1, NULL, eNetVarFromSGI, FALSE
},          /* LaserShutter, open/close 1/0 */
        /* 7*/ { "StepDelay", HDWR_TURRET_STEPDELAY,{0}, 2, NULL, eNetVarFromSGI,
FALSE },          /* TurretStepDelay */
20    /* 8*/ { "ObjStatus", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
ObjectiveSensorStat on/off =??? */
        /* 9*/ { "ClickStatus",1,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
clickSensorStatus on/off= ??? */
        /*10*/ { "HomeStatus", 1,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
25    HomeSensorStatus on/off= ??? */
        /*11*/ { "TurrPosStat",1,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
TurretPosStatus 1-5 */
        { NULL,      0,{0}, -1, NULL, FALSE, FALSE },
    } },
30    { "Cassette Switch",
    {
        /* 10, Cassette switch */
        /* 0*/ { "SoftWr Ver", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev, */
        /* 1*/ { "ResetStatus",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
        /* 2*/ { "ModeCmd", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* Mode, cmd */
35    /* 3*/ { "Status", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
        /* 4*/ { "CassPresent",HDWR_CASS_NOT_PRESENCE,{0}, 1, NULL, eNetVarToSGI, FALSE
},          /*O CassettePresent Yes/no=??? */

```

```

/* 5*/ { "SizeStatus", HDWR_CASS_SIZE_8IN,{0}, 1, NULL, eNetVarToSGI, FALSE },/*O
CassetteSizeStatus      0-3 */
/* 6*/ { "InterlockStat",0,{0}, 1, NULL, eNetVarToSGI, FALSE },      /* Mode, cmd */
      { NULL,          0,{0}, -1, NULL, FALSE, FALSE },
5  } },
  { "Solenoid Ctrl",
    {
      /* 11, Solenoid controller */
/* 0*/ { "SoftWr Ver",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev, */
/* 1*/ { "ResetStatus", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
10 /* 2*/ { "ModeCmd",    0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* Mode, cmd */
/* 3*/ { "Status",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
/* 4*/ { "DoorCmd",      HDWR_SOL_WFDRCMD_CLS, {0}, 1, NULL, eNetVarFromSGI,
FALSE }, /* CmdWaferDoor Open/close= ??? */
/* 5*/ { "ChuckVacCmd",HDWR_SOL_CHKVAC_NONE, {0}, 1, NULL, eNetVarFromSGI,
15 FALSE }, /* CmdWaferChuckSolenoid      0-3 */
/* 6*/ { "DoorStat",      HDWR_SOL_WFDRSTAT_CLS,{0}, 1, NULL, eNetVarToSGI,
FALSE }, /*O Stat WaferDoor open/close=??? */
/* 7*/ { "ChuckVacStat",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
StatChuckVac      on/off= ??? */
20 /* 8*/ { "HouseVacStat",0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O
StatHouseVac      on/off= ??? */
      { NULL,          0,{0}, -1, NULL, FALSE, FALSE },
    } },
  { "Aperture Ctrl",
25  {
      /* 11, Solenoid controller */
/* 0*/ { "SoftWr Ver",    0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O SW Rev, */
/* 1*/ { "ResetStatus", 0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O ResetStat, */
/* 2*/ { "ModeCmd",    0,{0}, 1, NULL, eNetVarFromSGI, FALSE }, /* Mode, cmd */
/* 3*/ { "Status",      0,{0}, 1, NULL, eNetVarToSGI, FALSE }, /*O Status, */
30 /* 4*/ { "Position Cmd",0,{0}, 1, NULL, eNetVarFromSGI, FALSE },
/* 5*/ { "Pos status", 0,{0}, 1, NULL, eNetVarToSGI, FALSE },
/* 6*/ { "Home Status", 0,{0}, 1, NULL, eNetVarToSGI, FALSE },
/* 7*/ { "Step Count", 0,{0}, 1, NULL, eNetVarToSGI, FALSE },
/* 8*/ { "Travel Lim", 0,{0}, 1, NULL, eNetVarFromSGI, FALSE },
35 { NULL,          0,{0}, -1, NULL, FALSE, FALSE },
    } },
  { NULL,

```

```

    {
        { NULL,    0, -1, NULL, FALSE, FALSE },
    } },
};
5

10.  #define AF_OBJ_PWR_NEED_INTEN 10
    /*
    * stage autofocus
    * Laser Intensity/pmt gain setting are from configuration file
    * Low power objectives: added extra travel
15  *****
    */

extern int lonui_StageAF(void)
{
    UV_CONFIG    *cfg    = cfg_GetLockRdConfig();
20  OPTICS_PARMS *pOptic    = da_GlobalLonPointer();
    int          OldLaserInten, AFLaserInten, OldPMTampAZero, iValue, iSpeed;
    BOOL          rt = FALSE;
    float         fZPos, fZLimUm, fCt2Um, fSpeed;
    int           iObj, iSpdMult, iOldYAmpl, iIntenOffset, iTimeCt, CurrFineZ, iObjPwr;
25

    lon_APIQuery( eLonAPITurretPosCmd,  &iObj,      0);
    iObj --;

    fZLimUm    = cfg->StageSpec.ZLimitUm[iObj] - pOptic->SampleThick;
30  iObjPwr = cfg->ObjSpec.Power[iObj];

    CurrFineZ = HDWR_FASTZ_MID_POS;
    lonui_SetFastZPos( CurrFineZ);

35  if ( iObjPwr < AF_OBJ_PWR_NEED_INTEN &&
        cfg->StageSpec.StgAFLoMagAlgor == eStgAFLoMagMoveToLim )
    {

```

```

        fZPos = fZLimUm - 150.0;
        n = lonuiMoveStageZUm( fZPos, 100, eStgSwitchAxis);
    }
else
5   {
        fCt2Um      = cfg-> StageSpec.ZMStep2um;
        OldLaserInten = UIImage.LaserPwr;
        AFLaserInten = lonui_GetIntensityFromLineObjSelected( pOptic, iObj, cfg);

10   lon_APIQuery( eLonAPIPageYAmpl,      &iOldYAmpl, 0);
        lon_APICmd ( eLonAPIPageYAmpl,      5);

        OldPMTampAZero = -1;
        if( lon_APIQuery( eLonAPIPMTampAZeroCmd, &iValue, 0) == TRUE)
15   OldPMTampAZero = iValue ;

        lonuiAZFunction( cfg-> StageSpec.StgAF_Zero);
        lonui_LaserPwrAndPMT( AFLaserInten);

20   if( iObjPwr < AF_OBJ_PWR_NEED_INTEN &&
            cfg-> StageSpec.StgAFLoMagAlgor == eStgAFLoMagIntArray)
        lonuiLoPwrStgAF( fCt2Um, fZLimUm, AFLaserInten, cfg-> ObjSpec.NA[iObj]);
    else
    {
25   n = stg_ReadStageZUm( &fZPos, TRUE);
        if( n == TRUE && (int)fZPos > BACKUP_BEFORE_AF_UM )
            fZPos = lonuiMoveStgDownBeforeAF( fZPos, fZLimUm, iObjPwr,
                eStgUseCurrentAxis);

30   iSpeed = (int)( (double)fZLimUm - fZPos);
        if( iSpeed < 1000)
            mot_SetCurrAxisSpeed( (int)( (double)iSpeed / fCt2Um));

        if( n == TRUE)
35   n = lonimg_OneAFCycle( fZLimUm);
        fprintf( stderr, "1n=%d,", n );
        if( n == TRUE)

```

```

        r = lonuiStgAF2And3( OldLaserInten, AFLaserInten, fZLimUm, fCt2Um, iSpdMult,
iObjPwr );
    }
    lonui_LaserPwrAndPMT( OldLaserInten);
5    if( OldPMTAmpAZero > 0)
        lon_APICmd( eLonAPIPMTAmpAZeroCmd, OldPMTAmpAZero);

    mot_SetDefaultSpeed( eAxisZ, FALSE);
    if( r == TRUE && iObjPwr >= OBJ_PWR_FAF_AFTER_STGAF )
10    CurrFineZ = lonuiLocalFineAF();
        lon_APICmd ( eLonAPIPageYAmpl, iOldYAmpl);
    }
    mot_EnableCurrentAxisJoystick();
    return CurrFineZ;
15 }

/*****
* do stage AF at low power obj
* steps:
20 * 1: go to fix position, LO_MAG_AF_RANGE below limit
* 2. Move to limit while watching for intensity
* 3. based on intensity read, set new laser intensity and
*    calc first gauss target
* 4. start from LO_MAG_2NDPASS_RANGE above and go to
25 *    LO_MAG_2NDPASS_RANGE below, while watching
*    for intensity.
* 5. based on intensity read, calc and goto final target
* 4a. if intensity read from 1st search pass is all 0
*    the second search range is LO_MAG_LO_LITE_RANGE
30 *****/
*/

#define LO_MAG_AF_RANGE    1000.0
#define LO_MAG_LO_LITE_RANGE    400.0
#define LO_MAG_2NDPASS_RANGE    120.0
35 #define LO_MAG_2NDPASS_2RANGE ( 2 * LO_MAG_2NDPASS_RANGE)

static BOOL lonuiLoPwrStgAF( float fCt2Um, float fZLimUm, int LastLaserI, float fNA)

```

```

{
float      fZPos, fMedPosUm, fTarget, fDelta;
int        iMaxInten, iRange, iPeak;
double     dHalfRange2ndPass;
5  BOOL     rt, bPMTClamped = FALSE;

dHalfRange2ndPass = LO_MAG_2NDPASS_RANGE;
if( fNA > 0.01 && fNA < 0.15)
{
10  dHalfRange2ndPass = 2.0 / ( (double)fNA * fNA);
    if( dHalfRange2ndPass > (LO_MAG_AF_RANGE / 2) )
        dHalfRange2ndPass = (LO_MAG_AF_RANGE / 2);
    }
fZPos = fZLimUm - LO_MAG_AF_RANGE;
15  rt = lonuiMoveStageZUm( fZPos, 100, eStgSwitchAxis);
    if( rt == TRUE)
        rt = mot_SetCurrAxisSpeed( (int)( LO_MAG_AF_RANGE/ fC2Um));
    if( rt == TRUE)
        rt = lonuiMoveStgAndReadInten( fZPos, fZLimUm, &iMaxInten, &fMedPosUm, FALSE);
20  if( rt == TRUE)
    {
        lon_APIQuery( eLonAPIPMTampClamped, &bPMTClamped, 0);
        if( bPMTClamped == TRUE)
            iMaxInten = 250;
25
            iPeak = 2 * iMaxInten ;
            if( iPeak > 400)
                iPeak = 400;
            if( iPeak < 20)
30  iPeak += 20;
            else
                if( iPeak < 50)
                    iPeak += ( 20 * (50 - iPeak)/30);

35  zr_SetTargetLaserIFrlmgInten( iPeak, LastLaserI);
    fTarget = (double)fZPos + fMedPosUm;
    rt = stg_ReadStageZUm( &fZPos, FALSE);

```

```

        if( iMaxInten < 1 || bPMTClamped == TRUE)
            fTarget = fZLimUm - LO_MAG_LO_LITE_RANGE;
        }
        if( r == TRUE && fZPos > (double)fTarget)
5         {
            if( iMaxInten > 0 && bPMTClamped != TRUE)
                {
                    fDelta = (double)fZPos - fTarget;
                    if( fDelta > dHalfRange2ndPass)
10                     {
                        fZPos = fZPos + dHalfRange2ndPass - fDelta;
                        r = lonuiMoveStageZUm( fZPos, 100, eStgUseCurrentAxis);
                    }
                    fTarget = fZPos - ( dHalfRange2ndPass + dHalfRange2ndPass);
15                 }
                if(r == TRUE)
                    {
                        iRange = (int)fZPos - (int)fTarget;
                        mot_SetCurrAxisSpeed( (int)( (double)iRange / fCt2Um));
20                 fprintf( stderr, "fr= %d,", (int)fZPos);
                        r = lonuiMoveStgAndReadInten( fZPos, fTarget, &iMaxInten, &fMedPosUm, FALSE);
                    }
                if( r == TRUE)
                    {
25                 fTarget = (double)fZPos - fMedPosUm;
                        r = lonuiMoveStageZUm( fTarget, 100, eStgUseCurrentAxis);
                    fprintf( stderr, "Target= %d\n", (int)fTarget);
                    }
                }
30         return r;
        }

/*****
* move stage down: ensure start from below focus
35 *****/
*/

```



```

static float lonuiMoveStgDownBeforeAF( float fZPos, float fZLimUm, int iObjPwr,
                                         STG_ISSUE_SW_AXIS_CMD AxisCmd)
{
    int iTimeCt;

5
    if( (int)fZPos > (int)fZLimUm)
        fZPos = fZLimUm;
    if( iObjPwr <= AF_OBJ_PWR_NEED_INTEN )
        fZPos -= WORKING_DIST_10X_UM;
10
    else
        fZPos -= WORKING_DIST_UM;
    lonuiMoveStageZUm( fZPos, 100, AxisCmd);
    return fZPos;
}

15

/*****
 * do step 2 and 3 of Stage AF
 *****/

20
*/

static BOOL lonuiStgAF2And3( int OldLaserInten, int AFLaserInten, float fZLimUm, float
                                         fCzUm, int iSpdMult, int iObjPwr )
{
    BOOL      r = FALSE;
25
    float      fZPos, fCurrLim;
    int        iOffset, iSpeed, iDiffUm;
    UV_CONFIG  *pCfg = cfg_GetLockRdConfig();

    r = stg_ReadStageZUm( &fZPos, FALSE);
30
    if( r == TRUE)
    {
        if( fabs( (double)fZPos - fZLimUm) > 50)
            iOffset = -1200;
        else
35
            iOffset = -200;
    }
    if( r == TRUE)

```

```

    {
        rt = lonuiAutoFocusCycle23( &AFLaserInten, iOffset, fZLimUm, fCt2Um, iSpdMult,
fZPos);
        fprintf( stderr, "2rt=%d,", rt);
5         if( rt == TRUE)
            rt = stg_ReadStageZUm( &fZPos, FALSE);
            if( rt == TRUE)
                rt = lonuiAutoFocusCycle23( &AFLaserInten, 140, fZLimUm, fCt2Um, iSpdMult,
fZPos);
10         }
        fprintf( stderr, "3rt=%d,", rt);
        return rt;
    }

15     /*
        if ( iObjPwr < AF_OBJ_PWR_NEED_INTEN &&
            pCfg->StageSpec.StgAFLoMagAlgor == eStgAFLoMagIntArray)
        {
            fprintf( stderr, "Pos=%d.Inten=%d,", (int)fZPos, AFLaserInten );
20             iOffset = -200;
            fCurrLim = (double)fZPos + iOffset;
            if( fCurrLim < 0.0)
                fCurrLim = 0;
            if( fCurrLim > (double)fZLimUm)
25                 fCurrLim = fZLimUm ;
            iDiffUm = abs( (int)( (double)fCurrLim - fZPos) );
            iSpeed = iDiffUm;
            if( iSpeed < 50)
                iSpeed = 50;
30             rt = lonuiAFSetLaserIntenAndSpeed( &AFLaserInten, (int)( (double)iSpeed / fCt2Um) );
            if( rt == TRUE)
                rt = lonuiLoObjPwrStgAF2And3(OldLaserInten, fZPos, fCurrLim, FALSE);
        }
    */

35     /*
        *****
        * run the 2nd, 3rd cycle of the autofocus
    */

```

```

*****
*/
static BOOL lonuiAutoFocusCycle23( int *AFLaserInten, int Offset, float fZLim, float fC2Um,
                                   int iSpdMult, float fZPos )
5   {
    BOOL r = FALSE;
    float fZCurrLim;
    int iDiffUm, iSpeed;

10   {
        fZCurrLim = (double)fZPos + Offset;
        if( fZCurrLim < 0.0)
            fZCurrLim = 0;
        if( fZCurrLim > (double)fZLim)
15         fZCurrLim = fZLim;
        iDiffUm = abs( (int)( (double)fZCurrLim - fZPos) );
        iSpeed = iDiffUm / 2;
        /* 2second to travell */
        if( iSpeed < 25)
            iSpeed = 25;
20         r = lonuiAFSetLaserIntenAndSpeed( AFLaserInten, (int)( (double)iSpeed / fC2Um) );
        fprintf( stderr, "Pos=%d,Inten=%d,", (int)fZPos, *AFLaserInten );
    }
    if( r == TRUE)
        r = lonimg_OneAFCycle( fZCurrLim);
25     return r;
    }

/*****
* stage AF step 2 and 3 for low power obj
30 * Step 2: move Z to target position, and at the same time read
    * intensity. At the end of the movement, get the readings
    * and determine, based on intensity array, where the peak is
    * Peak is determined to be at the median point
    * Step 3: base on the step 2 intensity array, in-focus
35 * stage position is determined, and stage is moved there.
    * Before moving, laser intensity, stage speed are restored.
    * Restore laser intensity to prevent PMT overload, restore

```

```

    * speed to improve throughput.
    *
    *****
    */
5  static BOOL lonuiLoObjPwrStgAF2And3( int OldLaserInten, float fCurrZUm, float fZLimUm,
                                         BOOL          bWakeAxis)
    {
        BOOL    rt, bMotoring;
        int      iMaxInten, iCt;
10     float     fMedPosUm, fNewTarget;

        rt = lonuiMoveStgAndReadInten( fCurrZUm, fZLimUm, &iMaxInten, &fMedPosUm,
        bWakeAxis);
        lonui_LaserPwrAndPMT( OldLaserInten);
15     if( rt == TRUE )
        {
            fNewTarget = (double)fCurrZUm - fMedPosUm;
            if( fNewTarget < 0.0)
                fNewTarget = 0;
20     fprintf( stderr, "target=%d:%d, ", (int)fMedPosUm, (int)fNewTarget );
            mot_SetDefaultSpeed( eAxisZ, FALSE);
            rt = stg_MoveStageZUm( fNewTarget, FALSE);
            for( iCt= 0, bMotoring = TRUE; iCt < 20 && bMotoring == TRUE; iCt++ )
            {
25         bMotoring = mot_IsMotorMoving();
                if( bMotoring == TRUE)
                    util_sginap( 2);
            }
            if( bMotoring == TRUE)
30         {
                mot_StopCurrentAxis();
                rt = FALSE;
            }
        }
35     return rt;
    }

```

```

typedef unsigned int  UINT;
#define STG_AF23_ARRAYS 5
#define AF23_SCALE      80    /* suppose to read 80, z stops at more/less 80 */
/* *****
5   * move stage and watch for intensity
   * check for cases where sum = 1, or only one nonzero
   *      number
   * *****
   */
10  static BOOL lonuiMoveStgAndReadInten( float fCurrZUm, float fZLimUm, int *iMaxInten,
                                         float *fMedPosUm, BOOL bWakeAxis )
    {
        BOOL      n = FALSE, bMotoring;
        int       iValue, iCt, iArrayCt;
15  BYTE         chArray [STG_AF23_ARRAYS][LON_ARRAY_DATA_SZ];
        int       *iArray;
        int       *iFiltered;

        int       iSum, iMedian, iMedCt, iNumArrays, iMaxI, iArraySize;
20  float        fZScale, fOffset, fNewTarget;

        iArraySize = STG_AF23_ARRAYS * LON_ARRAY_DATA_SZ;
        iArray  = (int *)XtMalloc( sizeof(int) * iArraySize);
        iFiltered = (int *)XtMalloc( sizeof(int) * iArraySize);
25

        iNumArrays = STG_AF23_ARRAYS;
        longmi_UpdateNV( FAST_Z_NODE_ADDR, FASTZ_OPCMD_ARR2_NV_ADDR, 6);
        longmi_UpdateNV( FAST_Z_NODE_ADDR, FASTZ_OPCMD_ARR1_NV_ADDR,
LON_ARRAY_DATA_SZ * iNumArrays);
30  fOffset = 5;
        fZScale = ( (double)fCurrZUm - fZLimUm) / AF23_SCALE;
        fZScale = fabs( fZScale);
        fprintf( stderr, "Target= %d.", (int)fZLimUm );
        n = stg_MoveStageZUm(fZLimUm, bWakeAxis);
35  n = lonimg_QueryNetVar( FAST_Z_NODE_ADDR, FASTZ_MODE_RD_INTENSITY,
                           &iValue, TRUE);

```

```

    bMotoring = mot_IsMotorMoving();
    if( bMotoring == TRUE)
        mot_StopCurrentAxis();
    fprintf( stderr, "motoring=%d,", bMotoring );
5      iSum = 0;
      iMaxI = 0;
      *iMaxInten = iMaxI;

    for( iCt = 0; iCt < iNumArrays && pi == TRUE; iCt++)
10      {
        if( lonimg_QueryNetVar16( FAST_Z_NODE_ADDR, FASTZ_INTEN_ARRAY0 + iCt,
                                &iValue, (char *)chArray[ iCt], TRUE) !=
        LON_ARRAY_DATA_SZ)
            pi = FALSE;
15      for( iArrayCt= 0; iArrayCt < LON_ARRAY_DATA_SZ; iArrayCt++)
        {
            if( iMaxI < (UINT)chArray[iCt][iArrayCt])
                iMaxI = chArray[iCt][iArrayCt];
            iArray[ iCt * LON_ARRAY_DATA_SZ + iArrayCt] = chArray[iCt][iArrayCt];
20        }
    }
    if( pi == TRUE)
    {
        *iMaxInten = iMaxI;
25      if( iMaxI < 3)
        memcpy( iFiltered, iArray, sizeof(int) * iArraySize);
    else
    {
        iFiltered[0] = (3 * iArray[0] + iArray[1]) /4;
30      for( iCt = 1; iCt < (iArraySize -1); iCt++)
        iFiltered[iCt] = (iArray[iCt -1] + (2 * iArray[iCt]) + iArray[iCt +1]) /4;
        iFiltered[iArraySize-1] = (3 * iArray[iArraySize -1] + iArray[iArraySize -2] )/4;

        iMaxI = 0;
35      for( iCt = 0; iCt < iArraySize; iCt++)
        {
            if( iMaxI < iFiltered[iCt])

```

```

        iMaxI = iFiltered[iCt];
    }
    *iMaxInten = iMaxI;
}
5   for( iCt =0; iCt < iArraySize; iCt++)
    iSum += iFiltered[ iCt];

    if( iSum > 1)
        iSum = iSum /2;
10
    iMedian = 0;
    for( iCt = 0; iCt < iArraySize && iMedian < iSum; )
    {
        iMedian += iFiltered[ iCt];
15     if( iMedian < iSum)
        iCt ++;
    }
    *fMedPosUm = fOffset + ( (double)fZScale * iCt);
}
20
fprintf( stderr, "\nPos= %d,", iCt);
for( iCt = 0; iCt < iArraySize; iCt++)
{
    fprintf( stderr, "%d ", iFiltered[ iCt]);
25     if( (iCt %16) == 15)
        fprintf( stderr, "\n");
}
fprintf( stderr, "\nMaxI = %d,", iMaxI);

30   XtFree( (char *)iArray);
    XtFree( (char *)iFiltered);
    return ( rt);
}

```

35

```

/*****

```

```

    * support functions for Autofocus
    * set laser intensity, Z speed
    * after each trial intensity is reduced
    *****

5   */
static BOOL lonuiAFSetLaserIntenAndSpeed( int *AFLaserInten, int SpeedCt)
{
    BOOL      rt;
    int      iValue;

10   rt = lon_APIQuery( eLonAPIPMTAmpClamped, &iValue, 0);
    if( rt == TRUE && iValue == PMT_CLAMPED)
        *AFLaserInten = *AFLaserInten * AF2_INTEN_DEC;
    else
15   *AFLaserInten = *AFLaserInten * AF2_INTEN_DEC_NORMAL;
    rt = lonui_LaserPwrAndPMT( *AFLaserInten);

    if( rt == TRUE)
        mot_SetCurrAxisSpeed( SpeedCt);
20   return rt;
}

/*****
    * doing an autozero function is very slow, so just set it
    *****/

25   */
static BOOL lonuiAZFunction( int iAZero)
{
    lon_APICmd( eLonAPIPMTAmpAZeroCmd, iAZero);

30   return TRUE;
}

/*****
    * move Z and wait till done
35   *****/

    */

```



```

static BOOL lonuiMoveStageZUm( float fZPos, int iMaxCt, STG_ISSUE_SW_AXIS_CMD
AxisCmd)
{
    BOOL    r ;
5    int     iTimeCt;

    r = stg_MoveStageZUm( fZPos, (AxisCmd == eStgUseCurrentAxis)? FALSE: TRUE );
    if( r == TRUE)
    {
10        iTimeCt = 0;
        while( mot_IsMotorMoving() == TRUE && iTimeCt < iMaxCt)
        {
            iTimeCt++;
            util_sginap(3);
15        }
        if( iTimeCt >= iMaxCt)
        {
            r = FALSE;
            mot_StopCurrentAxis();
20        fprintf( stderr, "time = %d.", iTimeCt);
        }
    }
    return r;
}

25
/*****
* one cycle auto focus
* 1. reset AF flipflop,
* 2. move Z,
30 * 3. wait till Ztarget done or AF threshold exceeded
*****/
*/

extern BOOL lonimg_OneAFCycle( float fZLimUm)
{
35    BOOL          r = FALSE, bMotoring;
    HARDWR_IMAGE_ENTRY *FastZimg, *PMTimg;
    int            iValue, iTimeCt;

```

```

FastZImg = &( HardWrlmgTable[ FAST_Z_NODE_ADDR].NVImg[
FASTZ_OPSTATUS_NV_ADDR]);
PMTImg = &( HardWrlmgTable[ PMT_AM_NODE_ADDR].NVImg[
5 PMTAM_CLAMPED_NV_ADDR));
    r = lonimg_QueryNetVar( FAST_Z_NODE_ADDR,
FASTZ_MODE_WAITFOCUS_NV_ADDR, &iValue, TRUE);
    if( r == FALSE)
        fprintf( stderr, "cannot send WaitForFoc.");
10    else
        if( iValue != C_AF_COMPARATOR_RDY)
            fprintf( stderr, "WaitForFoc=%d,", iValue);

        if( r == TRUE && iValue == C_AF_COMPARATOR_RDY)
15        {
            fprintf( stderr, "Tar=%d,", (int)fZLimUm );

            FastZImg->UpdatedFlg = FALSE;
            PMTImg->UpdatedFlg = FALSE;
20            r = stg_MoveStageZUm(fZLimUm, FALSE);
            if( r == FALSE)
                fprintf( stderr, "StageMoveErr.");
                if( r == TRUE)
                    {
25                        iTimeCt = 0;
                        for( bMotoring= TRUE; bMotoring == TRUE && FastZImg->UpdatedFlg == FALSE
&& PMTImg->UpdatedFlg == FALSE && iTimeCt < 100; )
                            {
                                util_sginap(1); /* dont put this at the end of the loop, too slow */
30                                LonTimeBlockReadOneAndDecodeMsg(AF_READ_TICK);
                                if( PMTImg->UpdatedFlg == FALSE && FastZImg->UpdatedFlg == FALSE)
                                    bMotoring = mot_IsMotorMoving();
                                fprintf( stderr, "move=%d,", bMotoring );
                                    iTimeCt++;
35                                }
                                if( bMotoring == TRUE)
                                    mot_StopCurrentAxis();

```

```

        if( PMTImg->UpdatedFlg == FALSE &&
            ( FastZImg->UpdatedFlg == FALSE || FastZImg->iNVValue !=
C_AF_COMP_TRIP) )
            rt = FALSE;
5
        if( PMTImg->UpdatedFlg == TRUE )
            fprintf( stderr, "PMT-OL=%d,", PMTImg->iNVValue);

        if( FastZImg->UpdatedFlg == FALSE)
10        fprintf( stderr, "NoResponseFrFastZ.");
        else
            if( FastZImg->iNVValue != C_AF_COMP_TRIP)
                fprintf( stderr, "FastZResp=%d,", FastZImg->iNVValue );

15        fprintf( stderr, "Time=%d,", iTimeCt );
            }
            if( FastZImg->UpdatedFlg == FALSE || FastZImg->iNVValue != C_AF_COMP_TRIP )

                lon_APICmd( cLonAPIFastZAbort, 0);
20        }
        return rt;
    }

    /*****
        * calc new laser intensity based on image intensity
25        * feedback
        * image Intensity here is heavily filtered, and is
        * usually lower than actual peak
        * if saturated, cut laser intensity by 5
        * if less than ideal, incr of 10( maybe in config) laser Inten = 100% img Inten
30        * for incr of less the 10 laser inten, each step = 10%( linear, not expotential)
        * return TRUE if
        *****/

        */

#define IDEAL_IMG_INTEN          195
35 #define HALF_OF_IDEAL          (IDEAL_IMG_INTEN /2)
#define MAX_TARGET_LASER_INTEN  75

```

```

extern BOOL zr_SetTargetLaserIfImgInten( int iImgI, int CurrLaserI)
{
    UV_CONFIG    *pCfg    = cfg_GetLockRdConfig();
    float        fNew, fCt;
5    BOOL        n = FALSE;
    int          iNewI;

    if( iImgI > IDEAL_IMG_INTEN)
        fNew = CurrLaserI - (iImgI - IDEAL_IMG_INTEN) /10;
10    else
    {
        fNew = CurrLaserI;
        if( iImgI < 15)
            iImgI = 15;
15    for( ; iImgI < HALF_OF_IDEAL; iImgI += iImgI)
        fNew += (double)pCfg->SysSpec.ZRangeLaserGain;
        fCt = (double)pCfg->SysSpec.ZRangeLaserGain * (IDEAL_IMG_INTEN - iImgI) / iImgI;

        fNew += (double)fCt;
20    if( (int)fNew > MAX_TARGET_LASER_INTEN)
        fNew = MAX_TARGET_LASER_INTEN;
        n = TRUE;
    }
    iNewI = fNew;
25    if(iNewI < 0)
        iNewI = 0;
    lomui_LaserPwrAndPMT( iNewI);
    fprintf( stderr, "newLI= %d, ", iNewI);
    return TRUE;
30 }

/*****
 * get laser intensity setting based on current laser line selected
 *****/
*/
35 static int lomui_GetIntensityFromLineObjSelected( OPTICS_PARMS *pOptic,
                                                    int iObj, UV_CONFIG *pCfg)
{

```

```

int      iIntenOffset, *pIOffset, iInten;

switch( UIImage.eLaserLine)
{
5   case eSrc458: iInten = pOptic->AFInten458; break;
    case eSrc488: iInten = pOptic->AFInten488; break;
    case eSrc515: iInten = pOptic->AFInten515; break;
    case eSrcWhite: iInten = pOptic->AFIntenAll; break;

10  default:      iInten = pOptic->AFInten488;
                    break;
}

pIOffset = Ionui_GetLaserIOffsetForLaserLine( &pCfg->ObjSpec, UIImage.eLaserLine);
iIntenOffset = pIOffset[iObj] - pIOffset[0];
15  iInten = iInten + iIntenOffset;
    if( iInten < 0)
        iInten = 0;
    if( iInten > 100)
        iInten = 100;
20  return iInten;
}

/*****
 * get the Laser I offset array based on laserLine specified
 *****/
25  */

extern int * Ionui_GetLaserIOffsetForLaserLine( UV_OBJ_SPEC *ObjSpec,
OPTICS_LASER_SRC LineFilter)
{
    int * pOffset;
30
    switch ( LineFilter)
    {
        case eSrc458:
            pOffset = ObjSpec->LaserI458Offset;
35            break;
        case eSrc488:
            pOffset = ObjSpec->LaserI488Offset;

```

```

        break ;
    case eSrc515:
        pOffset = ObjSpec->LaserI515Offset;
        break;
5   case eSrcWhite:
    default:
        pOffset = ObjSpec->LaserLAllOffset;
        break ;
    }
10  return pOffset;
    }

/*****
 * change the laser intensity ( PMT gain and laser power) only.
 * no laser line change
15  *
 * Specified LaserInten = 0-100
 * laser Power is set according to config, based on current laserline used
 *
 * if( LaserInten > LASER_I_PERCENT_CHG_PWR, laser power is increased upto
20  MAX_LASER_PWR
 *****/

*/
extern BOOL lonui_LaserPwrAndPMT( int LaserInten)
{
25  BOOL    n = FALSE;
    int     iValue, LaserPwrRaw, PMT, MinPMTGain, PMTRange;
    UV_CONFIG *cfg = cfg_GetLockRdConfig();
    OPTICS_LASER_SRC LaserLine;

30
    if( LaserInten > 100)
        LaserInten = 100;
    if( LaserInten < 0)
        LaserInten = 0;
35
    lon_APIQuery( eLonAPILaserLineFilterMove, (int *)&LaserLine, 0);
    switch( LaserLine)

```

```

    {
        case eSrc458: LaserPwrRaw = cfg->ConfocalSpec.LaserPwr458;
                        break;
        case eSrc488: LaserPwrRaw = cfg->ConfocalSpec.LaserPwr488;
5                        break;
        case eSrc515: LaserPwrRaw = cfg->ConfocalSpec.LaserPwr515;
                        break;
        case eSrcWhite:
        default:      LaserPwrRaw = cfg->ConfocalSpec.LaserPwrAll;
10                        break;
    }

    if( LaserInten > LASER_I_PERCENT_CHG_PWR && LaserPwrRaw < MAX_LASER_PWR
    )
    {
15        LaserPwrRaw += (double)( LaserInten - LASER_I_PERCENT_CHG_PWR) *
                        ( MAX_LASER_PWR - LaserPwrRaw) / ( 100 -
LASER_I_PERCENT_CHG_PWR );
    }

    n = lon_APIQuery( eLonAPILaserOnOff, &iValue, 0);
20    if( n == TRUE && iValue == HDWR_LASER_OFF)
        n = lon_APICmd( eLonAPILaserOnOff, HDWR_LASER_ON);
    if( n == TRUE)
        n = lon_APIQuery( eLonAPILaserRunIdle, &iValue, 0);
    if( n == TRUE && iValue == HDWR_LASER_IDLE)
25        n = lon_APICmd( eLonAPILaserRunIdle, HDWR_LASER_RUN);
    if( n == TRUE)
    {
        MinPMTGain = MIN_PMT_GAIN + cfg->ConfocalSpec.PMTMinGainOffset;
        PMT = MinPMTGain - ((LaserInten * cfg->ConfocalSpec.PMT_Range)/
30        LASER_PWR_SLIDER_RANGE);

        if( PMT < 50)      /* safety */
            PMT = 50;
        if( PMT > 255)
35            PMT = 255;
        if( LaserPwrRaw > MAX_LASER_PWR)
            LaserPwrRaw = MAX_LASER_PWR;
    }

```

```
    if( LaserPwrRaw < 0)
        LaserPwrRaw = 0;
    rt = lon_APICmd( eLonAPILaserPwr, LaserPwrRaw);
    if( rt == TRUE)
5      rt = lon_APICmd( eLonAPIPMTAmpVoltCmd, PMT);

    UIImage.LaserPwr    = LaserInten;
    }
    return rt;
10  }
    □
```



```

/*****
 * Ultrapointe Corp. copyright 1992,93,94
 * stage.c
 *
5  * UltraView1000 software
 * this module handles the high level stage interface:
 * stage coordinates are in microns, with XYoffset from different Objective
 * reference is based on 2nd obj: Obj1, ie: no correction is zero for obj1
 * If currently using other objective, the XOffsetCt/YOffsetCt from .ucvconfig
10 * is used to compensate
 *
 *****/

 */

#ifdef Irix5plus
15 #define NeedFunctionPrototypes 1
#else
#define FUNCPROTO
#endif

20 #include <stdio.h>
#include <stdlib.h>
#include <sys/ipc.h>
#include <sys/sem.h>
#include <sys/signal.h>
25 #include <math.h>
#include <unistd.h>
#include <limits.h>          /* sginap */
#include <X11/Intrinsic.h>
#include "uv_util.h"
30 #include "uvconfig.h"
#include "wfrfile.h"
#include "uv_data.h"
#include "motor.h"
#include "stage.h"
35 #include "deskew.h"
#include "stgchild.h"
#include "lon_msg.h"

```

```

#include "lon_img.h"
#include "lonuiimg.h"
#include "lon_gmi.h"

5  #define RD_DELAY_TIME          60    /* in ms */
    #define RD_DELAY_TICK          (RD_DELAY_TIME* CLK_TCK / 1000)
    #define HOME_XY_TIMEOUT_CT     10

    #define MAX_NO_MOVE_UM         1.0
10  #define MAX_NO_MOVE_ECT         2
    #define NUM_CLOSED_LOOPS       6
    #define MAX_ERROR_LIM_UM       2000.0

    #define OBJ_MOVE_CLOSE_ROBOT   0    /* first object as reference */
15  #define MIN_XY_SPEED_UM         100

    #define SENDER 1               /* direction for pipe */
    #define RECEIVER 0
    #define FILE_PTR int           /* file table index */
20

    #define BACKLASH_UM            40  /* amount of xy backlash */
                                     /* see EndPtPosCorrection() */

    #define INIT_BACKLASH_UP_UM    3000 /* back-and-fore stage init */
25  #define INIT_BACKLASH_DN_UM    1000 /* back-and-fore stage init */

    #define MAX_REL_MOVE           25000 /* max rel move is 20000um each time */
    #define STG_POS_END_POINT      120000.0 /* ends at +/- 120mm */
    #define STG_NEG_END_POINT      -120000.0
30  #define STG_XY_TRAVEL_LIM       110000 /* travel limit, with offset */
    #define STG_SQ_XY_TRAVEL_LIM   (110000.0 * 110000.0)

    #define STG_Z_TRAVEL_LIM        12000 /* um */

35  static STAGE_STATUS StageInfo;
    static BOOL stgBackgroundProc(STAGE_STATUS *pStgInfo, int iObj);
    static BOOL stgEndPtPosCorrection( MOT_AXIS Axis, double TargetUM, int iObj);

```

```

float EResol, float MResol);

static BOOL stg_MoveStageTillDone( STAGE_XY *StageXY, STG_MOVE_MODE MoveMode,
int iObj,

int XOffset, int YOffset);

5 static int stgUM2MotorCount( MOT_AXIS Axis, double um, int iObj);
static double stgEncoderCount2UM( MOT_AXIS Axis, double dCounts, int iObj);
static int stgUM2EncoderCount( MOT_AXIS Axis, double um, int iObj);
static int stgObjMOffsetCt( MOT_AXIS Axis, UV_STAGE_SPEC *StageSpec, int iObj);
static int stgObjEOffsetCt( MOT_AXIS Axis, UV_STAGE_SPEC *StageSpec, int iObj);
10 static BOOL stg_RdStgXYPosAndMotorEncoderOffset( STAGE_XY *StageXY,
STG_MOVE_MODE MoveMode,
int iObj,
int *XMotorCtOffset, int *YMotorCtOffset);
static BOOL stgSafeMoveXYStageRawAndNormal( STAGE_XY *StageXY, int iObj,
15 STG_MOVE_F_MODE eFMode, STG_MOVE_MODE MoveMode,
STG_MOVE_JOY_STATUS eMoveJoy );
static BOOL stgOrthoCorrAndDskToAppl( STAGE_XY *DestXY, STAGE_XY *SrcXY,
STG_MOVE_MODE MoveMode);
static BOOL stgOrthoCorrAndDskToStg ( STAGE_XY *DestXY, STAGE_XY *SrcXY,
20 STG_MOVE_MODE MoveMode);
static double stgIdealUMToMappedUM( MOT_AXIS Axis, UV_CONFIG *pCfg, double
dIdealUM);

static BOOL DebugPrtCh( char ch);
25

/*****
* stg_initialize( void)
* initialize XYZ stage
30 *****/
*/
extern BOOL stg_initialize( int hUldd)
{
UV_STAGE_SPEC StageSpec;
35
StageInfo.StgStatus = eStgIdle;
StageInfo.pidStage = -1;

```

```

    cfg_CopyStageConfig( &StageSpec);
    mot_Initialize( hUldd);
    dsk_ResetDeskewMatrix();
    StageInfo.CurrentXY.x = 0;
5   StageInfo.CurrentXY.y = 0;
    StageInfo.TargetXY.x = 0;
    StageInfo.TargetXY.y = 0;
    StageInfo.BacklashXY.x = 0;
    StageInfo.BacklashXY.y = 0;
10  StageInfo.CurrentZ      = 0;
    StageInfo.TargetZ      = 0;
    StageInfo.BacklashZ     = 0;
    StageInfo.bStage       = FALSE;
    StageInfo.XMResol       = StageSpec.XMStep2um;
15  StageInfo.YMResol       = StageSpec.YMStep2um;
    StageInfo.ZMResol       = StageSpec.ZMStep2um;
    StageInfo.XEResol       = StageSpec.XEStep2um;
    StageInfo.YEResol       = StageSpec.YEStep2um;
    StageInfo.ZEResol       = StageSpec.ZEStep2um;
20  StageInfo.BacklashCts = (double)BACKLASH_UM /StageSpec.XMStep2um;
    return TRUE;
}

/*****
25  * cleanup before program ends
*****/

*/

extern BOOL stg_TerminateStage( void)
{
30  mot_ClosePort();
    return TRUE;
}

extern BOOL stg_IsStageOnline( void)
35  {
    BOOL n = FALSE;
    n = StageInfo.bStage;

```

```

        return (rt);
    }

    /*****
5   * return TRUE if stage is put online
    * turret is also initialized here: home turret and then turn to low mag
    * set global update message
    *
    *****/

10  */

extern BOOL stg_PutStageOnline( void)
{
    int Pos, j, iJGain;
    UV_CONFIG    *pCfg    = cfg_GetLockRdConfig();
15  BOOL    rt;
    BYTE ch;

    if( StageInfo.StgStatus == eStgIdle )
    {
20      StageInfo.bStage = mot_PutStageOnline();

        rt = StageInfo.bStage;
        if( rt == TRUE)
            rt = mot_MoveToHome( eAxisZ);
25      longmi_UpdateNV( TURRET_NODE_ADDR, TURRET_POS_CMD_NV_ADDR,
        HDWR_TURRET_HOME_OBJ);
        longmi_UpdateNV( TURRET_NODE_ADDR, TURRET_POS_CMD_NV_ADDR,
        HDWR_TURRET_LO_OBJ);
        lonmi_SetLaserAzenToCurrObjAndLine();
30

        /* no gap!!! */
        if( rt == TRUE)
        {
            for( j = 0, rt = FALSE; rt == FALSE && j < HOME_XY_TIMEOUT_CT; j++)
35          {
                rt = mot_CurrAxisGetOneChar( &ch);
                if( rt == FALSE || (int)ch != 'F')

```

```

        rt = FALSE;
    }
    StageInfo.bStage = rt;
    if( j >= HOME_XY_TIMEOUT_CT )
5   fprintf( stderr, "zStillNotHome.");
    }
    if( rt == TRUE )
    {
        j = HDWR_TURRET_LO_OBJ -1;
10   if( j < 0 )
        j = 0;
        iJGain = pCfg->StageSpec.JoyStickGain[j];
        mot_SetStgJoyStickGain( eAxisX, iJGain);
        mot_SetStgJoyStickGain( eAxisY, iJGain);
15   iJGain = pCfg->StageSpec.JoyStickZGain[j];
        mot_SetStgJoyStickGain( eAxisZ, iJGain);
        stg_SetZJoystickLimitToCurrParms();          /* cfg/laserParms/obj set already */
    }
    dat_SendInfo_ObjChanged();
20
    }
    rt = StageInfo.bStage;
    return (rt );
}

25  /******
    * Put Stage Offline
    * return TRUE if online
    *   FALSE if offline
    *****/

30  */
extern BOOL stg_PutStageOffline( void)
{
    BOOL    rt = FALSE;

35   if( StageInfo.bStage == TRUE)
    {
        if( StageInfo.StgStatus == eStgIdle )

```

```

    {
        rt = mot_PutStageOffline();
        if( rt == TRUE)
            StageInfo.bStage = FALSE;
5      }
    }

    rt = StageInfo.bStage;
    return ( rt);
}

10

/*.....
 * home y stage and wait till done
 *.....
 */

15 extern BOOL stg_HomeXYAndWaitTillDone( void)
    {
        BOOL    rt, bMoving;
        BYTE     ch;
        int      iRetry;

20
        rt = mot_MoveToHome( eAxisX);
        if( rt == TRUE)
        {
            rt = FALSE;
25         for( iRetry = 0; rt == FALSE && iRetry < HOME_XY_TIMEOUT_CT; iRetry ++ )
            {
                rt = mot_CurrAxisGetOneChar( &ch);
                if( rt == FALSE || (int)ch != 'F')
                    rt = FALSE;
30            }
        }
        if( rt == TRUE )                /* wait till X done */
        {
            rt = mot_MoveToHome( eAxisY);
            if( rt == TRUE)
35            {
                stgWaitTillCurrAxisStop( 1);
                rt = mot_EnableJoystick();
            }
        }
    }

```

```

    }
    }
}

5   fprintf( stderr, "HomeCt: %d= %d.", iRetry, n );

    return n;
}

/*****
10  * return the current Z position in um
    * Send out WakeUnit if bAxisCmd == TRUE
    *****/

    */

extern BOOL stg_ReadStageZUm( float *Um, BOOL bAxisCmd)

15  {
    BOOL    n = FALSE;
    int     Pos;
    if( StageInfo.StgStatus == eStgIdle )
    {
20      if( bAxisCmd == TRUE)
          n = mot_ReportOnePos( eAxisZ, &Pos);
        else
          n = mot_CurrAxisReportMotorCt( &Pos);
        if( n == TRUE)
25          {
              *Um = Pos * StageInfo.ZMResol;
              if( *Um < 0.0 || (double)*Um > STG_Z_TRAVEL_LIM)
                  n = FALSE;
          }
30      }
    return n;
}

/*****
    * move Z stage to specified Location in um
35  * Send out WakeUnit if bAxisCmd == TRUE
    *****/

    */

```



```

extern BOOL stg_MoveStageZUM( float Um, BOOL bAxisCmd)
{
    BOOL    n = FALSE;
    int     Pos;
5   if( StageInfo.StgStatus == eStgIdle )
    {
        Pos = (double)Um / StageInfo.ZMResol;
        if( bAxisCmd == TRUE)
            n = mot_AbsMove( eAxisZ, Pos);
10        else
            n = mot_CurrAxisAbsMove( Pos);
    }
    return n;
}

15
/*****
 * read real current stage XY position
 * return TRUE if both XY valid
 *****/

20 */
extern BOOL stg_ReadStageXYPosition( STAGE_XY *StageXY, STG_MOVE_MODE
MoveMode, int iObj)
{
    BOOL    n = FALSE;
25    double  dEnc;
    int     RawEnc;
    STAGE_XY StageLocoXY;

    if( StageInfo.StgStatus == eStgIdle )
30    {
        n = mot_ReportOneEnc( eAxisX, &dEnc, &RawEnc);
        if( n == TRUE)
        {
            StageLocoXY.x = stgEncoderCount2UM( eAxisX, dEnc, iObj);
35            n = mot_ReportOneEnc( eAxisY, &dEnc, &RawEnc);
            if( n == TRUE)
            {

```

```

        StageLocoXY.y = stgEncoderCount2UM( eAxisY, dEnc, iObj);
        stgOrthoCorrAndDskToAppl( StageXY, &StageLocoXY, MoveMode );
    }
}
5    }
    return ( r);
}

/*****
* read real current stage XY position, and motor ct/encoder ct offset
10 * return TRUE if both XY valid
* offsets work this way:
*  ————0—————40000
*  —————————20030
*
15 * Moffset = -60 motor ct, assume 1 encoder ct = 2 motor cts.
* to go to target encoder ct = 20000, initial motor ct = 40000
* real motor ct to get to target = 40000 - 60
*
*****/
20 */
static BOOL stg_RdStgXYPosAndMotorEncoderOffset( STAGE_XY *StageXY,
STG_MOVE_MODE MoveMode,
                                int iObj,
                                int *XMotorCtOffset, int *YMotorCtOffset)
25 {
    BOOL    r = FALSE;
    double dEnc;
    int     RawEnc, MotCt, ConvertedMCt;
    STAGE_XY StageLocoXY;
30
    if( StageInfo.StgStatus == eStgIdle )
    {
        r = mot_ReportOneEnc( eAxisX, &dEnc, &RawEnc);
        if( r == TRUE)
35         r = mot_CurrAxisReportMotorCt( &MotCt);
        if( r == TRUE)
        {

```

```

        ConvertedMCt = RawEnc * StageInfo.XEResol/ StageInfo.XMResol;
        *XMotorCtOffset = MotCt - ConvertedMCt;
        StageLocoXY.x = stgEncoderCount2UM( eAxisX, dEnc, iObj);
        rt = mot_ReportOneEnc( eAxisY, &dEnc, &RawEnc);
5      }
      if( rt == TRUE)
        rt = mot_CurrAxisReportMotorCt( &MotCt);
      if( rt == TRUE)
      {
10      ConvertedMCt = RawEnc * StageInfo.YEResol/ StageInfo.YMResol;
        *YMotorCtOffset = MotCt - ConvertedMCt;
        StageLocoXY.y = stgEncoderCount2UM( eAxisY, dEnc, iObj);
        stgOrthoCorrAndDskToAppl( StageXY, &StageLocoXY, MoveMode );
      }
15    }
    return ( rt);
  }

  /*****
20  * Move stage with endpoint correction, or stopped by other cmd
  *****/

  /*
  static BOOL stg_MoveStageTillDone( STAGE_XY *StageXY, STG_MOVE_MODE MoveMode,
  int iObj,
25      int XMOffset, int YMOffset)
  {
    int    x, y;
    BOOL   rt= FALSE;
    int    PIDStat, Enc;
30    STAGE_STATUS *pStgInfo = &StageInfo;

    fprintf( stderr, "MoveTillDone=");

    if( StageInfo.bStage == TRUE)
35    {
      stgOrthoCorrAndDskToStg( &(StageInfo.TargetXY), StageXY, MoveMode );
      x = stgUM2MotorCount( eAxisX, StageInfo.TargetXY.x, iObj);

```

```

        y = stgUM2MotorCount( eAxisY, StageInfo.TargetXY.y, iObj);

        rt = mot_AbsXYMove( x + XOffset, y + YOffset);
        if( rt == TRUE)
5           rt = stgBackgroundProc( pStgInfo, iObj);
    }

    fprintf( stderr, "%d.", rt);
    return rt;
}

10
/*****
 * background process:
 * do backlash/inaccuracy correction
 *
15 * do just one command: wait till xy move done
 * !! assumes axis currently selected is Y
 *
 *****/

20 static BOOL stgBackgroundProc(STAGE_STATUS *pStgInfo, int iObj)
{
    char chBuf[30];
    BOOL rt = FALSE, bMoving;
    int iRetry;

25
    for( iRetry = 0; rt == FALSE && iRetry < 4; iRetry ++ )
    {
        rt = mot_CurrAxisGetOneChar( (BYTE *)chBuf);
        if( rt == TRUE)
30         if( (int)chBuf[0] != 'F')
            {
                rt = FALSE;
                DebugPrtCh( chBuf[0] );
            }
35     }

    if( rt == TRUE && (int)chBuf[0] == 'F') /* wait till Y done */
    {

```

```

    n = mot_SwitchAxis( eAxisX);
    stgWaitTillCurrAxisStop( 2);

    if( n == TRUE)
5      n = stgEndPtPosCorrection( eAxisX, pStgInfo->TargetXY.x, iObj,
                                   StageInfo.XEResol, StageInfo.XMResol);

    if( n == TRUE)
      n = mot_SwitchAxis( eAxisY);
10

    if( n == TRUE)
      n = stgEndPtPosCorrection( eAxisY, pStgInfo->TargetXY.y, iObj,
                                   StageInfo.YEResol, StageInfo.YMResol);
    }
15    return n;
    }

/*****
 * Assume the last move had moved stage to 2 * backlash pass
20 * the desirable position.
 * first move back to just one time backlash to fix the backlash
 * problem. Then move an amount base on encoder count to fix
 * the inaccuracy problem.
 *****/
25 */

static BOOL stgEndPtPosCorrection(MOT_AXIS Axis, double TargetUM, int iObj,
                                   float EResol, float MResol)
{
30   BOOL      n = FALSE, bDone;
   int      TargetMotCt, RawEnc, iLoopCt;
   double    dEnc, dTargetEnc, dDeltaEnc1;

   dTargetEnc = stgUM2EncoderCount( Axis, TargetUM, iObj);
35   bDone = FALSE;
   n = TRUE;

```

```

    util_sginap( 4);
    for( iLoopCt = 0; iLoopCt < NUM_CLOSED_LOOPS && bDone == FALSE && r ==
    TRUE; iLoopCt++)
    {
5      r = mot_CurrAxisReportEnc( &dEnc, &RawEnc);
      if( r == TRUE)
      {
          dDeltaEnc1 = dTargetEnc - dEnc;
          if( abs( dDeltaEnc1) > MAX_NO_MOVE_ECT)
10          {
              TargetMotCt = dDeltaEnc1 * EResol / MResol;
              r = mot_CurrAxisRelMove( TargetMotCt);
              if( r == TRUE)
              {
15                  stgWaitTillCurrAxisStop( 3);
                  util_sginap( 4);
              }
          }
          else
20          bDone = TRUE;
      }
    }
    return r;
}

25 /*****
    * MoveXYStage
    * send move command to xy-stage
    * but first call MoveZStage to move z-stage away from objective
    *****/
30 */
extern BOOL stg_SafeMoveRawXYStage( STAGE_XY *StageXY, int iObj,
    STG_MOVE_F_MODE eFMode)
    {
        BOOL    r = FALSE;
35
        r = stgSafeMoveXYStageRawAndNormal( StageXY, iObj, eFMode, cStgRawPos,
        cMoveEnableJoystick);
    }

```

```

        return rt;
    }
    /*****
    * MoveXYStage
    5  * send move command to xy-stage
    * but first call MoveZStage to move z-stage away from objective
    *****/
    */
extern BOOL stg_SafeMoveXYStage( STAGE_XY *StageXY, int iObj, STG_MOVE_F_MODE
10 eFMode )
    {
        BOOL rt = FALSE;

        fprintf( stderr, "StgMove,");
    15     rt = stgSafeMoveXYStageRawAndNormal( StageXY, iObj, eFMode, eStgNormalPos,
        eMoveEnableJoystick);
        return rt;
    }

    20 #define WAVE_LEN_UM      .45
    /*****
    * Safe move stage, with AutoFine Focus, normal/raw moves
    * if AF is specified: assume stage tilt is in .uvconfig
    *
    25 *          check depth-of-focus, do AF if XY movement exceed d-of-foc
    *****/
    */
static BOOL stgSafeMoveXYStageRawAndNormal( STAGE_XY *StageXY, int iObj,
        STG_MOVE_F_MODE eFMode, STG_MOVE_MODE MoveMode,
    30 STG_MOVE_JOY_STATUS eMoveJoy )
    {
        BOOL      rt;
        STAGE_XY  CurrXY;
        double    deltaX, deltaY;
    35     int      XMotorCtOffset, YMotorCtOffset;
        float     fNA2, fDepthFoc, fTilt;
        UV_CONFIG *cfg = cfg_GetLockRdConfig();
    
```

```

    n = stg_RdStgXYPosAndMotorEncoderOffset( &CurrXY, eStgRawPos, iObj,
                                              &XMotorCrOffset, &YMotorCrOffset);

    if( n == TRUE)
    {
5      deltaX = CurrXY.x - StageXY->x;
      deltaY = CurrXY.y - StageXY->y;
      if( ( fabs(deltaX) <= MAX_NO_MOVE_UM) && ( fabs(deltaY) <=
MAX_NO_MOVE_UM) )
        n = TRUE;
10     else
        {
            n = stg_MoveStageTillDone( StageXY, MoveMode, iObj, XMotorCrOffset,
YMotorCrOffset);
            if( n == TRUE)
15             {
                switch( eFMode )
                {
                    case eStgCoarseAF:
                        lomui_AFAfterXYMove( elomui_CoarseAF);
20                        dat_SendInfoToWindows( eFineZPosChanged);
                        break;
                    case eStgAutoFocus:
/*
                        fTilt = fTilt * cfg->StageSpec.LevelXOffsetZUm/200000 +
25                        (float)deltaY * cfg->StageSpec.LevelYOffsetZUm/200000 ;
*/
                        fTilt = sqrtf( ((double)deltaX * deltaX) + ((double)deltaY * deltaY) );
                        fTilt = (double)fTilt * cfg->StageSpec.LevelXOffsetZUm /200000;
                        fTilt = fabs(fTilt);
30                        fNA2 = 1;
                        if( cfg->ObjSpec.NA[iObj] != 0.0)
                            fNA2 = (double)cfg->ObjSpec.NA[iObj] * cfg->ObjSpec.NA[iObj];
                        fDepthFoc = (double)WAVE_LEN_UM / fNA2 /4;
                        if( (double)fTilt > fDepthFoc )
35                         {
                            lomui_AFAfterXYMove( elomui_FineAF);
                            dat_SendInfoToWindows( eFineZPosChanged);

```



```

        }
        break;
    }
    if( eMoveJoy != eMoveDisableJoystick)
5      mot_EnableXYJoystick();
    }
}

return rt;
10 }

/*.....
 * Turn on JoyStick
 *.....
*/

15 extern BOOL stg_EnableJoystick( void)
{
    return (mot_EnableJoystick() );
}

/*.....
20 * stg_relmoveXY() relatively move XY in small distance
 * distance clipped at 100um
 * otherwise no safety check
 * Used for interactive mouse/traceball moves
 *.....
25 */

extern BOOL stg_RelMoveXY( int Xum, int Yum, int iObj)
{
    STAGE_XY    CurrXY;
    BOOL        rt;
30 double      ChkLimX, ChkLimY;
    int         XStageDir, YStageDir;

    XStageDir = X_STAGE_DIR;
    YStageDir = Y_STAGE_DIR;
35

    rt = stg_ReadStageXYPosition( &CurrXY, eStgRawPos, iObj);
    if( rt == TRUE)

```

```

    {
        ChkLimX = CurrXY.x + Xum;
        ChkLimY = CurrXY.y + Yum;
        if( (ChkLimX * ChkLimX + ChkLimY * ChkLimY) > STG_SQ_XY_TRAVEL_LIM)
5         rt = FALSE;
        else
        {
            if( XStageDir < 0)
                Xum = -Xum;
10         if( YStageDir < 0)
                Yum = -Yum;

            rt = mot_RelXYMove( (double)Xum /StageInfo.XMResol,
                                (double)Yum /StageInfo.YMResol );

15         if( rt == TRUE)
            {
                stgWaitTillCurrAxisStop( 5);
                rt = mot_SwitchAxis( eAxisX);
                if( rt == TRUE)
20                 stgWaitTillCurrAxisStop( 6);
            }
        }
    }

    return rt;
25 }

/*
* relative move on current axis
*/

30 extern BOOL stg_RelMoveUm( MOT_AXIS Axis, int iUm, int iObj)
    {
        double    dEnc;
        int        RawEnc, iCurrUm;
        BOOL        rt = FALSE;
35         int        XStageDir, YStageDir;

        XStageDir = X_STAGE_DIR;

```

```

YStageDir = Y_STAGE_DIR;

rt = mot_ReportOneEnc( Axis, &dEnc, &RawEnc);
if( rt == TRUE )
5   {
      iCurrUm = stgEncoderCount2UM( Axis, dEnc, iObj);
      if( abs( iCurrUm + iUm) > STG_XY_TRAVEL_LIM)
          rt = FALSE;
      else
10   {
          if( Axis == eAxisX)
              {
                  if( XStageDir < 0)
                      iUm = -iUm;
15   rt = mot_CurrAxisRelMove( (double)iUm /StageInfo.XMResol);
              }
          else
              if( Axis == eAxisY)
                  {
20   if( YStageDir < 0)
                      iUm = -iUm;
                      rt = mot_CurrAxisRelMove( (double)iUm /StageInfo.YMResol);
                  }
          else
25   rt = FALSE;
              }
          }
      if( rt == TRUE)
          stgWaitTillCurrAxisStop( 7);
30  return rt;
    }

35  /*****
      * set XY axis to default speed
      *****/

```

```

    */
extern BOOL stg_SetXYToDefaultSpeeds( void)
{
    return (mot_SetXYToDefaultSpeeds() );
5   }

/*****
 * set XY to speed specified
 *****/

10  */
extern BOOL stg_SetXYSpeedsUm( int iInitSpdUm, int iFinalSpdUm)
{
    BOOL      rt ;
    int      iInitSpdCt, iFinalSpdCt;

15  if( iFinalSpdUm < MIN_XY_SPEED_UM)
        iFinalSpdUm = MIN_XY_SPEED_UM;
    if( iInitSpdUm < 0)
        iInitSpdUm = 0;

20  iInitSpdCt = iInitSpdUm / StageInfo.XMResol;
    iFinalSpdCt = iFinalSpdUm / StageInfo.XMResol;

    rt = mot_SetInitAndFinalSpeedsCt( eAxisX, iInitSpdCt, iFinalSpdCt);
25  if( rt == TRUE)
    {
        iInitSpdCt = iInitSpdUm / StageInfo.YMResol;
        iFinalSpdCt = iFinalSpdUm / StageInfo.YMResol;
        rt = mot_SetInitAndFinalSpeedsCt( eAxisY, iInitSpdCt, iFinalSpdCt);
30  }
    return rt;
}

/*****
35  * convert um to MOTOR count
 *****/

*/

```

```

static int stgUM2MotorCount( MOT_AXIS Axis, double um, int iObj)
{
    int      Counts, iObjOffset;
    double    fOffset, fRes;
5    UV_CONFIG    *cfg = cfg_GetLockRdConfig();
    UV_STAGE_SPEC *StageSpec;
    GLOBAL_DATA    *pUV = dat_ObtainGlobDataPtr();
    int      XStageDir, YStageDir, iMapCt;

10    XStageDir = X_STAGE_DIR;
    YStageDir = Y_STAGE_DIR;

    StageSpec = &(cfg->StageSpec);
    if( um > STG_XY_TRAVEL_LIM)
15    um = STG_XY_TRAVEL_LIM;
    else
        if( um < -STG_XY_TRAVEL_LIM )
            um = -STG_XY_TRAVEL_LIM;

20    if( Axis == cAxisY)
    {
        fOffset = YMOT_OFFSET + pUV->UVSysParms.YCenterOffsetUm;
        fRes    = StageInfo.YMResol;
        um = stgIdealUMToMappedUM( Axis, cfg, um);
25    if( YStageDir < 0)
        um = -um;
    }
    else
    {
30    fOffset = XMOT_OFFSET + pUV->UVSysParms.XCenterOffsetUm;
        fRes    = StageInfo.XMResol;
        um = stgIdealUMToMappedUM( Axis, cfg, um);
        if( XStageDir < 0)
            um = -um;
35    }

    Counts = (int) ((um + fOffset)/ fRes);

```

```

    iObjOffset = stgObjMOffsetCt( Axis, StageSpec, iObj);
    Counts += iObjOffset;
    return (Counts);
}

5
/*****
 * convert specified( ideal) position to mapped( real stage) position
 *****/

10 static double stgIdealUMToMappedUM( MOT_AXIS Axis, UV_CONFIG *pCfg, double
dIdealUM)
{
    STAGE_MAP_PARMS *pStgMap = &(pCfg->StageMappingParms);
    double    dConverted, dGridLo, dGridHi, dIdealLower;
15    int      iCt;

    dConverted = dIdealUM;
    if( Axis == eAxisY)
    {
20        if( pStgMap->iYPosCt > 1)
        {
            if( dIdealUM > pStgMap->YGridTable[0].RawPosUm &&
                dIdealUM < pStgMap->YGridTable[pStgMap->iYPosCt -1].RawPosUm )
            {
25                iCt = (dIdealUM - pStgMap->YGridTable[0].RawPosUm) /
pStgMap->fYGridSizeUm;
                dIdealLower = pStgMap->YGridTable[ 0 ].RawPosUm + ((double)iCt *
pStgMap->fYGridSizeUm);
                dGridLo    = pStgMap->YGridTable[ iCt  ].RawPosUm;
30                dGridHi    = pStgMap->YGridTable[ iCt + 1].RawPosUm;
                dConverted = dGridLo +
                    ( (dIdealUM - dIdealLower) * (dGridHi -dGridLo)/
pStgMap->fYGridSizeUm);
            }
35        }
    }
    else

```

```

    {
        if( pStgMap->iXPosCt > 1)
        {
            if( dIdealUM > pStgMap->XGridTable[0].RawPosUm &&
5             dIdealUM < pStgMap->XGridTable[pStgMap->iXPosCt -1].RawPosUm )
            {
                iCt = (dIdealUM - pStgMap->XGridTable[0].RawPosUm) /
                pStgMap->fXGridSizeUm;
                dIdealLower = pStgMap->XGridTable[ 0 ].RawPosUm + ((double)iCt *
10 pStgMap->fXGridSizeUm);
                dGridLo = pStgMap->XGridTable[ iCt].RawPosUm;
                dGridHi = pStgMap->XGridTable[ iCt + 1].RawPosUm;
                dConverted = dGridLo +
                    ( (dIdealUM - dIdealLower) * (dGridHi -dGridLo)/
15 pStgMap->fXGridSizeUm);
            }
        }
        return dConverted;
20 }

/*****
 * find the XY offset in motor ct due to Object changes
 *****/

25 */
static int stgObjMOffsetCt( MOT_AXIS Axis, UV_STAGE_SPEC *StageSpec, int iObj)
{
    int iOffset = 0;
    if( iObj >= 0 && iObj < CONFIG_NUM_OBJ)
30 {
        if( Axis == eAxisX)
            iOffset = StageSpec->XOffsetMCt[iObj];
        else
            iOffset = StageSpec->YOffsetMCt[iObj];
35 }
    return iOffset;
}

```

```

/*****
 * find the XY offset in encoder ct due to Object changes
 *****/
*/
5 static int stgObjEOffsetCt( MOT_AXIS Axis, UV_STAGE_SPEC *StageSpec, int iObj)
{
    int iOffset = 0;
    if( iObj >= 0 && iObj < CONFIG_NUM_OBJ)
    {
10         if( Axis == eAxisX)
            iOffset = StageSpec->XOffsetECt[ iObj];
        else
            iOffset = StageSpec->YOffsetECt[ iObj];
    }
15     return iOffset;
}

/*****
 * convert ENCODER count to um
 *****/
20
*/
static double stgEncoderCount2UM( MOT_AXIS Axis, double dCounts, int iObj)
{
    double    UM, fOffset, fRes;
25     int      Dir = 1, iObjOffset;
    UV_CONFIG *cfg = cfg_GetLockRdConfig();
    UV_STAGE_SPEC *StageSpec;
    GLOBAL_DATA *pUV = dat_ObtainGlobDataPtr();

30     StageSpec = &(cfg->StageSpec);
    iObjOffset = stgObjEOffsetCt( Axis, StageSpec, iObj);
    dCounts -= iObjOffset;
    if( Axis == eAxisY)
    {
35         fOffset = YMOT_OFFSET + pUV->UVSysParms.YCenterOffsetUm;
        fRes = StageInfo.YEResol;
        Dir = Y_STAGE_DIR;
    }
}

```



```

    }
    else
    {
        fOffset = XMOT_OFFSET + pUV->UVSysParms.XCenterOffsetUm;
5       fRes   = StageInfo.XEResol;
        Dir    = X_STAGE_DIR;
    }
    UM = (dCounts * fRes) - fOffset;
    if( Dir < 0)
10    UM = -UM;
    return UM;
}

/*****
 * convert um to Encoder count
15 *****/

static int stgUM2EncoderCount( MOT_AXIS Axis, double um, int iObj)
{
    int      Counts, iObjOffset;
20    double  fOffset, fRes;
    UV_CONFIG *cfg = cfg_GetLockRdConfig();
    UV_STAGE_SPEC *StageSpec;
    GLOBAL_DATA      *pUV = dat_ObtainGlobDataPtr();
    int      XStageDir, YStageDir;
25
    XStageDir = X_STAGE_DIR;
    YStageDir = Y_STAGE_DIR;

    StageSpec = &(cfg->StageSpec);
30    if( um > STG_XY_TRAVEL_LIM)
        um = STG_XY_TRAVEL_LIM;
    else
        if( um < -STG_XY_TRAVEL_LIM )
            um = -STG_XY_TRAVEL_LIM;
35
    if( Axis == eAxisY)
    {

```

```

        fOffset = YMOT_OFFSET + pUV->UVSysParms.YCenterOffsetUm;
        fRes = StageInfo.YEResol;
        if( YStageDir < 0)
            um = -um;
5      }
    else
    {
        fOffset = XMOT_OFFSET + pUV->UVSysParms.XCenterOffsetUm;
        fRes = StageInfo.XEResol;
10     if( XStageDir < 0)
        {
            um = -um;
        }

        Counts = (int) ((um + fOffset)/ fRes);
15     iObjOffset = stgObjEOffsetCt( Axis, StageSpec, iObj);
        Counts += iObjOffset;
        return (Counts);
    }

20
    /*
    * move XY relative ,
    * to compensate for objective changes
    */
25
extern BOOL stg_MoveXYToCompensateObjChanges( int iOldObj, int iNewObj)
{
    int iOldCt, iNewCt, iXTargetCt, iYTargetCt;
    UV_CONFIG *cfg = cfg_GetLockRdConfig();
30     BOOL    n = FALSE;

    iXTargetCt = 0;
    iYTargetCt = 0;
    if( iOldObj >= 0 && iNewObj >= 0 && iOldObj < CONFIG_NUM_OBJ && iNewObj <
35     CONFIG_NUM_OBJ)
    {
        iOldCt = cfg->StageSpec.XOffsetMCt[ iOldObj];

```

```

        iNewCt = cfg->StageSpec.XOffsetMCt[ iNewObj];
        iXTargetCt = iNewCt - iOldCt;
        iOldCt = cfg->StageSpec.YOffsetMCt[ iOldObj];
        iNewCt = cfg->StageSpec.YOffsetMCt[ iNewObj];
5       iYTargetCt = iNewCt - iOldCt;
        rt = mot_RelXYMove( iXTargetCt, iYTargetCt);
        }
        return rt;
    }
10

/*****
 * move stage within the robot's range
 * z should be home as part of the safe procedure of known good Z pos
15  *
 * but some control takes extra long( after any actual motion) to home
 * it is takenout for now,
 * but then the motor slip, dangerous, put back in.
 *****/
20  */
extern BOOL stg_MoveCloseForRobot( STG_MOVE_JOY_STATUS eMoveJoy)
{
    STAGE_XY      stagexy_data;
    BOOL          rt = FALSE;
25  int           iRetry;
    BYTE          ch;
    int           XStageDir, YStageDir;
    GLOBAL_DATA   *pUV = dat_ObtainGlobDataPtr();

30  XStageDir = X_STAGE_DIR;
    YStageDir = Y_STAGE_DIR;

    /* home z*/
    if( XStageDir >= 0)
35  stagexy_data.x =  ROBOT_LOAD_XPOS - pUV->UVSysParms.XCenterOffsetUm;
    else
        stagexy_data.x =  ROBOT_LOAD_XPOS + pUV->UVSysParms.XCenterOffsetUm;

```

```

    if( YStageDir >= 0)
        stagexy_data.y =    ROBOT_LOAD_YPOS - pUV->UVSysParms.YCenterOffsetUm;
    else
        stagexy_data.y =    ROBOT_LOAD_YPOS + pUV->UVSysParms.YCenterOffsetUm;
5
    rt = stg_HomeZAndWaitTillDone();
    if( rt == TRUE)
    {
        mot_EnableCurrentAxisJoystick();
10    rt = stgSafeMoveXYStageRawAndNormal( &stagexy_data, OBJ_MOVE_CLOSE_ROBOT,
                                           eStgNoAF, eStgRawPos, eMoveJoy);

    }
    return rt;
15 }
    /*
    *****
    * move z to home and wait till done, or timeout
    *****
    */

20 extern BOOL stg_HomeZAndWaitTillDone( void)
    {
        BOOL    rt = FALSE;
        int      iRetry;
        BYTE      ch;
25
        if( mot_MoveToHome( eAxisZ) == TRUE)
        {
            for( iRetry = 0; rt == FALSE && iRetry < HOME_XY_TIMEOUT_CT; iRetry ++ )
            {
30                rt = mot_CurrAxisGetOneChar( &ch);
                if( rt == FALSE || (int)ch != 'F')
                    rt = FALSE;
            }
        }
35 if( iRetry > 1)
    fprintf( stderr, "homeZ.NumRetry = %d,", iRetry);

```

```

        return rt;
    }

    .....

5    * move stage within the robot's range
    .....

    */

extern BOOL stg_IsStageCloseForRobot( void)
{
10    int rt = FALSE;
    return( TRUE);
}

    .....

    * relative move Z in um
15    * only move within limit
    .....

    */

extern BOOL stg_SafeRelMoveStageZUm( int iUm )
{
20    UV_CONFIG    *pCfg = cfg_GetLockRdConfig();
    OPTICS_PARMS *pOptic = da_GlobalLonPointer();
    BOOL    rt = FALSE;
    float    fZPos, fZLimUm, fTarget;
    int    iObj;

25    lon_APIQuery( eLonAPITurretPosCmd,  &iObj,    0);
    iObj --;
    if( iObj < 0)
        iObj = 0;
30    if( iObj >= CONFIG_NUM_OBJ)
        iObj = CONFIG_NUM_OBJ -1;

    fZLimUm = pCfg->StageSpec.ZLimitUm[iObj] - pOptic->SampleThick;
    if( stg_ReadStageZUm( &fZPos, TRUE) == TRUE)
35    {
        fTarget = fZPos + (double)iUm;
        if( fTarget > 0.0 && (double)fTarget < fZLimUm)

```

```

    {
        rt = stg_MoveStageZUm( fTarget, TRUE);
        mot_EnableCurrentAxisJoystick();
    }
5   }

    return rt;
}

/*****
 * set z joystick travel limit
10  * based on current obj, sample thickness, cfg limits
 *****/

*/

extern BOOL stg_SetZJoystickLimitToCurrParms( void )
{
15  UV_CONFIG    *pCfg = cfg_GetLockRdConfig();
    OPTICS_PARMS *pOptic = da_GlobalLonPointer();
    int          iObj;
    float        fZLimUm;
    BOOL          rt;

20
    lon_APIQuery( eLonAPITurretPosCmd, &iObj, 0);
    iObj--;
    fZLimUm = pCfg->StageSpec.ZLimitUm[iObj] - pOptic->SampleThick;
    rt = mot_SetStgPositiveLim( eAxisZ, (int)( (double)fZLimUm /
25  pCfg->StageSpec.ZMStep2um)) ;
    return rt;
}

30

static BOOL DebugPrtCh( char ch)
{
    int iCh;
35  iCh = ch;
    if( iCh > ' ' && iCh < 0x7f)
        fprintf( stderr, "%c", ch);
}

```

```

        else
            fprintf( stderr, "< %x> ", iCh);
        return TRUE;
    }

5
:

/*****
 * check move
 * code is used only to identify caller, each call
10  * should have an unique number
 *****/

extern BOOL stgWaitTillCurrAxisStop( int Code)
{
15  int iRetry;
    BOOL    bMoving;

    for( bMoving = TRUE, iRetry = 0; iRetry < 200 && bMoving == TRUE; iRetry++)
    {
20      bMoving = mot_IsMotorMoving();
        if( bMoving == TRUE)
            util_sginap( RD_DELAY_TICK );
        fprintf( stderr, "ChkMove%d= %d,", Code, bMoving );
    }
25  return TRUE;
}

/*****
 * correct for stage orthogonality
 * from encoder to application
30  *****/

static BOOL stgOrthoCorrAndDskToAppl( STAGE_XY *DestXY, STAGE_XY *SrcXY,
STG_MOVE_MODE MoveMode)
{
35  STAGE_XY LocXY;
    UV_CONFIG    *pCfg    = cfg_GetLockRdConfig();

```

```

LocXY.y = SrcXY->y;
LocXY.x = SrcXY->x - (SrcXY->y * pCfg->StageSpec.stgOrthoTangent);

if( MoveMode == eStgNormalPos)
5   dsk_InvDeskew( DestXY, &LocXY );
else
    memcpy(DestXY, &LocXY, sizeof( STAGE_XY ));
return TRUE;
}

10
/*****
* correct for stage orthogonality
* and deskew from appl to motor
*****/

15 */
static BOOL stgOrthoCorrAndDskToStg( STAGE_XY *DestXY, STAGE_XY *SrcXY,
STG_MOVE_MODE MoveMode )
{
    STAGE_XY    LocXY;
20   UV_CONFIG   *pCfg  = cfg_GetLockRdConfig();

    LocXY.y = SrcXY->y;
    LocXY.x = SrcXY->x + (SrcXY->y * pCfg->StageSpec.stgOrthoTangent);

25
    if( MoveMode == eStgNormalPos)
        dsk_Deskew( &LocXY, DestXY);
    else
        memcpy( DestXY, &LocXY, sizeof( STAGE_XY ));
30   return TRUE;
}

```

□



```

/*****
 * (c) Copyright 1992,1993, Ultrapointe Corp.
 * All rights reserved
 *
5  * motor.c
 *
 * UltraView1000 software
 * this module handles the low level motor interface:
 *
10 *****/

#ifdef Irix5plus
#define NeedFunctionPrototypes 1
#else
15 #define FUNCPROTO
#endif

#include <stdio.h>
#include <stdlib.h>
20 #include <X11/Intrinsic.h>
#include "uv_util.h"
#include "uvconfig.h"
#include "motor.h"
#include "rs232.h"
25 #include "lon_msg.h"
#include "lon_img.h"
#include "options.h"

#define TIMEOUT_LINE 5 /*timeout: 5second for line */
30 #define TIMEOUT_MOV_DONE 4
#define TIMEOUT_CHAR 2

#define MAX_MOTOR_SPEED 150000
#define DEFAULT_ACC 240000 /* 7um/sec/sec */
35 #define START_MTR_SPEED 2500 /* um/sec */

#define MAX_ZMOTOR_SPEED 3000

```

```

#define START_ZMTR_SPEED 300      /* this is important for AutoFocus */
#define DEFAULT_ZACC 20000 /* um/sec/sec */

#define XY_POSITIVE_LIM_UM 230000
5  #define Z_POSITIVE_LIM_UM 8000

#define CH_BUFFER_SIZE 80
static HANDLE      hdStageX = -1;
static HANDLE      hdStageY = -1;
10 static HANDLE      hdStageZ = -1;
static RS232_PORT_SELECT SGIPortX; /* GMI_RS2_0, 1, SGI */
static RS232_PORT_SELECT SGIPortY;
static RS232_PORT_SELECT SGIPortZ;

15 static MOT_AXIS   CurrAxis = eAxisNone;
static RS232_PORT_SELECT CurrSGIPort;
static HANDLE      CurrhdStage;

20 static BOOL      bXOnline = FALSE;
static BOOL      bYOnline = FALSE;
static BOOL      bZOnline = FALSE;

#define CMD_ABS_MOVE      "MA%d\r"
25 #define CMD_REL_MOVE      "MR%d\r"
#define CMD_MOVE_HOME      "MH\r"
#define CMD_REP_POS      "RP\r\n"
#define CMD_REP_ENC      "RE\r\n"
#define CMD_SET_SPEED      "VF%d\r"
30 #define CMD_SET_INTT_SPEED "VI%d\r"
#define CMD_SET_ACC      "AC%d\r"
#define CMD_SET_POSITION  "SP%d\r"
#define CMD_SET_ENCODER   "SE%d\r"
#define CMD_WAKE_X      "WU1\r"
35 #define CMD_WAKE_Y      "WU2\r" /* some delay */
#define CMD_WAKE_Z      "WU3\r"
#define CMD_WAKE_NOBODY   "WU9\r"

```

```

#define CMD_SHORT_CMD      "MED\r"
#define CMD_STOP           "ST\r"
#define CMD_LO_IDLE       "ICD\r"
#define CMD_ENABLE_JOY     "AE\r"
5  #define CMD_DISABLE_JOY  "\33\r"
#define CMD_JOYSTICK_GAIN  "JG%d\r"
#define CMD_JOY_DEADBAND   "JD15\r"
#define CMD_POSITIVE_LIM   "YP%d\r"
#define CMD_NEGATIVE_LIM   "YN%d\r"
10
#define CMD_MOV_DONE_ACK   "MFE\r"
#define CMD_MOV_DONE_NAK   "MFD\r"
#define CMD_REPORT_STATUS  "RS\r"
#define CMD_NOTHING        "\r" /* delay, space filler */
15 #define EOL               0xd

static BOOL motProcessReportMessage( int *Pos);
static BOOL motOnlineOneAxisInit( HANDLE hdStage, RS232_PORT_SELECT SGIPort,
                                int InitSpd, int FinalSpd, int Acc, int PosLimUm, float MStep2um);
20 static BOOL motInitAccSpeed( int InitSpd, int FinalSpd, int Acc, float MStep2um);
static BOOL motWakeupAxisAndCmd(MOT_AXIS Axis, char *pCmd);
static BOOL motWakeupAxisFast(MOT_AXIS Axis);
static BOOL motCurrAxisSendStringAndGetLine( char *Cmd, char *chBuff, int iSize);
static BOOL motNoCheckWakeUnit( HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
25 *WakeCmd);
static void DebugDisplayOneChar( BYTE chBuf);
static BOOL motSendAndWaitStar( HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
*pCmd);
static BOOL motCurrAxisSendCmdToReadInt( char *pQuery, int *iValue);
30 static double motMapEncCt( int iNumPos, double dGridSizeEnc,
                            STAGE_MAP_GRID_PAIR *GridTable, int iRawCt);
static BOOL motCurrAxisRdMappedEncCt( double *dEncCt, int *iRawCt);

35
/***** these function is for NEAT compatibility *****/
static BOOL command_noecho(HANDLE hdStage, RS232_PORT_SELECT SGIPort, char

```

```

*cmdstr);
static BOOL command_noresp(HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
*cmdstr);
static RS2_CMD_STATUS send_echostr(HANDLE hdStage, RS232_PORT_SELECT SGIPort,
5 char *str);
static BOOL wake_unit_funct_fast(HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
*WakeCmd);

/*****/
10

/*****
* initialization once per program
*****

15 */
extern BOOL mot_Initialize( int hUldd) /* program start */
{
    hdStageZ = hUldd;
    hdStageX = hUldd;
20 hdStageY = hUldd;
    return TRUE;
}
extern BOOL mot_ClosePort( void) /* program end */
{
25 mot_DisableJoystick();
    rs2_ClosePort(hdStageZ, SGIPortZ );
    if( hdStageZ != hdStageX)
    {
        rs2_ClosePort(hdStageX, SGIPortX );
30 rs2_ClosePort(hdStageY, SGIPortY );
    }
    hdStageX = -1;
    hdStageY = -1;
    hdStageZ = -1;
35 return TRUE;
}
/*****

```

```

    * user turn on stage
    .....

    */

extern BOOL mot_PutStageOnline( void)

5      {
        UV_CONFIG      *pCfg = cfg_GetLockRdConfig();
        UV_STAGE_SPEC *StageSpec;
        BOOL            r = FALSE;
        char            chBuf[CH_BUFFER_SIZE + 10];
10     volatile int j;

        if( opt_HasOptionHardware() != TRUE)
            return FALSE;

15     StageSpec = &(pCfg-> StageSpec);

        if( StageSpec->PortNum >= 0 && StageSpec->PortNum < 6 )
        {
            switch( StageSpec->PortNum)
20         {
                case 0:
                case 1: SGIPortX = eSGI_RS2_Port;
                       SGIPortY = eSGI_RS2_Port;
                       SGIPortZ = eSGI_RS2_Port;
25                 break;

                case 2: SGIPortX = eGMI_RS2_0;
                       SGIPortY = eGMI_RS2_0;
                       SGIPortZ = eGMI_RS2_0;
30                 break;

                case 4:
                case 5:
                       SGIPortX = eSGI_RS2_Port;
35                 SGIPortY = eSGI_RS2_Port;
                       SGIPortZ = eGMI_RS2_1;
                       break;
            }
        }
    }

```

```

        case 3:
        default:SGIPortX = eGMI_RS2_1;
                SGIPortY = eGMI_RS2_1;
                SGIPortZ = eGMI_RS2_1;
5         break;
    }
    if( StageSpec->PortNum == 4 || StageSpec->PortNum == 5)
    {
        hdStageZ = rs2_initialize( hdStageZ, 3, BAUDRATE96);
10     if( StageSpec->PortNum == 4)
        {
            hdStageX = rs2_initialize( hdStageZ, 0, BAUDRATE96);
            hdStageY = rs2_initialize( hdStageZ, 1, BAUDRATE96);
        }
15     else
        {
            hdStageX = rs2_initialize( hdStageZ, 10, BAUDRATE96);
            hdStageY = rs2_initialize( hdStageZ, 11, BAUDRATE96);
        }
20     }
    else
    {
        hdStageZ = rs2_initialize( hdStageZ, StageSpec->PortNum, BAUDRATE96);
        hdStageX = hdStageZ;
25     hdStageY = hdStageZ;
    }
    fprintf( stderr, "initialized axisxyz = %x, %x, %x, ". hdStageX, hdStageY, hdStageZ );

    /*****
30     * all these extra wakeups seem neccessary to ensure
    * reliable init, especially at controller powe-up
    *****/

    /*
        bXOnline = bYOnline = bZOnline = FALSE;
35     if( StageSpec->bZOnline == TRUE)
        mot_SendWakeUnitString( eAxisZ);    /* z axis is very difficult to talk to */

```

```

    if( StageSpec->bXOnline == TRUE)
    {
        rt = motNoCheckWakeUnit( hdStageX, SGIPortX, CMD_WAKE_X);
        if( rt == TRUE)
5           rt = motOnlineOneAxisInit( hdStageX, SGIPortX,
                                         StageSpec->iXYInitSpdUm, StageSpec->iXYFinalSpdUm,
                                         StageSpec->iXYAccUm,
                                         XY_POSITIVE_LIM_UM, StageSpec->XMStep2um);

        bXOnline = rt;
10    }

    if( StageSpec->bYOnline == TRUE && rt == TRUE)
    {
        rt = motNoCheckWakeUnit( hdStageY, SGIPortY, CMD_WAKE_Y);
        if( rt == TRUE)
15           rt = motOnlineOneAxisInit( hdStageY, SGIPortY,
                                         StageSpec->iXYInitSpdUm, StageSpec->iXYFinalSpdUm,
                                         StageSpec->iXYAccUm,
                                         XY_POSITIVE_LIM_UM, StageSpec->YMStep2um);

        bYOnline = rt;
20    }

    if( StageSpec->bZOnline == TRUE && rt == TRUE )
    {
        rt = motNoCheckWakeUnit( hdStageZ, SGIPortZ, CMD_WAKE_Z);
        if( rt == TRUE)
25           rt = motOnlineOneAxisInit( hdStageZ, SGIPortZ,
                                         START_ZMTR_SPEED ,MAX_ZMOTOR_SPEED,
                                         DEFAULT_ZACC,
                                         StageSpec->ZLimitUm[0], StageSpec->ZMStep2um);

        bZOnline = rt;
30    }

    fprintf( stderr, "XYZOnline=%d,%d,%d,", bXOnline, bYOnline, bZOnline);
    if( rt == TRUE)
        rt = (bXOnline || bYOnline || bZOnline);
    }
35    return rt;
}

```

```

/*****
* initialize one axis
*****/

5 static BOOL motOnlineOneAxisInit( HANDLE hdStage, RS232_PORT_SELECT SGIPort,
                                   int InitSpd, int FinalSpd, int Acc, int PosLimUm, float MStep2um)
{
    BOOL    rt;
    char chBuf[100];

10    rt = motSendAndWaitStar( hdStage, SGIPort, CMD_SHORT_CMD);
    if( rt == TRUE)
        rt = motSendAndWaitStar( hdStage, SGIPort, CMD_STOP);
    if( rt == TRUE)
15    {
        sprintf( chBuf, CMD_SET_INTT_SPEED, (int)( (double)InitSpd / MStep2um));
        rt = motSendAndWaitStar(hdStage, SGIPort, chBuf );
    }
    if( rt == TRUE)
20    {
        sprintf( chBuf, CMD_SET_SPEED, (int)( (double)FinalSpd / MStep2um));
        rt = motSendAndWaitStar(hdStage, SGIPort, chBuf );
    }
    if( rt == TRUE)
25    {
        sprintf( chBuf, CMD_SET_ACC, (int)( (double)Acc / MStep2um));
        rt = motSendAndWaitStar( hdStage, SGIPort, chBuf);
    }
    if( rt == TRUE)
30    rt = motSendAndWaitStar( hdStage, SGIPort, CMD_LO_IDLE); /* idle curr disabled */
    if( rt == TRUE)
        rt = motSendAndWaitStar( hdStage, SGIPort, CMD_JOY_DEADBAND);
    if( rt == TRUE)
        rt = motSendAndWaitStar( hdStage, SGIPort, CMD_MOV_DONE_ACK);
35    if( rt == TRUE)
        rt = mot_SetCurrentStgNegativeLim( 0);
    if( rt == TRUE)

```



```

    r1 = mot_SetCurrentStgPositiveLim( (double)PosLimUm / MStep2um );
    if( r1 == TRUE)
        r1 = mot_SetCurrentStgJoyStickGain( 1);
    return r1;
5   }
    /*****
    * set XY axis to default speed
    *****/
    */
10  extern BOOL mot_SetXYToDefaultSpeeds( void )
    {
        char          chBuf[100];
        UV_CONFIG      *cfg = cfg_GetLockRdConfig();
        BOOL           r1 = FALSE;
15  float             XMStep2Um, YMStep2Um;

        XMStep2Um = cfg->StageSpec.XMStep2um;
        YMStep2Um = cfg->StageSpec.YMStep2um;

20  r1 = motWakeupAxisAndCmd( eAxisX, NULL);
    if( r1 == TRUE)
        {
            sprintf( chBuf, CMD_SET_INIT_SPEED, (int)( (double)cfg->StageSpec.iXYInitSpdUm/
            XMStep2Um) );
25  r1 = motSendAndWaitStar(CurrhdStage, CurrSGIPort, chBuf);
        }
    if( r1 == TRUE)
        {
            sprintf( chBuf, CMD_SET_SPEED, (int)( (double)cfg->StageSpec.iXYFinalSpdUm/
30  XMStep2Um) );
            r1 = motSendAndWaitStar( CurrhdStage, CurrSGIPort, chBuf);
        }
    if( r1 == TRUE)
        r1 = motWakeupAxisAndCmd( eAxisY, NULL);
35  if( r1 == TRUE)
        {
            sprintf( chBuf, CMD_SET_INIT_SPEED, (int)( (double)cfg->StageSpec.iXYInitSpdUm/

```

```

YMStep2Um) );
    rt = motSendAndWaitStar(CurrhdStage, CurrSGIPort, chBuf);
    }
    if( rt == TRUE)
5      {
        sprintf( chBuf, CMD_SET_SPEED, (int)( (double)cfg->StageSpec.iXYFinalSpdUm/
YMStep2Um) );
        rt = motSendAndWaitStar( CurrhdStage, CurrSGIPort, chBuf);
        }
10    return rt;
    }

extern BOOL mot_PutStageOffline( void)
15    {
        mot_DisableJoystick();
        if( bZOnline)
            rs2_ClosePort(hdStageZ, SGIPortZ );
        if( hdStageZ != hdStageX && bXOnline && bYOnline)
20        {
            rs2_ClosePort(hdStageX, SGIPortX );
            rs2_ClosePort(hdStageY, SGIPortY );
        }
        CurrhdStage = -1;
25        bXOnline = FALSE;
        bYOnline = FALSE;
        bZOnline = FALSE;
        return TRUE;
    }
30

/*****
* high level functions
*****/

*/
35 extern BOOL mot_RelMove(MOT_AXIS Axis, int Pos)
    {
        char chBuf[80];

```

```

        BOOL    r = FALSE;

        sprintf( chBuf, CMD_REL_MOVE, Pos);
        r = motWakeupAxisAndCmd( Axis, chBuf);
5      return r;
    }

extern BOOL mot_CurrAxisRelMove( int Pos)
    {
        char chBuf[80];
10      BOOL    r = FALSE;

        sprintf( chBuf, CMD_REL_MOVE, Pos);
        r = motSendAndWaitStar( CurrhdStage, CurrSGIPort, chBuf);
        return r;
15    }

extern BOOL mot_AbsMove(MOT_AXIS Axis, int Pos)
    {
20      char chBuf[80];
        BOOL    r = FALSE ;

        sprintf( chBuf, CMD_ABS_MOVE, Pos);
        r = motWakeupAxisAndCmd( Axis, chBuf);
25      return r;
    }

/******
 * move both XY
*****
30  */

extern BOOL mot_AbsXYMove( int XPos, int YPos)
    {
        char chBuf[80];
        BOOL    r = FALSE;
35

        sprintf( chBuf, CMD_ABS_MOVE, XPos);
        r = motWakeupAxisAndCmd( eAxisX, chBuf);

```

```

    if( rt == TRUE)
    {
        sprintf( chBuf, CMD_ABS_MOVE, YPos);
        rt = motWakeupAxisAndCmd( eAxisY, chBuf);
5      }
    return rt;
}

/*****
10  * this function do relative move,
    * will not wait for stop
    *****/

    /*
extern BOOL mot_RelXYMove( int XPos, int YPos)
15  {
    char chBuf[80];
    BOOL    rt = FALSE ;

    sprintf( chBuf, CMD_REL_MOVE, XPos);
20    rt = motWakeupAxisAndCmd( eAxisX, chBuf);
    if( rt == TRUE)
    {
        sprintf( chBuf, CMD_REL_MOVE, YPos);
        rt = motWakeupAxisAndCmd( eAxisY, chBuf);
25    }
    return rt;
}

/*****
    * enable XYZ joystick
30  *****/

    /*
extern BOOL mot_EnableJoystick( void)
    {
        volatile int j;
35    BOOL    rt = FALSE;

    rt = motWakeupAxisAndCmd( eAxisX, CMD_ENABLE_JOY);

```

```

    if( r1 == TRUE)
        r1 = motWakeupAxisAndCmd( eAxisY, CMD_ENABLE_JOY);
    if( r1 == TRUE)
        r1 = motWakeupAxisAndCmd( eAxisZ, CMD_ENABLE_JOY);
5   return r1;
    }

/*****
* enable XYZ joystick and set xy joystick gain
10 *****/
*/
extern BOOL mot_EnableJoystickAndSetXYGain( int iJGain)
{
    char      chBuf[80];
15   BOOL      r1 = FALSE;

    if( iJGain < 1) iJGain = 1;
    if( iJGain > 60) iJGain = 60;

20   sprintf( chBuf, CMD_JOYSTICK_GAIN, iJGain);
    r1 = motWakeupAxisAndCmd( eAxisX, chBuf);
    if( r1 == TRUE)
        r1 = motSendAndWaitStar(CurrhdStage, CurrSGIPort, CMD_ENABLE_JOY);
    if( r1 == TRUE)
25     r1 = motWakeupAxisAndCmd( eAxisY, chBuf);
    if( r1 == TRUE)
        r1 = motSendAndWaitStar(CurrhdStage, CurrSGIPort, CMD_ENABLE_JOY);
    if( r1 == TRUE)
        r1 = motWakeupAxisAndCmd( eAxisZ, CMD_ENABLE_JOY);
30   return r1;
    }

/*****
* enable XY joystick
35 *****/
*/
extern BOOL mot_EnableXYJoystick( void)

```

```

    {
        volatile int j;
        BOOL    r = FALSE;

5       r = motWakeupAxisAndCmd( eAxisX, CMD_ENABLE_JOY);
        if( r == TRUE)
            r = motWakeupAxisAndCmd( eAxisY, CMD_ENABLE_JOY);
        return r;
    }
10    /*****
        * disable joystick
        *****/

    /*
extern BOOL mot_DisableJoystick( void)
15    {
        BOOL    r = FALSE;
        if( bXOnline == TRUE)
            {
                if( hdStageZ == hdStageX)
20                 r = motSendAndWaitStar( hdStageX, SGIPortX, CMD_WAKE_X);
                else
                    r = TRUE;
                if( r == TRUE)
                    r = motSendAndWaitStar( hdStageX, SGIPortX, CMD_DISABLE_JOY);
25                if( r == TRUE)
                    r = rs2_WaitForResp( hdStageX, SGIPortX, '*', TIMEOUT_CHAR);
            }
        if( r == TRUE && bYOnline == TRUE)
            {
30                if( hdStageZ == hdStageX)
                    r = motSendAndWaitStar( hdStageY, SGIPortY, CMD_WAKE_Y);
                if( r == TRUE)
                    {
                        r = motSendAndWaitStar( hdStageY, SGIPortY, CMD_DISABLE_JOY);
35                        if( r == TRUE)
                            r = rs2_WaitForResp( hdStageY, SGIPortY, '*', TIMEOUT_CHAR);
                    }
            }
    }

```

```

    }
    if( r == TRUE && bZOnline == TRUE)
    {
        if( hdStageZ == hdStageX)
5         r = motSendAndWaitStar( hdStageZ, SGIPortZ, CMD_WAKE_Z);
        if( r == TRUE)
        {
            r = motSendAndWaitStar( hdStageZ, SGIPortZ, CMD_DISABLE_JOY);
            if( r == TRUE)
10             r = rs2_WaitForResp( hdStageZ, SGIPortZ, '*', TIMEOUT_CHAR);
        }
    }
    return r;
}
15  /*****
    * enable joystick, no checking
    *****/
    */
extern BOOL mot_EnableCurrentAxisJoystick( void)
20  {
    return( motSendAndWaitStar( CurrhdStage, CurrSGIPort, CMD_ENABLE_JOY) );
}
/*****
    * stop current axis
25  * need to send without normal char-by-char echo check
    * because we may get the 'F' status in the middle of sending
    * the stop cmd
    *****/
    */
30  extern BOOL mot_StopCurrentAxis( void)
    {
        BOOL r;

        r = motSendAndWaitStar( CurrhdStage, CurrSGIPort, CMD_STOP);
35  return r;
    }

```

```

/*****
* send command to home stage
*****/

5  extern BOOL mot_MoveToHome(MOT_AXIS Axis)
    {
        BOOL    r = FALSE;

        r = motWakeupAxisAndCmd( Axis, CMD_STOP);
10  if( r == TRUE)
        r = motSendAndWaitStar( CurrhdStage, CurrSGIPort, CMD_MOVE_HOME);
        return r;
    }
/*****
15  * send command to home stage
*****/

    extern BOOL mot_CurrAxisMoveHome( void )
    {
20  BOOL    r = FALSE;

        r = motSendAndWaitStar(CurrhdStage, CurrSGIPort, CMD_MOVE_HOME);
        return r;
    }
25
/*****
* Send Stop Command to the specified axis
*****/

30  extern BOOL mot_StopOneAxis(MOT_AXIS Axis)
    {
        BOOL    r = FALSE ;

        r = motWakeupAxisAndCmd( Axis, CMD_STOP);
35  return r;
    }

extern BOOL mot_SetPos(MOT_AXIS Axis, int Pos)

```



```

    {
        char chBuf[80];
        BOOL    r = FALSE;

5       sprintf( chBuf, CMD_SET_POSITION, Pos);
        r = motWakeupAxisAndCmd( Axis, chBuf);

        return r;
    }
10    extern BOOL mot_SetAcc(MOT_AXIS Axis, int Acc)
    {
        char chBuf[80];
        BOOL    r = FALSE;

15       sprintf( chBuf, CMD_SET_ACC, Acc);
        r = motWakeupAxisAndCmd( Axis, chBuf);
        return r;
    }
    extern BOOL mot_SetEnc(MOT_AXIS Axis, int Pos)
20    {
        char chBuf[80];
        BOOL    r = FALSE;

        sprintf( chBuf, CMD_SET_ENCODER, Pos);
25       r = motWakeupAxisAndCmd( Axis, chBuf);

        return r;
    }

30    /*****
        * query position/encoder
        *****/

    /*
    extern BOOL mot_ReportOnePos(MOT_AXIS Axis, int *Pos)
35    {
        BOOL r = FALSE;

```

```

    *Pos = 0;
    if( motWakeupAxisAndCmd(Axis, NULL) == TRUE)
        r = motCurrAxisSendCmdToReadInt( CMD_REP_POS, Pos);
    return r;
5   }

/*****
 * query current axis motor counts
 *****/

*/

10  extern BOOL mot_CurrAxisReportMotorCt( int *MotCt)
    {
        BOOL    r = FALSE;

        r = motCurrAxisSendCmdToReadInt( CMD_REP_POS, MotCt);
15  return (r);
    }

/*****
 * report encoder count
 *****/

20  */

extern BOOL mot_ReportOneEnc(MOT_AXIS Axis, double *Pos, int *RawPos)
    {
        BOOL r = FALSE;

25  if( motWakeupAxisFast( Axis) == TRUE)
        r = motCurrAxisRdMappedEncCt( Pos, RawPos);
    return r;
    }

/*****
 * report current axis encode count
 *****/

*/

30  extern BOOL mot_CurrAxisReportEnc( double *Pos, int *RawPos)
    {
35  BOOL r = FALSE;

        r = motCurrAxisRdMappedEncCt( Pos, RawPos);

```

```

    return rt;
}

/*****
* apply stage mapping to encoder read
5  * encoder count read from stage is matched to
*      RawCt in uvconfig,
* mapped encode cts are from Computed Counts
*
* Computed  ---0---1-----n-1---
10 * raw Ct   ---0---1---n-1---
* case 1      *
* case 2          *
* case 3              *
*
15 *
* No correction done to case 1,3
*****/

*/

static BOOL motCurrAxisRdMappedEncCt( double *dEncCt, int *iRawCt)
20 {
    UV_CONFIG    *cfg = cfg_GetLockRdConfig();
    BOOL          rt = FALSE, bFound;
    int           iNumPos, iPosCt;
    double        dGridSizeEnc;
25  STAGE_MAP_GRID_PAIR    *GridTable;

    *dEncCt = 0;
    *iRawCt = 0;

    rt = motCurrAxisSendCmdToReadInt( CMD_REP_ENC, iRawCt);
30  if( rt == TRUE && (CurrAxis == eAxisX || CurrAxis == eAxisY ) )
    {
        *dEncCt = *iRawCt;
        if( CurrAxis == eAxisX)
        {
35          iNumPos    = cfg->StageMappingParms.iXPosCt;
          dGridSizeEnc    = cfg->StageMappingParms.dXGridSizeEncCt;
          GridTable    = cfg->StageMappingParms.XGridTable;

```

```

    }
    else
    {
        iNumPos    = cfg->StageMappingParms.iYPosCt;
5       dGridSizeEnc    = cfg->StageMappingParms.dYGridSizeEncCt;
        GridTable   = cfg->StageMappingParms.YGridTable;
    }
    *dEncCt = motMapEncCt( iNumPos, dGridSizeEnc, GridTable, *iRawCt);
    }
10    return rt;
    }

/*****
 * convert raw encoder ct to Mapped Enc ct
15  * depending on stage axis direction
 * the array of Counts can be either increasing/
 * or decreasing
 *
 *****/
20  */
static double motMapEncCt( int iNumPos, double dGridSizeEnc,
                          STAGE_MAP_GRID_PAIR *GridTable, int iRawCt)
{
    double    dMappedEnc, dComputedCtLo;
25    BOOL      bFound;
    int       iPosCt, iRawCtLo, iRawCtHi;

    dMappedEnc = iRawCt;
    if( iNumPos > 1)
30    {
        if( GridTable[0].RawEncCt < GridTable[1].RawEncCt) /* normal.. number increasing */
        {
            if( iRawCt > GridTable[0].RawEncCt && iRawCt < GridTable[iNumPos -1].RawEncCt)
            {
35                bFound = FALSE;
                for( iPosCt= 1; iPosCt < (iNumPos -1) && bFound == FALSE; )
                {

```

```

        if( iRawCt < GridTable[ iPosCt].RawEncCt)
            bFound = TRUE;
        else
            iPosCt++;
5      }

        dComputedCtLo    = GridTable[0].RawEncCt + ((iPosCt-1) * dGridSizeEnc);
        iRawCtLo          = GridTable[ iPosCt -1].RawEncCt;
        iRawCtHi          = GridTable[ iPosCt  ].RawEncCt;
10      dMappedEnc = dComputedCtLo + ((double)( iRawCt - iRawCtLo) * dGridSizeEnc /
                                   (iRawCtHi - iRawCtLo));
    }
}
else
15  {
        /* decreasing */
        if( iRawCt < GridTable[0].RawEncCt && iRawCt > GridTable[iNumPos -1].RawEncCt)
        {
            bFound = FALSE;
            for( iPosCt= 1; iPosCt < (iNumPos -1) && bFound == FALSE; )
20          {
                if( iRawCt > GridTable[ iPosCt].RawEncCt)
                    bFound = TRUE;
                else
                    iPosCt++;
25          }

            dComputedCtLo    = GridTable[0].RawEncCt - ((iPosCt-1) * dGridSizeEnc);
            iRawCtLo          = GridTable[ iPosCt -1].RawEncCt;
            iRawCtHi          = GridTable[ iPosCt  ].RawEncCt;
30      dMappedEnc = dComputedCtLo - ((double)( iRawCt - iRawCtLo ) * dGridSizeEnc /
                                   (iRawCtHi - iRawCtLo));
    }

    /*
    fprintf( stderr, "Ct=%d, Raw/Computed Ct=%d,%d \n", iPosCt, iRawCt, iMappedEnc);
35  */
}
}

```

```

        return dMappedEnc;
    }

5  /*****
   * send cmd(must be a query cmd), and expect to get an integer
   *****/
   */

static motCurrAxisSendCmdToReadInt( char *pQuery, int *iValue)
10  {
    BOOL rt = FALSE;
    if( hdStageZ != hdStageX)
    {
        if ( rs2_FlushQueSendString( CurrhdStage, CurrSGIPort, pQuery) == TRUE)
15      rt = motProcessReportMessage( iValue);
    }
    else
    {
        if( rs2_SendStringWithPerCharHandShake(CurrhdStage, CurrSGIPort, pQuery) ==
20      eRS2StatusOK)
            rt = motProcessReportMessage( iValue);
    }
    return rt;
}

25  /*****
   * the report message has the command mixed in with the result
   * assuming the real reports are all in numbers
   *****/

30  */

static BOOL motProcessReportMessage( int *Pos)
{
    BOOL rt;
    char chBuf[CH_BUFFER_SIZE + 2];
35  int j, iLen, Sign= 1;

    rt = rs2_GetOneLine(CurrhdStage, CurrSGIPort, chBuf, CH_BUFFER_SIZE,

```

```

TIMEOUT_LINE);
    if( r == TRUE)
        r = rs2_WaitForResp( CurrhdStage, CurrSGIPort, '*', 2);
    iLen = strlen( chBuf);
5    if( r == TRUE && iLen > 0 )
        {
            for( j=0; j < iLen && r == TRUE && (int)chBuf[j] != EOL; j++)
                if( !isdigit( chBuf[j] ) && !isalpha( chBuf[j] ) && chBuf[j] != (unsigned)'\'' )
                    r = FALSE;
10
            for( j=0; j < iLen && !isdigit( chBuf[j] ) && chBuf[j] != (unsigned)'\'' ; j++)
                {};
            if( j < iLen && r == TRUE )
                {
15                 *Pos = atoi( &chBuf[j]);
                    r = TRUE;
                }
            else
                r = FALSE;
20    if( r == FALSE)
        fprintf( stderr, "str=%s,\n", chBuf);
        }
        return r;
    }
25    /*****
        * send a string to neat. and get a response
        *****/

    /*
static BOOL motCurrAxisSendStringAndGetLine( char *Cmd, char *chBuff, int iSize)
30    {
        BOOL r = FALSE;
        if( hdStageZ != hdStageX)
            {
                if( rs2_FlushQueSendString( CurrhdStage, CurrSGIPort, Cmd) == TRUE)
35                 r = rs2_GetOneLine( CurrhdStage, CurrSGIPort, chBuff, iSize, TIMEOUT_LINE);
            }
        else

```

```

    {
        if( rs2_SendStringWithPerCharHandShake( CurrhdStage, CurrSGIPort, Cmd) ==
eRS2StatusOK)
            rt = rs2_GetOneLine( CurrhdStage, CurrSGIPort, chBuff, iSize, TIMEOUT_LINE);
5        }
        if( rt == TRUE)
            rt = rs2_WaitForResp( CurrhdStage, CurrSGIPort, '*', 2);
        return rt;
    }
10

extern BOOL mot_SetSpeed( MOT_AXIS Axis, int Speed)
{
    char chBuf[80];
15    BOOL    rt = FALSE;

    sprintf( chBuf, CMD_SET_SPEED, Speed);
    rt = motWakeupAxisAndCmd( Axis, chBuf);
    return rt;
20    }

/*****
 * set init and final speed to specified axis
 *****/

25    */
extern BOOL mot_SetInitAndFinalSpeedsCt( MOT_AXIS Axis, int iInitSpdCt, int iFinalSpdCt)
{
    char chBuf[80];
    BOOL    rt;
30

    sprintf( chBuf, CMD_SET_INIT_SPEED, iInitSpdCt);
    rt = motWakeupAxisAndCmd( Axis, chBuf);
    if( rt == TRUE)
    {
        35    sprintf( chBuf, CMD_SET_SPEED, iFinalSpdCt );
        rt = motSendAndWaitStar(CurrhdStage, CurrSGIPort, chBuf );
    }

```



```

        return r;
    }

extern BOOL mot_SetCurrAxisSpeed( int Speed)
5   {
    char chBuf[80];
    BOOL    r = FALSE;

    sprintf( chBuf, CMD_SET_SPEED, Speed);
10   r = motSendAndWaitStar( CurrhdStage, CurrSGIPort, chBuf);
    return r;
    }

/*****
 * reset speed to default
15  * if( bAxisCmd == TRUE, send WakeUnit command
 *****/

*/
extern BOOL mot_SetDefaultSpeed( MOT_AXIS Axis, BOOL bAxisCmd)
    {
20   UV_CONFIG    *cfg    = cfg_GetLockRdConfig();
    char chBuf[80];
    BOOL    r = FALSE;
    int Speed;

25   switch( Axis)
    {
        case eAxisX: Speed = (double)cfg->StageSpec.iXYFinalSpdUm/
cfg->StageSpec.XMStep2um;
            break;
30   case eAxisY: Speed = (double)cfg->StageSpec.iXYFinalSpdUm/
cfg->StageSpec.YMStep2um;
            break;

        case eAxisZ:
35   default: Speed = (double)MAX_ZMOTOR_SPEED / cfg->StageSpec.ZMStep2um;
            break;
    }
}

```

```

    sprintf( chBuf, CMD_SET_SPEED, Speed);
    if( bAxisCmd == TRUE)
        rt = motWakeupAxisAndCmd( Axis, chBuf);
    else
5        rt = motSendAndWaitStar(CurrhdStage, CurrSGIPort, chBuf);

    return rt;
}

/*****
10  * che f the specified axis is moving
    * return TRUE if yes
    *****/

*/
extern BOOL mot_IsAxisMoving(MOT_AXIS Axis)
15  {
    BOOL    rt = FALSE;

    if( motWakeupAxisAndCmd( Axis, NULL) == TRUE)
        rt = mot_IsMotorMoving();
20  return (rt);
}

/*****
    * check if the currently selected motor is moving
25  * return FALSE if no motor connected
    * NOTE:
    * !!!!
    * this is control specific. Check here if need to change
    * controller
30  * !!!!
    * while Report Status does not return 'F', wait and retry
    * after max retry, still invalid response, then assume stopped
    *****/

*/
35  extern BOOL mot_IsMotorMoving( void)
    {
        char    chBuf[CH_BUFFER_SIZE + 2];

```

```

int      count, j;
BOOL      r = FALSE, bSure = FALSE;

if( bXOnline == TRUE || bYOnline == TRUE || bZOnline == TRUE )
5      {
        chBuf[0] = 0;
        chBuf[1] = 0;
        chBuf[2] = 0;
        r = motCurrAxisSendStringAndGetLine( CMD_REPORT_STATUS, chBuf,
10      CH_BUFFER_SIZE);
        if( r == TRUE && ( (int)chBuf[0] == 'F' || (int)chBuf[0] == 'M' ||
                        (int)chBuf[2] == 'F' || (int)chBuf[2] == 'M' ) )
            {
                for( j=0; j < strlen( chBuf); j++)
15      DebugDisplayOneChar( chBuf[j]);
                r = ( (int)chBuf[0] == 'M' || (int)chBuf[2] == 'M')? TRUE: FALSE;
            }
        else
            {
20      for( count = 0; count < 30 && bSure == FALSE; count++)
                {
                    chBuf[0] = 0;
                    chBuf[1] = 0;
                    chBuf[2] = 0;
25      fprintf( stderr, "%d,", count);
                    r = motCurrAxisSendStringAndGetLine( CMD_REPORT_STATUS, chBuf,
                        CH_BUFFER_SIZE);
                    if( r == TRUE && ( (int)chBuf[0] == 'F' || (int)chBuf[0] == 'M' ||
                        (int)chBuf[2] == 'F' || (int)chBuf[2] == 'M' ) )
30      bSure = TRUE;
                }
                r = ( (int)chBuf[0] == 'M' || (int)chBuf[2] == 'M')? TRUE: FALSE;
            }
        }
35      return( r);
    }

```

```

/*****
 * switch to specified axis
 *****/

*/
5  extern BOOL mot_SwitchAxis( MOT_AXIS Axis)
    {
        BOOL    r = FALSE;

        r = motWakeupAxisAndCmd( Axis, NULL);
10  return r;
    }

#define WAKEUNIT_BUF_SIZE 70
extern BOOL mot_SendWakeUnitString( MOT_AXIS Axis)
15  {
        BOOL    r = FALSE;
        char chBuf[WAKEUNIT_BUF_SIZE + 2];
        char     *chPtr;

20  CurrAxis = Axis;
        switch( Axis)
        {
            case eAxisX: CurrSGIPort = SGIPortX;
                        CurrhdStage = hdStageX;
25  chPtr = CMD_WAKE_X;
                        break;

            case eAxisY: CurrSGIPort = SGIPortY;
                        CurrhdStage = hdStageY;
30  chPtr = CMD_WAKE_Y;
                        break;

            case eAxisZ: CurrSGIPort = SGIPortZ;
                        CurrhdStage = hdStageZ;
35  chPtr = CMD_WAKE_Z;
                        break;
        }
        if( chPtr)
            r = rs2_FlushQueSendString( CurrhdStage, CurrSGIPort, chPtr);

```

```

    if( rt == TRUE)
        rt = rs2_GetOneLine(CurrhdStage, CurrSGIPort, chBuf, WAKEUNIT_BUF_SIZE, 2);
    if( rt == TRUE)
    {
5      fprintf( stderr, "str= %s\n", chBuf);
        rt = rs2_FlushQueSendString( CurrhdStage, CurrSGIPort, CMD_SHORT_CMD);
    }
    return rt;
}

10  /******
    * move motor to specified position, on currently selected
    * axis
    *****/

    /*
15  extern BOOL mot_CurrAxisAbsMove( int TargetPos)
    {
        BOOL    rt = FALSE;
        char chBuf[30];

20      sprintf( chBuf, CMD_ABS_MOVE, TargetPos);
        rt = motSendAndWaitStar( CurrhdStage, CurrSGIPort, chBuf);
        return rt;
    }

    /******
25  * get one char from current axis
    *****/

    /*
    extern BOOL mot_CurrAxisGetOneChar( BYTE * ch)
    {
30      BOOL rt;
        rt = rs2_GetOneChar( CurrhdStage, CurrSGIPort, ch, TIMEOUT_CHAR);
        return( rt );
    }

35  /******
    * send command to wakeup the axis

```

```

* and if command is specified, send it also
*****

*/

static BOOL motWakeupAxisAndCmd(MOT_AXIS Axis, char *pCmd)
5   {
    BOOL    r = FALSE;
    char     *pWake = NULL;

    if( (Axis == eAxisX && bXOnline == TRUE) ||
10    (Axis == eAxisY && bYOnline == TRUE) || (Axis == eAxisZ && bZOnline ==
TRUE) )
    {
        CurrAxis = Axis;
        switch( Axis)
15    {
            case eAxisX: CurrhdStage = hdStageX;
                        CurrSGIPort = SGIPortX;
                        pWake = CMD_WAKE_X;
                        break;
20    case eAxisY: CurrhdStage = hdStageY;
                        CurrSGIPort = SGIPortY;
                        pWake = CMD_WAKE_Y;
                        break;
25    case eAxisZ: CurrhdStage = hdStageZ;
                        CurrSGIPort = SGIPortZ;
                        pWake = CMD_WAKE_Z;
                        break;
        }
        if( hdStageZ != hdStageX)
30    r = TRUE;
        else
        {
            if( pWake)
35    r = motNoCheckWakeUnit( CurrhdStage, CurrSGIPort, pWake);
        }
        if( r == TRUE && pCmd)
            r = motSendAndWaitStar(CurrhdStage, CurrSGIPort, pCmd);
    }
}

```

```

    }
    if( rt == FALSE)
    fprintf( stderr, "wake %d = %d,", Axis, rt);
    return rt;
5    }

/*****
* send command to wakeup the axis
10 *****/
*/

static BOOL motWakeupAxisFast(MOT_AXIS Axis)
{
    BOOL    rt = FALSE;
15    char    *pWake;

    if( (Axis == eAxisX && bXOnline == TRUE) ||
        (Axis == eAxisY && bYOnline == TRUE) || (Axis == eAxisZ && bZOnline ==
TRUE) )
20    {
        CurrAxis = Axis;
        switch( Axis)
        {
            case eAxisX: CurrhdStage = hdStageX;
25                CurrSGIPort = SGIPortX;
                pWake = CMD_WAKE_X;
                break;

            case eAxisY: CurrhdStage = hdStageY;
                CurrSGIPort = SGIPortY;
30                pWake = CMD_WAKE_Y;
                break;

            case eAxisZ: CurrhdStage = hdStageZ;
                CurrSGIPort = SGIPortZ;
35                pWake = CMD_WAKE_Z;
                break;

        }

        if( hdStageZ != hdStageX)

```

```

        rt = TRUE;
    else
    {
        if( pWake)
5         rt = wake_unit_funct_fast( CurrhdStage, CurrSGIPort, pWake);
    }
}

if( rt == FALSE)
fprintf( stderr, "wake %d= %d,", Axis, rt);
10

return rt;
}

/*****
 * send a command to neat and wait for a *
15 *****/

*/

static BOOL motSendAndWaitStar( HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
*pCmd)
{
20     BOOL    rt = FALSE;

    if( hdStageZ != hdStageX)
    {
        if( rs2_FlushQueSendString(hdStage, SGIPort, pCmd) == TRUE)
25         rt = rs2_WaitForResp(hdStage, SGIPort, '*', TIMEOUT_CHAR);
    }
    else
    {
        if( rs2_SendStringWithPerCharHandShake(hdStage, SGIPort, pCmd) == eRS2StatusOK)
30         rt = rs2_WaitForResp(hdStage, SGIPort, '*', TIMEOUT_CHAR);
    }

    return rt;
}

35

/*****
 * This function is copied from NEAT 3/25/93

```



```

*
* wake_unit_func: wakes 310 unit given axis index, returns TRUE if OK, *
* FALSE if bad axis index value or I/O error *
*****
5  */
static BOOL motNoCheckWakeUnit( HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
*WakeCmd)
{
    int      count,dval;
10  BOOL      r = FALSE;

    command_noecho(hdStage, SGIPort, "");    /* send out CR to clear and delay */
    command_noecho(hdStage, SGIPort, CMD_WAKE_NOBODY);
    command_noecho(hdStage, SGIPort, "");    /* send out CR to clear and delay */
15  command_noecho(hdStage, SGIPort, WakeCmd);    /* send out command, don't expect echo
*/
    if(command_noresp(hdStage, SGIPort, "") == TRUE)    /* expect echo */
        r = TRUE;
20  else
    {    /* unit not responding--retry */
        for(count=0; count<5 && r == FALSE; count++)
        {
            command_noecho(hdStage, SGIPort, WakeCmd);    /* resend wake command--don't
25  check echo */
            if(command_noresp(hdStage, SGIPort, WakeCmd) == TRUE ) /* resend again checking
echo */
                r = TRUE;
        }
30  if( count > 1)
        fprintf( stderr, "WakeCt=%d,", count);
    }
    return r;
}
35
/*****
* simplified wakeup, good for

```

```

    * getting XY stage position
    *****

    */

static BOOL wake_unit_func_fast(HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
5  *WakeCmd)
    {
        int          count,dval;
        BOOL          rt = FALSE;

10     for(count=0; count<5 && rt == FALSE; count++)
        {
            command_noecho(hdStage, SGIPort, WakeCmd);          /* send out command, don't
            expect echo */
            if(command_noresp(hdStage, SGIPort, WakeCmd) == TRUE ) /* resend again checking
15     echo */
                rt = TRUE;
        }

        if( count > 1)

20     fprintf( stderr, "WakeCt= %d,", count);

        return rt;
    }

25

    /*****
    *****/

    */

30     #define NO_ECHO_LOOP_CT          5
        #define NO_ECHO_WAIT_MS        5

static BOOL command_noecho(HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
    *cmdstr)
35     {
        BOOL          rt = FALSE;
        BYTE          ch;

```

```

    int      TimeoutCt;

    rt = rs2_FlushQueSendString( hdStage, SGIPort, cmdstr);
    if( rt == TRUE)
5      {
        TimeoutCt = NO_ECHO_LOOP_CT;
        do
        {
            if(rs2_NonblockGetOneChar( hdStage, SGIPort, &ch) == TRUE)
10          TimeoutCt = NO_ECHO_LOOP_CT;
            util_WaitMs( NO_ECHO_WAIT_MS);
            }while(--TimeoutCt > 0);
        }
    return( rt );
15 }

/*****
 * This function is for NEAT compatibility
 *****/

20 */

static BOOL command_noresp(HANDLE hdStage, RS232_PORT_SELECT SGIPort, char
*cmdstr)
{
25  BOOL  rt;

    rt = motSendAndWaitStar( hdStage, SGIPort, cmdstr);
    return rt;
}

30 /*****
 * This function is here for NEAT compatibility
 * send_echostr: sends given string over serial port, confirming an echo *
 *               on each character sent, and appends a CR onto the      *
 *               string and looks for no more characters, returns TRUE  *
35  *               if no I/O errors and all characters in string were    *
 *               echoed, returns FALSE if I/O error or bad or no echo,  *
 *               stops sending characters after the 1st bad one          *
 *****/

```

```

*****
*/
static RS2_CMD_STATUS send_echostr(HANDLE hdStage, RS232_PORT_SELECT SGIPort,
char *str)
5   {
    RS2_CMD_STATUS rtStat ;
    rtStat = rs2_SendStringWithPerCharHandShake( hdStage, SGIPort, str);
    return rtStat;
  }
10

/*****
* testing gmi board:
15 * PORT a and B need to looped back
* LaserPwrControl node must be present
*
*****
*/
20 #define RANDOM_SIZE 500
#define READ_LOOP (RANDOM_SIZE / 20)
#define PORTA 2
#define PORTB 3

25 extern BOOL mot_UltrapGMITest( int iLoops)
{
    BYTE *InputA , *OutputA, *InputB, *OutputB;
    int iCt, IAIndex, IBIndex, Vers;
    int iChkCt, iRdCt, iChCt;
30    BOOL rt = FALSE, bHold, bOK;

#ifdef DO_GMI_TEST

    InputA = UvMalloc( RANDOM_SIZE +10 );
35    InputB = UvMalloc( RANDOM_SIZE +10 );
    OutputA = UvMalloc( RANDOM_SIZE +10 );
    OutputB = UvMalloc( RANDOM_SIZE +10 );

```

```

    fprintf( stderr, "\nGMI rs232 ports need to be connected to each other.\n");
    fprintf( stderr, "Laser power control node must be present.\n");
    fprintf( stderr, "This test sends random char between the 2 port and read sw\n");
5   fprintf( stderr, "version number from laser power control node.\n");
    fprintf( stderr, "Will do %d times.\n", iLoops );

    for( iCt= 0; iCt < RANDOM_SIZE; iCt++)
    {
10      OutputA [iCt] = rand() & 0xff;
      OutputB [iCt] = rand() & 0xff;
    }
    if( rs2_initialize( hdStage, PORTA, BAUDRATE96) == hdStage )
      if( rs2_initialize( hdStage, PORTB, BAUDRATE96) == hdStage)
15      {
        bOK = TRUE;
        rt = TRUE;
        for( iCt = 0; iCt < iLoops && rt == TRUE; iCt++)
        {
20      fprintf( stderr, "%d.", iCt);
          IAIndex = 0;
          IBIndex = 0;          /* swap */
          memcpy( InputA, OutputA, RANDOM_SIZE);
          memcpy( OutputA, OutputB, RANDOM_SIZE);
25      memcpy( OutputB, InputA, RANDOM_SIZE);

          rs2_SendOneLine( hdStage, eGMI_RS2_0, (char *)OutputA, RANDOM_SIZE);
          rs2_SendOneLine( hdStage, eGMI_RS2_1, (char *)OutputB, RANDOM_SIZE);
          for( iRdCt = 0; iRdCt < READ_LOOP &&
30              ( IAIndex < RANDOM_SIZE || IBIndex < RANDOM_SIZE ) &&
              bOK == TRUE; iRdCt++ )
          {
            for( iChCt = 0; iChCt < 50 && IAIndex < RANDOM_SIZE; iChCt++)
            {
35              bHold = rs2_NonblockGetOneChar( hdStage, eGMI_RS2_0, &InputA[IAIndex]);
              if( bHold == TRUE)
                IAIndex++;
            }
          }
        }
      }

```

```

    }
    for( iChCt=0; iChCt < 50 && IBIndex < RANDOM_SIZE; iChCt++)
    {
        bHold = rs2_NonblockGetOneChar( hdStage, eGMI_RS2_1, &InputB[IBIndex]);
5      if( bHold == TRUE)
        IBIndex++;
    }
    bOK = lonimg_QueryNetVar( LASER_PWR_NODE_ADDR,
    LASER_SW_REV_NV_ADDR, &Vers, TRUE);
10    }
    if( iRdCt >= READ_LOOP)
        fprintf( stderr, "timeout, no responds rx: A=%d,B=%d,\n", IAIndex, IBIndex );

    for( iChkCt= 0; iChkCt < RANDOM_SIZE && r == TRUE; iChkCt++)
15    {
        if( (int)InputA[iChkCt] != (int)OutputB[iChkCt] )
        {
            fprintf( stderr, "rs232 error,rx from PORTA at %d,", iChkCt);
            r = FALSE;
20        }
        if( (int)InputB[iChkCt] != (int)OutputA[iChkCt] )
        {
            fprintf( stderr, "rs232 error,rx from PORTB at %d,", iChkCt);
            r = FALSE;
25        }
    }
    fprintf( stderr, "firstchar=%X,%X,%X,%X,\n", InputA[0], OutputB[0],
        InputB[0], OutputA[0] );

    if( bOK != TRUE)
30    {
        fprintf( stderr, "lon_error,");
        r = FALSE;
    }
}
35 }

UvFree( InputA);

```

```

        UvFree( InputB);
        UvFree( OutputA);
        UvFree( OutputB);
    #else
5      fprintf( stderr, "gmi test not implemented for now. \n");
    #endif

        return rt;
    }
10  /******
    * set current stagelimit to specified counts
    *****/

    /*
    extern BOOL mot_SetCurrentStgPositiveLim( int StgCt)
15    {
        char chBuf[80];
        BOOL    rt;

        sprintf( chBuf, CMD_POSITIVE_LIM, StgCt);
20    rt = motSendAndWaitStar(CurrhdStage, CurrSGIPort, chBuf );
        return rt;
    }

    extern BOOL mot_SetStgPositiveLim(MOT_AXIS Axis, int iCt)
    {
25    char chBuf[80];
        BOOL    rt;

        sprintf( chBuf, CMD_POSITIVE_LIM, iCt);
        rt = motWakeupAxisAndCmd( Axis, chBuf);
30    return rt;
    }

    extern BOOL mot_SetCurrentStgNegativeLim( int StgCt)
    {
35    char chBuf[80];
        BOOL    rt;

```

```

    sprintf( chBuf, CMD_NEGATIVE_LIM, StgCt);
    rt = motSendAndWaitStar(CurrhdStage, CurrSGIPort, chBuf );
    return rt;
}

5  extern BOOL mot_SetCurrentStgJoyStickGain( int iJGain)
    {
        char chBuf[80];
        BOOL    rt;

10     sprintf( chBuf, CMD_JOYSTICK_GAIN, iJGain);
        rt = motSendAndWaitStar( CurrhdStage, CurrSGIPort, chBuf );
        return rt;
    }

extern BOOL mot_SetStgJoyStickGain( MOT_AXIS Axis, int iJGain)
15     {
        BOOL    rt = FALSE;

        if( motWakeupAxisAndCmd( Axis, NULL) == TRUE)
            rt = mot_SetCurrentStgJoyStickGain( iJGain);

20     return rt;
    }

    /******
    * debug
    *****/

25     */

static void DebugDisplayOneChar( BYTE chBuf)
    {
        if( ( (int)chBuf > ' ') && ((int)chBuf < 0x7f) )
            fprintf( stderr, "%c", chBuf);

30     else
        fprintf( stderr, "< %x>", chBuf);
    }

```

□



```

/*****
* Ultrapointe Corp. copyright 1992,93,94
*
* Interface to rs232
5  * for info: see man page serial/termio:
*   man serial
*   man termio
* and page 75-89 of Advance UNIX programming by Marc J. Rochkind
*
10  * Note:
*   4 ports supported in UV:
*   PortNum = 0,1,2,3 = SGI_PORT1, SGI_PORT2, ULTRA_PORT0, ULTRA_PORT1
*   port assignment:
*   stage: PortNum 0,1,2
15  * robot: PortNum 0,1,3
*
* SGI_PORT1 and SGI_PORT2 may be used as terminal by the
* system. Make sure serial port used here is turned OFF by the
* system, and uv has r/w access to /dev/ttyd1/2
20  *
*****
*/

/*
25  #define MOTOR_DEBUG
*/

#ifdef Irix5plus
#define NeedFunctionPrototypes 1
#else
30  #define FUNCPROTO
#endif

#include <stdio.h>
#include <termio.h>
35  #include <errno.h>
#include <setjmp.h>
#include <signal.h>

```

```

rs232.c

#include <unistd.h>
#include <limits.h>          /* sginap */
#include <sys/stat.h>
#include <fcntl.h>
5  #include <unistd.h>
    #include <X11/Intrinsic.h>

#include "uv_util.h"
#include "rs232.h"
10  #include "uldd.h"

/* cannot find the prototype include file for ioctl..... workaround and defined here....
*/
extern int ioctl (int fildes, int request, ...);
15

#define MAX_LINE 100

#define PORT1 "/dev/ttyd1"
20  #define PORT2 "/dev/ttyd2"
    #define PORT1X "/dev/ttyd05%d"

#define TIMEOUT_CHAR      1    /* something that should have no delay */
#define TIMEOUT_LINE      10   /* something that is slow */
25  #define TIMEOUT_RESPONSE 10

#define NEAT_INTER_CHAR_WAIT 1    /* wait 1ms per char */

static jmp_buf env_alarm;
30

static BOOL rsSetRaw(HANDLE hdRS232, int BaudRate);
static void sig_alarm( int signo);
static void DebugDisplayOneChar( BYTE chBuf);
static BOOL rs2SendOneChar( HANDLE hdRS232, RS232_PORT_SELECT SGIPort, char ch);
35

/*****

```

rs232.c

```

* initialization:
*   PortNum: 0= PORT1,1=PORT2
*   BaudRate: 0= 9600, 1=19200, 2=38400
* see: man page on termio
5  *****

*/
extern int rs2_initialize( int hUldd, int PortNum, int BaudRate )
{
  char  *pPort, chBuf[80];
10  int  hdRS232 = -1;
  CMD_RS2_CONTROL RstCtl;

  switch (PortNum)
  {
15  case 0:
    case 1:
    case 10:
    case 11:
    case 12:
20  case 13:
    case 14:
    case 15:
    case 16:
    case 17:
25  if( PortNum >= 10 && PortNum <= 17 )
    {
      sprintf( chBuf, PORT1X, PortNum -10);
      pPort = chBuf;
    }
30  else
    pPort = ( PortNum == 1)? PORT2 :PORT1;

    hdRS232 = open( pPort, O_RDWR | O_NOCTTY);
    if( hdRS232 >= 0)
35  {
    if( rsSetRaw( hdRS232, BaudRate) == FALSE)
    {

```

rs232.c

```

        close( hdRS232);
        hdRS232 = -1;
    }
}
5      break;
case 2:
case 3:
    RstCtl.CtlCode = ( PortNum == 2)? eCmdRS2InitPort0 : eCmdRS2InitPort1;
    switch( BaudRate)
10      {
        case BAUDRATE192: RstCtl.BaudRate = GMI_BAUD_19200; break;
        case BAUDRATE384: RstCtl.BaudRate = GMI_BAUD_38400; break;
        case BAUDRATE96:
        default:          RstCtl.BaudRate = GMI_BAUD_9600; break;
15      }

    if( ioctl( hUldd, eCmdRS2Control, &RstCtl ) >= 0)
    {
        hdRS232 = hUldd;
20      }
    else
        fprintf( stderr, "ioctlerr=%X, ", errno );

        break;
25      }
    if( hdRS232 < 0)
        fprintf( stderr, "cannot open, error=%d,", errno);
    return hdRS232;
}
30

extern BOOL rs2_ClosePort( HANDLE hdRS232, RS232_PORT_SELECT SGIPort)
{
    if(SGIPort == eSGI_RS2_Port && hdRS232 >=0)
35      close(hdRS232);
    return TRUE;
}

```

```

rs232.c

/*****
 * send one string to rs232.. For the NEAT controller
 * check for echo of each char
 * up to 2 different char recieved per one sent are tolerated
5  * NEAT's reponse may have 'F' mixed in the char stream
 *****/

*/

#define NUM_DIFF_TOLERATED      3
extern RS2_CMD_STATUS rs2_SendStringWithPerCharHandShake( HANDLE hdRS232,
10 RS232_PORT_SELECT SGIPort, char * chPtr)
{
    RS2_CMD_STATUS rtStat= eRS2StatusBad;
    BYTE chTerm = '\r', chBuf;
    int      iDiffCt;
15    BOOL      rt;

    if( hdRS232 >=0 && strlen(chPtr) > 0 )
    {
        rs2_FlushInputQueue(hdRS232, SGIPort);
20        rtStat = eRS2StatusOK;
        do
        {
            rs2SendOneChar( hdRS232, SGIPort, *chPtr);
            rtStat = eRS2StatusBad;
25        for( iDiffCt = 0; iDiffCt < NUM_DIFF_TOLERATED && rtStat == eRS2StatusBad;
            iDiffCt++)
        {
            rt = rs2_GetOneChar( hdRS232, SGIPort, &chBuf, TIMEOUT_CHAR);
            if ( rt == FALSE)
30                rtStat = eRS2StatusTimeout;
            else
            {
                if( ((int)chBuf & 0x5f) == ((int)*chPtr & 0x5f) )
                    rtStat = eRS2StatusOK;
35                else
                {
                    fprintf( stderr, "chrx= %X,", chBuf);

```

```

        util_WaitMs( NEAT_INTER_CHAR_WAIT);
    }
}

5      chPtr++;
    } while( rtStat == eRS2StatusOK && (int)*chPtr != '\r' && (int)*chPtr != '\0' );
    rs2SendOneChar( hdRS232, SGIPort, chTerm) ;
}
return ( rtStat);
10 }

/*****
 * SGIPort = eSGI_RS2_Port: from standard indigo port
 *      = eGMI_RS2_0/1: GMI board, hdRS232=0: channel A
 *****/

15 */
static BOOL rs2SendOneChar( HANDLE hdRS232, RS232_PORT_SELECT SGIPort, char ch)
{
    CMD_RS2_WR_PACKET WrPack;
    int iWr= 0;
20    BOOL    r = FALSE;

    if( SGIPort == eSGI_RS2_Port )
        iWr = write( hdRS232, &ch, 1 ) ;
    else
25    {
        WrPack.Line    = (SGIPort == eGMI_RS2_0)? 0: 1;
        WrPack.iMaxChar    = 1;
        WrPack.iWr    = 0;
        WrPack.chBuf    = &ch;
30    if( ioctl( hdRS232, eCmdRS2BlockWrite, &WrPack ) >= 0)
        iWr = WrPack.iWr;
    }
    return (r = ( iWr> 0)? TRUE: FALSE );
}

35 /*****
 * bSGIPort = TRUE: from standard indigo port
 *      = FLASE: GMI board, hdRS232=0: channel A
 *****/

```

```

***** rs232.c

*/
extern int rs2_SendOneLine( HANDLE hdRS232, RS232_PORT_SELECT SGIPort, char *chPtr,
int NumChar)
5   {
    CMD_RS2_WR_PACKET WrPack;
    int iWr = 0;

    if( SGIPort == eSGI_RS2_Port)
10      iWr = write( hdRS232, chPtr, NumChar ) ;
    else
    {
        WrPack.Line    = (SGIPort == eGMI_RS2_0)? 0: 1;
        WrPack.iMaxChar = NumChar;
15      WrPack.iWr      = 0;
        WrPack.chBuf    = chPtr;
        if( ioctl( hdRS232, eCmdRS2BlockWrite, &WrPack ) >= 0)
            iWr = WrPack.iWr;
    }
20  return iWr;
    }

/*****
* send a string to rs232, dont wait for response
25  *****/
*/
extern BOOL rs2_FlushQueSendString( HANDLE hdRS232, RS232_PORT_SELECT SGIPort,
char * chPtr)
{
30  BOOL    r1 = FALSE;

    if( hdRS232 >= 0 && strlen(chPtr) > 0 )
    {
        rs2_FlushInputQueue( hdRS232, SGIPort);
35      if( rs2_SendOneLine( hdRS232, SGIPort, chPtr, strlen( chPtr ) > 0)
            r1 = TRUE;
    }
}

```

rs232.c

```

    return rt;
}
/*****
 * send a string to rs232, no flush, no wait for response
5 *****/
*/
extern BOOL rs2_SendStringNoWaitNoFlush( HANDLE hdRS232, RS232_PORT_SELECT
SGIPort, char * chPtr)
{
10     BOOL    rt = FALSE;

    if( hdRS232 >=0 && strlen(chPtr) > 0 )
    {
        if( rs2_SendOneLine( hdRS232, SGIPort, chPtr, strlen( chPtr ) > 0)
15         rt = TRUE;
    }
    return rt;
}
/*****
20 * wait for the correct respond
 * return TRUE if received correct responds
 *     FALSE if timed out
 *****/
*/
25 extern BOOL rs2_WaitForResp(HANDLE hdRS232, RS232_PORT_SELECT SGIPort, char ch,
int Seconds)
{
    BYTE chBuf[10];
    BOOL    rt= FALSE;
30
    if( hdRS232 >=0)
    {
        do
        {
35         rt = rs2_GetOneChar( hdRS232, SGIPort, chBuf, Seconds);
        } while ( rt == TRUE && (int)chBuf[0] != (int)ch );
    }
}

```



rs232.c

```

    return r;
}

/*****
 * wait to get one char
5  * return 1 if got a char,
 *    0 if not: file closed, or timeout
 *****/

*/

extern BOOL rs2_GetOneChar( HANDLE hdRS232, RS232_PORT_SELECT SGIPort, BYTE
10 *ch, int Seconds)
{
    BOOL r = FALSE;
    CMD_RS2_RD_PACKET RdPack;

15    if( hdRS232 >=0)
    {
        if( SGIPort == eSGI_RS2_Port)
        {
            signal( SIGALRM, sig_alarm);
20            if( setjmp( env_alarm) != 0)
                return FALSE;          /* return 0 if timeout */
            else
            {
                alarm( Seconds);
25                if( read( hdRS232, ch, 1) > 0)
                    r = TRUE;
                alarm( 0);
            }
        }
30    else
    {
        RdPack.Line      = ( SGIPort == eGMI_RS2_0)? 0: 1;
        RdPack.iMaxChar = 1;
        RdPack.iRd       = 0;
35        RdPack.chBuf    = (char *)ch;
        RdPack.iTimeout = Seconds * CLK_TCK;
        if( ioctl( hdRS232, eCmdRS2BlockRead, &RdPack ) >= 0)

```

rs232.c

```

        {
            if( RdPack.iRd == 1)
                r = TRUE;
        }
5      }

    }

    return r;
}

10  /*****
    * get one char if available, assume rtyport ( if used) set to non-block
    * return TRUE  if got a char,
    *      FALSE if not: empty, file closed, or timeout
    *****/

15  */

extern BOOL rs2_NonblockGetOneChar( HANDLE hdRS232, RS232_PORT_SELECT SGIPort,
BYTE *ch)
{
    BOOL r = FALSE;
20  CMD_RS2_RD_PACKET RdPack;

    if( hdRS232 >=0)
    {
        if( SGIPort == eSGI_RS2_Port)
25      {
            if( rs2_SetToNonBlockedRead( hdRS232) == TRUE)
            {
                if ( read( hdRS232, ch, 1) > 0)
                    r = TRUE;
30      }
            rs2_SetToBlockedRead( hdRS232);
        }
    }
    else
    {
35  RdPack.Line      = ( SGIPort == eGMI_RS2_0)? 0: 1;
        RdPack.iMaxChar = 1;
        RdPack.iRd      = 0;
    }
}

```

rs232.c

```

    RdPack.chBuf      = (char *)ch;
    if( ioctl( hdRS232, eCmdRS2NonBlockRd, &RdPack ) >= 0)
    {
        if( RdPack.iRd > 0)
5         rt = TRUE;
    }
}
return rt;
10 }

15
/*****
 * software interrupt
 * from alarm timer
 *****/

20 */
static void sig_alarm( int signo)
{
    longjmp( env_alarm, 1);
}

25
extern BOOL rs2_FlushInputQueue( HANDLE hdRS232, RS232_PORT_SELECT SGIPort)
{
    BOOL      rt = FALSE;
    char chBuf[MAX_LINE];
30    CMD_RS2_RD_PACKET RdPack;

    if( SGIPort == eSGI_RS2_Port)
    {
        if( rs2_SetToNonBlockedRead( hdRS232) == TRUE)
35        {
            while( read( hdRS232, chBuf, 1) > 0)
                DebugDisplayOneChar( chBuf[0]);
        }
    }
}

```

```

        rt = rs2_SetToBlockedRead( hdRS232);
    }
}
else
5   {
    rt = TRUE;
    do
    {
        RdPack.Line      = (SGIPort == eGMI_RS2_0)? 0: 1;
10   RdPack.iMaxChar = MAX_LINE;
        RdPack.iRd  = 0;
        RdPack.chBuf  = chBuf;
        if( ioctl( hdRS232, eCmdRS2NonBlockRd, &RdPack ) < 0)
            rt = FALSE;
15   } while ( rt == TRUE && RdPack.iRd > 0);
    }
    return rt;
}

/*****
20  * get one line, terminated by LF or CR
   * return FALSE if TIMEOUT, or max line length exceeded
   *****/

*/

extern BOOL rs2_GetOneLine( HANDLE hdRS232, RS232_PORT_SELECT SGIPort, char
25  *chBuf, int iMax,
                               int Seconds)
{
    int iCt;
    BOOL    rt = FALSE, bEnd = FALSE;
30
    if( hdRS232 >=0 && iMax > 0)
    {
        rt = TRUE;
        for( iCt=0; iCt < iMax && bEnd == FALSE && rt == TRUE; iCt++)
35   {
            rt = rs2_GetOneChar( hdRS232, SGIPort, (BYTE *)&chBuf[iCt], Seconds);
            if( rt == TRUE)

```

rs232.c

```

    {
        DebugDisplayOneChar( chBuf[iCt]);
        if( (int)chBuf[iCt] == 0xd || (int)chBuf[iCt] == 0xa)
        {
5           bEnd = TRUE;
            chBuf[iCt] = '\0';
        }
    }
10    }

    if( rt == FALSE )
        chBuf[0] = '\0';
    return( rt);
}

15  /*****
    * set to raw mode
    * BaudRate: 0= 9600;
    *           1= 19200
    *****/

20  */

static BOOL rsSetRaw(HANDLE hdRS232, int BaudRate)
{
    struct termio tbuf;
    BOOL      rt= FALSE;
25  char ibuf[MAX_LINE];

    if( hdRS232 >=0)
    {
        if( ioctl( hdRS232, TCGETA, &tbuf) == -1)
30         fprintf(stderr, "ioctl get:error ");
        else
        {
            tbuf.c_iflag &= ~(INLCR | ICRNL | IUCLC | ISTRIP | IXON | BRKINT);
            tbuf.c_oflag &= ~OPOST;
35         tbuf.c_lflag &= ~(ICANON | ISIG | ECHO | XCASE);
            tbuf.c_cflag &= ~(CBAUD | CSIZE | PARENB ); /* baud, 8bit parity */
            switch ( BaudRate )

```

rs232.c

```

    {
        case BAUDRATE192:
            tbuf.c_cflag |= ( B19200 | CS8 | CREAD | CLOCAL);
            break;
5       case BAUDRATE384:
            tbuf.c_cflag |= ( B38400 | CS8 | CREAD | CLOCAL);
            break;
        case BAUDRATE96:
        default:
10         tbuf.c_cflag |= ( B9600 | CS8 | CREAD | CLOCAL);
            break;
    }

    tbuf.c_cc[0] = 0xff;
    tbuf.c_cc[1] = 0xff;
15    tbuf.c_cc[2] = 0xff;
    tbuf.c_cc[3] = 0xff;
    tbuf.c_cc[VMIN] = 1; /* min */
    tbuf.c_cc[VTIME] = 1; /* time */
    if( ioctl(hdRS232, TCSETAF, &tbuf) == -1)
20     fprintf(stderr, "ioctl set: error ");
    else
    {
        rt = rs2_FlushInputQueue( hdRS232, eSGI_RS2_Port );
    }
25 }

return rt;
}

/*****
30  * set to BLOCKed read
  * see man page termio. and fcntl
  *****/

extern BOOL rs2_SetToBlockedRead( HANDLE Receiver)
35 {
    int FFlag;
    FFlag = fcntl( Receiver, F_GETFL, 0);

```

rs232.c

```

    if( fcntl( Receiver, F_SETFL, FFlag &( ~FNDELAY) ) == 0)
        return TRUE;
    else
        return FALSE;
5   }
    /*****
    * set to non-BLOCKed read
    * see man page termio, and fcntl
    *****/
10  */
    extern BOOL rs2_SetToNonBlockedRead( HANDLE Receiver )
    {
        int  FFlag;
        FFlag = fcntl( Receiver, F_GETFL, 0);
15  if( fcntl( Receiver, F_SETFL, FFlag | FNDELAY) == 0)
        return TRUE;
        else
            return FALSE;
    }
20  /*****
    * debug
    *****/
    */
    static void DebugDisplayOneChar( BYTE  chBuf)
25  {
    #ifdef MOTOR_DEBUG
        if( chBuf > ' ' )
            fprintf( stderr, "%c", chBuf);
        else
30  . fprintf( stderr, "< %x>", chBuf);
    #endif
    }

```

5

**A METHOD AND APPARATUS FOR PERFORMING  
AN AUTOMATIC FOCUS OPERATION**

10

**Christopher R. Fairley  
Timothy V. Thompson  
Ken K. Lee**

15

20

**APPENDIX F**

25

**M-2464-1P US**



1/6/95

/\*\*\*\*\*

\* (c) Copyright 1994

\* Ultrapointe Corp.

\*

5 \* page scan fine Z autofocus

\* host computer support functions

\*

\* functions used but not listed:

\*

10 \* lon\_APICmd()

\* set the specified subsystem parameter to specified value

\*

\* lon\_APIQuery()

\* query the status of the the specified subsystem parameter

15 \*

\* cfg\_GetLockRdConfig

\* Get pointer to system configuration data structure

\*

\*\*\*\*\*

20 \*/

/\* system hardware image \*/

static UI\_LON\_IMAGE UIImage;

/\*\*\*\*\*

25 \* controlling the page scanning hardware

\*

\*

\*

\*

30 \*

\*

\*

\*

\*

\*

35 \*

\*

\*

\*

| 1 |

2

| 3 |

1

| 2

regions

1/6/95

\*

\* region 2 gives the visible data

\*

\* scanner f= 7.875 kHz, duty cycle for y=256 is 92%, y=512 is 96%

5

\*

\* example: y=256 (region 2) = 256 cycles

\* retrace( region 3+1)= 5ms = 39.4 cycles

\*

\*

10

\*  $\frac{\text{region}(1+2)}{\text{region } 1+2+3} = \frac{\text{region}(1+2)}{295.4} = .92 \rightarrow \text{region3} = 23.6$ 

\* for a 92% duty cycle = ----- = .92 = ----- -&gt; region3 = 23.6

\*  $\frac{\text{region } 1+2+3}{295.4} \text{ then, region1} = 39.4 -$ 

\*

23.6=15.7

\* converting these into RAMPCLK, 8times

15

\* ie: y=256: VSYNC = 15.7 \* 8 = 125.6 = 126

\* TRACE = 256 \* 8 = 2048

\* RETRACE = 23.6 \* 8 = 188.8 = 189

\* peak-to-peak-ampl = Incr \* ( VSYNC + TRACE + 1) = Decr \* ( RETRACE + 2)

\* for peak-to-peak = about 12288 -&gt; 6 \* (126 + 2048 + 1) = 13050

20

\* and adjust DECR to fit,

\* after adjusting( duty cycle, peak-to-peak)

\* for y=128: VSYNC = 131 duty cycle = 90.4 %

\* TRACE = 1024

\* RETRACE = 202

25

\* INCR = 12

\* DECR = 68

\* Total Ampl = (1024+131+1) \* 12 = 13872 -&gt; half height 6936 -&gt; start =

\*

1256

\* for y=256: VSYNC = 127 92%

30

\* TRACE = 2048

\* RETRACE = 190

\* INCR = 6

\* DECR = 68

\* Total Ampl = (2048+127+1) \* 6 = 13056 -&gt; half height 6528 -&gt; start = 1664

35

\*

1/6/95

```

* for y=512: VSYNC = 135                      96%
*           TRACE = 4096
*           RETRACE = 182
*           INCR = 3
5 *           DECR = 69
*           Total AmpI = (4096+135+1) * 3 = 12696 -> half height 6348 -> start =
*                                           1844
*****

```

\*/

10

typedef struct PAGE\_SCAN\_PARAMS\_tag

```

{
    int    iStartValue;
    int    iIncr;
15  int    iDecr;
    int    iVSync;
    int    iTrace;
    int    iRetrace;
    }PAGE_SCAN_PARAMS;

```

20

```

static PAGE_SCAN_PARAMS PageScanParms128 = { 1256, 12, 68, 131, 1024, 202 };
static PAGE_SCAN_PARAMS PageScanParms256 = { 1664, 6, 68, 127, 2048, 190 };
static PAGE_SCAN_PARAMS PageScanParms512 = { 1844, 3, 69, 135, 4096, 182 };
static PAGE_SCAN_PARAMS PageScanSuperAF = { 0, 32, 32, 9, 246, 254 };

```

25

/\*\*\*\*\*

\* super fine Autofocus

\* with page scan

\*\*\*\*\*/

30

\*/

```

static BOOL IonuiSuperFineAF( UV_CONFIG *pCfg, OPTICS_PARMS *pLaser, int AFCode,
int *PeakInten)

```

{

BOOL r = FALSE;

35

int iZIndex, YAmp, YStart;

1/6/95

```

/* get laser current zoom factor*/
iZIndex = lonui_GetIndexFrPixZoomEnum( pLaser->PixZoom);
/* get YStart value, and YAmpl from configuration */
if( pLaser->YDim == 512)
5   {
      YStart = pCfg->StageSpec.FineAFPScanYStart512[iZIndex];
      YAmpl  = pCfg->StageSpec.FineAFPScanYAmpl512 [iZIndex];
    }
else
10  {
      YStart = pCfg->StageSpec.FineAFPScanYStart256[iZIndex];
      YAmpl  = pCfg->StageSpec.FineAFPScanYAmpl256 [iZIndex];
    }

15  n = lon_APICmd( eLonAPIPageYAmpl, 0);
    if( n == TRUE)
        n = lon_APICmd ( eLonAPIPageControlReg, PAGE_SCAN_CTRL_SCAN_ON |
                           PAGE_SCAN_CTRL_STOP_DAC);

    if( n == TRUE)
20    n = lon_APICmd( eLonAPIPageVSync,    PageScanSuperAF.iVSync);
    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageTraceTime, PageScanSuperAF.iTrace);
    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageRetrace,  PageScanSuperAF.iRetrace);
25  if( n == TRUE)
        n = lon_APICmd( eLonAPIPageStartAmpl, YStart);
    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageIncrAmpl, PageScanSuperAF.iIncr);
    if( n == TRUE)
30    n = lon_APICmd( eLonAPIPageDecrAmpl, PageScanSuperAF.iDecr);
    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageYAmpl,    YAmpl);
    if( n == TRUE)
        n = lon_APICmd ( eLonAPIPageControlReg, PAGE_SCAN_CTRL_SCAN_ON );
35

```

1/6/95

```

if( rt == TRUE)

```

```

    rt = lonimg_QueryNetVar( FAST_Z_NODE_ADDR, AFCode, PeakInten, TRUE);

```

```

    lonui_PageScan( pLaser->YDim, pLaser->PixZoom);

```

```

5    return( rt);

```

```

    }

```

```

/*****

```

```

    * convert from enum to index

```

```

10  *****/

```

```

    */

```

```

extern int lonui_GetIndexFrPixZoomEnum( LASER_PIX_ZOOM ePixZoom)

```

```

{

```

```

    int ZoomIndex;

```

```

15

```

```

    switch( ePixZoom)

```

```

    {

```

```

        case ePixZoom12nm: ZoomIndex = 3;      break;

```

```

        case ePixZoom25nm: ZoomIndex = 2;      break;

```

```

20

```

```

        case ePixZoom50nm: ZoomIndex = 1;      break;

```

```

        case ePixZoom100nm:

```

```

            default:      ZoomIndex = 0;      break;

```

```

    }

```

```

    return ZoomIndex;

```

```

25

```

```

}

```

```

/*****

```

```

    * initialize /serup page scanner

```

```

    * page scanner is now alway on

```

```

    *****/

```

```

30

```

```

    */

```

```

extern BOOL lonui_PageScan( int iYDim, LASER_PIX_ZOOM ePixZoom)

```

```

{

```

```

    PAGE_SCAN_PARAMS      *PScanParms;

```

```

    BOOL                  rt = FALSE;

```

```

35

```

```

    int                  XAmp, YAmp, YIndex, ZoomIndex, TurPos, PageCtrlReg;

```

```

    UV_CONFIG            *pCfg    = cfg_GetLockRdConfig();

```

1/6/95

```

    if( lon_APIQuery( eLonAPITurretPosCmd, &TurPos, 0) == TRUE)
        TurPos--;
    else
        TurPos = 0;

    if( TurPos < 0 || TurPos >= (HDWR_TURRET_HOME_OBJ -1) )
        TurPos = 0;

    switch ( iYDim)
    {
    case 256: PScanParms = &PageScanParms256;
                break;
    case 512: PScanParms = &PageScanParms512;
                break;
    default:   PScanParms = NULL;
                break;
    }

    ZoomIndex = lonui_GetIndexFrPixZoomEnum( ePixZoom);

    if( PScanParms)
    {
        UImage.iYCl      = iYDim;
        UImage.XYZoom    = ePixZoom;

        n = lon_APICmd( eLonAPIPageYAmpl, 0);
        if( n == TRUE)
            n = lon_APICmd ( eLonAPIPageControlReg, PAGE_SCAN_CTRL_SCAN_ON |
                            PAGE_SCAN_CTRL_STOP_DAC);

        if( n == TRUE)
            n = lon_APICmd( eLonAPIPageVSync, PScanParms->iVSync);
        if( n == TRUE)
            n = lon_APICmd( eLonAPIPageTraceTime, PScanParms->iTrace );
        if( n == TRUE)
            n = lon_APICmd( eLonAPIPageRetrace, PScanParms->iRetrace);
    }

```

1/6/95

```

    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageStartAmpl, PScanParms->iStartValue);
    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageIncrAmpl, PScanParms->iIncr);
5    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageDecrAmpl, PScanParms->iDecr);
    if( n == TRUE)
        n = lon_APICmd( eLonAPIPageXAmpl,
                        pCfg->ConfocalSpec.XScanAmpl[TurPos][ZoomIndex] );
10    if( n == TRUE)
        {
            YAmpl = pCfg->ConfocalSpec.YScanAmpl[TurPos][ZoomIndex];
            if( iYDim == 256 )
                YAmpl = YAmpl/2;
15
            n = lon_APICmd( eLonAPIPageYAmpl, YAmpl);
        }
    n = lon_APICmd ( eLonAPIPageControlReg, PAGE_SCAN_CTRL_SCAN_ON );
    }
20    return n;
    }

```

5

**A METHOD AND APPARATUS FOR PERFORMING  
AN AUTOMATIC FOCUS OPERATION**

10

15

**Christopher R. Fairley  
Timothy V. Thompson  
Ken K. Lee**

20

**APPENDIX G**

25

**M-2464-1P US**



5

**Ultrapointe Model 1000  
Laser Imaging System  
System User's Manual**

10

15

Copyright 1993, 1994  
Ultrapointe Corporation  
163 Baypointe Parkway  
San Jose, CA 95134

(408) 894-7080

20

P/N 000621

Rev. 1.3

**DRAFT**

Revised: 3 June, 1994

5

### Foreword

10

Revision 1.3 of the UV1000 software incorporates all changes and features provided in various unreleased versions through 1.2T(u), and issued before 18 May, 1994. This manual revision 1.3 documents all those features. In addition, customer preventative maintenance procedures have been incorporated as Appendix E. Finally, while the manual refers to the LIS 1000 hardware, the manual applies also to the model LIS 500. The manual assumes the LIS 500 has been ordered with full 3-D capability, but if this is not the case, then those features not ordered will not be available.

20

## SAFETY

All user interaction with the instrument is via the front panel switches, the mass storage devices in the drive bay, or the keyboard and display console. The user must be familiar with all external system components, per the Safety and Operating sections of this manual. Servicing is to be performed by factory-authorized personnel only.

The Ultrapointe Model 1000 has been designed to minimize the hazards presented to the user during normal operation of the instrument. All servicing is to be performed by factory-authorized personnel only.

### Laser Safety

In accordance with the Regulations for the Administration and Enforcement of the Radiation Control for Health and Safety Act of 1968 (applicable to laser products), a Laser Product Report for this instrument has been forwarded to the Center for Devices and Radiological Health (CDRH). This instrument falls under Class I definitions because there is no user access to laser radiation.

**CAUTION:** Use of controls or adjustments, or performance of procedures other than those specified herein may result in hazardous radiation exposure.

The protective housing of the instrument should be kept intact during normal operation. Internal to the instrument is a Class IIIb argon-ion laser rated at 25 milliwatts with a maximum possible output of 500 milliwatts. Removal of any part of the upper protective housing can expose the user unnecessarily to hazardous laser radiation. Danger

labels reproduced below indicate covers that provide access to hazardous laser light exposure.

5

**DANGER - Laser radiation when  
open and interlock defeated.  
AVOID DIRECT EXPOSURE TO BEAM**

10

**DANGER - Laser radiation when  
open. AVOID DIRECT EXPOSURE  
TO BEAM**

## **Safety, Continued**

### **High Voltage**

Internal to the instrument are components using voltages up to 1000 VDC. Access to and servicing of all internal components is by factory-authorized personnel only.

5

### **Moving Parts**

As part of routine instrument operations, various parts of the machine make movements automatically and when commanded by the user. The user should take care to keep body parts and clothing away from the robot area. The user should also ensure that the current robot operation is complete before adding or removing a cassette from the cassette locator plate.

10

Much of the instrument is attached to an internal frame that is mounted to the external frame and protective covers via pneumatic vibration isolators. The "floating" internal frame can therefore move freely with respect to the covers and creates possible pinch points between floating components and the covers or external frame. Externally accessible floating components include the cassette platform, the robot, and the flat finder.

15

Inside the instrument, are many moving parts which could create possible hazards, including motors, shutters, and an X/Y/Z stage. Access to components inside the machine with system power on is restricted to factory-authorized service personnel.

20

Moving parts on the outside of the system include the robot handler, pre-aligner, and wafer shutter door, as well as the floating optics platform. The user should take care to keep body parts and clothing clear of all moving parts at all time.

25

### Table of Contents

	1. Introduction . . . . .	307
	2. Principles of operation . . . . .	311
5	3. Safety . . . . .	315
	3.1. OPERATING SAFETY: . . . . .	315
	3.1.1. Laser Safety: . . . . .	315
	3.1.2. High Voltage: . . . . .	316
	3.1.3. Moving Parts: . . . . .	316
10	4. Software Installation . . . . .	318
	5. System Usage . . . . .	319
	5.1. Introduction . . . . .	319
	5.1.1. Basic Controls . . . . .	319
	5.1.1.1. Mouse . . . . .	319
15	5.1.1.2. Menu Bars . . . . .	320
	5.1.1.3. Option Menus . . . . .	320
	5.1.1.4. Buttons . . . . .	320
	5.1.1.5. Slider Bars . . . . .	321
	5.1.1.6 Entry Field . . . . .	322
20	5.2. Data files and organization . . . . .	322
	5.2.1. Volume and Surface . . . . .	323
	5.2.2. Wafer map File . . . . .	324
	5.2.2.1 Defect information . . . . .	324
	5.2.2.2 Wafer information . . . . .	325
25	5.2.2.3 Identification . . . . .	325
	5.2.3. Recipe . . . . .	325
	4.2.3.1 Screen layout . . . . .	325
	4.2.3.2 Laser scan parameters . . . . .	326
	4.2.3.3 Defect Sort Criterion . . . . .	326
30	5.2.4. Bit map files: . . . . .	326

	5.2.5. Defect Code Table files: . . . . .	326
	5.3. Top Menu Bar . . . . .	327
	5.3.1. File: . . . . .	327
5	5.3.1.1. Image: . . . . .	327
	5.3.1.2. Recipe: . . . . .	327
	5.3.1.3. Wafer Map: . . . . .	328
	5.3.1.4. Defect Code: . . . . .	329
	5.3.1.5. Exit: . . . . .	330
	5.3.2. Screen: . . . . .	330
10	5.3.3. Setup: . . . . .	330
	5.3.3.1. LaserScan Parms: . . . . .	330
	5.3.3.2. Deskew Options . . . . .	331
	5.3.3.3 Edit Wafer Map: . . . . .	331
	5.3.3.4. Defect List Sort: . . . . .	333
15	5.3.3.5 Rearranging Wafers . . . . .	336
	5.3.4. User: . . . . .	336
	5.3.4.1. Login: . . . . .	337
	5.3.5. Config: . . . . .	337
	5.3.5.1. SECS: . . . . .	337
20	5.3.5.2. Printer: . . . . .	337
	5.3.5.3 Misc. . . . .	337
	5.3.6. Maintenance: . . . . .	340
	5.3.6.1. Manual stage contr . . . . .	340
	5.3.6.2. Manual robot control: . . . . .	340
25	5.3.6.3. File Transfer: . . . . .	340
	5.3.6.4. UserAcct: . . . . .	342
	5.3.6.5. LON Diagnostics: . . . . .	344
	5.3.6.6. Calibration: . . . . .	344
	5.3.7. PrtScreen: . . . . .	344
30	5.4. Window and Screen . . . . .	344
	5.4.1. Common window attributes . . . . .	345

	5.4.1.1. Window sizes. . . . .	346
	5.4.1.2. Swapable to other window sites on-screen. . . . .	346
	5.4.1.3. Global and local data. . . . .	347
5	5.4.1.4. Inter-window data exchange. . . . .	347
	5.4.2. Window types . . . . .	348
	5.4.2.1. Cassette Window . . . . .	348
	5.4.2.2. Wafer Window . . . . .	351
	5.4.2.3. Laser Window . . . . .	353
10	5.4.2.4. Microscope window . . . . .	361
	5.4.2.5. 2D Image Window . . . . .	361
	5.4.2.6. 3D Surface Image Window . . . . .	363/2
	5.4.2.7. Library Image Window . . . . .	363/6
	5.4.2.8. Text and status window . . . . .	363/7
15	5.4.3. Screen . . . . .	363/13
	5.4.3.1. Switching to a different screen . . . . .	363/13
	5.5. User privilege . . . . .	363/14
	5.6. Operator Mode . . . . .	363/14
	5.6.1. System Startup . . . . .	363/14
20	5.6.1.1. Laser imaging system Login: . . . . .	363/15
	5.6.1.2. Hardware initialization: . . . . .	363/15
	5.6.2. Wafer loading, cassette window . . . . .	363/16
	5.6.3. Deskew, Text Status window . . . . .	363/16
	5.6.4. Defect location, Wafer window . . . . .	363/17
25	5.6.5. Image Acquisition, Laser window . . . . .	363/18
	5.6.5.1. White light only. . . . .	363/18
	5.6.5.2. White light and laser . . . . .	363/18
	5.6.6. Screen Printing . . . . .	363/20
	5.6.7. Image storage . . . . .	363/20
30	5.6.8 Defect Classification . . . . .	363/20



	5.7. Engineering Mode . . . . .	363/21
	5.7.1. User type and Operator access control . . . . .	363/21
	5.7.2. Recipe preparation . . . . .	363/22
	5.7.2.1. Screen setup and sequencing . . . . .	363/22
5	5.7.2.2. Laser parameter control . . . . .	363/23
	5.7.2.3 Sort criteria . . . . .	363/23
	5.7.3. Manual Robot control . . . . .	363/23
	5.7.4. Manual stage control . . . . .	363/24
	5.7.4.1. Enable/Disable joystick, all axes . . . . .	363/24
10	5.7.4.2. Report Position/encoder count . . . . .	363/24
	5.7.4.3. Report sample thickness . . . . .	363/24
	5.7.4.4. XYstage burn-in . . . . .	363/24
	5.7.4.5. Stop . . . . .	363/24
	5.7.4.6. XYTo Load/unload . . . . .	363/25
15	5.7.4.7. Move, absolute or relative . . . . .	363/25
	5.7.4.8. Set speed . . . . .	363/25
	5.7.4.9. Set acceleration . . . . .	363/25
	5.7.4.10. Move to Home . . . . .	363/25
	5.7.5. Hardware diagnostics . . . . .	363/25
20	5.7.6. System Calibration . . . . .	363/26
	5.7.6.1. Laser scan XY dimension calibration . . . . .	363/26
	5.7.7 System configuration . . . . .	363/26
	5.7.7.1 Screen/Laser saver timer interval . . . . .	363/26
25	5.7.7.2 File DeskewDie Offset . . . . .	363/26
	5.3.5.3.3 Flat angle offset( deg.) . . . . .	363/27
	5.3.5.3.4 White light mode . . . . .	363/27
	5.3.5.3.5 Next die position . . . . .	363/28
30	5.7.8. System Startup . . . . .	363/28

	5.7.8.1. Operating system and windowing system . . . . .	363/28
	5.7.8.2. IRIX Computer user autologin and autostart: . . . . .	363/29
5	6. Troubleshooting . . . . .	363/30
	6.1. Robot operation . . . . .	363/30
	6.2. Stage Operation . . . . .	363/30
	6.3. Power fail Recovery . . . . .	363/30
	6.3.1. Wafer recovery . . . . .	363/31
10	6.3.2. System component reset . . . . .	363/31
	6.4. LON diagnostics dialog box . . . . .	363/31
	6.4.1 Wafer door interlock. . . . .	363/31
	6.4.2 Critical parameters . . . . .	363/31
	6.5. System Service . . . . .	363/31
15	6.5.1. Software configuration . . . . .	363/31
	6.6. System Messages and Response . . . . .	363/31
	6.6.1. Caution/Information . . . . .	363/32
	6.6.2. Error Condition and Recovery . . . . .	363/32
	6.6.2.1 Floppy disk not responding . . . . .	363/32
20	6.6.2.2 Printer not printing . . . . .	363/32
	7. Applications Information . . . . .	363/34
	7.1. Wavelength Selection . . . . .	363/34
	7.2. Image interpretation . . . . .	363/34
	7.2.1. Surface reflectivity . . . . .	363/34
25	7.2.2. Limits of performance . . . . .	363/35
	7.3. Complex surfaces . . . . .	363/36
	7.3.1. Transparent coatings . . . . .	363/36
	7.3.2. Multilayer structures . . . . .	363/36
	7.3.3. High contrast variations . . . . .	363/36
30	7.3.4. Vertical Surfaces . . . . .	363/37

	8. System Maintenance . . . . .	363/38
	8.1. System software maintenance . . . . .	363/38
	8.1.1 Setting Date and Time . . . . .	363/38
	8.1.2. System shutdown . . . . .	363/38
5	8.2. Maintenance and Laser Safety . . . . .	363/38
	8.3. Customer Preventative Maintenance Procedures	363/39
	8.4. Extraordinary maintenance . . . . .	363/39
	8.4.1. Wafer removal . . . . .	363/39
	8.4.2. Other . . . . .	363/39
10	9. Appendices . . . . .	363/40
	Appendix A . . . . .	363/40
	A.1 To Install the Indigo System Software . . . . .	363/40
	A.2 To Install the Ultrapointe Application Software .	363/41
	A.3 InstallApplic Script . . . . .	363/42
15	A.3.1 Important Files . . . . .	363/43
	A.4 To Install the Peripherals . . . . .	363/44
	A.4.1 Tape . . . . .	363/45
	A.4.2 Printer (Kodak XL7700) . . . . .	363/45
	A.4.3 Printer (Codonics NP600) . . . . .	363/46
20	A.4.4 Floppy . . . . .	363/47
	A.4.5 Network to Other Hosts . . . . .	363/47
	A.4.6 Serial Ports . . . . .	363/47
	A.5 To Install Incremental Software Upgrades . . . . .	363/47
	A.6 To Archive Data Files onto QIC-150 1/4" Tapes	363/48
25	A.7 To Prepare an InstallSystem Tape . . . . .	363/48
	A.8 PrepareSystemTape Script . . . . .	363/49
	A.9 To Prepare an InstallApplication Tape . . . . .	363/49
	A.10 PrepareApplicTape Script . . . . .	363/49
	A.11 To Format the HardDisk . . . . .	363/49
30	A.12 To Build a New Kernel (Operating System) . .	363/50

	A.13 BuildNewKernel Script . . . . .	363/50
	Appendix B . . . . .	363/52
	B.1. Introduction: . . . . .	363/52
	B.2. Screen spec: . . . . .	363/52
5	B.3. User_spec: . . . . .	363/53
	B.4. Confocal_spec: . . . . .	363/54
	B.5. Robot spec: . . . . .	363/57
	B.6. Stage_spec . . . . .	363/57
	B.6.1. Procedure to determine stage_af_zero .	363/61
10	B.6.2. Procedure to determine stage_af_inten .	363/61
	B.6.3 Procedure to determine	
	stg_orthogonality_deg . . . . .	363/62
	B.6.4 Procedure to determine	
	stage_angle_tangent . . . . .	363/63
15	B.7. Objective Specs . . . . .	363/63
	B.8. Lon Specs . . . . .	363/67
	B.9. System Specs . . . . .	363/67
	Appendix C: IRIS INDIGO BASIC COMMANDS . . . . .	363/70
	C.1. Introduction . . . . .	363/71
20	C.2. Login in . . . . .	363/71
	C.3. Help . . . . .	363/71
	C.4. File manipulation . . . . .	363/71
	C.5. File information . . . . .	363/72
	C.6. Directory manipulation . . . . .	363/72
25	C.7. Find and search . . . . .	363/73
	C.8. System information . . . . .	363/73
	C.9. User information . . . . .	363/73
	C.10. Command information . . . . .	363/74
	C.11. Process control . . . . .	363/74
30	C.12. Time/Date . . . . .	363/74
	C.13. File conversion . . . . .	363/75

	C.14. Piping commands . . . . .	363/75
	C.15. Sourcing and Sinking . . . . .	363/75
	C.16. Miscellaneous Commands . . . . .	363/75
	C.17. System shutdown procedure . . . . .	363/76
5	C.18. Text editing using Jot . . . . .	363/76
	C.19. DOS Floppy Disk . . . . .	363/77
	C.20. Tapes . . . . .	363/77
	C.21. Adding users setup . . . . .	363/78
	C.22. Networking setup . . . . .	363/78
10	C.23 IMPORT directory setup . . . . .	363/78
	C.23-1 Network setup on the remote machine .	363/78
	C.23-2 Add the remote machine to LIS Network	
	setup . . . . .	363/79
	C.23-3, Make the IMPORT directory . . . . .	363/79
15	C.23-4 Make the IMPORT directory accessible	
	by the remote machine . . . . .	363/79
	C.23-5, Tell the system start using	
	this new table: . . . . .	363/79
	C.24 REMOTE/REMOTE2 directory setup . . . . .	363/79
20	C.24-1 Network setup on remote machine . . .	363/79
	C.24-2 Add the remote machine to LIS Network	
	setup . . . . .	363/79
	C.24-3 Make the LIS REMOTE directory . . . . .	363/80
	C.24-4, Enable network connection . . . . .	363/80
25	C.25. File transfer . . . . .	363/80
	C.26. Logging on remotely . . . . .	363/81
	C.27. Home Directory . . . . .	363/81
	Appendix D: Support for KLA 20xx Defect File . . . . .	363/82
	Appendix E: Customer Preventative Maintenance	
30	Procedures . . . . .	363/83
	E.1. Scope: . . . . .	363/83

	E.2. Safety Precautions: . . . . .	363/83
	E.2.1 Laser Safety: . . . . .	363/83
	E.2.2 High Voltage: . . . . .	363/86
	E.2.3 Moving Parts: . . . . .	363/86
5	E.3. PM Schedule: . . . . .	363/87
	E.3.1 System Startup . . . . .	363/87
	E.3.2 System Shutdown . . . . .	363/88
	E.3.3 Lockup Recovery . . . . .	363/91
	E.3.3.1 Restarting the Ultrapointe	
10	Application . . . . .	363/92
	E.3.3.2 Executing a Five Finger Reset . . . . .	363/93
	E.3.3.3 Rebooting the Computer . . . . .	363/93
	E.3.3.4 Cycling Power to the Instrument . . . . .	363/94
	E.3.4 Wafer Removal . . . . .	363/94
15	E.3.4.1 Wafer Removal Using the	
	Application Software . . . . .	363/95
	E.3.4.2 Manual Wafer Removal . . . . .	363/96
	E.3.5 System Wipe Down . . . . .	363/97
	E.3.6.2 System Restore . . . . .	363/98
20	E.3.7 Chuck Cleaning . . . . .	363/98
	E.3.8 Laser Lifetime and Interlock Inspection . . . . .	363/99
	E.3.9 Robot Alignment Check . . . . .	363/101
	E.3.9.1 Robot to Cassette	
	Alignment Check . . . . .	363/102
25	E.3.9.2 Robot to Pre-Aligner	
	Alignment Check . . . . .	363/102
	E.3.9.3 Robot to XYZ Stage	
	Alignment Check . . . . .	363/103
	E.3.10 Changing the Halogen Bulb . . . . .	363/103
30	E.3.11 Ultrapointe Factory PM . . . . .	363/106
	E.4. Preventative Maintenance Log Sheet . . . . .	363/107
	E.5. Problem Report Worksheet . . . . .	363/108

Ultrapointe Model 1000Laser Imaging SystemSystem User's ManualRev. 1.3

5

1. Introduction

The Ultrapointe Model 1000 Laser Imaging System is designed to replace and out perform optical microscope review stations now  
10 utilized in the Semiconductor Fab environment to examine optical anomalies previously detected by Patterned Wafer Defect Scanning Instruments and to perform a variety of microscopic inspection functions. Significantly, the Ultrapointe system is the first defect review tool whose optics and functionality, have been designed  
15 explicitly for efficient performance of the dedicated production review and inspection tasks.

The Ultrapointe Model 1000 produces a three dimensional image, which may be rotated or tilted or shaded with correct perspective  
20 maintained with simple image rendering techniques without ever moving or disturbing the sample and which may be stored and recalled for later viewing. Also, the unit provides quantitative dimensional information (x, y, and z). The Ultrapointe Laser Imaging System has an ability the SEM cannot match, subsurface viewing of defects lying  
25 beneath dielectric layers. Combined with 3-D analysis software, the user is able to examine cross sections of the defect and surrounding material, and to assess the impact on the circuit layers.

The operator views the image on a CRT, with comfortable ergonomics,  
30 and the wafers are not be exposed to operator contamination or airflow. Unlike SEMs, the laser imaging system operates in air with

class 1 cleanroom compatibility. While possessing metrological capabilities, the Laser Imaging System is not intended to compete with metrology systems now on the market, which also use confocal technologies, but are not designed for defect evaluation.

5

Cassettes of wafers (typically ranging from 3 inches (75mm) to 200mm diameter) are presented to the machine. At the same time, a file of data from the wafer defect scanner, specifying the coordinates of the detected defects is transferred to the review station computer, either by diskette or other media, or by communication via a link or network such as Ethernet, RS232, etc.

10

A precision, high reliability robotic wafer handler removes wafers from the cassette and performs a prealignment step, sensing the notch and/or flat(s) on the wafer. The wafer is then loaded onto the optical unit's XYZ stage and is secured on the vacuum chuck. The operator accomplishes fine alignment ("de-skew") by lining up the visible light microscope field of view with either etched fiducial marks or other pre-specified structures on the wafer surface. Then, the unit can accurately translate any specified location on the wafer into the field of view with an accuracy of a few microns.

15

20

The operator examines each defect image as the tool presents it. The operator has three modes from which to select viewing. This examination can be either with the white-light conventional microscope optics, or in a real time laser scanning mode, or simultaneously with both laser and white light through the same optics. In white light, the operator can select from one of several objective lenses, varying effective magnification of the image. (The laser image scales simultaneously with the white-light image.) The white-light illumination may be polarized, and a cross polarized filter inserted in

25

30



the observed light , causing optically active objects (such as quartz) to appear brighter. Optional bright-field-dark field objectives and illumination may be used to help locate particles and other defects on the wafer surface.

5

If the defect is not in view, or the operator wishes to examine a larger or different area, a joystick controller allows the operator to "cruise" the wafer.

10 Automatic operation features include autofocus, auto gain, and auto ranging of the vertical scanning. Engineering mode allows recipe setup to pre-specify operating parameters for use by operators working on a specific process level and product. Utilities are available in pop-up dialog boxes allowing manual control of the stage, polling of stage  
15 variables, plus robot manual control (e.g. allowing movement of a wafer from the flat finder or stage after a power failure). Diagnostics are called up via a pop-up window displaying all LON nodes and system variables, plus providing direct control of system functions.

20 By utilization of "recipes", the system may be pre-programmed for automatic operation. on a specific type of wafer, and at a specific process step for that wafer. The laser wavelength, and power, to be used may be selected. The number of slices of data and their spacing (in nanometers) may be specified. Finally, autofocus may be specified,  
25 and the offset in z (vertical) from the autofocus position to the ideal viewing position may also be preset.

A second major use of conventional microscope stations has been the more general inspection function, where pre-selected sites on a wafer  
30 are inspected for efficacy of a previous process step. The Ultrapointe Laser Imaging System is directly usable in this application, which faces

Ultrapointe LIS 1000 User's Manual, Rev. 1.3

Revised 3 June, 1994

exactly the same dilemma of defect imaging: Decreasing size of objects of interest, a lack of resolution, and the need for 3-D imaging.

- 5 Viewing through transparent layers allows determination of the vertical site of defects, or provide some specialized inspection. For example, stress voids in metal layers below dielectric layers may be viewed. The depth of metal plugs in glass insulators may also be seen.

## 2. Principles of operation

The system functional diagram (Figure 23) illustrates the general principles of operation. The Laser Imaging System uses the basic principles of confocal microscopy, where the illuminating light passes through a spatial filter (pinhole), and the image of this pinhole is then cast on the sample to be viewed. The light scattering from the sample returns through the optics to the pinhole, but only light from the focal plane of the imaging lens returns through the pinhole.

The light source is a multi-line laser (1) which produces polarized light at several discrete wavelengths which passes through a selectable filter (2) to isolate a single laser line or more than one line. A polarizing beam splitter (3) preferentially reflects light only of the proper polarization and directs it to the spatial filter (4). The spatial filter consists of optics which image the laser light on a pinhole whose diameter has been selected to re-image through the downstream optics and objective lens (5) producing a diffraction-limited spot on the sample (6). The x-y beam scanner (7) consists of two mirrors which can oscillate their angle with respect to the beam in a manner to scan a raster pattern in space, when operated in conjunction with the scan lens (8). The raster scanned image of the pinhole is produced between the scan lens and beam splitter cube (9). In our first embodiment of this design, the slow scan is varied at 13 to 26 hertz, while the fast unit operates at 8 kHz. A raster scan of 256 or 512 lines is produced at approximately 26 or 13 frames per second, and is imaged at the back focal plane of the tube lens (11).

A quarter wave plate (10) is positioned to convert the linearly polarized laser light to circularly polarized light. The tube lens (11) works with the objective (5) to de magnify the image of the raster scanned pinhole

image and project it on the sample. A portion of this light is scattered back into the objective and returns through the optic train. As it passes through the quarter wave plate (10), is converted to light linearly polarized at right angles to the light emitted by the laser. Thus  
5 when the light reaches the polarizing beam splitter (3), it passes through and reaches the photo detector (12).

As the raster is scanned, the surface data processor (13) keeps track of the returning intensity, and a two-dimensional intensity map of the  
10 sample is made. This image passes through the system computer (14) and is produced on the video monitor (15) as a digital image.

To obtain a three dimensional image, the height of the sample is varied by stepping the sample height small distances vertically using the fine  
15 z stage (16) between raster scans and keeping track of the vertical position for each horizontal point on the scan which produces the maximum intensity. This results in a 3-D map of the surface.

White light (conventional) microscope images are produced  
20 simultaneously with the live laser image by a video camera (19) which views the sample in white light emitted by the illuminator (20), an inserted in the optical path by the beam splitter (21). A filter which blocks the laser line in use but passes broad bands of light otherwise prevents the laser from saturating the image with reflected light.

25 A turret (23) allows objectives of different magnifying power to be used.

The optics unit must work in conjunction with a precision x-y stage  
30 (18) and "coarse" z stage (17). This module must accommodate from 3 inch (75 mm) to 200 mm diameter semiconductor wafers, working

also with the robotic handler to load and unload the wafers with high reliability (< 1 wafer drop per million transfers), and class 1 compatible cleanliness.

- 5      The electronics package (Figure 24-1, 24-2, 24-3) contains all analog and digital electronics, plus power supplies, for complete operation of the tool. The unit operates on 220 v.(200-240 v. nominal), 50/60 Hz single phase electric power (or the European and Japanese equivalents).

10

The Surface Data Processor (SDP) (1) is a proprietary Ultrapointe design which interfaces with the photo detector and is synchronized with the scanner electronics, and fine z stage control, producing an x, y, & z map of intensity which can either be stored directly in the  
15      computer memory, or processed to extract a surface image on the fly. This interfaces through the SDP Interface (2) to the system computer (3).

20

The system computer (3) is a high speed RISC graphics workstation, capable of handling concurrent tasks of robot functions, stage motion, operator interface, and optics control, while also performing image processing functions. In addition, it functions with a windowing user interface and high resolution color graphics.

25

The X, Y, and coarse Z stage controllers (4) communicate with the system computer via an RS-232C interface (5), as do the robot and pre-aligner controllers (6).

30

The balance of the system electronic functions communicate through a Local Operating Network (LON) Interface (7) built on the same interface slot as the RS-232C interface. The LON itself (8) is a pair of

wires that plug into each node serially around the system. Each node contains a local processor and firmware for LON communications, self diagnosis, and local operation of certain functions.

5 All user interface is via the system monitor, mouse/trackball, joystick/controller, and system keyboard. Through these controls and the windowing software, the operator can set up, program and operate the complete system, wafer selection and handling, defect editing and selection, automatic and/or manual wafer loading, defect classification,  
10 etc. In addition, he can utilize the tool as a microscope, selecting views, translating and focusing on details, changing magnification and zooming in on a sample. The operator console may optionally be remotely mounted. In addition, image processing and analysis functions may be performed from the console.

15

The overall system software is designed for operation on the review station in both operator and engineering mode, optionally using defect files supplied by various defect detection tools in a variety of formats. Both modes are password protected.

### **3. Safety**

5 All user interaction with the instrument is via the front panel switches, the mass storage devices in the drive bay, or the keyboard and display console (see Figure 22). Servicing is to be performed by factory-authorized personnel only.

#### **3.1. OPERATING SAFETY:**

10 The Ultrapointe Model 1000 has been designed to minimize the hazards presented to the user during normal operation of the instrument. All servicing is to be performed by factory-authorized personnel only.

##### **3.1.1. Laser Safety:**

15 In accordance with the Regulations for the Administration and Enforcement of the Radiation Control for Health and Safety Act of 1968 (applicable to laser products), a Laser Product Report for this  
20 instrument has been forwarded to the Center for Devices and Radiological Health (CDRH). This instrument falls under Class I definitions because there is no user access to laser radiation.

25 **CAUTION:** Use of controls or adjustments, or performance of procedures other than those specified herein may result in hazardous radiation exposure.

30 The protective housing of the instrument should be kept intact during normal operation. Internal to the instrument is a Class IIIb argon-ion laser rated at 25 milliwatts with a maximum possible output of 500

milliwatts. Removal of any part of the upper protective housing can expose the user unnecessarily to hazardous laser radiation. Danger labels reproduced below indicate covers that provide access to hazardous laser light exposure.

5

**DANGER - Laser radiation when  
open and interlock defeated.  
AVOID DIRECT EXPOSURE TO BEAM**

10

**DANGER - Laser radiation when  
open. AVOID DIRECT EXPOSURE  
TO BEAM**

15

### 3.1.2. High Voltage:

Internal to the instrument are components using voltages up to 1000 VDC. Access to and servicing of all internal components is by factory-authorized personnel only.

20

### 3.1.3. Moving Parts:

As part of routine instrument operations, various parts of the machine make movements automatically and when commanded by the user.

25

The user should take care to keep body parts and clothing away from the robot area. The user should also ensure that the current robot operation is complete before adding or removing a cassette from the cassette locator plate.

30

Much of the instrument is attached to an internal frame that is mounted to the external frame and protective covers via pneumatic



vibration isolators. The "floating" internal frame can therefore move freely with respect to the covers and creates possible pinch points between floating components and the covers or external frame. Externally accessible floating components include the cassette platform, the robot, and the flat finder.

Inside the instrument, are many moving parts which could create possible hazards, including motors, shutters, and an X/Y/Z stage. Access to components inside the machine with system power on is restricted to factory-authorized service personnel.

Moving parts on the outside of the system include the robot handler, pre-aligner, and wafer shutter door, as well as the floating optics platform. The user should take care to keep body parts and clothing clear of all moving parts at all time.

#### 4. Software Installation

Each Ultrapointe Model 1000 is pre-configured at the factory with all necessary system hardware and software. The system is ready to run  
5 after on-site installation by Ultrapointe Corp. No software installation is required.

However, for the purposes of system archiving and upgrades, each  
10 Ultrapointe Model 1000 is supplied with two 1/4 inch QIC-150 formatted tapes. One tape contains the generic computer operating system software, and the other contains the basic LIS software with some sample image files. See Appendix A for system backup, restore, and hardware upgrade procedures.

## 5. System Usage

### 5.1. Introduction

5 The Ultrapointe Model 1000 has a graphical user interface, and all functions can be accessed via "point-and-click" with the mouse. Within the LIS software, a screen is the collective look of the whole monitor screen, and it consists of a title bar, a top-menu bar and a collection of windows. The windows provide all primary functions, and  
10 the top menu bar provides access to the rest of the functions that are used infrequently. Windows are movable and sizable rectangles in various screens. A window encapsulates a major primary function and provides the necessary user interface. User can specify the arrangement of the windows in a given screen. Each item on the top  
15 menu bar pulls down a menu, and each item in the menu performs a specific task. The Model 1000 supports three levels of users: Normal user, administrator (engineering level), and service engineer. Normal users have no access to diagnostics and calibration functions.

#### 20 5.1.1. Basic Controls

⇒ See Figure 25 for Basic Controls

##### 5.1.1.1. Mouse

25 The red arrow on the screen will follow the movement of the mouse. By moving the mouse, one may position the red arrow over a widget such as a push-button, a menu bar item, an option menu, a radio button, a toggle button, or a slider  
30 bar.

Then, by clicking (pressing and releasing) the left button of the mouse, one activates the action associated with that particular widget.

5      5.1.1.2. Menu Bars

By clicking on a menu bar item, an associated pull down menu will appear. Then the user can activate one of the pull down menu items by clicking on it. Alternatively, one may press  
10      down on a menu bar item and not release the mouse button until the screen arrow is moved onto the desired pull down menu item.

15      5.1.1.3. Option Menus

An option menu is similar to a menu bar item, but it doesn't have to appear along the top of the screen, and it's the last selected item that will remain displayed. By clicking on an option menu, an associated pull down menu will appear.  
20      Then, the user can activate one of the pull down menu items by clicking on it. Alternatively, one could press down on an option menu and not release the mouse button until the screen arrow is moved onto the desired pull down menu item. The selected item will then remain visible as the other items in the  
25      pull down menu disappear.

#### 5.1.1.4. Buttons

##### 5.1.1.4.1. push-button

5 By clicking on a push-button, one would activate the action associated with the push-button.

##### 5.1.1.4.2. Toggle

10 By clicking on a toggle button, one would turn the switch on, if it was originally off, or turn the switch off, if it was originally on.

##### 5.1.1.4.3. Radio

15 Radio buttons are a group of toggle buttons where only one of the buttons can be 'on' at any one time. Therefore, by clicking on one button, the other buttons will all be turned off.

##### 20 5.1.1.5. Slider Bars

One may move the slider bar within a scale by 2 methods:

25 1. By clicking on and dragging the slider bar to the desired value.

2. By clicking on the area in between the slider bar and arrow button. This will increment/decrement the slider bar by a value of 10. (Most scales have 10 as a default)

30

#### 5.1.1.6 Entry Field

5 An **Entry Field** is a rectangular area where user can enter text and numbers. An Entry field accepts keyboard input only if it has the input focus, i.e.: the border of the entry field is a thick dark frame. Clicking within the entry field make the entry field to have input focus.

### 5.2. Data files and organization

10 The Model 1000 stores data( surface, volume, wafer map, recipe, and bit map image) in one of the following directories:

15	System	Accessible by all users
	User	Normal user can access his own user directory
	FLOPPY	DOS formatted 3.5" 1.44 megabyte floppy disk
	OPTICAL	3.5" 128 megabyte read-write optical disk
	IMPORT	Files to be imported to LIS in this directory. This directory is accessible from another machine. Linked via NFS( Ethernet, TCP/IP)
	REMOTE	Physically located in another machine. Linked via NFS
	REMOTE2	Physically located in another machine. Linked via NFS

20 A directory is created for each user, and this is the default directory where all data is stored for that user. A normal user can access only his own user directory, while an administrator can access all user directories. The **System**, **FLOPPY**, **OPTICAL**, **IMPORT**, **REMOTE**, and **REMOTE2** directories can be accessed by all users. All file names are

- keyed with a file extension. For example: all surface file names end with **.sur** and all volume files have names ended with **.vol**. In the LIS, a directory listing shows only those file names with extension that matches the file type selected. The **FLOPPY** directory allows reading from and writing to MS-DOS formatted 3-1/2 inch floppies. Normally, the **FLOPPY** directory lists only those file names that match file type selected. If wafer map type is selected, however, all file names in the **FLOPPY** directory are listed.
- 10 The **REMOTE**, **REMOTE2** and **IMPORT** directories are designed for interfacing to a data base host computer, for example, the Tencor Swift Station or the KLA 2550/2551. Some host computers require two directories for file transfer: The LIS receives a map file in the **IMPORT** directory, and the LIS sends the reviewed map file back to the
- 15 host computer in the **REMOTE** directory. Other host computers need only one directory for both sending and receiving, and this directory is usually located in the host computer and accessed from LIS as the **REMOTE** directory. Generally, the **IMPORT** and **REMOTE** directories list only those file names that match the file type selected. If wafer
- 20 map type is selected, however, all file names in the **IMPORT**, **REMOTE** or **REMOTE2** directory are listed. As part of the installation, the **REMOTE**, **REMOTE2** and **IMPORT** directories are configured to the correct network addresses. If the remote directories are not configured, or if the remote computer is not available, the **REMOTE** or
- 25 **REMOTE2** choices are not displayed. The **IMPORT** directory is made accessible to another machine.

### 5.2.1. Volume and Surface

- 30 The basic image from laser scanning is a two-dimensional image at the focus plane. A volume data set is a series of 2D

images taken at regular vertical intervals or "slices". From the volume data set, the system extracts the surface. The Model 1000 allows the user to view, manipulate, and store the surface as well as the volume data sets. Volume data set size is generally large (512 by 512 by 64 slices means 16 megabytes), and the corresponding surface data is generally smaller (512 by 512 by 2 bytes = 512 kilobytes.) Both volume and surface files are stored with their history records: date and time of image, wafer map used, operator, product ID, wafer slot, and defect number. If the image file is loaded from the top menu bar, the named wafer map is also loaded. Loading images from individual 2D or 3D windows will not update the wafer map.

#### 5.2.2. Wafer map File

A wafer map file contains information relevant to the whole lot of wafers. Generally a wafer map file can handle a lot of two cassettes, and each cassette can have up to 25 wafers. The Ultrapointe 1000 supports many standard wafer map formats. Including Tencor<sup>tm</sup>, KLA<sup>tm</sup>, and Inspec<sup>tm</sup>, among others. Typical contents of wafer map files are:

##### 5.2.2.1 Defect information

Each defect is specified in defect size, intensity, location, and classification code. Defect information is displayed in a list in the defect list window, and graphically in the wafer window.



#### 5.2.2.2 Wafer information

5 Die size, flat orientation and wafer size are displayed in the wafer window. Deskew locations are displayed in the deskew window.

#### 5.2.2.3 Identification

10 Product ID, lot ID and process ID are displayed in the status window.

#### **5.2.3. Recipe**

15 Recipes are used to customize and automate many of the instrument's functions. By loading a properly specified recipe, the user can setup the LIS with the exact screen layout and scan parameters, as well as defect sort criteria, without any user input. When a user logs in, the system will search the user's directory for a recipe with the user's name. If the  
20 recipe is found, it is loaded automatically. A recipe stores:

#### 4.2.3.1 Screen layout

25 Different steps in the defect review process require different screen layouts for optimal results. Instead of manually arranging the screens, the user can design a number of screens in a recipe and select those screens as needed. Each recipe's screen is given a descriptive name, and the list of recipe's screen is accessible from an option menu at the lower  
30 right corner of the monitor. In addition to specifying the

layout of screen's windows, the user can also specify the options within some of the windows, e.g.: 3D 2D and status.

#### 4.2.3.2 Laser scan parameters

5

Most of the optics parameters can be specified in the recipe: Camera intensity, laser intensity, number of slices, scan dimensions, bright and dark field selection, camera filter, laser scan xy resolution, laser line, sample thickness, laser display color option, as well as stage autofocus intensities and stage-center-offsets. These parameters are loaded into the LIS when the recipe is loaded.

10

#### 4.2.3.3 Defect Sort Criterion

15

Defects from a wafer map file can be filtered with a set of sort criterion so that only the desired defects are reviewed. The sort criterion can be specified in a recipe.

20

#### **5.2.4. Bit map files:**

25

Bit map images can be loaded and displayed in the library window. These are intended as reference images for defect classification. Bit map images are in Silicon Graphics RGB format.

#### **5.2.5. Defect Code Table files:**

30

The definitions of defect codes are stored in **Defect Code Table files**. User can edit, create, save and load these files.

### 5.3. Top Menu Bar

⇒ See Figure 25: Top menu bar

- 5 Clicking on an item on the top menu bar opens a menu of choices within that item. The Top menu bar provides the following groups of functions:

#### 5.3.1. File:

10

The *File* item on the top menu bar provides access to load, save, and/or edit a variety of Ultrapointe 1000 files via a series of dialog boxes.

15

##### 5.3.1.1. Image:

20

Popup the *Image Open* dialog box for loading previously stored volume and surface image files, and for saving current volume images. Note that surface images are saved from within a 3D window.

##### 5.3.1.2. Recipe:

25

Popup the *Recipe Edit* dialog box for editing, loading and saving recipes.

##### 5.3.1.2.1 Screen editing functions:

###### 5.3.1.2.1.1 Add

###### 5.3.1.2.1.2 Cut

###### 5.3.1.2.1.3 Rename

30

###### 5.3.1.2.1.4 Paste

###### 5.3.1.2.1.5 Window swapping

##### 5.3.1.2.2 laser Params

Popup optics parameters editing dialog box: Note that this dialog box is for editing the recipe's optics parameters, and does not control the parameters while it is being edited.

5                   5.3.1.2.3 Sort Params

Popup the sort criterion editing dialog box.

10                   5.3.1.2.4 File

Popup the file recipe file access dialog box. One can load and save recipe files.

15                   5.3.1.3. Wafer Map:

Popup the Wafer Map **FileSelection** dialog box. In this dialog box, the user can:

20                   5.3.1.3.1. Defect map file list

Display the list of available defect map files.

= See Figure 26

25                   The user can view the list of available maps from different directories. The System, Floppy and Optical Disk directories contain files that are available to all users. The **User** directory is the user's own directory. Each user generally can access his own directory, unless he is logged in as administrator, and  
30                   in this case, he can select any user's directory.

#### 5.3.1.3.2. Load the selected map.

5 Highlight the desired map from the list box and clicking on the **load** button loads the selected map. Any previously loaded map will be erased. If there is a wafer on the stage, the user will be prompted to unload wafer first.

#### 5.3.1.3.3. Save currently loaded map.

10 Click on the **selection** text entry field and type in a file name to save the current map. The defect classification code is saved with the map.

#### 5.3.1.3.4. Edit currently loaded map.

15 Some of the fields in a defect map can be edited: **Wafer number, Lot ID, and File Description.**

#### 5.3.1.4.5. Clear Map:

20 **Clear Map** clears the current wafer map, and creates a new map. The user will be prompted for the desired Wafer Size (4, 5, 6, 8 in.), the Flat Orientation Angle (0, 90, 180, 270 deg.), and the Type of Map.

#### 25 5.3.1.4. Defect Code:

30 Popup the **Defect Code Edit** dialog box for editing, loading, and saving the defect code set.

⇒ See Figure 27

#### 5.3.1.5. Exit:

Popup the ***Confirm Before Quitting*** dialog box. Clicking the ***OK*** button allows the user to exit from the Model 1000 software. The ***Quit*** button cancels the command. Exit is available only to users logged in as Administrator or as Service Engineer.

#### 5.3.2. Screen:

***Screen*** Shows currently configured screen and allows the user to select from the list.

Clicking on a screen item brings up that screen.

#### 5.3.3. Setup:

##### 5.3.3.1. LaserScan Parms:

Popup the ***LaserScan Parms*** dialog box for controlling laser scanning. Most of the functions in this dialog box are the same as in a ***Laser*** window, except that Sample thickness, live laser color, stage autofocus parameters and Stage center offset can only be changed here.

##### 5.3.3.1.1 Stage Center Offset

Set this parameter to make LIS1000 and scanner stage origin match. Setting this parameter correctly will ensure that the deskew points are within field-of-view of a 5X objective, thereby greatly ease the deskew process.

### 5.3.3.2. Deskew Options

5        Popup the *Deskew Options* dialog box . In this dialog box user  
can select the mode of deskew: Use Map Deskew Info, or  
Defects as deskew points.

#### 5.3.3.2.1 Use map deskew info

10       In this mode, the deskew process requires user manually  
matching three predefined deskew points . These deskew  
points' stage coordinates are specified in the defect map file.  
This is the normal deskew mode for patterned wafers.

#### 5.3.3.2.2 Defects as deskew pts

15       Bare wafers with defect maps can be reviewed by the  
LIS1000. On a bare wafer, there is no prominent pattern that  
can be used for deskew. Instead, the defects themselves can  
be used for deskew.

20

### 5.3.3.3 Edit Wafer Map:

25       Popup the *Edit Wafer Map* dialog box for changing current  
wafer map file. Changes made here are those that affect the  
whole file: wafer size, flat orientation, die size, deskew points  
and wafer IDs. Individual defect editing can be done from a  
*Textstatus* window.

#### 5.3.3.3.1 Widgets:

30       5.3.3.3.1.1 FileName: shows the currently loaded wafer map file name.

5.3.3.3.1.2 Lot ID: allows editing of the LotId

5.3.3.3.1.3 Description: file description.

5.3.3.3.1.4 WaferSize option menu: allows choices of 100, 125, 150, 200 mm

5.3.3.3.1.5 Flat Orientation option menu:

allows choices of Down, Left, Up, and Right.

5 5.3.3.3.1.6 Die Pitch:

X and Y die size in microns; not editable for Tencor formatted files.

• 5.3.3.3.1.7 Deskew1

10 First Deskew location, in microns; this field is not editable for KLA formatted files.

5.3.3.3.1.8 Deskew2

15 Second deskew location, in microns; this field is not editable for KLA formatted files.

5.3.3.3.1.9 Registration:

Displays the registration location. The registration location should NOT be changed unless it is known to be incorrect.

20 Die patterns and defect locations are all relative to the registration point. Changing the location of the registration will shift the overall wafer map.

5.3.3.3.1.10 Wafer ID test entry field:

25 Allows user to enter wafer ID to specified slot. A valid wafer file must have at least one slot with a WaferID.

5.3.3.3.1.11 WaferID, Enter To List button:

30 Click on this button to accept the contents in the WaferID text entry field.



5.3.3.3.1.12 WaferID list box:

Displays WaferIDs for all slots. For a slot that has no map, "[No Map]" is displayed.

5

5.3.3.3.1.13 OK:

Accept changes as displayed and close the dialog box. All affected windows will update according to the new parameters.

10

5.3.3.3.1.14 Quit:

Ignore any changes made, and close the dialog box.

5.3.3.4. Defect List Sort:

15

Popup the *Defect selection* dialog box for editing defect sorting criteria. Sort criteria can be cascaded, and the Criteria list box displays the whole sort criteria.

5.3.3.4.1 Widget:

20

5.3.3.4.1.1 Criteria list box

Displays the whole list of sort criteria.

5.3.3.4.1.2 Criteria choice

25

Displays and allows selections of provided sort criteria.

5.3.3.4.1.3 Include and Exclude

30

Displays and allows selections of qualifier to the criteria. For the criterion **ALL**, defects are always included and the choices of include and exclude do not apply.

5.3.3.4.1.4 Minimum

Inputs the minimum requirement for the criteria. For *All*, *First*, and *Random* criterion, the minimum requirement does not apply.

5

5.3.3.4.1.5 Maximum

Inputs the maximum requirement for the criteria.

5.3.3.4.1.6 Cut

10 Clicking on this button cuts away the selected criteria list item and copies it to the widget displays

5.3.3.4.1.7 Paste Up

15 Clicking on this button pastes criterion as displayed by the widgets to the criteria list box. New list box entries appear above the currently selected list box item.

5.3.3.4.1.8 Paste Down

20 Clicking on this button pastes criterion as displayed by the widgets to the criteria list box. New list box entries appear below the currently selected list box item.

5.3.3.4.2 Supported criterion :5.3.3.4.2.1 All

25 Selects all defects from the wafer map.

5.3.3.4.2.2 First n.

Selects first n defects from the wafer map.

30

5.3.3.4.2.3 Random n.

Selects n randomly selected defects from the wafer map.

5.3.3.4.2.4 X size.

Selects defects with x size within the minimum and maximum from the wafer map.

5

5.3.3.4.2.5 Y size.

Selects defects with y size within the minimum and maximum from the wafer map.

10

5.3.3.4.2.6 Intensity.

Selects defects with intensity within the minimum and maximum from the wafer map.

15

5.3.3.4.2.7 Defect code.

Selects defects with defect code within the minimum and maximum from the wafer map.

20

5.3.3.4.2.8 X location.

Selects defects with x location coordinate within the minimum and maximum from the wafer map.

5.3.3.4.2.9 Y location.

Selects defects with y location coordinate within the minimum and maximum from the wafer map.

25

5.3.3.4.2.10 Every n

Select only one out of every n defects, i.e. for n=5, defect number 5,10,15,20.... are selected

30

5.3.3.4.3 Qualifier.

The criteria can be qualified with:

5.3.3.4.3.1 Include

The sorted list includes the defects that passed the sort criterion.

5

5.3.3.4.3.2 Exclude

The sorted list excludes the defects that passed the sort criterion.

5.3.3.5 Rearranging Wafers

10

This dialog box is provided to rearrange the slot number of wafer maps in a wafer map file, in case that there is a mismatch between actual wafer placement and map file contents. Rearrangement is made only to the map file, no actual wafer movement is involved. Rearrangement is unnecessary if the wafer slot location matches those specified in the map file. In this dialog box, 25 buttons are provided corresponding to the 25 slots in a cassette. The slot with a wafer map is displayed in GREEN, slot with no map is colored WHITE. To move the wafer map from one slot to another: Click on the button for the first slot; hold down the mouse button while moving the mouse to the button of the second slot; release mouse button. The color of the button will change to reflect file changes. This Click-n-drag operation checks for the follow error conditions: Source slot has no map, destination slot already has a map, Click or release button on space other than the 25 slot buttons.

15

20

25

30

**5.3.4. User:**

5.3.4.1. Login:

⇒ See Figure 28

5        Popup the **Login** dialog box. This dialog box allows user to  
logout, so that others will have to log in to use the system, as  
well as log in as another user.

5.3.5. Config:

10

The **Config** menu is used for option configuration.

5.3.5.1. SECS:

15        **SECS** is planned to be used for the SECS II option, not  
currently installed.

5.3.5.2. Printer:

20

⇒ See Figure 29

**Printer** is used to setup the print options: Printer type, and  
number of copies.

25

5.3.5.3 Misc.

System Configuration dialog box allows Administrators or  
service engineers to edit the following system parameters:

#### 5.3.5.3.1 Screen/laser saver timer, in minutes

5 This item specifies the number of minutes of mouse or key board inaction before the system goes to screen/laser saver mode. In the screen/laser saver mode, the Laser is switched to idle to prolong laser life, and the Laser window is turned to dark to save the monitor screen.

#### 10 5.3.5.3.2 Deskew die Offset

15 For a defect file that contains no deskew information (i.e. KLA formatted files.) The LIS1000 generates the three deskew points at the lower-left-corner of three key dies. The left-most middle die is selected as the first deskew die, the right-most middle die is selected as the second deskew die, while the center, 00 die is used as the registration die. Sometimes, a wafer may have no pattern on the outer-most die locations, and the use of a non-existent die for deskew will cause problem. A *Deskew Die Offset* of 1 makes the dies next to the  
20 outer-most dies as the deskew die, an offset of 2 selects the dies 2 die space away from the left and right die.

25 The die offset adjustments applies even to those files that do have deskew information. The deskew points in a Tencor file, for example, sometimes point to an incomplete die, or a test die. In this case, a LIS1000 user can use the *Deskew die offset* to move the deskew points to another die.

#### 30 5.3.5.3.3 Flat angle offset (deg)

If the wafer pattern is slightly skewed with respect to the wafer flat, or if there is an error in the robot teach angle, the wafer pattern on the monitor screen will appear to be tilted.

A stage movement along the wafer's horizontal street will cause the wafer pattern to move vertically. This XY tiling can be removed by adjusting the *Flat angle offset*.

This offset parameter allows fine tuning of the angle in which a wafer is loaded onto the LIS1000 stage. The adjustment resolution is .125 degree( 1/8 of a degree) and typical value is within +/- 0.5 degrees.

#### 5.3.5.3.4 White light mode

Two white-light-only options are available: white light with a beam cube *in* or with it *out*. The beam splitter is required in the laser scan mode to get simultaneous white light and laser viewing, as well as in autofocus. However, with the cube in, the white light intensity is cut by 50%( although image quality remains unchanged, and perceived image brightness is down by marginal amount.) For very dark samples, one therefore would want to switch out the cube in white-light-only mode. But, since the cube has to be *in* in the laser mode, the beam splitter cube will be switching in and out every time the laser and white-light-only mode is changed. The beam splitter switching is slow and can be a throughput bottleneck. To improve throughput, user has the option of keeping the beam splitter *in* in the white-light-only mode.

Cube Out: brightest white light image, but switching between white-light and laser can be slow.

Cube in: Slightly dimmer image for the same camera intensity, but otherwise no changes in image quality. Significantly better throughput and reliability.

#### 5.3.5.3.5 Next die position

When viewing a defect, sometimes It is useful to compare the image with those at the same die location, at the next die. The LIS1000 has a button allowing user to go to the same location on the next die. User uses the *Next Die position* to specify where the next die is located: Bottom/Left/Top/Right.

#### 5.3.6. Maintenance:

##### 5.3.6.1. Manual stage control:

⇒ See Figure 30

Popup the *Manual Stage Control* dialog box. This dialog box is used to initialize the XYZ stage, as well as to perform stage diagnostics.

##### 5.3.6.2. Manual robot control:

⇒ See Figure 31

Popup the *Manual Robot Control* dialog box. In this dialog box, users can initialize the robot, and command elementary robot functions. This is used, for example, to remove a wafer left on the stage because of a power failure.

##### 5.3.6.3. File Transfer:

⇒ See Figure 32



Popup the *File transfer* dialog box to *move, copy, rename* and *delete* files.

#### 5.3.6.3.1 File lists:

5 Files from two directories, selectable using the *From* and the *To* Option menu, are listed in this dialog box. File types listed are selectable from the *FileType* option menu.

#### 5.3.6.3.2 File Type

10 The file type of a file is generally distinguishable from its file name extension, the file name's last three characters. The File transfer dialog box lists only one file type at one time, and user can choose surface, volume, wafer map, bit map, recipe, or defect code to list the specific file type. An *All* file type  
15 choice is provided for listing all files in any directory. When the Wafer Map file type is selected along with the *FLOPPY, IMPORT, REMOTE* or *REMOTE2* directory, then all file names are listed.

#### 5.3.6.3.3 Target Name

20 Clicking on a file from the *From* directory selects that file for later delete, move, copy or rename. The name of the selected file is also copied to the *TargetName* entry field. The name in the entry field is used as the new file name in *Rename*, and as  
25 the target name in copy and move. Users can edit the name in the entry field for the different name desired.

#### 5.3.6.3.4 Delete

30 The selected file from the *From* directory is deleted.

**5.3.6.3.5 Move**

The selected file from the *From* directory is moved to the *To* directory with a file name from the *TargetName* entry field. The selected file in the *From* directory will be deleted.

5

**5.3.6.3.6 Copy**

The selected file from the *From* directory is copied to the *To* directory with a file name from the *TargetName* entry field.

10

**5.3.6.3.7 Rename**

The selected file from the *From* directory is renamed to be the name from the *TargetName* entry field.

**5.3.6.3.8 File selection**

15

User selects from the *From* file name list box by clicking on the desired file name. Multiple files can be selected: click on the first one, hold down the mouse button and drag the mouse to the last one. A warning message is displayed for deleting, coping, moving multiple files. Only one file can be selected for Rename.

20

**5.3.6.4. UserAcct:**

⇒ See Figure 33

25

Popup the *User Accounting* dialog box. Each model 1000 user is assigned a name, password (can be blank), and a privilege level (operator or administrator). Users can be added and deleted in this dialog box.

30

**5.3.6.4.1 User Names list box**

This list box displays the list of currently configured users.

**5.3.6.4.2 Copy**

5 Copies the user info of the high-lighted user in the user names list box. The **Name**, **Password** and **privilege level** will be updated.

**5.3.6.4.3 Cut**

10 Copies the user info of the high-lighted user in the user names list box. The **Name**, **Password** and **privilege level** will be updated. The high-lighted user will be deleted from the user list.

**5.3.6.4.4 Paste**

15 Create a new user and put into the user list, from the **Name**, **Password** and **Privilege** displayed.

**5.3.6.4.5 Name, Password, Privilege level**

20 Displays and Allows editing of user information. After editing a new user, the new user info needs to be pasted into the user list, and the user list must be saved.

**5.3.6.4.8 Save**

25 Save all the changes made, and the LIS software start using the new user list.

**5.3.6.4.9 Quit**

30 Ignore all changes to the user list, and quit the User Accounting dialog box.

#### 5.3.6.5. LON Diagnostics:

⇒ See Figure 34

5        Popup the ***LON Diagnostics*** dialog box. This dialog box provides diagnostics and service control for most of the system hardware.

#### 5.3.6.6. Calibration:

10

⇒ See Figure 35

      Popup the ***Laser Scan Calibration*** dialog box. This dialog box is used to input the XY laser scan calibrations.

15

#### **5.3.7. PrtScreen:**

***PrtScreen*** prints the whole screen on a digital printer.

20

      It takes up to 10 seconds to capture the screen image. The user is asked to confirm printing when screen is captured, but before the image is send to the printer. If more than one digital printer is installed, the user is asked to specify the destination. The system will run slower than normal with printing in the background.

25

#### **5.4. Window and Screen**

30

      There are eight window types in the Model 1000 graphical user interface: ***3D, 2D, Cassette, Wafer map, Status, Laser Scan, Image Library, and uscope*** (microscope). These windows are sizable and movable and they may be configured

by the user in a non-overlapping fashion. A specific screen is a particular arrangement of windows. The following screens have been pre-configured:

5                    *FourView*  
                      *BigSurface*  
                      *TwoSliceSurface*  
                      *3View*

1                    10                    Different screens can be installed by the user.

#### 5.4.1. Common window attributes

1                    15                    ⇒ See Figure 36

2                    20                    Each major functional module is encapsulated in a window. For some window types, one can have multiple windows of that type in a screen. For others, there must be one and only one of that window type in a screen. 2D, 3D, library and text status windows support multiple windows. There must, however, be one and only one of the following window type in a screen: Laser, microscope, wafer, and cassette. When switching screens, the contents of screen's windows carry over to the next if there is a corresponding window type in the new screen. It is more convenient therefore, to make all screens to have the same number and types of windows. For example, when one switches from a screen with three 2D windows to a screen with two 2D windows and back again, one of the three 2D window will lose its image. Generally, one would want to have three 2D windows to show the top ("XY"), side ("YZ"), and front ("XZ") views of the sample, two

25                    25                   

30                    30

text status windows and two 3D windows. Although the functionality is different, all window types share the following attributes:

5      5.4.1.1. Window sizes.

All window types have three distinctive looks, depending on the size of the window. Size can be:

10      5.4.1.1.1. Large:

When the size is large enough to show the full functionality of the window. Typically a large window is wider than a third of the screen. All controls for a large window are directly accessible

15      5.4.1.1.2. Small:

When a window is large enough to show its major functions and a few controls but not large enough for full functionality. A small window allows access to its full range of functions indirectly through a control dialog box.

20      5.4.1.1.3. Icon:

When a window is less than 95 pixels wide. An iconized window displays only a button with the window name on it. To use the window, it must be swapped to a larger window location.

5.4.1.2. Swapable to other window sites on-screen.

30      All windows uses the click-drag method to swap between different windows.

Click and hold the mouse on the swap button in one window and drag it to another window. The two windows will interchange positions. The swap button is labeled "L".

5      5.4.1.3. Global and local data.

Operations originated from the top menu bar are global in nature, i.e.: they affect all windows. Operations originated from within a particular window affect that window only. For example, loading a surface file from the top menu puts the same surface on all 3D and 2D windows, as well as loading the wafer map file for that image. Loading a surface from within a 3D or 2D window updates the surface only in that window and will not update wafer map.

15      Volume data sets are generally large, and there is a limit of only one volume data set (file) loaded in the system memory at any one time. Loading a volume file from within a 2D window therefore also updates all 2D windows.

20      5.4.1.4. Inter-window data exchange.

There are some user interface features that require the cooperation of multiple windows; these features are implemented with the technique of inter-window data exchange. For example, the displaying of 2D XY cross sectional views using the top file load, allows the "*cut*" cursor in the XY display to control the slices viewed in the XZ and YZ displays in the other 2D windows.

30

## 5.4.2. Window types

### 5.4.2.1. Cassette Window

5           ⇒ See Figure 37

10           The **Cassette** window provides the user interface for loading and unloading wafers in a cassette. It gives a graphical display of the current status of the wafer transport subsystem.

#### 5.4.2.1.1. Wafer status:

15           There are three icons on the left of the cassette window, indicating the location status of a wafer among the robot arm, stage, and the pre aligner.

            This window also displays the current cassette size, and the name of current wafer map file.

#### 20           5.4.2.1.2. Cassette status and wafer selection

            In the middle of the window, there are 25 rows of buttons and color bars corresponding to the 25 wafers in a cassette.

25           Buttons are for selecting the wafer, to be used for loading and unloading. In addition, prior to actually loading any wafer, one can also **preview** the different maps in the file, by selecting its corresponding slot. If the selected slot does not have a map, the user will be prompted if an empty map for that slot is to be created.

30           Color bars indicate the status of the cassette slot:



**Present:** A wafer is considered to be present in a slot on the cassette when that a wafer has been successfully loaded or unloaded in that slot.

5           **Inspected:** A wafer is considered to be inspected when it has been successfully loaded on the stage and returned (unloaded) back to the cassette.

10           **File Present:** Valid defect information for the slot is available from the currently loaded defect map file.

**5.4.2.1.3. Command Buttons:**

5.4.2.1.3.1 Lower25/Upper25

15           For a wafer map that has wafer information for 2 cassettes, this button Toggles between using the first/second cassette.

5.4.2.1.3.2. LoadWafer:

20           **LoadWafer** loads the wafer from the selected slot onto the stage. The load wafer sequence consists of:

Cassette presence check, chuck vacuum check, opening the wafer door, and moving the stage to the load/unload position.

25           The robot first takes the wafer from the specified cassette slot with its vacuum actuated end effector, and then transfers it to the flat finder (pre aligner).

30           The flat finder orients the wafer according to the specified flat position.

The robot retrieves the wafer from the flat finder, and puts it onto the stage chuck, where it is secured and sensed by the robot vacuum system.

5      5.4.2.1.3.3. Unload:

**Unload** removes the wafer from the stage and returns it to the selected slot. The unload sequence involves:

10      Robot picks up wafer from chuck, which has moved to the unload position.

Check for cassette status: Displays a warning message if the cassette has been removed.

15      Robot puts wafer into cassette slot.

5.4.2.1.3.4. Stop:

20      **STOP** stops current robot operation. A warning message is displayed and the user can **reset** the robot, or **resume** the original operation.

5.4.2.1.3.5. Initialize:

**Reset** initializes and reset robot to its reset position.

25      5.4.2.1.3.6. UsingMap/NotUsingMap:

A wafer map file can have multiple wafer maps: one for each slot. Two modes are provided for the different ways wafer map files are used.

30      In the **UsingMap** mode, loading a wafer from a slot with defect map information will also load that wafer map. If the map for that slot is not found, a new one is created. This

new map will contain the same die information as the other maps in the file. In the *NotUsingMap* mode, wafer loading has no effect on current wafer map.

5        5.4.2.1.3.7. Flat orientation

This button tells the robot to orient wafer to one of the four directions: 0/90/180/270 degrees. "0" is with the primary flat/notch at the bottom, and the subsequent positions proceed clock-wise. The robot will use this new orientation to place the next wafer. Note: Most wafer maps will specify flat/notch orientation, and user normally needs not to change the wafer orientation from those specified.

10       5.4.2.2. Wafer Window

15       ⇒ See Figure 38, Wafer Map Window

This window displays the name of the current wafer map file, wafer size, slot number and the XY coordinates of the current stage position, and selected target stage position.

20       There are two graphical display sections for the defects. The one on the left displays the overall view of the wafer. The one on the right displays a zoomed view of the section selected by the zoom box. The Zoomed view tracks the zoom box and update in real time.

25       5.4.2.2.1. Rubber band zoom box :

30       5.4.2.2.1.1. Zoom Box Size:

Click and hold on the lower right hand corner of the zoom box and then drag the mouse to change the size of the zoom box.

5.4.2.2.1.2. Zoom Box Positioning:

Click and hold on the upper left hand corner of the zoom box and then drag the mouse to move the zoom box.

5      **5.4.2.2.2. Selecting a defect location**

The wafer window is used to select an XY location to visit. In this window, there are two ways to select a defect location:

10      5.4.2.2.2.1. Graphical Selection.

Clicking on the desired position (defect sites are shown as colored dots) on the overview defect display or the zoomed defect display selects the location. Due to screen pixel size limitations, selecting a defect from the overview section will only be precise to several hundreds of microns. From a greatly zoomed display, however, one can generally reduce placement uncertainties to a few tens of microns.

Clicking on the **GoTo** button moves the stage to this selected position.

5.4.2.2.2.2. Type-in the desired location.

Type in the XY location on the XY text entry field, then click on the **GoTo** button to move the stage to this location

25      **5.4.2.2.3. Operation without a wafer map**

Without a wafer map, the wafermap window shows only the wafer outline without defect nor die grids. One can still use the wafer window, in this case, to move the stage and to read current stage positions.

#### 5.4.2.2.4 Current XY stage position

Current XY stage coordinates are displayed in 2 ways: as numbers in the top left portion of the window, and as a cross in the graphic display. Current stage positions are updated at the end of all programmed stage moves. Also, they are updated when user clicks anywhere inside the overview defect display or the zoomed defect display .

#### 5.4.2.3. Laser Window

⇒ See Figure 39: Laser Scan Window

This window provides the user interface for most of the white-light imaging control as well as the laser scanning functionality.

##### 5.4.2.3.1 White light image

The live white light image is normally displayed on a separate monitor. The user has the option of selecting **LaserScan** or **whiteLight**. In LaserScan mode, both the white light image and the laser image are displayed. The white light image is yellow filtered and thus appears yellowish. The WhiteLight mode gives the best white light image, but with only white light image displaying, the laser image is blocked,

##### 5.4.2.3.2. Live laser imaging

The live laser scan image is displayed here. The user can use the joystick to move the stage in XYZ direction and watch for their effects on this window. When the laser window is large enough, a vertical scale showing the scan dimension is displayed on the right side of the Live Laser display.

#### 5.4.2.3.3. Autofocus

There are two Z stages: the coarse and the fast.

5 The slower, Coarse Z has a travel of six millimeters, and is used to move the wafer surface into general focus from the load position of the stage (full-down).

10 The faster, fine Z has a travel of about 50 microns and is used for precision focus and to step the wafer through different Z positions during the volume and surface image acquisitions.

Two means of autofocus are provided: stage autofocus and fine autofocus.

##### 15 5.4.2.3.3.1. Stage Autofocus

The **Stage AF** button initiates this autofocus: The fine Z stage is first moved to the middle of its travel, the coarse Z stage is moved three times at successively slower speed to find the best focus. A stage autofocus generally takes about 4  
20 seconds. Coarse autofocus is generally used once, after the wafer is loaded, but may be required occasionally during inspection if the range of the fine focus is exceeded. Parameters for coarse auto focus can be pre-calibrated as part of the recipe, and should work for all samples, without regular  
25 user adjustments.

##### 5.4.2.3.3.2. Fine Autofocus

30 The **AutoFocus** button initiates fine autofocus. The fine Z stage is first moved to its lowest position, then moves up until focus is within reach. A more detailed search of the best focus is then performed. A fine Z autofocus generally takes less than a second. The user may need to set the laser

intensity for samples that are particularly dark or particularly bright. If fine autofocus fails, the Current Z slider goes to the top of its range ( 50 micron) or bottom of the range( 0 micron).

5

#### 5.4.2.3.4. Objective Selection

Click on the **Objectives** option menu to select desired objective. Several things happen when objective is changed: XYZ stage moves to effect parfocality and parcentricity, Laser intensity adjustments, camera attenuation for 5X objective and fine-Z-autofocus. All operations for objective changes are transparent to user, and is configurable by system administrator or service engineers.

10

15

#### 5.4.2.3.5 Camera control parameters

##### 5.4.2.3.5.1 Camera intensity

Adjusts between the nominal camera intensity between 0 to 255 for different samples.

20

##### 5.4.2.3.5.2 Camera mode

Provides three options: normal( bright field), dark field and polarizer. Normal and dark field are applicable to both white-light-only and laser scan. Polarizer only works in white-light-only. If currently in laser scan mode, selecting polarizer camera mode switches system to white-light-only.

25

#### 5.4.2.3.6. Laser scan parameters

These controls affect the operation of the laser and its scanner.

30

5.4.2.3.6.1. (Laser) Scan area

XY Laser Scan area can be selected by changing objective, XY scan resolution, and data size. The **XY Res** option menu displays the current laser scan XY resolution, and the horizontal and vertical range scales display the corresponding scan range.

A 20X objective gives e.g., 5 times the scan length and width as 100X objective. Click on the **Objective** button to select objective size.

Four XY Scan zoom settings, (corresponding to pixel sizes in laser display) are available: 100, 50, 25 and 12.5 nanometers (at 100X). The pixel size scales with objective magnification. Click on the **XY Res** button to select a pixel size/zoom factor.

Two data sizes are selectable: 512 and 256. Using the smaller data set requires less data storage and less processing and gives a faster response. Click on the **Dimensions** button to select.

5.4.2.3.6.2. Wave length selection

The Ultrapointe Model 1000 uses a multi-line argon laser for its laser source. User can select one of the three dominant wavelengths (458, 488, 515 nm), or all lines from the laser. (Additional argon lines are available as factory-installed options.) The **LaserLine** button allows wavelength selection.

In general, shorter wavelength gives better resolution, but there are several other reasons for selecting a particular wavelength. For example, sample spectral reflectivity,



spectral opacity, and spectral absorption, all may play a role in selecting the wavelength.

#### 5.4.2.3.6.3. Image intensity

Image intensity is affected by many factors: laser power, laser sensor gain, and sample reflectivity for the wavelength. The ***LaserIntensity*** slider control combines laser power and sensor gain in a normalized scale of 1 to 100. The user can use this control independently with the wavelength selection. The laser image sensor has a built-in protection mechanism that will shut itself down when intensity becomes too bright. The laser scan display becomes dark when the sensor shuts down, and the color of the ***LaserIntensity*** slider label changes to red. Clicking on the ***LaserIntensity*** slider bar at a lower value resets the sensor and the label color.

#### 5.4.2.3.6.4. Focus change

Laser scanning has a very short depth of focus. By adjusting the ***Current Z*** slider control, the user can obtain a view of the sample at the desired depth. The fine Z stage control allows user to position the sample in steps of about 12 nm over a range of about 50 microns.

#### 5.4.2.3.7. Image Positioning

##### 5.4.2.3.7.1. Joystick

There is a toggle switch on the joystick, when this switch is in the up or "Z" position, moving the joystick up and down moves the coarse Z stage up (+Z) or down (-Z). In this position, the left-right motion of the joystick has no effect. The Z control with the joystick is intended for engineering use only. All routine focusing operation should be performed with

the coarse and fine autofocus functions, and the fine z slider control. Crash protection firmware prevents accidental "crashes" of the objective against the wafer.

5 When the toggle switch is in the down or "XY" position, moving the joystick up and down moves the stage in the +Y and -Y directions and moving the joystick left and right moves the stage in the +X and -X directions. Joystick Z travel limits are set based on configurable sample thickness.

10 5.4.2.3.7.2. Point and go\*

(Not yet implemented). Click on the screen commands the stage to move that point to the center of display.

15 5.4.2.3.8. Surface Imaging Mode

5.4.2.3.8.1. Principles

20 The underlying principle of scanning laser imaging is the short depth of focus of the confocal optics. By taking images at regular intervals of Z, the sample's XYZ intensity information is obtained. From this set of 3D volume data, the Model 1000 extracts the sample's 3D surface. In a 3D window, the sample surface can be viewed with different tilt and rotation angles, with zoom, different illumination effects as well as  
25 with different color schemes

A conventional microscope can only look at the top view of the sample, without Z information.

30 5.4.2.3.8.2. Slices

The number of slices selected is the number of data points in the Z direction. The Model 1000 has a scan speed of about

13 slices per second at 512 by 512 pixel image size and a maximum number of slices of 128. (At 256 by 256 pixels, the scan speed is 26 frames per second.)

5      5.4.2.3.8.3. Range Selection

The Z distance between slices is selectable to be between 1 to 30 basic steps, where a single step is approximately 12 nm. Use the **Slices** slider to select the number of slices desired. The number of slices and step size together  
10      determines the range of the Z motion in an image acquisition.

As a convenient feature, the user can use the **S**(tart) and **F**(inish) buttons to define the start and finish of a scan, and the system adjusts scan step size. The user utilizes the  
15      **Current Z** slider to move the laser focus above the region of interest. Clicking the **S** button sets the current z position as the starting point. Next, the user moves the slider down so that the laser focus is now below the region of interest. Clicking the **F** button sets the finishing point. Scan step size  
20      is then automatically adjusted to match this selected scan range, and the focus returns to the mid point between start and finish positions. Note that the **S** and **F** button are used for calculating the step size only. Actual data acquisition  
25      starts at half of the range above current Z position and steps at step-size for the number of steps specified.

5.4.2.3.8.4 SetZ

Instead of manually setting the Z scan parameters, the user can click on the **SetZ** button to automatically set the Z scan  
30      start position, range, and laser intensity. When user clicks on **SetZ**, the system initiates a surface acquisition around the current Z over a 20 micron range. From the resultant surface,

the system calculates and sets the best starting point, range and the optimal laser intensity. It also sets the fine Z stage to approximately mid point between the start and stop positions. The **SetZ** function takes about 3 seconds to complete for the 20 micron range. The actual range and time is configurable.

*5.4.2.3.8.4.1 Limitations of the SetZ function:*

Although the **SetZ** function has been optimized for most typical use, there are cases where the Z scan parameters should be set manually:

The **SetZ** function searches for the sample's top and bottom within a range of about 20 microns. It can give erroneous results for defects taller than 15 microns.

The minimum range set by **SetZ** is about 2.5 microns (0.04 microns step size for 64 slices), and the resolution for this range may be too coarse for some flat defects.

5.4.2.3.8.5. Image Acquisition

Clicking the **Surf Acq** button starts the acquisition of a surface image at the half scan range above current Z position, using the pre-selected scan range. In this mode a volume scan is performed, and the surface data processor automatically produces the 3D surface.

5.4.2.3.9. Advanced features

5.4.2.3.9.1. Volume Data Acquisition

For a compositional complex sample, such as semi-transparent, multilayer structures, subsurface defects, etc., there may be a need to look at the basic volume data set directly, rather than the reduced surface. Clicking the **VolAcq** starts a volume acquisition scan without the surface extraction stage.

5.4.2.3.9.2. Multi surface analysis

Complicated samples with multiple layers can be viewed in the Model 1000. The following features have been implemented for extracting a particular layer surface:

5

4.4.2.3.9.2.1 Different view direction

Views of the volume data in 2D windows from different viewing directions.

10

4.4.2.3.9.2.2 New surface generation

From a 2D window displaying a XZ or YZ cut, generate a new surface with restricted start and finish slices.

5.4.2.4. Microscope window

15

⇒ See Figure 40

This window shows the white light video image in real time. It is not yet implemented.

20

5.4.2.4.1. Simultaneous Viewing5.4.2.4.2. Objective Selection5.4.2.4.3. Illumination Control5.4.2.4.4. Image Filtering\*

25

5.4.2.4.5. Bright field/Dark field Operation5.4.2.4.6. Image Storage\*5.4.2.4.7. High-resolution Mode (no laser image)5.4.2.5. 2D Image Window

30 ⇒ See Figure 41

The 2D window can display surface images and volume images in any one of the 3 cutting planes: XY, XZ, YZ.

#### 5.4.2.5.1. Multiple views within window

Within each 2D window, a sequence of 1, 4, 9, or 16 views can be displayed. Clicking the **Sections** button cycles through 1, 4, 9, and 16 views displays. The 1st Slice and Spacing slider bar determines the position of the first view and the interval between views.

Volume data size tends to be large in general, and there is a limit of one volume file in the system at a time. Because of this limit, loading a volume file, or acquiring a new volume data updates all 2D windows. No such limit exists for surface data.

#### 5.4.2.5.2. Cursors

Clicking on anywhere within the view area brings up the cursor. There are two cursor types: sizing and cutting. Size and Cut buttons toggle between the two choices.

The sizing cursor is a "rubber band rectangle". Clicking and dragging on the upper left corner moves the cursor, clicking and dragging the lower right corner size the cursor. The vertical and horizontal dimensions of the cursor are displayed in the large screen window.

The cutting cursor is a "rubber band cross hair". Clicking and dragging on the horizontal or the vertical line moves the respective cursor line. The slices (cutting planes) addressed are displayed in the large screen window.

#### 5.4.2.5.2.1. Multi-window operation, the cut cursor

For a complete 2D view of a 3D object, one probably wants to look from three different directions: top, side, and front. Having an

indication of the positions of the cutting planes for the side and front view may also be useful. Rather than having a dedicated window having these complex functionality, the Model 1000 uses the concept of multi-window operation to achieve the same goal. For a complete  
5 2D view, a Model 1000 user can use 3 2D windows, one for each cutting plane and selects the cutting cursor for the window displaying the XY view. The cursor information of this window is dynamically broadcast to other 2D windows. One window will therefore display the XY plane with the lines showing the cutting planes for the other  
10 two 2D windows.

#### 5.4.2.5.3. Image appearance functions

##### 5.4.2.5.3.1. Smoothing

15 Clicking the Smooth button smoothes the displayed data. If global surface data is displayed, displays in other windows may share this effect.

##### 5.4.2.5.3.2. Color

20 Clicking the Color button toggles the display between false color and gray scale. (This function is different in the 3D window).

##### 5.4.2.5.4. Viewing directions: Surface

##### 5.4.2.5.4.1. Top view (XY)

30 In top view, the 2D window displays the surface in top view, in false color, or in gray scale. Since there is only one surface, displaying more than one view means display the same surface multiple times.

#### 5.4.2.5.4.2. Vertical Slices (XZ, YZ)

A surface profile is displayed in XZ or YZ plane.

#### 5 5.4.2.5.4.3. Multi-surface

More than one surface can be displayed in a 2D window, in XZ and YZ plane.

10 These surfaces may be generated from the same volume data using different surface extraction methods. Or, these surfaces may be similar but from different samples. Since each surface file holds only one surface, for multi-surface, all surfaces files must be loaded locally, one at a time, with the Add/Replace option set to Add.

15

#### 5.4.2.5.5. Viewing directions: Volume File

##### 5.4.2.5.5.1. XY, Top View

20 In top view, **XY** the 2D window displays the intensity volume data in projected top view, at the selected Z slice number.

##### 5.4.2.5.5.2. XZ, front view

25 In front view, **XZ**, the 2D window displays the intensity volume data as viewed from the front, at the selected Y slice number.

##### 5.4.2.5.5.3. YZ, side view

30

In side view **YZ**, the 2D window displays the intensity volume data as viewed from the left side, at the selected X slice number.



#### 5.4.2.5.5.4. XZPeaks, YZPeaks

5

**XZPeaks, YZPeaks** are views from the same direction and at the same slice number as XZ and YZ. However, the extracted peaks, indicating possible surfaces, are displayed. The user may need to set the **Threshold** slider bar to cut out unwanted noise.

#### 5.4.2.5.5.5. Z Profile

10

In the **Z Profile** view, the intensity profile along the Z axis at the selected XY location is plotted. This profile is useful as a diagnostics tool, as well as in multilayer viewing.

15

#### **5.4.2.5.6. Image History**

20

Surface and volume image are saved with their history records. Clicking the History button pops up the **Image File History** dialog box. The image file name, comments, date and time of acquisition, wafer map used, operator name, wafer slot, defect number, and Laser line used are displayed in the **Image File History** dialog box.

25

#### **5.4.2.6. 3D Surface Image Window**

⇒ See Figure 42

30

This window displays the surface in 3D. Various image manipulation tools are provided.

#### **5.4.2.6.1. Image manipulation**

## 5.4.2.6.1.1. Rotate

The **Rot-** button rotates the image counter clock wise, the + button rotates clockwise.

5

## 5.4.2.6.1.2. Tilt

The **Tilt-** button tilts the image toward top view, the + button towards the front view.

10

## 5.4.2.6.1.3. Zoom Function

Each click of the **zoom+** button brings the image closer toward viewer, so to make the image larger. The **zoom-** button reverse the effects of the previous zoom+ action. Click on **Reset** to undo zoom.

15

Detail control of Zoom is accomplished with the zoom cursor box, activated by clicking anywhere on the 3D window image. This is a special "rubber band box" whose size is controlled by depressing the mouse button, and sliding the mouse left or right. To position the re-sized zoom cursor, release the mouse button, slide the arrow cursor to the desired location and click once momentarily. The zoom cursor will move to the new location, centered on the arrow cursor position. Finally, click on **Zoom+** to implement the zoom.

20

25

## 5.4.2.6.1.4. Translation

The ← → and Up Down buttons move the image vertically and horizontally.

30

## 5.4.2.6.1.5. Z enhancement

Each click on the Z enhancement button enlarges the Z scale.

5

## 5.4.2.6.2. Image appearance

10

Surface data, height and intensity are displayed in the 3D windows in a photo-realistic algorithm. The image is rendered with illumination from a stimulated light source. The direction of the light source can be adjusted to appear increasingly from the left or right by clicking the **Xdir-** and **+** buttons, respectively.

15

## 5.4.2.6.2.1. Brightness and contrast

Clicking the **Bri-** and **+** button adjusts the overall image brightness.

20

Clicking the **Contr-** and **+** button adjusts the contrast.

Clicking the **Color** button repeatedly cycles through the four available color schemes:

25

- Height and intensity in false color,
  - Height and intensity in gray scale,
  - Height only, with height used as intensity, in false color
  - Height only, no color or intensity information.
- Changing the lighting direction may be used here to

30

highlight small step changes, emulating Differential Interference Contrast

#### 5.4.2.6.2.2. Smooth

5

Clicking the Smooth button smoothes the surface. Clicking Reset will not undo the smoothing effect.

#### 5.4.2.6.3. Sectioning

10

Moving this **Sectioning** slider selects the ending point of the display, thus gives a sectioning effect at that point.

#### 5.4.2.6.4. Reset

15

Clicking on the **Reset** button resets all the image manipulation steps to the default condition. Note however, the effect of smoothing cannot be reset.

#### 5.4.2.6.5. Image History

20

Surface and volume image are saved with their history records. Clicking the History button pops up the **Image File History** dialog box. The image file name, comments, date and time of acquisition, wafer map used, operator name, wafer slot, defect number, and Laser line used are displayed in the **Image File History** dialog box.

25

#### 5.4.2.6.6. Image Storage

30

A surface image can be saved by clicking on the **Open** button to popup the **FileSelection** dialog box. Click on the **Selection** text entry field, entering a file name, and then click on the

**Save Surface** button to save the displayed surface under that file name. The bit map of the currently displayed 3D image can also be saved by clicking the **Save Bit map** button.

5      5.4.2.7. Library Image Window

⇒ See Figure 43

10      The user can load and display previously saved bit map images in a **Library** window.

5.4.2.7.1. Library image selection

15      Click on the **File** button opens the **File Selection** dialog box. Highlight the desired bit map file name with the arrow cursor and click on the **Load** button to load the image.

5.4.2.7.2. Defect classification

20      In this mode, the general defect map information is displayed. The user can enter a description for the map, as well as a defect code for each defect visited. Clicking the **List** button pops up the defect list dialog box. When a defect map is saved, the entered description and the defect code is also  
25      saved.

5.4.2.7.2.1. Automatic use of defect library Function

30      The definition of the defect code is determined by the defect code file used. A defect code file can be created, edited, loaded and saved by clicking on the **File** menu on the top menu bar, and selecting the **Defect Code** item.

#### 5.4.2.8. Text and status window

⇒ See Figure 44

5 This window displays text based status information. This window can be in one of four display modes: ***Defect status, deskew, defect selection, and wafer map edit.***

##### **5.4.2.8.1. Defect status:**

10

This mode is optimized for routine defect reviews. Defect code entry, defect selection and deskew can all be initiated here. The status window displays name of wafer map, product ID, Lot ID, process ID, number of defects displayed, number of total defects, operator name, defect site number, defect die number, defect XY coordinates.

15

##### 5.4.2.8.1.1 Routine defect classification

20

Entering a defect code to the defect entry field tags the defect code to current defect and move the stage to next defect in the defect list. If another status window is visible showing the defect list, the defect code will be updated and the next defect will be high-lighted. Defect code entries can be made using only the numeric key pad, in ***NumLock*** mode. Pressing "/" key in the numeric key pad high-lights the defect code entry field.

25

##### 5.4.2.8.1.2 EditDefect

30

The user can **append** a new defect or **remove** an existing defect. Bring up the Edit Defect Dialog Box by clicking on the

EditDefect push-button. To add a defect, type in its XY coordinates, its XY size, and its Intensity, and then click on the **AppendDefect** button. This defect will be appended to the Defect List. To remove an existing defect, select it in the defect list, and then click on the **RemoveDefect** button. This will change its Event Code to negative 1, which indicated this defect is not to be reviewed. Defects with their defect codes set to -1, will not be displayed if the Sorting routines are set to exclude defects with negative defect codes.

5.4.2.8.1.3 DeskewOff

5.4.2.8.1.4 CalcMapOffset

5.4.2.8.1.5 StartDeskew/ StartDefDeskew

See 5.4.2.8.2 for more detail

5.4.2.8.1.5 Prev Def

Click on this button to move the stage to the previous defect in the wafer defect list.

5.4.2.8.1.6 Next Def

Click on this button to move the stage to the next defect in the wafer map defect list.

**5.4.2.8.2. Deskew**

Deskew is the process of determining the amount of correction needed for pattern and wafer placement errors, so

that the XY coordinates in the wafer map correspond precisely to the LIS1000 stage.

5 The deskew procedure consists of matching three reference points in the wafer map: the two deskew points and the one event origin point. Optionally, the user may use the offset feature to position defect to a particular part of the screen. After the deskew process, all stage coordinates are corrected to match those in the wafer map.

10 The LIS1000 system supports two deskew modes: patterned wafer deskew and bare wafer deskew. In a patterned wafer deskew, the deskew information from the wafer map is used to deskew. In a bare wafer deskew, the deskew information  
15 from the wafer map is not used, instead, user selected defects are used as deskew points.

#### 5.4.2.8.2.1. DeskewOff

20 The deskew procedure involves many steps, any time within the procedure, one can click on the **DeskewOff** button to cancel the effect of deskew, so one can start the process over again.

#### 5.4.2.8.2.2. Start Deskew( patterned wafer deskew)

25 Click on **StartDeskew** to start the normal deskew process. The XY Stage will be moved to the first deskew point, usually at the lower left hand corner of the first die from the left, on  
30 the middle row of the wafer. After autofocus, manually move the stage, with the joystick, until this lower left corner is at the center of the screen.



5.4.2.8.2.3. Accept Deskew1( patterned wafer deskew)

5 Click on the ***Accept Deskew1*** button to confirm the first  
deskew point and to move the stage to deskew point 2,  
usually at the lower left hand corner of the last die on the  
middle row. Manually move the stage with the joystick until  
this lower left corner is at the center of the screen. Generally,  
the amount of adjustment in the X (horizontal) direction will be  
small, perhaps up to few ten's of microns (within a quarter of  
10 the screen at 5X). An adjustment much larger than this  
indicates error, perhaps using the wrong die. The Y (vertical)  
direction will be larger, because wafer rotation is being  
corrected.

15 5.4.2.8.2.4. Accept Deskew2( patterned wafer deskew)

Click on the ***Accept Deskew2*** button to move stage to the  
reference event origin. This reference point is usually the  
middle of a street intersection, and is defined by the patterned  
20 wafer scanner used. Using the joystick, move this origin to  
the center of the stage.

5.4.2.8.2.5. Accept Deskew 3( patterned wafer deskew)

25 Click on the ***Accept Deskew3*** button completes the deskew  
process.

5.4.2.8.2.6 StartDefDeskew( bare wafer deskew)

30 Click on this button to start the process of bare wafer deskew.  
There is no stage movement, but the button's label changes to  
AcceptDsk1. The User is then expected to move the stage to  
a large (one micron or larger) defect from the defect list.

5.4.2.8.2.7 AcceptDsk1( bare wafer deskew)

After the user is sure that the defect displayed on screen is the one specified in the defect list, he can click on the ***AcceptDsk1*** button to accept that for the first deskew. After clicking the ***AcceptDsk1*** button, the stage will not move, but the label changes to ***AcceptDsk2***. The user is then expected to move the stage to another large (one micron or larger) defect from the defect list, preferably one that is a large distance to the right of the first.

5.4.2.8.2.8 AcceptDsk2( bare wafer deskew)

After the user is sure that the defect displayed on screen is the one specified in the defect list, he can click on the ***AcceptDsk2*** button to accept that for the second deskew. The deskew process is completed after the second deskew is accepted, and the button label changes to ***StartDefDskew***

5.4.2.8.2.9 SetDeskew1

For non-KLA formatted maps, the user can set the deskew points to current stage coordinates. This feature is useful in creating a defect map, setting the deskew points and then adding the defects.

Click on the SetDeskew1 button to define or reset the first deskew point to current stage position. A confirmation dialog box will pop up and the user needs to confirm to actually set the deskew1 position.

5.4.2.8.2.10 SetDeskew2

Click on the SetDeskew2 button to define or reset the second deskew point to current stage position. A confirmation dialog box will pop up and user needs to confirm to actually set the

deskew2 position. Note that once the deskew1 and deskew2 are defined from actual stage position, the third deskew point (the registration point) is not needed for deskew. Although the deskew process still involves going through the three deskew points, the user should accept deskew3 without moving the stage.

#### 5.4.2.8.3. Defect list

Display the **Defect List**. The user can select a defect from this list and command the stage to view the selection, as well as view the next one on the list.

#### 5.4.2.8.4. Wafer Map Edit.

Some of the fields in the currently loaded wafer map which can be edited are: **Wafer number**, **Lot ID**, and **File Description**.

#### 5.4.2.8.5 Search

A method of systematic spiraling search movement is implemented through this option. The user can specify the step size in each move: One die distance or 80% of the current field of view of the laser scan. The total area covered is through the spiral moves of a 5 by 5 move pattern. A set of 25 buttons are assigned for the corresponding search sections. Clicking on one of the 25 button moves stage relatively to that point. The button corresponding to the current search position is colored green, the rest white. Clicking on the **StartSearch** button advances to the next step in the search sequence. Clicking on the **StartSearch** button and holding down the button starts the search and continues the search, until the button is lifted. Clicking on the

**ResetCenter** button sets the current stage position as the center, the starting point of the search pattern, without any actual stage movements. All search movements are relative to current position, i.e. the user can independently use the joystick to move the stage.

#### 5.4.3. Screen

A screen is a particular layout of windows: window type, size, and positions. A set of default screens are provided and can be selected from the **Screen** menu item on the top menu bar. When switching screens, the window contents are also transferred. No data reloading is required if the old and new screen have the same number and type of windows. All pre configured screens have the same number and type of windows:

- Three 2D windows
- Two 3D windows
- One cassette window
- One wafer map window
- One laser window
- One microscope window
- One image library window
- Two text and status windows

##### 5.4.3.1. Switching to a different screen

There are two ways to switch to a different screen:

##### 5.4.3.1.1. Select a screen from the *Screen* pull-down menu

#### 5.4.3.1.2. Select a screen from the current recipe:

If a recipe has been selected, three buttons will appear at the bottom of the screen. These buttons select screens from the screen list in current recipe. *Next* selects next screen in the list. *Last* selects screen in front of current screen in the list. The middle button allows user to select from the list of available screens.

### 10 5.5. User privilege

The Ultrapointe Model 1000 software supports multiple user with different privileges and password. Each user has his own file directory where other normal user can not access. Three user privilege levels are supported: normal, administrator, and service engineer. Normal users are blocked from accessing functions such as calibration, and diagnostics. Administrators can access these functions, but are blocked from user account changes. Service engineer can access all functions. In the *User Accounting* dialog box, one can add, delete, modify users.

### 5.6. Operator Mode

25 The Ultrapointe Model 1000 should be started by a user with administrator or service privilege.

#### 5.6.1. System Startup

30 The Model 1000 is configured so that on power up, the Model 1000 software is automatically started and ready for user

Login. Users are not expected to have to deal with the underlining operating system.

5.6.1.1. Laser imaging system Login:

5

⇒ See Figure 45: Copyright Screen

This copyright notice screen appears first and the user needs to click on the *Continue* button to go onto next step.

10

⇒ See Figure 46: Login dialog box

Click on one of the user names in the *user name list box* to select a user. If password is used, enter the password by clicking on the *password entry field* and then type in the password. Click on the *login* button to login as the selected user.

15

As shipped from the factory, only one user, "service" is installed with no password assigned.

20

5.6.1.2. Hardware initialization:

⇒ See Figure 47

25

Press *OK* to initialize Hardware dialog box.

After checking that the system is safe to start, click on the *OK* button. System hardware initialization takes several seconds. After hardware initialization, the manual stage control dialog box appears.

30

#### 5.6.1.2.1. Stage initialization

5 Click on the **InitStage** button to initialize the XYZ stage. Stage initialization can take up to 20 seconds. After the stage is initialized, the user must approve homing of the X and Y axes. It can take up to 20 seconds to home all axes.

#### 5.6.1.2.2. robot initialization

10 Swap in a **Cassette** window, if necessary. Click on the **Initialize** button in the cassette window to initialize the robot.

#### 5.6.2. Wafer loading, cassette window

15 If a wafer map is available, it must be loaded before the wafer itself is loaded on the stage. To load a wafer map: click on the **WaferMap** item on the **File** menu. This pops up the **wafer map file selection** dialog box. Select the desired wafer slot in the **Cassette** window, and then click on the **LoadWafer**  
20 button. The robot will fetch the selected wafer from cassette, pre-align, and place it on the stage. The wafer map corresponding to the selected wafer slot is also loaded and displayed in the **Wafer** window.

#### 25 5.6.3. Deskew, Text Status window

30 If a wafer map is used, the map's defect locations must first be corrected by the deskew process. The deskew process involves user visually matching 3 predefined wafer locations using the joystick. The first two are generally the lower left corners of the middle left and right die. The third deskew point

is generally the lower left corner of the center die. The *textstatus* window provides the user-interface for the deskew process.

5           After the wafer is loaded, the system starts the deskew process automatically: The stage moves to the first deskew point, and autofocuses. The deskew button on the textstatus window displays ***AcceptDeskew1***. After manually moving the stage to precisely the first deskew point, the user clicks on  
10           the ***AcceptDeskew1*** button. The system then moves to the second deskew point and the deskew button displays ***AcceptDeskew2***. The user again moves the stage to match the second deskew point and then clicks on ***AcceptDeskew2***.  
15           The deskew process is completed when the third deskew point is accepted. After the deskew process is completed, the system moves to the first defect on the wafer map.

#### 5.6.4. Defect location, Wafer window

20           The user can easily select a defect location to review with the Model 1000. He can point-click and GoTo the defect in the wafer window's two wafer map views. Or, he can use the defect map list in the text status window to select a site in the defect list. In addition, entering a defect code to the  
25           current defect also moves the stage to the next defect. Once the site is selected and the stage has moved to that site, the system does an autofocus automatically. All imaging control functions, laser and white light, are provided in the ***Laser*** window. In the ***Laser*** window, there are two autofocus  
30           buttons: The fine autofocus and the stage autofocus. One starts with the ***stage autofocus***, and then uses ***fine autofocus***



to bring sample into precise focus. The workings of the stage autofocus are pre configured, and no normal user adjustment is required. Fine autofocus has a range of about 50 microns and works over a wide range of sample intensity. However, if the sample is very dark or very bright, the user may need to adjust the laser intensity in the *Laser* window.

#### 5.6.5. Image Acquisition, Laser window

Once the desired sample is brought into focus, it can be laser imaged. The Laser window can be in one of two modes: *White Light Only* to get the best white light camera image, and *Laser and White Light* for both laser scan and camera image.

##### 5.6.5.1. White light only.

Select this mode when the best white light microscope image is desired. The laser shutter is closed, and only the white-light illumination falls on the sample. No filters are inserted except when selected by the user. Use:

5.6.5.1.1. *CamIntensity* slider to select camera image intensity.

5.6.5.1.2. *CameraMode* to select normal, dark field and *polarizer*.

##### 5.6.5.2. White light and laser

In the *White light and laser* mode both the white light and laser image are available, and one can access all laser scan functions. The yellow filter is inserted in the white-light light path to prevent the laser light from saturating the video camera.

NOTE: When the image becomes much brighter than saturation, the laser image sensor may shut-down, causing the laser image to become dark green. Lowering the *LaserIntensity* slider setting will reset the sensor.

5

#### 5.6.5.2.1. Setting scan range and scan start, manual method

Move the *Current Z* slider up and down and observe that the image intensity changes. The intensity displayed is adjustable with the *LaserIntensity* slider. Bright white color indicates saturation and dark green means no light sensed. On the brightest slice, the image intensity should be adjusted to be as bright as possible, without being saturated. Laser scan always starts down from the current fine Z position. Use the *Slices* and *StepSize* sliders to set the scan range and step size. Or, move the *CurrentZ* slider to the top of scan and click on the *S* button, then move the *CurrentZ* to the bottom of the scan and click the *F* button. The *F* button sets the *StepSize* slider for the selected slices to cover the whole range.

20

#### 5.6.5.2.2 Setting scan parameters, automatic method

Clicking on the *SetZ* button starts the automatic scan parameter setting procedure. The system performs a search 20 microns around current Z, sets scan start, scan range, and optimal laser intensity. The *SetZ* command takes about 3 seconds.

25

#### 5.6.5.2.3. Image acquisition

Clicking the *StartSurface* button starts the surface acquisition and clicking the *VolAcq* button starts the volume acquisition. In volume and surface acquisition, slices of 2D images are

30

taken starting at current fine Z. position and stepping down at *StepSize* interval for *Slices* number of slices. In volume acquisition, only 2D windows are updated with the new volume data, surface image is not generated. In surface acquisition, a surface is generated and displayed on all 2D and all 3D windows, but the volume itself is discarded.

#### 5.6.6. Screen Printing

Clicking on the ***PrtScreen*** menu on the top menu bar prints the whole screen. After a screen is captured, taking about 10 seconds, a confirmation dialog box is opened automatically. The user can then click on the ***OK*** button to confirm the print or ***Quit*** to abort.

The number of copies and the type of printer can be adjusted from the ***Printer*** dialog box, opened by clicking on the ***Printer*** item from the ***Config*** Top menu bar.

#### 5.6.7. Image storage

Both surface and volume data can be stored on the hard disk. Since the size of a volume data is generally large, one may want to limit the number of volume data files saved on disk at any time. Current volume data can be save by popping up the ***Image*** data file selection dialog box from the top menu. Surface data can be saved by popping up the image data file selection dialog box from a 3D window.

#### 5.6.8 Defect Classification

A defect classification code is entered into the text status window's defect code entry field. After the code is entered,

the stage moves to the next defect in the defect list. Windows displaying the defect list are updated with the entered defect code. An **Entry field** accepts keyboard input only if it has the input focus, i.e., the border of the entry field is a thick dark frame. Clicking within the entry field make the entry field to have input focus. Most of the defect code entry can be done using only the numeric key pad portion of the keyboard. Pressing the "/" key on the numeric key pad at any time forces input focus to the defect code entry field.

### 5.7. Engineering Mode

Engineering mode operations are those that deal with system configuration, recipe editing, calibration and maintenance. Most of these functions are restricted to administrators and service engineers.

#### 5.7.1. User type and Operator access control

The Model 1000 defines three classes of users: regular user, administrator, and service engineer. Regular users has no access to other user's directory. The Administrator has access to all directories.

Regular users also are restricted from using most of the engineering functions.

#### 4.7.1.1 User accounting

Clicking on the **UserAcct** item from the User menu on the top menu bar pops up the user accounting dialog box. From this dialog box, the Service Engineer can add or delete a user. A

user's password and privilege status can also be changed here.

Changes will only take effect after being saved.

5

In case of error conditions, or for diagnostics, a user logged in as administrator can popup the manual control dialog boxes: the robot, and the stage manual control dialog box.

10

### 5.7.2. Recipe preparation

15

Many Model 1000's functions can be customized for the user and stored in a recipe file. A recipe with a file name the same as the user is automatically searched for in the user's directory when the user logs in. If found, this recipe is loaded and the system configure itself according to this recipe. Any Recipe can be edited by popping up the ***Recipe Edit*** dialog box.

20

#### 5.7.2.1. Screen setup and sequencing

25

In the recipe dialog box, one can create a sequence of screens. The screen layouts are selected from a set of pre-installed screens. The user can swap the window in the screen. Clicking on the window brings up an ***Edit Window parameter*** dialog box where one can set the default conditions for that window. Current recipe screens are accessible via three buttons in the bottom screen bar

#### 5.7.2.2. Laser parameter control

Laser scan parameters can also be customized in the recipe. Clicking on the **LaserParams** button brings up the **Laser and Imaging control** dialog box. One can set the desired laser scan parameters here so that the Model 1000 will come up this way when user logs in. A warning message: **EditingOnly** is displayed indicating no actual hardware is being controlled.

#### 5.7.2.3 Sort criteria

All defects from the wafer map are filtered using the sort parameters. A number of sort criterion are provided, and they can be edited into the recipe. Clicking on the **SortParms** button brings up the **SortParametersEdit** dialog box. In this dialog box, one can edit the sort parameters. See **DefectList Sort** menu item from the top Menu bar for more information on editing sort parameters.

#### 5.7.3. Manual Robot control

The manual robot control dialog box is provided for single stepping the robot operation. The user can manually instruct the robot to get or put a wafer. Wafer source and destination are selected from the list box: stage, pre-aligner (flat-finder), and the cassette slots.

During robot movement, clicking the **Stop** button stops the robot. After the robot is stopped, a dialog box pops up so that user can continue the interrupted move, or reset the robot.

#### 5.7.4. Manual stage control

5        *Manual stage control* provides a list of elementary stage functions. All operations in this dialog box are in motor or encoder counts, except where noted. The user can use this dialog box to:

##### 5.7.4.1. Enable/Disable joystick, all axes

##### 10        5.7.4.2. Report Position/encoder count

##### 5.7.4.3. Report sample thickness

15        Manually focus on the sample, then select the ***Z Dist Fr Config Limit*** function to display sample thickness. Since Z config Limit is the Z stage position at which the system is focused on the stage itself, Z Distance Fr Config Limit (when sample is focused) equals wafer thickness, in microns. This function assumes the maximum z travel has been set properly.

##### 20        5.7.4.4. XYstage burn-in

25        The ***XYstage burn-in*** function commands the XY stage to move to random positions for a number of times. User specifies the number of times in the text entry.

##### 5.7.4.5. Stop

30        This function stops all axes.

#### 5.7.4.6. XYTo Load/unload

This function moves the XYZ stage to the loading/unloading position.

5

#### 5.7.4.7. Move, absolute or relative

Move current axis to the specified count number, or to the specified number relative to current position.

10

#### 5.7.4.8. Set speed

#### 5.7.4.9. Set acceleration

15

#### 5.7.4.10. Move to Home

Move current axis to its home position. Home for Z is at the z axis's full-down position.

20

Looking at the machine from the front, home for X is at the back, home for Y is at the right.

25

### **5.7.5. Hardware diagnostics**

The *LonDiagnostics* dialog box is provided for subsystem diagnostics. In this dialog box, most of the common functions are directly accessible through slider or option menus. All functions are accessible through the Node and NetVariables list box.

30



### 5.7.6. System Calibration

5       The Laser scan dimensions can be calibrated by popping up the **Calib** dialog box. Clicking on the **Calib** item from the Maintenance menu on the top menu bar pops up the dialog box.

#### 5.7.6.1. Laser scan XY dimension calibration

10       The model 1000 maintains a calibration table for the five objective's four scanner resolutions. This table can be calibrated by using the **LaserScan Calibration** dialog box.

### 5.7.7 System configuration

15       Some of the items in system configuration can be accessed from within the LIS1000 operating software. Currently, the following items can be edited in the Misc dialog box:

20               Laser saver timer  
                Deskew die offset  
                Pre-aligner (Flat finder) angle offset  
                White light mode  
                Next die position

Click on the **Misc** menu item from the **Config** top menu pops up the System Configuration dialog box.

#### 25       5.7.7.1 Screen/Laser saver timer interval

Set the timer to the desired interval in minutes,

#### 5.7.7.2 File DeskewDie Offset

30       Moves deskew points away default, The number entered determines the number of die deskew points are to be moved toward the wafer center.

#### 5.3.5.3.3 Flat angle offset( deg.)

If the wafer pattern is slightly skewed with respect to the wafer flat, or if there is an error in the robot teach angle, the wafer pattern on the monitor screen will appear to be tilted.

5 A stage movement along the wafer's horizontal street will cause the wafer pattern to move vertically as well. This XY tiling can be removed by adjusting the *Flat angle offset*.

10 This offset parameter allows fine tuning of the angle in which a wafer is loaded onto the LIS1000 stage. The adjustment resolution is .125 degree( 1/8 of a degree) and typical value is within +/- 0.5 degrees.

#### 5.3.5.3.4 White light mode

15 Two white-light-only options are available: white light with a beam splitter cube *in* or with it *out*. The beam splitter is required in the laser scan mode to get simultaneous white light and laser viewing, as well as in autofocus. However, with the cube in, the white light intensity is cut by 50% (although image quality remains unchanged, and perceived image  
20 brightness is down by marginal amount.) For very dark samples, one therefore would want to switch out the cube in white-light-only mode. But, since the cube has to be *in* in the laser mode, the beam splitter cube will be switching in and out every time the laser and white-light-only mode is changed.  
25 The beam splitter switching is slow and can be a throughput bottleneck. To improve throughput, user has the option of keeping the beam splitter *in* in the white-light-only mode.

30 Cube Out: Brightest white light image, but switching between white-light and laser can be slow.

Cube in: Slightly dimmer image for the same camera intensity, but otherwise no changes in image quality. Significantly better throughput and reliability.

5      5.3.5.3.5 Next die position

When viewing a defect, sometimes it is useful to compare the image with those at the same die location, at the next die. The LIS1000 has a button allowing user to go to the same location on the next die. The user uses the *Next Die position* to specify where the next die is located: Bottom/Left/Top/Right.

15      **5.7.8. System Startup**

The Ultrapointe Model 1000 is configured so that on power up, the Model 1000 software is automatically started.

20      In the case where user may need to get out of the Model 1000 software for computer system usage, the following is provided as a guide to manually starting the model 1000 software.

25      5.7.8.1. Operating system and windowing system

The Ultrapointe Model 1000 software runs on the Silicon Graphics workstation over the IRIX operating system, SGI's version of UNIX. User interface is implemented using the X-Window windowing system, and Motif tool kit. The Model 30      1000 software is designed to provide all necessary functions without user having to interact with the underlining operating

system or its user interface. However, IRIX also provides a graphical user interface, the WorkSpace, for system access, Please refer to the IRIS WorkSpace User's Guide for more information on the **WorkSpace** and **X-Windows** windowing system.

5

#### 5.7.8.2. IRIX Computer user autologin and program autostart:

IRIX is a multi-user operating system, and on power up, it will ask for user login. One can by pass this computer login step by assigning a user with Autologin option. A user "uv" has been installed, without password, but with autologin at the factory. This manual is written assuming you have logged in as "uv". Additional users and passwords can be setup using the **System Manager** tool from the **Toolchest** system menu on screen. After logging in, the **X Window** windowing system is started, and the **toolchest**, **console icon**, and the **workspace window** appear. The model 1000 is also configured at the factory to Auto-start the Model 1000 software. With autologin and auto start, the Model 1000 software runs after power up without user intervention. The file **.sgisession** in the user's home directory specifies the program to run when user logs in. Alternatively, especially after exiting from the Model 1000 software, double-clicking the mouse's left button on the **uv1000** icon in the WorkSpace window starts the Ultrapointe model 1000 system software.

10

15

20

25

## 6. Troubleshooting

### 6.1. Robot operation

5           The *manual robot control* dialog box is provided for single stepping the robot operation. The user can manually instruct the robot to get or put a wafer. *Source* and *destination* are selected from the list box: stage, flat finder, and the cassette slots.

10

During robot movement, clicking the *Stop* button stops the robot. After the robot is stopped, a dialog box pops up so that user can continue the interrupted move, or reset the robot. The Model 1000 communicates with the robot through a RS232. port. Robot can be disabled and assigned to different port by editing the */usr/uv.uvconfig* file. The positions of the robot arm are programmable and are stored in the robot controller. Please refer to the MECS robot user manual for detail.

15

20

### 6.2. Stage Operation

A list of elementary stage functions are provided in the manual stage control dialog box. All operations in this dialog box are in motor or encoder counts, except where noted. The Model 1000 communicates with the stage through the RS232 port. The stage can be disabled and configured by editing the */usr/uv.uvconfig* file.

25

30

### 6.3. Power fail Recovery

### 6.3.1. Wafer recovery

### 6.3.2. System component reset

## 5 6.4. LON diagnostics dialog box

All system functions can be accessed through this dialog box.  
Please use this dialog box with caution, as no safety check is  
made when accessing the system in this way.

10

### 6.4.1 Wafer door interlock.

Wafer door and laser shutter are wired together, so that  
when the wafer door is open, the laser is blocked.

15

### 6.4.2 Critical parameters

Some parameters work in groups, arbitrary resetting one  
without also changing the others in the group may cause  
damage to the system. The page scanner controller's Vertical  
Sync, TraceTiming, Retrace timing, StartAmp offset, Incr and  
20 Decr are all working together as a group.

## 6.5. System Service

### 6.5.1. Software configuration

25

All configurable software parameters are stored in  
*/usr/uv/.uvconfig*.

See appendix B for details.

30

## 6.6. System Messages and Response

### 6.6.1. Caution/Information

### 6.6.2. Error Condition and Recovery

#### 5      6.6.2.1 Floppy disk not responding

The Model 1000 assumes the floppy drive is enabled, or "**Mounted**". If the floppy is unmounted, it will be inaccessible. Use the system manager in the Toolchest to check if floppy is mounted. See Appendix A.4.4 for floppy drive installation.

#### 10     6.6.2.2 Printer not printing

##### 6.6.2.2.1 SCSI connected printer, i.e.: Kodak XL7700:

The power-up sequence of the printer and the Model 1000 is important: the printer must be ready before the Model 1000. That is: make sure the printer is turned on first, wait till the front panel on the printer says it is ready, then turn on the Model 1000. The SCSI cable between the Model 1000 and the printer can be connected or disconnected only when both are turned off. If the printer has been off or the Model 1000 has been off, the order of power-up must be observed.

20     This printer is managed by the IRIX **Print Manager**, and printer status is available from the print manager. Printing is spooled by the IRIX in the hard disk. Do not queue up more than a few prints, as doing so may severely slow down the system, and may cause printer failure. Click on the system item on ToolChest and select **Print Manager** for status.

25     Trouble shooting: If after checking everything and the printer still is not printing, cycle power to the printer, wait till printer is ready, shutdown the Model 1000 and start again. The prints that are queued can be deleted in the **Print Manage** dialog box.

30

##### 6.6.2.2.2 Network connected printer, i.e.: Codonics NP600

A network printer needs to have the correct address configured before it can accept print files. See Appendix A.4.3 for network printer installation. Once installed, the network print should need no special power-up procedure.

5

Trouble shooting: If after checking everything and printer still is not printing, cycle power to the printer, wait for 3 minutes and print again. If the printer still is not printing, reset the print queue by re-starting the UV1000 and try again.



## **7. Applications Information**

### **7.1. Wavelength Selection**

5 The Model 1000 allows selection of different laser wavelengths in order to optimize viewing of the wafer surface. Selection depends on the spectral reflectivity, absorption, and transparency of the materials in the field of view, and the purpose of the observation.

10 The Model 1000 provides for selection of three primary argon ion wavelengths, 458 nm (violet), 488 nm (blue), and 515 nm (green). (Optional extra lines are also available). All give good resolution because they are at the shorter wavelength end of the visible spectrum, but the response of various materials and film thickness is quite different for each line. For example, polysilicon reflects the 515  
15 nm line strongly, but not the 488 or 458.

The selection of wavelength also depends on whether the top surface is to be viewed, or whether a sub-surface feature is of interest. One value of wavelength almost invariably is most penetrating in a given situation, another is most reflected. Since materials and processes  
20 vary so much, experimentation is required in most cases to select the optimum. This information may then be entered into a recipe assigned to the product and process step of interest, and recalled by the operator during routine review and inspection operations.

### **7.2. Image interpretation**

#### **7.2.1. Surface reflectivity**

25 The Model 1000 maps gray scale intensity into false color to enhance contrast and visibility of small differences in reflectivity.

### 7.2.2. Limits of performance

XY resolution and vertical resolution are a trade-off for very small structures. While the system can "see" individual  $0.1\ \mu\text{m}$  particles and resolve slightly larger structures horizontally, it cannot simultaneously resolve in the z (vertical) dimension until the x-y structure is on the order of  $0.3$  to  $0.4\ \mu\text{m}$ . When the x-y structure is of sufficient area  $>0.5\ \mu\text{m}^2$ , then the high resolution of the system vertically (approx.  $25\ \text{nm}$ ).

The system in its simplest mode works on opaque, diffusely scattering surfaces.

Vertical surfaces are not well resolved.

Calibration latex spheres are a problem (as with all imaging optical instruments), because they act as tiny lenses, refracting the light and causing false images. Small, rounded transparent defects can have a similar effect..

Too small an object is a scatterer only, and shows up as an extinction object against a flat background.

Imaging down small diameter, high aspect ratio cylindrical holes of round, square, rectangular, or other cross sections can suffer when multiple reflections create false focuses.

Light scattering and interference effects can cause distorted images of sharp edges or high contrast change boundaries.

### 7.3. Complex surfaces

#### 7.3.1. Transparent coatings

5 Many types of photoresist are transparent as are dielectric materials. These and other transparent media can cause confusion with the simple surface extraction process of the SDP. In the surface acquisition mode, the SDP selects the z coordinate at each xy coordinate which returns the brightest signal, and interprets it as a  
10 point on the surface. All surface points are connected to produce the display. If a (semi)transparent material covers a dark area, the brightest signal will come from the transparent surface. If it covers a highly reflective material, then the lower surface will be brightest and represent a surface point. If contrast varies on the lower surface, the  
15 SDP will produce an interpretation of the surface which may be an amalgam of regions of the two physical surfaces. In these cases, some of the advanced processing features may be needed to extract a single surface.

#### 20 7.3.2. Multilayer structures

In the case of two or more layers to be viewed, the data should be taken with the *VolAcq* control, and viewed with the *2D XZ* or *YZ* views. The *New Surface* feature may be used to extract a single surface, or the *XZPeaks* or *YZPeaks* feature may be used to view the  
25 multiple surface profiles.

#### 7.3.3. High contrast variations

30 Some surfaces, such as rough metal lines, will have a very large range of reflectivity in the field of view, which exceeds the dynamic range of the photo detector electronics. In order to see some of the lower

contrast area, it may be necessary to increase *Laser Intensity*. This, in turn may cause the highest reflectivity portions of the image to saturate. If the saturation occurs over more than one step in Z, the SDP will be unable to select the brightest surface point at that XY location, and will treat the surface as lying at the first point where saturation occurred. This may distort the surface shape in the saturated region, which will appear white in color in the false color intensity display.

#### 7.3.4. Vertical Surfaces

Virtually no reflected light is returned by vertical or near-vertical surfaces. While the surface algorithm shows a vertical surface, it is dark green to black in color, indicating no light was detected. In very dark areas, background noise may be detected by the surface algorithm, depicting a very rough, non-physical surface. Again, this is shown as dark green to black in color, indicating no light was detected.

## **8. System Maintenance**

### **8.1. System software maintenance**

#### **5        8.1.1 Setting Date and Time**

Exit Model 1000 software, and open a *Winterm* for a command prompt by clicking on the *Tool* item on the *Toolchest* menu and select *shell*. From the command prompt enter:

10                date mmddhhmm < enter>

example

date 06091415     < enter>     setting date/time to  
June 9,2:15

15                Close the Winterm by double clicking the icon on the top left corner of the window.

#### **8.1.2. System shutdown**

User must follow this system shutdown procedure before turning power off. Exit Model 1000 software, Click on the System item on the Toolchest and select System shutdown.

20                A confirmation dialog box will popup for system shutdown confirmation: click Yes to shutdown or No to cancel request. The message " Okay to power off the system now" appears when the system is ready for power down.

### **8.2. Maintenance and Laser Safety**

25

All maintenance procedures are to be performed with the laser power off. The system is equipped with interlocks to prevent accidental power-on status of the laser with the system covers removed. Never attempt to defeat the interlocks during normal use or maintenance.

30

### **8.3. Customer Preventative Maintenance Procedures**

See Appendix E

### **5 8.4. Extraordinary maintenance**

These procedures are to be followed in unusual circumstances.

#### **8.4.1. Wafer removal**

10

#### **8.4.2. Other**

## 9. Appendices

### Appendix A

- 5 The *InstallSystem* Tape is for system software installation and the *InstallApplication* Tape is for application software installation.

One would install the SystemTape, then the ApplicationTape, and then the Peripherals.

10

#### A.1 To Install the Indigo System Software

1. If the Hard Disk needs to be Formatted, please see Section A.11

2. If the Hard Disk already has the operating system installed,

15 Skip to Section A.2, otherwise continue.

3. The system normally has the TapeDrive already on the system's SCSI bus with a default SCSI ID of 7.

Load the InstallSystem Tape into the TapeDrive.

4. Bring up the System Maintenance Menu by:

20 Powering on the system.

Pressing the Esc key, when prompted.

5. Load in the installation program by:

Selecting Option 2 (*Install System Software*) at the System Maintenance Menu.

25 6. If the Hard Disk has just been formatted, the installation program will ask if it's OK to make New File Systems, Answer Yes, each time.

7. When the Installation Menu appears:

Type 'install eoe2.man.bsdldr'

30 8. Type 'install eoe2.sw.bsdldr'

9. Type 'go'

10. After the 'Inst> prompt ' reappears, Type 'go' once again.

This is to install the maintenance software.

11. After the 'Inst> prompt ' reappears, Type 'quit'

12. Answer **Yes** to the Restart prompt.

5

## **A.2 To Install the Ultrapointe Application Software**

A.2.1. Add a user named 'uv' by:

Logging in as '**root**'

10 Select **System**, in the **ToolChestWindow** at the TopLeftHandCorner  
of the screen.

Select the System Manager

Select the **Users** icon.

Select **Add**.

Type in 0000 for **UserID**. (This gives 'uv' superuser privileges).

15 Select **Other** for Shell

Type in '/bin/ksh' for the **ShellName**

A.2.2. Bring up a **ShellWindow** by:

Selecting **Tools**, in the **ToolChestWindow** at the TopLeftHandCorner  
of the screen.

20 Select **Shell**

A.2.3. Login as 'uv' by Typing 'login uv'

A.2.4. Load the InstallApplication Tape into the TapeDrive.

A.2.5. Type 'tar xv'. (This program will extract files from the Tape  
onto the HardDisk)

25 A.2.6. Login once again by Typing 'login uv'

A.2.7. Execute the command file (script) 'InstallApplic' by Typing  
'InstallApplic'

A reboot will occur as the last step of the command file.

A.2.8. Turn on Autologin and turn on Workspace by:

30 Select **System**, in the **ToolChestWindow** at the TopLeftHandCorner  
of the screen.



Select the System Manager

Select the *Users* icon.

Select *uv*

Select 'On' for Workspace

5 Select 'On' for Autologin

Select Accept

A.2.9. To start the Ultrapointe Application program via the autologin,  
just type 'reboot'

### A.3 InstallApplic Script

10

```
# The following script assumes the printer's SCSI ID has
# been set to the default value of 6, and
# we are logged in as 'uv'
```

```
chmod a+rw /dev/ttyd1
```

15

```
chmod a+rw /dev/ttyd2
```

```
cp /usr/people/uv/Install/KodakXL7700.info /usr/lib/print
```

```
cp /usr/people/uv/Install/k77print /usr/lib/print
```

```
cp /usr/people/uv/Install/centface /usr/spool/lp/model
```

```
cp /usr/people/uv/Install/mkcentpr /usr/sbin
```

20

```
chown lp /usr/lib/print/KodakXL7700.info
```

```
chown lp /usr/lib/print/k77print
```

```
chown lp /usr/spool/lp/model/centface
```

```
chown root /usr/sbin/mkcentpr
```

```
chgrp lp /usr/lib/print/KodakXL7700.info
```

25

```
chgrp lp /usr/lib/print/k77print
```

```
chgrp lp /usr/spool/lp/model/centface
```

```
chgrp sys /usr/sbin/mkcentpr
```

```
chown lp /dev/scsi/sc0d6l0
```

```
chgrp lp /dev/scsi/sc0d6l0
```

30

```
# cp /usr/people/uv/Install/hosts /etc
```

```
cp /usr/people/uv/Install/printcap /etc
```

```

cp /usr/people/uv/Install/bsdipr /etc/init.d
mkdir /usr/spool/np1
cp /usr/people/uv/Install/kernel /usr/sysgen/master.d
mknod /dev/uldd c 70 0
5 mknod /dev/ulsd c 71 0
cat /usr/people/uv/Install/systemAddOn > >
    /usr/sysgen/system
cp /usr/people/uv/Install/uldd /usr/sysgen/master.d
cp /usr/people/uv/Install/ulsd /usr/sysgen/master.d
10 cp /usr/people/uv/Install/uldd.a /usr/sysgen/boot
cp /usr/people/uv/Install/ulsd.a /usr/sysgen/boot
rm /etc/inittab.O
rm /etc/ioctl.syscon.O
rm /etc/passwd.N
15 rm /usr/sysgen/master.d/kernel.O
rm /usr/sysgen/system.O
cp /usr/people/uv/Install/.Xdefaults /usr/people/uv
cp /usr/people/uv/Install/.sgisession /usr/people/uv
cp /usr/people/uv/Install/.profile /usr/people/uv
20 # cp /usr/people/uv/Install/uv /usr/people/uv
# cp /usr/people/uv/Install/.uvconfig /usr/uv
autoconfig -f
reboot

```

### A.3.1 Important Files

```

25 /usr/people/uv/.Xdefaults      (X Window Config File)
   /usr/people/uv/.sgisession    (Auto Login File)
   /usr/people/uv/.cshrc         (User Config File)
   /usr/people/uv/.login         (User Config File)
   /usr/people/uv/.profile       (User Config File)
30 /usr/people/uv/uv              (Application Program)
   /usr/people/uv/Install/InstallApplic (InstApplicSoftwrScript)

```

	/usr/people/uv/Install/InstallUpgrd	(InstUpgrdSoftwrScript)
	/usr/people/uv/Install/PrepareSystemTape	(PrepSystemTapeScript)
	/usr/people/uv/Install/PrepareApplicTape	(PrepApplicTape Script)
	/usr/people/uv/Install/KodakXL7700.info	(Printer Install File)
5	/usr/people/uv/Install/k77print	(Printer Install File)
	/usr/people/uv/Install/centface	(Printer Install File)
	/usr/people/uv/Install/mkcentpr	(Printer Install File)
	/usr/people/uv/Install/hosts	(Printer Install File)
	/usr/people/uv/Install/printcap	(Printer Install File)
10	/usr/people/uv/Install/bsdldr	(Printer Install File)
	/usr/people/uv/Install/kernel	(Kernel Build File)
	/usr/people/uv/Install/system	(Kernel Build File)
	/usr/people/uv/Install/uldd	(Kernel Build File)
	/usr/people/uv/Install/ulsdp	(Kernel Build File)
15	/usr/people/uv/Install/uldd.a	(Kernel Build File)
	/usr/people/uv/Install/ulsdp.a	(Kernel Build File)
	/usr/uv/.uvconfig	(Application Pgm Config File)
	/usr/uv/data/*.bmp	(Bitmap Files)
	/usr/uv/bruce/data/*.sur	(Surface Files)
20	/usr/uv/bruce/data/*.vol	(Volume Files)
	/usr/uv/bruce/data/*.map	(Wafer Map Files)
	/usr/uv/bruce/data/*.rcp	(Recipe Files)

#### A.4 To Install the Peripherals

25

All peripherals can be installed after the System Tape has been installed, except for the printer, which should be installed after both the System and Application Tapes have been installed.

#### A.4.1 Tape

A.4.1.1. The system normally has the TapeDrive already on the system's SCSI bus with the default SCSI ID of 7.

5 A.4.1.2. The TapeDrive is automatically installed when it's Detected on SCSI Bus.

#### A.4.2 Printer (Kodak XL7700)

1. Please note that the printer should be installed after the Application Tape has been installed.

10 2. Set the Printer's SCSI ID to 6 (Default)

3. Connect the printer to the system's external SCSI port. Power up the printer, 'Ready' should be displayed in the printer's LCD Window.

15 4. Reboot the system. Typing 'hinv' should output a line similar to

Printer: unit 6 on SCSI controller 0

5. Select System, in the ToolChestWindow at the TopLeftHandCorner of the screen.

6, Select the System Manager

20 7, Select the Printers Icon

8, Select Add

9, Type in a PrinterName (e.g. KodakColor)

10, If the printer is on the SCSI Bus, Select Parallel

11, Select the correct PrinterType, Select Accept

25 12, If the printer is on a remote Host, Select Network

Type in the RemoteHostName and the RemotePrinterName,

13, Select Accept

14, Select the Printer that has been just added.

15, Select Yes for 'Accepting Requests'

30 16, Select Yes for 'Printing Requests'

**A.4.3 Printer (Codonics NP600)**

1. Please note that the printer should be installed after the Application Tape has been installed.

2. Install Files from Codonics Installation Diskettes onto the Printer by:

Turning off the printer using the power switch at the BACK.

Inserting Diskette #1 into the disk drive.

Turning on the power. The printer will boot up and load in the software.

Waiting until there are no more beeps (5 minutes).

Ejecting the floppy and inserting the next diskette.

Continuing this until all the diskettes have been installed.

3. Enter its IP address (or network address) by:

Obtaining an IBM PC/AT keyboard.

Turning off the printer using the power switch at the BACK.

Plugging the keyboard into Port D at the back.

Turning on the power using the power switch at the BACK.

Turning on the power switch at the FRONT.

Typing Control-C when the front display is flashing the word 'NETWORK'.

Pressing ESC on the keyboard and waiting for 3 beeps. Redo this step if a typing mistake is made. Typing

'i' <ENTER>

Typing 127.0.2.200 <ENTER>

Connecting the ethernet network cable to the COMM PORT.

The word 'NETWORK' should be displayed (Not Flashing) on the front panel.

Turning off the power using the power switch at the BACK

Unplugging the keyboard.

Turning on the the power using the power switch at the BACK.

**A.4.4 Floppy**

1. Select System, in the ToolChestWindow at the  
TopLeftHandCorner of the screen.
- 5 2. Select the System Manager
3. Select the Disks&Files Icon
4. Select Floppy
5. Type in '/FLOPPY' for the mount directory
6. Select Mounted
- 10 7. Select ShowDetailedInfo
8. Select Yes for 'Mount At Startup'
9. Select Accept

**A.4.5 Network to Other Hosts**

- 15 1. Select System, in the ToolChestWindow in  
TopLeftHandCorner of the screen
2. Select the SystemManager
3. Select the Networking Icon
4. Select Add,
- 20 5. Type in the New HostName
6. Type in the IP Address
7. Select Accept

**A.4.6 Serial Ports**

- 25 1. Select System, in the ToolChestWindow in  
TopLeftHandCorner of the screen
2. Select the SystemManager
3. Select Serial Port1, Set Port1 to Off
4. Select Serial Port2, Set Port2 to Off

**A.5 To Install Incremental Software Upgrades**

- 30 1. Load the InstallUpgrade Tape into the TapeDrive.

2. Execute the script 'InstallUpgrd' by:

Typing 'InstallUpgrd'

#### A.5.1 InstallUpgrd Script

# The following script assumes we are logged in as 'uv'

tar xv

### A.6 To Archive Data Files onto QIC-150 1/4" Tapes

To archive or store Data Files on QIC-150 1/4" Tapes, please use the 'tar' command as described in Appendix C20.

### A.7 To Prepare an InstallSystem Tape

1. Using a ReadyToRun Indigo XS24 (R3000), with a CDROM Drive and a TapeDrive attached:

Load a CDROM containing the operating system IRIX 4.0.5F. into the CDROM Drive,

Load a blank tape into the TapeDrive.

2. Mount the CDROM onto the directory /CDROM by:

Selecting System, in the ToolChestWindow at the TopLeftHandCorner of the screen.

Select the System Manager

Select the Disks&Files Icon

Select CDROM

Type in '/CDROM' for the mount directory

Select Mounted

3. Copy Files From the CDROM to a Tape via a Local Directory (/usr/TAPE).

Make sure have enough Disk Space under /usr by:

Typing: 'df'

Execute the script 'PrepareSystemTape' by:

Typing 'PrepareSystemTape'

### A.8 PrepareSystemTape Script

```
# The following script assumes we are logged in as 'uv'
mkdir /usr/TAPE
5      cp /CDROM/RELEASE.info /usr/TAPE
      cp /CDROM/dist/eoe*      /usr/TAPE
      cp /CDROM/dist/motif*    /usr/TAPE
      cp /CDROM/dist/sa        /usr/TAPE
      cp /CDROM/dist/mr        /usr/TAPE
10     distcp -vw /usr/TAPE     /dev/nrtape
      rm /usr/TAPE/*
      rmdir /usr/TAPE
```

### A.9 To Prepare an InstallApplication Tape

- 15 1. Load a blank tape into the TapeDrive.
2. Execute the script 'PrepareApplicTape' by:  
Typing 'PrepareApplicTape'

### A.10 PrepareApplicTape Script

```
20  # The following script assumes we are logged in as 'uv'
      tar cv /usr/people/uv /usr/uv
```

### A.11 To Format the HardDisk

0. One would normally not have to do this.
- 25 1. Bring up System Maintenance Menu by:  
Powering on the system.  
Pressing the Esc key, when prompted.
2. Enter the Command Mode by:  
Selecting Option 5 (CommandMonitor) at the System
- 30 Maintenance Menu.



3. Load in the DiskManagement program from the tape by:  
Typing "boot -f tpssc(0,7)(fx.IP12)"  
where the number '7' is the default SCSI ID of the  
TapeDrive.
- 5 4. When asked to be in ExpertMode, answer Yes.
5. When the DiskManagement(fx) Menu appears:  
Setup the Disk automatically by:  
Typing 'auto'. This will take approximately 30 minutes.
6. Exit from the fxMenu to the System Maintenance Menu  
10 by:  
Typing 'exit'
7. Then continue with Section A.1, Step 5.

#### A.12 To Build a New Kernel (Operating System)

- 15 0. Normally one would not have to build a new Kernel,  
unless there are new Device Drivers to be added in.  
Ultrapointe's 2 Device Drivers are already included in  
the current Kernel.  
They are for the GMI Board and the SDP Board.
- 20 The GMI Board resides on the GIO Bus and has LON  
Network, and Serial Ports functionalities.  
The SDP Board has FrameGrabbing and  
SurfaceDataProcessing functionalities.
1. Execute the script 'BuildNewKernel' by:  
25 Typing 'BuildNewKernel'

#### A.13 BuildNewKernel Script

# The following script assumes we are logged in as 'uv'  
# The following 6 commands have already been executed  
30 # by the InstallAppl Script.  
# They are included here for information purposes.

```
# cp /usr/people/uv/Install/kernel /usr/sysgen/master.d
# cp /usr/people/uv/Install/system /usr/sysgen
# cp /usr/people/uv/Install/uldd /usr/sysgen/master.d
# cp /usr/people/uv/Install/ulsdp /usr/sysgen/master.d
5 # cp /usr/people/uv/Install/uldd.a /usr/sysgen/boot
# cp /usr/people/uv/Install/ulsdp.a /usr/sysgen/boot
#
# Other Vendors commands should go here.
cd /
10 autoconfig -f
reboot
```

## Appendix B

### **B.1. Introduction:**

5 The file: */usr/uv/.uvconfig* contains all software configuration information. It is an ASCII file that can be edited by using the *jot* editor. Currently, the *.uvconfig* contains the following parameters:

B.1.1 screen spec

10 B.1.2 user spec

B.1.3 confocal spec

B.1.4 stage spec

B.1.5 robot spec

B.1.6 objectivesspec

15 B.1.7 turret step size

B.1.8 system specs

The general format of *.uvconfig* is:

20 B.1.2.1 Each line is at most 80 character long, characters after this limit are ignored.

B.1.2.2 Lines that start with # are treated as comments.

B.1.2.3 Empty line is also treated as comment.

25 B.1.2.4 Every non-comment line is expected to be in the form of:

key\_word = values

### **B.2. Screen spec:**

30 line 0: screen\_spec = nn  
nn = number of lines that follows.

line 1:     screen\_name = name  
           name = descriptive name of the screen..  
 line 2...   window\_spec = WindowTypeKey xStart yStart width  
   height

5           WindowTypeKey is one of the following:

          slices  
           surface  
           cassette  
           rt\_confocal <the laser window>  
 10          referencelmng  
           defect\_map  
           text\_status  
           microscope

          xStart, yStart, width, and height are in screen pixels, in the  
 15          range of x: 0 to 1262, y: 0 to 913. These are less than  
           1280x1024, because of the border and top/bottom displays.

Example:

20          screen\_spec = 11  
           screen\_name = 2Surfaces  
           window\_spec = cassette     5   3 90 40  
           window\_spec = surface     100 3 90 40  
           window\_spec = microscope 5 46 90 40  
           window\_spec = text\_status 100 46 90 40  
 25          window\_spec = slices       5 89 90 40  
           window\_spec = rt\_confocal 5 132 90 40  
           window\_spec = defect\_map 240 4 411 266  
           window\_spec = surface     4 273 647 640  
           window\_spec = slices     654   4 604   453  
 30          window\_type = slices     654 460 604 453

### B.3. User\_spec:

line 0: user\_spec = nn  
 35

nn = number of lines that follows

line 1+: name directory password privilege

name is the name of the user

directory is the name of the user's directory from /usr/uv

5 password can be alphabet or number, or a single '\*' for none.

privilege is either 's' for administrator, or nothing for normal user.

Example:

```
10      user_spec = 2
      bruce bruce abc s
      joe   joe   *
```

User specs can also be modified from inside of UV: click on the UserAcct item from the Maintenance menu.

15

#### B.4. Confocal\_spec:

line 0: confocal\_spec = nn

20 nn = number of lines that follows.

line 1: confocal\_xy\_res = aa bb cc dd

numbers are with 0.01 resolutions for the four scan resolutions aa,bb,cc,dd are nominal resolutions, after the machine is

25 calibrated.

line 2: confocal\_z\_res = ee

actual fineZ step size, in nanometer

line 3: confocal\_x\_amp\_obj0 = m n o p

line 4: confocal\_x\_amp\_obj1 = m n o p

30 line 5: confocal\_x\_amp\_obj2 = m n o p

line 6: confocal\_x\_amp\_obj3 = m n o p

line 7: confocal\_x\_amp\_obj4 = m n o p

line 8: confocal\_y\_amp\_obj0 = m n o p

line 9: confocal\_y\_amp\_obj1 = m n o p

line 10: confocal\_y\_amp\_obj2 = m n o p

line 11: confocal\_y\_amp\_obj3 = m n o p

5 line 12: confocal\_y\_amp\_obj4 = m n o p

for all objectives, and page-scan XY amplitudes.

numbers are in integer, for the 4 scan resolutions, These are the actual page-scan x amplitude used to get the nominal XY resolution.

10 line 13: laser\_pwr\_at\_line = xx yy zz nn

For each of the 4 laser line choices, the laser power lever is specified here.

xx = laser power setting for the 458 line, yy zz nn are for 488, 515, and all. range is 0 to 125 for 0 to 25mw.

15 line 14: laser\_atten\_at\_458line = xx yy zz mm nn

line 15: laser\_atten\_at\_488line = xx yy zz mm nn

line 16: laser\_atten\_at\_515line = xx yy zz mm nn

line 17: laser\_atten\_at\_Allline = xx yy zz mm nn

For each of the 4 laser line and 5 objective choices, the laser power filter selection is specified here.

xx = laser power attenuation filter position setting for the first objective position, yy zz nn are for second, third fourth and fifth. range is 1 to 6.

1 = 0 attenuation, 2 = minimum, 3 = 0, 4 = medium, 5 = 0, 6 = high attenuation.

25 line 18: laser\_inten\_dark\_field\_offset = xx

Compensate for the differences in laser intensity when switching between bright and dark field.

xx = (laser intensity at bright field ) -

30 ( laser intensity setting to get equivalent intensity at dark field)

line19: pmtamp\_min\_offset\_n\_range = xx yy

xx is the minimum gain, can be -10 to 10. This is the gain sent to hardware if user selects 0 in LaserIntensity slider.

yy is the range of LaserIntensity slider. The LaserIntensity slider always display a range of 100, internally, 100 is linearly map into yy. yy can be between 80 to 120.

line 20: white\_light\_dark\_field\_power = nn

the number nn is the camera intensity used in dark field.

line 21: laser\_white\_switch\_cube\_1\_0 = n

n = 0: the beam splitter cube is always in, no switching between laser/white light mode chnages.

n = 1: beam splitter is in for laser, out for white light. This setting degrade system throughput and reliablitiy.

#### Example:

```

confocal_spec = 21
confocal_xy_res = 100.00 50.00 25.00 12.50
confocal_z_res = 13.35
confocal_x_amp_obj0 = 1426 713 356 280
confocal_x_amp_obj1 = 1426 713 356 280
confocal_x_amp_obj2 = 1426 713 356 280
confocal_x_amp_obj3 = 1426 713 356 280
confocal_x_amp_obj4 = 1426 713 356 280
confocal_y_amp_obj0 = 1426 713 356 280
confocal_y_amp_obj1 = 1426 713 356 280
confocal_y_amp_obj2 = 1426 713 356 280
confocal_y_amp_obj3 = 1426 713 356 280
confocal_y_amp_obj4 = 1426 713 356 280
laser_pwr_at_line = 125 30 100 0
laser_atten_at_458line = 1 1 1 1 1
laser_atten_at_488line = 1 1 1 1 1
laser_atten_at_515line = 1 1 1 1 1
laser_atten_at_Allline = 1 1 1 1 1
laser_inten_dark_field_offset = -8
pmtamp_min_offset_n_range = 0 90

```

```
white_light_dark_field_power = 200
laser_white_sw_cube = 0
```

Confocal xy amps of all objectives can be calibrated from  
within UV:

Click on the calibration item from the maintenance menu.  
Confocal z resolution must be calculated manually and entered  
into .uvconfig with *jot*.

#### B.5. Robot spec:

```
line 0: robot_spec = 2
```

```
line 1: port_number = a
```

a = -1, 0, 1, 2, 3   -1 = no robot, 0, 1, 2, 3 = SGI port 1, 2, 3, 4

```
line 2: 456inch_wafer_flat_size_deg = xx yy zz
```

The minimum wafer flat angle size, in degrees for the 4, 5 and  
6 inch wafers. If one uses single flat wafer exclusively,  
xx,yy,zz can all be set to 0.5 degrees. If, however, two flat (  
major and minor) wafers are expected, then xx yy and zz  
should be set to angle that is half way between the major and  
minor flats.

Example:

```
robot_spec = 3
port_number = 2
456inch_wafer_flat_size_deg = 29.000 34.500 34.500
```

#### B.6. Stage\_spec



line 0: stage\_spec = 24  
line 1: port\_number = n  
          n = -1           no stage.  
          n = 0,1,2, 3    single rs232 port on SGI port 1,2,3, or 4.  
5           n = 4           use 3 rs232 to communicate to the stage,  
                          SGI port 1/2/4 connected to X/Y/Z axis.  
line 2: x\_online = a  
line 3: y\_online = b  
line 4: z\_online = c  
10           a,b,c = 0,1; 0 = offline, 1 = online  
line 5: x\_resolution = xx  
line 6: y\_resolution = yy  
line 7: z\_resolution = zz  
line 8: x\_encoder\_res = xe  
15   line 9: y\_encoder\_res = ye  
line 10: z\_encoder\_res = ze  
          xx,yy,zz are micron per motor step, xe,ye,ze are micron per  
          encoder step.  
line11: max\_z\_um = aa bb cc dd ee  
20           aa bb cc dd ee are the maximum stage autofocus z travel  
          limits, in micron, for each objective.  
line12: stage\_af\_zero = zz  
25           xx is the dark level set during coarse autofocus  
line13: stage\_af\_inten = aa bb cc dd  
          intensity values to autofocus for 458/488/515/All laser lines  
line14: xy\_200mm\_level\_zoffset = xx yy  
30           xx yy are the Z position difference for a travel of 200mm, in  
          micron.

xx = focus position( z) at Right of stage xy position - Z at Left  
 yy = focus position( z) at Top of stage xy position - Z at  
 bottom.

line15: xy\_center\_offset\_um = xx yy

5

If the center of wafer is not exactly at stage coordinate 0,0,  
 user can change this parameter to make it center.

line 16: joystick\_gain = aa bb cc dd ee

XYJoystick gain settings for all objectives are specified here.

10

Range = 1 - 5

line 17: joystick\_z\_gain = aa bb cc dd ee

Joystick gain of Z axis for all objectives. Range = 1 - 5

line 18: stage\_angle\_tangent = xx

15

Wafers are placed by the robot and the position of the wafer  
 are programmed into the robot. The robot position is coarse but  
 is repeatable. Placement skew can be fine tuned with this  
 parameter. The wafer rotation, in tangent is calculated: For  
 the left and right points on a horizontal line( Pt1, pt2),  $xx =$   
 $(Y2 - Y1) / (X2 - X1)$

20

line 19: stg\_orthogonality\_deg = xx

The angle between the stage's X and Y axes is machined into  
 the stage and can not be changed. Adjustment of this angle,  
 however, can be made in software to compensate for  
 differences in tolerance. The degrees of XY orthogonality,  
 expressed in degrees from 90 is specified here.

25

line 20 fine\_af\_offset = n

Diagnostics purposes, should always set to 0.

line 21: fine\_af\_algorithm = n

Diagnostics purposes, should always set to 0.

30

line22: stg\_xy\_init\_final\_spd\_acc\_um = xx yy zz

xx is the initial xy speed, yy is the final speed, and zz is the acceleration of the xy axis for any programmed moves. all in microns.

line 23: stg\_af\_lo\_mag\_algorithm = x

5           x = 0: use first low-mag stage autofocus algorithm, x = 1:  
          use second lo-mag stage autofocus. Second lo-mag stage  
          autofocus requires fastZ firmware version 17(hex)

line 23: stg\_af\_lo\_mag\_algorithm = n

n can be either 0,1, or 2.

10          n = 0 old stage autofocus algorithm, to be used only for  
          systems with fast Z board version 1.6 or older.

n = 1 new stage autofocus algorithm, this is a more reliable  
algorithm, but requires fast Z board version 1.7 or newer.

15          n = 2: no stage autofocus, instead, move to Z limit, less 150  
          microns.

#### Example:

```

20           stage_spec = 24
          port_number = 3
          x_online = 1
          y_online = 1
          z_online = 1
          x_resolution = 0.2500
25           y_resolution = 0.2500
          z_resolution = 0.1016
          x_encoder_res = 0.5000
          y_encoder_res = 0.5000
          z_encoder_res = 0.1016
30           max_z_um = 6950 6900 6950 6990 6500
          stage_af_zero = 110
          stage_af_parms = 58 40 40 40 40
          xy_200mm_level_zoffset = 10 20
          xy_center_offset_um = 100 -312
35           joystick_gain = 4 4 2 1 4
          joystick_z_gain = 1 1 1 1 1

```

stage\_angle\_tangent = 0.0000000  
stg\_orthogonality\_deg = 0.0000000  
fine\_af\_offset = 0  
fine\_af\_algorithm = 0  
5 stg\_xy\_init\_final\_spd\_acc\_um = 5000 150000 250000  
stg\_af\_lo\_mag\_algorithm = 1

#### B.6.1. Procedure to determine stage\_af\_zero

10 B.6.1.1 Manually focus on a dark sample, with 100X  
objective at 458 line

B.6.1.2 Set laser intensity to just before PMT overload

B.6.1.3 Manually move sample away from focus.

15 B.6.1.4 Go to Ion diagnostic and select PMTAmp dialog box.  
Adjust AZOffsetCmd to a number that causes the laser image  
to be BLACK.

B.6.1.5 Back to laser window, manually focus sample again,  
adjust laser intensity to just before PMT overload. Redo step  
6.1.3 and 6.1.4

20 B.6.1.6 Use the AZOffsetCmd setting for stage\_af\_zero, and  
the laser intensity setting for first number in the  
stage\_af\_inten.

B.6.1.7 Exit uv, enter the stage\_af\_zero, and 458  
stage\_af\_inten into .uvconfig

#### B.6.2. Procedure to determine stage\_af\_inten

25 Note that the recipe also holds the stage\_af\_inten, and loading  
a recipe will overwrite current stage\_af\_inten the user login  
procedure includes a loading of the recipe with the user's  
name. The stage\_af\_zero is loaded only during program  
30 startup.

6.2.1 Start uv again, go to the LaserScanParm dialog box to  
select stage\_af\_inten for different laser line. Experiment with

different settings to try to autofocus on dark and bright samples.

5 6.2.2 bring up the recipe's LaserScanParm dialog box from within the recipe dialog box, and enter the autofocus intensities.

6.2.3 Save the recipe.

10 **B.6.3 Procedure to determine stg\_orthogonality\_deg**

Set stg\_orthogonality\_deg to 0.000000 before start.

6.3.1 Obtain a test wafer with die patterns that are known to be orthogonal.

15 6.3.2 Move stage to the lower left corner of the middle left die, record actual position, Leftx = \_\_\_\_\_, Lefty = \_\_\_\_\_ um

6.3.3 Move stage to the lower left corner of the middle right die, record actual position, Rightx = \_\_\_\_\_, Righty = \_\_\_\_\_ um

20 6.3.4 Move stage to the lower left corner of the middle top die, record actual position, Topx = \_\_\_\_\_, Topy = \_\_\_\_\_ um

25 6.3.5 Move stage to the lower left corner of the middle bottom die, record actual position, bottomx = \_\_\_\_\_, bottomy = \_\_\_\_\_ um

RightY - LeftY

6.3.6  $\tan(\text{Horizontal}) = \frac{\text{RightY} - \text{LeftY}}{\text{RightX} - \text{LeftX}}$ , Horizontal angle = \_\_\_\_\_ deg

RightX - LeftX

30 TopX - BottomX

6.3.7  $\tan(\text{Vertical}) = \frac{\text{TopX} - \text{BottomX}}{\text{RightX} - \text{LeftX}}$ , Vertical angle

= \_\_\_\_\_ deg

TopY - BottomY

6.3.8 Orthogonality = vertical + horizontal, in degrees. Use this number to set stg\_orthogonality\_deg in /usr/uv/.uvconfig.

5

#### B.6.4 Procedure to determine stage\_angle\_tangent

This parameter is used to fine tune robotic wafer placement.

6.4.1 Load a test patterned wafer, move to the middle left side of the wafer.

10 6.4.2 Focus on the lower left corner of a die with 100X objective, write down the stage coordinates: LeftX = \_\_\_\_\_, LeftY = \_\_\_\_\_ um.

6.4.3 Move to the lower left corner of a die in the right side of the wafer, staying on the same street. Stage coordinates at 100X, RightX = \_\_\_\_\_, RightY = \_\_\_\_\_.

15

RightY - LeftY

stage\_angle\_tangent = -----

RightX - LeftX

20

#### B.7. Objective Specs

line 0: objective\_spec = aa bb cc dd ee

aa,bb,cc,dd,ee are text labels for the name of the objectives

line 1: objective\_power = p1,p2,p3,p4,p5

25 nominal power of the objectives

line 2: objective\_na = n1,n2,n3,n4,n5

nominal NA of the objectives

line 3: objective\_offset\_nm = o1,o2,o3,o4,o5

30 moved after each objective change. The first objective is used as reference, so that o1 should be zero(0), Unused objectives should also be zero.

line 4: objective\_xoffset\_um = ox1 ox2 ox3 ox4 ox5

line 5: objective\_yoffset\_um = oy1 oy2 oy3 oy4 oy5

Stage XY moved after each objective change. The first objective is used as reference, so that o1 should be zero(0),

5 Unused objectives should also be zero.

line 6: obj\_Laser\_int\_458\_offset = a458 b458 c458 d458 e458

line 7: obj\_Laser\_int\_488\_offset = a488 b488 c488 d488 e488

line 8: obj\_Laser\_int\_515\_offset = a515 b515 c515 d515 e515

10 line 9: obj\_Laser\_int\_ALL\_offset = aall ball call dall eall

The amount of relative laser intensity adjustment for objective changes at 458/488/515/ALL laser line. a458 should be zero(0).

15 Procedure to setting obj\_Laser\_int\_xxx\_offset, assuming objectives are arranged as follows: 5X, none, 50X, 100X, none:

20 B7.9.1 Load a typical sample, and focus on a typical feature in laser-scan mode with 5X objective at 458 laser line.

B7.9.2 Adjust laser intensity slider to get optimum illumination. Intensity = \_\_\_\_\_

25 B7.9.3 Change to different objectives setting and adjust intensity slider to get the same level of intensity. Record intenisty settings:

4585 = \_\_\_\_\_, 45850 = \_\_\_\_\_, 458100 = \_\_\_\_\_

B7.9.4 Change to 488 line and repeat intensity adjustments for all objectives.

4885 = \_\_\_\_\_, 48850 = \_\_\_\_\_, 488100 = \_\_\_\_\_

30 B7.9.5 Change to 515 line and repeat intensity adjustments for all objectives.

5155 = \_\_\_\_\_, 51550 = \_\_\_\_\_, 515100 = \_\_\_\_\_

B7.9.6 Change to All line and repeat intensity adjustments for all objectives.

All5 = \_\_\_\_\_, All50 = \_\_\_\_\_, All100 = \_\_\_\_\_

5. B7.9.7 set:

obj\_laser\_int\_458\_offset = 0 0 (45850- 4585) (458100- 4585) 0

obj\_laser\_int\_488\_offset =  
(4885- 4585) 0 ( 4885- 4585) ( 488100- 4585) 0

10 obj\_laser\_int\_515\_offset =  
(5155-4585) 0 (51550- 4585) (515100- 4585) 0

obj\_laser\_int\_all\_offset = (All5- 4585) 0 ( All50- 4585) (All100- 4585) 0

15 line 10: objective\_camera\_int\_all\_offset = aall ball call dall

line 11: objective\_camera\_int\_pol\_offset = apol bpol cpol dpol

line 12: objective\_camera\_int\_yel\_offset = ayel byel cyel dyel

line 13: objective\_camera\_int\_dim\_offset = adim bdim cdim ddim

The amount of relative camera intensity adjustment for objective changes. aall should be zero( 0)

20

Procedure to setting objective\_camera\_int\_xxx\_offset:

25 B7.13.1 Load a typical sample, and focus on a typical feature in white-light-only mode with 5X objective.

B7.13.2 Adjust camera intensity slider to get optimum illumination. Intensity = \_\_\_\_\_

B7.13.3 Change to different objectives setting and adjust intensity slider to get the same level of intensity. Record intenisty settings:

30

All5 = \_\_\_\_\_, All50 = \_\_\_\_\_, All100 = \_\_\_\_\_



B7.13.4 Change to polarizer and repeat intensity adjustments  
for all objectives.

Pol5 = \_\_\_\_\_, Pol50 = \_\_\_\_\_, Pol100 = \_\_\_\_\_

B7.13.5 Change to Yellow and repeat intensity adjustments  
for all objectives.

Yel5 = \_\_\_\_\_, Yel50 = \_\_\_\_\_, Yel100 = \_\_\_\_\_

B7.13.6 Change to Dim and repeat intensity adjustments for  
all objectives.

Dim5 = \_\_\_\_\_, Dim50 = \_\_\_\_\_, Dim100 = \_\_\_\_\_

B7.13.7 set:

objective\_camera\_int\_all\_offset = 0 0 (All50 - All5) (All100-  
All5) 0

objective\_camera\_int\_pol\_offset =

(Pol5- All5) 0 ( Pol5- All5) ( Pol100- All5) 0

objective\_camera\_int\_yel\_offset =

(Yel5-All5) 0 (Yel50-All5) (Yel100- All5) 0

objective\_camera\_int\_dim\_offset =

(Dim5- All5) 0 ( Dim50- All5) (Dim100- All5) 0

line 14: cam\_mode\_norm\_use\_dim = aa bb cc dd ee

where aa,bb,cc,dd,ee are either 0 or 1, for each objective.

0 means using all-pass camera filter in normal camera mode,

1 means using dim camera filter in normal camera mode.

Example: :

```

25      objective_spec = 5      00      50      100      00
      objective_power = 20 20 50 100 20
      objective_na = 0.30 0.55 0.85  0.95  0.95
      objective_offset_nm = 0 0 48400 48400 0
30      objective_xoffset_um = 0  0  -200 234 0
      objective_yoffset_um = 0  0   120 99  0
      obj_Laser_int_458_offset =  0 0 5 20 0
      obj_Laser_int_488_offset =  0 0 5 20 0
      obj_Laser_int_515_offset =  0 0 5 20 0
      obj_Laser_int_all_offset =  0 0 5 20 0

```

```
objective_camera_int_all_offset = 0 0 20 30 0
objective_camera_int_pol_offset = 70 0 90 100 0
objective_camera_int_yel_offset = 20 0 40 50 0
objective_camera_int_dim_offset = 50 0 70 80 0
5 cam_mode_norm_use_dim = 1 0 0 0 0
```

### B.8. Lon Specs

miscellaneous configuration for the LON network.

10 line 0: lon\_spec

line 1: turret\_step = xx

xx = turret zstep delay, set to 50

line2: filter\_wheel\_delay\_laser\_line = 600

line3: filter\_wheel\_delay\_laser\_atten = 600

15 line4: filter\_wheel\_delay\_cam\_filter = 600

Not currently used, set to 600

line5: cam\_filter\_ap\_pos\_5\_or\_6 = 6

use filter wheel position 5 or 6 as the ALL-PASS position.

20 Using position 5 improves the camera filter wheel response  
time, but require new camera filter wheel firmware changes.

Example:

```
lon_spec
25 turret_step = 50
filter_wheel_delay_laser_line = 600
filter_wheel_delay_laser_atten = 600
filter_wheel_delay_cam_filter = 600
cam_filter_ap_pos_5_or_6 = 6
30
```

### B.9. System Specs

line 0: system spec

line 1: laser\_saver\_minutes = xx

xx = the minutes of no mouse/key click before putting laser to idle, white light intensity to 0 and PMT gain to minimum. Any mouse click/key will restore previous settings

5      line2:    kla\_deskew\_die\_offset = xx

The two deskew points, deskew1 and deskew2 are moved in toward the wafer center by XX number of die.

10                    line3: cfg\_next\_die\_position = n

n = 0: next die is the one below the current die.

n = 1: to the left of current die

15                    n = 2: to the top

n = 3: to the right

20                    line 3: enable\_bin\_image\_tx = n

n = 0: disable the saving of images in binary form.

n = 1: Enable the saving of images in binary form, in tiff format. Images in this form are in 8 bit gray scale image in top view.

25                    line 4: zrange\_search\_num\_steps = nn

nn = the number of steps in the **SetZ** procedure. default is 44 steps. The step size is fixed at 30 micro steps( approximately 0.5 microns per step.) Scan speed is about 13 steps per second.

30                    line 5: zrange\_search\_laser\_gain = x

where x is in the range of 7.0 to 13.0.

This is the fine tune adjustment mapping between laser intensity setting and the resultant image intensity. If this value is too low or too high, user may need to perform the SetZ function to get a reliable optimal laser intensity setting.

5

To check for zrange\_search\_laser\_gain setting:

10

Set laser intensity to some low value so that image appears light green( about 10 for typical sample) at focus. click on SetZ, record the laser intensity setting. Click on SetZ again for 2 more times. If laser intensity is increasing after the second and third SetZ, then zrange\_search\_laser\_gain is too low, increment by steps of 0.5

15

Example:

20

```
system spec
laser_saver_minutes = 5
kla_deskew_die_offset = 1
cfg_next_die_position = 3
enable_bin_image_tx = 0
zrange_search_num_steps = 44
zrange_search_laser_gain = 9.0
```

**Appendix C: IRIS INDIGO BASIC COMMANDS****Contents:**

	C1. Introduction
5	C2. Logging in
	C3. Help
	C4. File Manipulation
	C5. File Information
	C6. Directory manipulation
10	C7. Find and search
	C8. System information
	C9. User information
	C10. Command information
	C11. Process control
15	C12. Time/Date
	C13. File conversion
	C14. Piping commands
	C15. Sourcing and Sinking
	C16. Miscellaneous Commands
20	C17. System shutdown procedure
	C18. Text editing using Jot
	C19. DOS Floppy Disk
	C20. Tapes
	C21. Adding users setup
25	C22. Networking setup
	C23. IMPORT directory setup
	C24. REMOTE/REMOTE2 directory setup
	C25. File transfer
	C26. Logging on remotely
30	C27. Home Directory

### C.1. Introduction

The basic operating system, IRIX 4.0.5F, is already preloaded onto the Indigo. IRIX 4.0.5F is Silicon Graphic's version of Unix.

A list of basic command follows.

5 For information on such topics as Floppy Disk Setup,  
please refer to the Ultrapointe Model 1000 User's Manual,  
Appendix A.

There are also very useful information in the Online Manual Pages.

### 10 C.2. Login in

login "login\_name" login to the Unix system

passwd change one's passwd

### C.3. Help

15 man "command" invokes online manual pages for a particular  
command

man -k "keyword" search manual pages for a particular keyword

### C.4. File manipulation

20 cp "srcfile" "destfile"

copy the file named "srcfile" to a file named "destfile"

mv "srcfile" "destfile"

move or rename a file named "srcfile" to a file named  
"destfile"

25 rm "filename" remove or delete a file.

can be very destructive, no deleted file recovery in Unix.

chmod "flags" "filename"

change permissions of a file

r = read, w = write, x = executable

30 u = user, g = group, o = other, a = all

chmod a+rw "filename"

give read and write permissions to everyone  
chmod o-x "filename"  
take away execute permission from others  
chown "NewOwner" "filename"  
5 In "srcfile" "destfile"  
instead of copying files, one can link files, and save space  
and redundancy.  
"destfile" is becomes a pointer to "srcfile", so whenever  
one accesses "destfile", one is actually accessing  
10 "srcfile".

### C.5. File information

more "filename" same as MORE in DOS world, stop at each page  
cat "filename" same as TYPE in DOS world  
15 head "filename" show the first few lines in a file  
tail "filename" show the last few lines in a file  
diff "file1" "file2"  
show differences between files, line by line  
cmp "file1" "file2"  
20 compare two files, byte by byte  
wc "filename" show number of characters, words, and lines in a  
file

### C.6. Directory manipulation

25 ls list the contents of current directory  
"ls -al" to show dotfiles, permissions, ownership, dates,  
filesizes.  
pwd print working directory, (show current directory)  
cd "directory" same as CD in DOS world, change to another  
30 directory  
cd .. change to parent directory

cd /      change to the root directory  
cd ~      change to your home directory  
mkdir "directory"    same as MD in DOS world, make a new directory  
rmdir "directory"    same as RD in DOS world, remove directory

5

### C.7. Find and search

grep -i "keyword" \*.c  
                search for the "keyword" in all files with extension  
                ".c" while being case insensitive.  
10    find "directory\_to\_start\_search\_from" -name "filename" -print  
                search for "filename" starting at this "directory"

### C.8. System information

hostname    display the name of the this computer, eg. "irisdemo"  
15    df -k            show amount of disk space free, (in kilobytes)  
      du -k            show amount of disk space used, (in kilobytes)  
      env             show your environment settings, such as PATH, SHELL,  
                      TERM, HOME, prompts, etc.  
      hinv            show hardware configuration, such as CPU type,  
20                    memory size, serial ports, ethernet, disk drives,  
                      graphics board options, etc.  
      versions        show versions of installed software  
      uptime          show how long the system has been up, and the number  
                      of users logged on

25

### C.9. User information

whoami      show whose account this is  
who          shows who is currently logged on  
finger "user\_name"  
30            show real name of "user\_name", etc.



**C.10. Command information**

which "command" show exactly which executable will be invoked  
when "command" is typed.

whereis "command" show where the location of this "command" is,  
this includes its executable, source code, and manual  
page.

w show which commands are being executed, and by which  
users, etc

**10 C.11. Process control**

reboot Shutdown system and then Boot up again. No musical  
tune is played. To power down, choose the "System  
Shutdown" option on the "System" popup menu at the  
top, left-hand corner of screen.

15 Ctrl-C kill current process

ps show the processes running currently, along with their  
PID, (process identification number)  
"ps -e" to show every process on system.

kill "process\_id\_number"

20 kill the specified process identified by its PID,  
(process\_id\_number)

"command" & execute "command in the background

fg bring background process back to the foreground.

**25 C.12. Time/Date**

date show the current date and time

date 06091445 change date and time to June 9, 2:45PM

time "command" show the amount of time it takes to execute  
"command"

30 cal show current month's calendar

**C.13. File conversion**

file "filename" show what type file this is,

ascii text, C program code, binary data, executable,  
image, archive, Postscript, packed, compressed, etc

5 to\_dos "unixfile" "dosfile"

convert Unix file to DOS file, these are ascii text files

convert <CR> <LF> to <LF>, append Ctrl-Z  
(EOF\_Marker)

to\_unix "dosfile" "unixfile"

10 convert DOS file to Unix file, these are ascii text files

convert <LF> to <CR> <LF>, strip off Ctrl-Z  
(EOF\_Marker)

**C.14. Piping commands**

15 "command1" | "command2"

use the output of "command1" as input to "command2"

diff "file1" "file2" | more

show the differences between "file1" and "file2", page by  
page.

20

**C.15. Sourcing and Sinking**

"command" > "filename"

put the output from "command" into "filename"

"command" < "filename"

25 use contents of "filename" as input to "command"

**C.16. Miscellaneous Commands**

su become the super-user

Ctrl-D quit being the super-user,

30 can also be used to destroy a Terminal(BlueScreen).

RightMouseButton press the right mouse button when in a

BlueScreen to bring up a popup menu.

Options in this Popup Menu will let the user change the font, or create another Terminal(BlueScreen).

5           Windows   All windows have the same basic functionalities, which  
                          are Move, Size, Minimize, Maximize, Raise, Lower,  
                          Quit. These functionalities can be activated either  
                          with the WindowManager's Popup Menu in the  
                          TopLeftHand corner or graphically with the mouse.

10           On the Ultrapointe Model 1000, we can also use the following  
                  commands:

          ll           list the contents of the current directory with detailed info  
          hh           show the previous few commands  
          Ctrl-P       retype the previous command automatically

15           **C.17. System shutdown procedure**

          0. Do not just turn off the power switch on the SGI Indigo computer.

          1. Select the "System" menu in the top left hand corner.

          2. Select the "System Shutdown" option.

          Do not choose the "Log Out" option.

20           3. Wait for screen message saying "OK to power off the system"  
                  appears, then turn off machine.

**C.18. Text editing using Jot**

25           Jot is an X Window application, easy to use, can cut and paste,  
                  1 level of undo, search for a string, pull in a file.

          jot "filename"   start the jot text editor.

                  Use the left mouse button for selecting regions of text.

                  Use the right mouse button for bringing up the popup  
                  menu.

30

**C.19. DOS Floppy Disk**

cp "filename" /FLOPPY     copy a file from current directory to the  
                             Floppy Drive  
cp /FLOPPY/"filename" .   copy a file from the Floppy Drive to current  
5                               directory  
eject                       eject the Floppy Diskette from the Drive

**C.20. Tapes**

Tar is most universal type of tape drive format.

10     One usually copies files from the specified directories, onto a tape  
(archiving), or extracts files from the tape, back into their original  
directories.

tar tv                    list the contents of this tape

tar cv "Filename1" "Filename2"

15                        will format the tape to the tar format, then write the  
specified files onto the tape.

tar cv "Directory1" "Directory2"

20                        will format the tape to the tar format, then write all files  
from the specified directories and their subdirectories  
onto the tape.

tar cv .                will write all the files from the current directory and its  
subdirectories onto the tape.

tar xv "Filename"     extracts the file "Filename" from the tape back  
into its original directory.

25     tar xv "Directory1" "Directory2"

extract all the files from the specified directories and  
their subdirectories from the tape.

tar xv                will extract all the files from the tape, back into their  
original directories, and subdirectories.

30

**C.21. Adding users setup**

1. Select the "System" menu in the top left hand corner.
2. Select the "System Manager" option.
3. Select the "Users" icon.
- 5 4. Select the "Add" pushbutton.
5. Enter login name
6. Fill in the detailed information such as password, full name,  
desired Unix shell. For Korn shell, select other, and enter /bin/ksh
7. Select Accept.

10

**C.22. Networking setup**

1. Select the "System" menu in the top left hand corner.
2. Select the "System Manager" option.
3. Select the "Networking" icon.
- 15 4. Select Add.
5. Type in the HostName of the workstation to be added.
6. Type in its IP address.
7. Turn on TCP/IP, turn off NFS and NIS.
8. Select Accept.

20

**C.23 IMPORT directory setup**

The IMPORT directory is a directory located physically in LIS but can be accessed by another machine as if it is its own directory. It is intended for another machine to IMPORT files into the LIS.

25 Configuration:

**C.23-1 Network setup on the remote machine**

Obtain from network administrator the remote machine's name and its IP address, the LIS' name and IP address. Setup the remote machine so that it can access LIS' IMPORT directory. The LIS' IMPORT directory is

30 /usr/uv/IMPORT.

**C.23-2 Add the remote machine to LIS Network setup**

see C22 Network setup for adding the remote machine

**C.23-3, Make the IMPORT directory**

make the /usr/uv/IMPORT directory with read/write access by:

```
mkdir /usr/uv/IMPORT  
chmod +rw /usr/uv/IMPORT
```

**C.23-4 Make the IMPORT directory accessible by the remote machine**

Add the IMPORT directory to the export table, so it is accessible to other machine. Edit the /etc/exports file:

```
jot /etc/exports  
add the line:  
/usr/uv/IMPORT
```

save the file

quit

**C.23-5, Tell the system start using this new table:**

```
exportfs -a
```

**C.24 REMOTE/REMOTE2 directory setup**

The REMOTE and REMOTE2 are two directories that LIS users can access in all file access dialog box. These two directories are intended to be mapped into directories that are physically located in different machines via TCP/IP and NFS. These directories are hidden until they are configured, Configuration:

**C.24-1 Network setup on remote machine**

Obtain from system administrator the remote machine's name and its IP address, i.e: host1 and 192.0.2.10, its directory LIS is expected to access. As an example, it is in host1, directory /usr/xxx/yyy/zzz. Make sure that this directory is exported to LIS with both read write access.

**C.24-2 Add the remote machine to LIS Network setup**

see C22 Network setup for adding the remote machine

**C.24-3 Make the LIS REMOTE directory**

make the /REMOTE directory with read/write access by:

```
mkdir /REMOTE
chmod +rw /REMOTE
```

**5 C.24-4, Enable network connection**

Edit the file /etc/fstab:

```
jot /etc/fstab
add the line:
```

```
Host : /usr/xxx/yyy/zzz /REMOTE nfs
10 rw,soft,intr,bg,0 0
```

save the file

quit

mount and use the directory:

```
mount /REMOTE
```

15

**C.25. File transfer**

There are 2 methods for transferring files to a remote system: rcp and ftp.

1. rcp hostname1:/usr/people/user1/filename1

20 hostname2:/usr/people/user2/filename2

note: if one is copying to the remote system, one must have permission to write into the remote directory.

2. ftp "hostname"

25 help type "help" for a list of available commands

dir list directory contents, can also use "ls"

get "filename" receive a file

put "filename" send a file

mget "\*.c" receive all files with the extension .c

mput "\*.h" send all files with the extension .h

30 cd "directory" change to another directory on remote system

lcd "directory" change to another directory on local system

pwd show current directory on remote system

bye disconnect from remote system

**C.26. Logging on remotely**

1. rlogin "hostname"
2. after being prompted for a login name and a password,  
it would seem as though one is logged on locally.

5

**C.27. Home Directory**

There are some dotfiles which need to be configured to one's taste.

The home directory is the directory one is placed in upon login in.

It is usually /usr/people/your\_login\_name

10 Dotfile are configuration files which are executed when one is  
login in.

Dotfiles are similar to autoexec.bat and config.sys in the  
DOS world.

15 Shells are powerful command interpreters. They are similar to  
command.com in the DOS world, but much more powerful.  
Can modify one's prompt, set aliases for long commands, recall  
previous commands, edit the command line, etc.

20 .login for C shell, configuration file when logging in  
.cshrc for C shell, configuration file when new C shell is started  
.profile for Korn shell, configuration file when logging in  
.Xdefaults configuration file for X Window applications  
.rhosts permission file for remote logging in or remote file  
25 transfer

For example, if one wishes to let others perform a remote  
login (rlogin) or perform a remote file transfer (rcp),  
append the following lines to the file .rhosts

30 hostname1 username1  
hostname1 username2  
hostname2 username3



**Appendix D: Support for KLA 20xx Defect File**

The Ultrapointe Model 1000 supports the KLA 20xx series defect scanners with its ability to read and write 20xx formatted wafer files. Due to the 20xx series unique features, the model 1000 implemented the following user interface to automate defect review using 20xx series wafer maps:

D1 Wafer map contains no wafer size information, wafer size is set to the cassette size currently used.

D2. Toggling between test die and reference die is provided through a button in the textstatus window. This button moves stage alternately between defect locations calculated with test die and reference die.

D3. Actual defect location, in test die or reference die, is updated when user enters a defect code for the defect. When a defect code is entered in the textstatus window, the defect list window, and the wafer map window will display the changes.

**Appendix E: Customer Preventative Maintenance Procedures****E.1. Scope:**

These Preventative Maintenance (PM) procedures describe the maintenance operations on the Ultrapointe Laser Imaging System (LIS) that can be routinely performed by on-site customer personnel (referred to as customer service personnel). The operations are low-risk in the sense that they do not involve delicate adjustments or potential contamination of critical optical components. Furthermore, there is no requirement to bypass or defeat any of the safety interlocks of the machine. At all times, however, the service person is required to observe common safety practices and to avoid altering any adjustments inside the machine that do not pertain to the procedure at hand. These operations are intended to be supplemented by Ultrapointe factory PM operations performed on a less frequent basis.

**E.2. Safety Precautions:**

The Ultrapointe LIS Model 1000 has been designed to minimize the hazards presented to the user during normal operation of the instrument. The customer service personnel will not be exposed to hazards if the following procedures are followed carefully, but customer service personnel should exercise due caution and common safety practices when working on the machine. Refer to figure 1 for identification of the major components of the instrument.

**E.2.1 Laser Safety:**

In accordance with the Regulations for the Administration and Enforcement of the Radiation Control for Health and Safety

Act of 1968 (applicable to laser products), a Laser Product Report for this instrument has been forwarded to the Center for Devices and Radiological Health (CDRH). This instrument falls under Class I definitions because there is no user access to laser radiation.

5

**CAUTION: Use of controls or adjustments, or performance of procedures other than those specified herein may result in hazardous radiation exposure.**

10

The protective housing of the instrument should be kept intact during normal operation. Internal to the instrument is a Class IIIb argon-ion laser rated at 25 milliwatts with a maximum possible output of 500 milliwatts. Removal of any part of the upper protective housing can expose the user unnecessarily to hazardous laser radiation. Danger labels reproduced below indicate covers that provide access to hazardous laser light exposure.

15

5

**DANGER - Laser radiation when  
open and interlock defeated.  
AVOID DIRECT EXPOSURE TO BEAM**

**DANGER - Laser radiation when  
open. AVOID DIRECT EXPOSURE  
TO BEAM**

-363/85-

RECTIFIED SHEET (RULE 91)

Customer service personnel will be required to remove the top protective cover in order to perform some of the procedures described here. The cover is interlocked with the laser power supply so that laser radiation is not present after the cover is removed. **The top cover interlocks are defeatable but should not, under any circumstances, be defeated by customer service personnel.**

#### **E.2.2 High Voltage:**

Internal to the instrument are components using voltages up to 1000 VDC and 240VAC. Where the procedures in this document permit the customer service personnel to have access to components that may be at hazardous voltages, there will either be warnings or instructions to shutdown, turn off, and lock out the AC input to the machine. In this manner, the customer service personnel will not be exposed to hazardous voltages within the instrument.

#### **E.2.3 Moving Parts:**

As part of routine instrument operations, various parts of the machine make movements automatically and when commanded by the user. The user should take care to keep body parts and clothing away from the robot area. The user should also ensure that the current robot operation is complete before adding or removing a cassette from the cassette locator plate.

Much of the instrument is attached to an internal frame that is mounted to the external frame and protective covers via pneumatic vibration isolators. The "floating" internal frame can therefore move freely with respect to the covers and creates possible pinch points between floating components

and the covers or external frame. Externally accessible floating components include the cassette platform, the robot, and the flat finder.

5 Inside the instrument, are many moving parts which could create possible hazards, including motors, shutters, and an X/Y/Z stage. Where the procedures in this document permit the customer service personnel to have access to components that might involve mechanical hazards or pinch points, there will be instructions to shutdown, turn off, and lock out the AC  
10 input to the machine. In this manner, the customer service personnel will not be exposed to mechanical hazards or pinch points.

Moving parts on the outside of the system include the robot handler, pre-aligner, and wafer shutter door, as well as the  
15 floating optics platform. The user should take care to keep body parts and clothing clear of all moving parts while the instrument is under power.

### E.3. PM Schedule:

20 Figure E2 represents the recommended PM intervals for the Ultrapointe Model 1000 and the responsibility for each procedure.

#### E.3.1 System Startup

25 The system is started from a power off condition by pressing the green power-on button at the front of the machine. The large red emergency off (EMO) mushroom button must be released before power will be applied; release the emergency off by rotating it counter-clockwise.

30 The computer operating system will go through its normal initialization and boot process, and finally the Ultrapointe

application will appear on the screen. A series of initialization windows will appear; follow the instructions in the windows as they appear. Two components in the system have to be explicitly initialized: the XYZ stage and the robot. The last window to pop up as part of the initialization process is a Manual Stage Control window; click on the 'InitStage' button to initialize the stage. The robot must be initialized via the cassette window: click on 'Initialize' to initialize the robot. Neither the stage nor the robot will be operable until they are initialized.

### E.3.2 System Shutdown

The system may be powered off at any time using the red 'off' button or the large EMO switch on the front panel. However, the UNIX operating system used internally on the computer platform will incur small amounts of 'damage' to its file system if the system is turned off without shutting down the computer in an orderly manner. This damage is routinely repaired as the computer system boots up at power on. If the system is repeatedly turned off without an orderly shutdown, the accumulated damage may be enough so that the automatic repair process fails. Before turning off power to the system, a system shutdown process is therefore recommended.

To perform a system shutdown, lower the Ultrapointe application: click and hold on the small square at the extreme top left of the screen. A pull-down menu appears; select 'Lower', to lower the application. After the Ultrapointe application is lowered, a set of buttons appears at the top left of the screen:

Procedure Name	Reference Paragraph	Description	Interval	Responsibility
System Startup	3.1	How to start the Ultrapointe application program	N/A	Customer
System Shutdown	3.2	Shutdown procedure for the LIS in preparation for power off	N/A	Customer
Lockup Recovery	3.3	Procedures for recovery from software lockups	N/A	Customer
Wafer Removal	3.4	Wafer removal procedures in the event of a lockup or failure	N/A	Customer
Wipe Down	3.5	Wipe down of the exterior of the machine	Daily	Customer
Backup Data	3.6	Procedure to backup critical data to tape	Weekly	Customer
Chuck Cleaning	3.7	Cleaning procedures for the chuck and robot end effector	Weekly	Customer
Laser Life and Inter-lock Check	3.8	Test procedure for the laser system to determine whether to schedule laser service, and to verify interlock integrity	Monthly	Customer

Ultrapointe LIS 1000 User's Manual, Rev. 1.3

Revised 3 June, 1994

-363/89-

RECTIFIED SHEET (RULE 91)



5

Robot Alignment Check	3.9	Test to verify that the robot handles wafers and interfaces appropriately to the cassette, pre-aligner and XYZ stage	Quarterly	Customer
Bulb Change	3.10	Install new halogen bulb for white light illumination	6 months	Customer
Bi-annual PM	3.11	Optics cleaning and alignment, system PM	6 months	Ultrapointe

**Figure E2. Recommended PM Intervals.**

10

Click on 'System' and another pulldown menu appears. Click on 'System Shutdown'. A window will pop up asking if you really want to shut the system down; click on 'Yes'. The computer system will now begin its shutdown process. After a minute or so, a message appears indicating that shutdown is complete and that it is OK to turn off the machine. At this point, turn off the machine by pressing the red 'Off' button or the EMO on the front panel.

### E.3.3 Lockup Recovery

It is possible for the system to lock up to varying degrees. There have been reports of cases where the system slows down noticeably, of cases where none of the on-screen buttons in the Ultrapointe application respond to the cursor, and of cases where the mouse and keyboard don't seem to work at all. There is a special case of lockup that can occur during a volume acquisition. These incidents should be rare occurrences and Ultrapointe is very interested in knowing about and fixing lockups that have any frequency. The possible recovery actions include: re-starting the Ultrapointe application program, re-booting the computer system, re-starting the XWindows server (five finger reset), and cycling power of the entire instrument. Whatever the recovery method, the problem and the corrective action should be recorded for maintenance purposes (see sections 4.0 and 5.0). The following table lists lockup symptoms and recovery actions recommended. For the sake of simplicity, the operator can attempt to recover from all lockups using the reboot procedure - this takes slightly longer than just re-starting the Ultrapointe application.

Lockup Symptom	Recovery Action	Paragraph #
System slows down	Restart the application	3.3.1
Application doesn't respond	Restart the application	3.3.1
Mouse doesn't respond	Five finger reset	3.3.2
Lockup during volume acquire	Reboot the computer	3.3.3
Cannot lower the application	Five finger reset	3.3.2
Nothing works	Cycle power	3.3.4

#### **E.3.3.1 Restarting the Ultrapointe Application**

To restart the Ultrapointe application, lower the Ultrapointe application: click and hold the square at the extreme top left of the screen. A pull-down menu appears; select 'Lower' to lower the application. After the Ultrapointe application is lowered, a set of buttons appears at the top left of the screen; click on 'Tools' and another pulldown menu appears. Click on 'Shell' and after a few seconds, the outline of a UNIX shell window will appear. Position the window appropriately and click the mouse to create the window. Position the cursor anywhere in the UNIX shell window and type 'uv1000' to restart the Ultrapointe application.

If the application cannot be lowered, attempt the five finger reset as described in section 3.3.2.

### **E.3.3.2      Executing a Five Finger Reset**

If the application cannot be lowered because the box at the extreme top left of the screen does not respond to the mouse, the operator can try to kill and restart the XWindows server using the keyboard. This involves pressing five keys simultaneously, and is referred to as the 'five finger reset'. To execute a five finger reset, press the 'Ctrl', 'Shift', 'Esc', 'F12', and '/' keys simultaneously, where the '/' key is on the numeric keypad. This should cause a dramatic change in the display, followed by a recovery to a logon screen known as the 'Pandora' screen. The Pandora screen has several icons with user names under them. Double click on the icon labeled 'uv' and the system will restart normally.

If the lockup occurred during a volume acquisition, the operator must still reboot the system even though the application has restarted. The reboot procedure is described in section 3.3.3.

If the five finger reset does not work, cycle power as described in section 3.3.4.

### **E.3.3.3      Rebooting the Computer**

To reboot the workstation, lower the Ultrapointe application: click and hold the square at the extreme top left of the screen. A pull-down menu appears; select 'Lower' to lower the application. After the Ultrapointe application is lowered, a set of buttons appears at the top left of the screen; click on 'Tools' and another pulldown menu appears. Click on 'Shell' and after a few seconds, the outline of a UNIX shell window will appear. Position the window appropriately and click the mouse to create the window. Position the cursor anywhere in the window and type 'reboot' to reboot the computer.

If the application cannot be lowered, attempt the five finger reset as described in section 3.3.2, and then execute the reboot procedure.

5      **E.3.3.4      Cycling Power to the Instrument**

10      If none of the above measures is successful in returning the application to normal, the only alternative is to cycle power to the instrument. Needless to say, this will be done without a proper shutdown process. To cycle power, shut off power using the red power 'off' button on the front panel, wait 30 seconds, then power on by pressing the green 'on' button. The system should come up as described in section 3.1. Be sure that the EMO switch is released by turning counter-clockwise, or the 'on' switch will have no effect.

15

**E.3.4      Wafer Removal**

20      There are two conditions under which a wafer can be left in the instrument. First, the system may lose AC power or the application software may lockup while a wafer is inside the optics chamber or in the pre-aligner; the operator will therefore have to remove the wafer (using the application software) immediately after re-starting the Ultrapointe application (see section 3.4.1). Second, the instrument may become non-functional (due to a component failure) while a wafer is in the optics chamber or in the pre-aligner, and the operator may want to remove the wafer manually because the application software cannot be re-started, even by cycling power to the instrument (see section 3.4.2).

25

#### E.3.4.1 Wafer Removal Using the Application Software

Wafers can be removed from the optics chamber or the pre-aligner using the Manual Robot Control window which appears as a menu selection under the 'Maintenance' button on the menu bar across the top of the screen. This menu selection is available only to users with service or administrator permissions. However, if the Ultrapointe application is started with a wafer on the XYZ stage (in the optics chamber) or on the pre-aligner, then the Manual Robot Control menu will appear automatically when the operator attempts to initialize the robot.

With the Manual Robot Control menu on the screen, the operator clicks on the a wafer source or destination in the list at the left of the window, then clicks on 'Robot Get' or 'Robot Put' to command the robot to either retrieve or deposit a wafer at that location. **This window should be used with great care, because the usual safeguards in the applications software are not active in this window.** For example, the Manual Robot Control menu will allow you to attempt to put two wafers in the same slot, or replace a wafer into a cassette of the wrong size.

To remove a wafer from the XYZ stage, click on the 'Stage' selection in the source/destination list at the left of the window, then click on 'Robot Get'. The robot will enter the optics chamber and retrieve the wafer. To put the wafer back into the cassette, first inspect the cassette to find a slot that is empty, then click on the appropriate slot number in the source/destination list at the left of the window. Click on 'Robot Put' to deposit the wafer in that slot. **Be sure that the cassette size matches the wafer, that the slot in question is empty, and that the correct slot has been highlighted in the**

list at the left of the window before clicking on 'Robot Put'.

The procedure for removing a wafer from the pre-aligner is similar, except that a 'Robot Get' operation is performed with the 'Flat Finder' selection highlighted in the source/destination list.

5

#### **E.3.4.2 Manual Wafer Removal**

If the Ultrapointe application software is inoperable, the user may remove a wafer from any of the major components in the system manually. The operator should be aware that he is interacting with mechanical components that move under computer control, so it is highly recommended that the system be shut down using the procedure of section 3.2 before attempting manual wafer removal.

10

If a wafer is on the robot, it can be easily removed with a vacuum wand. When the system is powered off, the operator can manually move the robot arm to a more appropriate position to remove the wafer. The default state of the vacuum hold-downs in the system is 'on', so it is highly recommended that the operator disconnect the vacuum supply at the rear of the instrument before attempting to remove a wafer from the end effector or the chucks. The same procedure applies to the pre-aligner module - if a wafer is left on its chuck, it can be removed manually with a vacuum wand.

15

20

25

If a wafer is to be removed from the XYZ stage in the optics chamber and the system is inoperable, the first step is to remove the facia plate that surrounds the wafer door. The facia plate is fastened with two screws near its base, remove these two screws and gently tilt out the bottom of the plate to remove it. Behind the facia plate is the wafer door; this

30

can be raised manually. With the wafer door raised, reach into the system with a vacuum wand and remove the wafer from the XYZ chuck.

#### 5 E.3.5 System Wipe Down

System wipe down consists of cleaning the exterior of the machine with a cleanroom-compatible cleaning cloth and a dilute solution of isopropyl alcohol or plain water. The user should check with the Ultrapointe factory if other solvents or  
10 cleaning solutions are to be used.

#### E.3.6 System Backup and Restore

To avoid loss of data due to a hard disk failure or other disasters, the user should periodically backup critical data.  
15 The following procedures describe backup and restore operations.

##### **E.3.6.1 System Backup**

To use the backup facilities, place a 150 megabyte quarter  
20 inch tape in the tape drive within the drive bay and close the door of the tape drive. Lower the Ultrapointe application program: click and hold on the small square at the extreme top left of the screen. A pull-down menu appears; click on 'Lower', to lower the application. A set of buttons appears at  
25 the top left corner of the screen, click on 'System', then 'System Manager'. A window appears with several icons and names above them. Double click on 'Backup&Restore'. Another window appears for setting up the backup operation; be sure that 'QIC150#7' is selected as the tape drive, and  
30 'Backup' as the operation. Click on 'Accept' to accept the setup. Yet another window appears for specifying the files to



5 backup. Select 'by Directory' as the files to backup, type in  
'/usr/uv' as the directories and files to be backed up, then  
click on 'Start Backup'. **Note: this is a typical backup  
operation, more or less data may be backed up depending on  
the user's preference - contact the Ultrapointe factory for  
recommendations.** The backup now begins, with the 'Backup  
is' indicating 'Active'. The backup proceeds with the list of  
files being backed up scrolling through the window until the  
10 backup is complete. At this point the 'Backup is' status  
changes to 'Inactive'.

#### E.3.6.2 System Restore

15 The restore procedure is very similar to that of section 3.6.1.  
Instead of selecting 'Backup' as the operation in the first  
window, select 'Restore'. When the setup is accepted, a  
window will pop up asking the operator to wait while the tape  
is being read. A window then appears to specify what files to  
restore. Be sure to specify that files should be restored to the  
same directory they came from, then click on 'Start Restore'  
20 to begin the restoration process. Just as before a status  
message indicates whether the restore process is 'Active' or  
'Inactive' (complete).

#### E.3.7 Chuck Cleaning

25 The backside of wafers come into contact with three pieces  
of equipment in the LIS: the robot end effector, the XYZ stage  
chuck, and the pre-aligner chuck. These three surfaces must  
be cleaned periodically. Cleaning the robot end effector is  
straightforward using a cleanroom compatible clean wipe and  
30 a cleaning solution of either isopropyl alcohol or water. A  
special tool is supplied by Ultrapointe for cleaning the two

chucks. It should be fitted with a fresh cleaning sleeve and used with a alcohol or water to wipe off the surface of both chucks. Do not use acetone with the chuck cleaning tool.

Cleaning the pre-aligner chuck is also straightforward. If the robot is initialized, it will move down as far as possible so that it does not interfere with the cleaning process. The chuck is visible through the slot through which the wafer enters the pre-aligner.

To clean the XYZ stage chuck, first initialize the robot to get it out of the way using the 'Initialize' button in the cassette window, select a turret position containing no objective ('00' selection from 'Objective' button at top right of LaserScan window), then open the wafer door using the LON Diagnostics window. The LON Diagnostics window is accessed by clicking on 'Maintenance' on the top menu bar, and then on 'LONDiagnostics'. A button called 'WaferDoor' appears roughly in the middle of the window; click on it to open or close the wafer door. With the wafer door open, the operator may reach in with the chuck cleaning tool to clean the surface of the chuck.

The operator should be aware that he is interacting with mechanical components that move under computer control, so it is highly recommended that the system be shut down using the procedure of section 3.2 before attempting chuck cleaning.

Close the wafer door using the LON Diagnostics window when the cleaning operation is complete.

#### E.3.8 Laser Lifetime and Interlock Inspection

The normal aging process of the laser tube can be monitored by measuring the current required to obtain a given power

output. The procedure involves setting the laser power to maximum (25 milliwatts) and measuring the required tube current. With a new laser this is typically 6 Amps. The laser power supply is capable of supplying a maximum of 10.5 amperes. When the laser tube requires 9.5 amps to achieve 25 milliwatts of power it is close to the end of its life (est. 1000 hours remaining) and a replacement should be scheduled.

10 The procedure is as follows:

- 1) Click on the 'Maintenance' selection on the menu bar at the top of the screen. Select 'LONDiagnostics' from the resultant pull-down menu.
- 15 2) On the right hand side of the LON Diagnostics screen are two columnar lists: the left column (labeled 'Nodes') lists the I/O boards in the instrument, the right column (labeled NetVariables) lists the available network variables for the I/O board selected in the left column. Select the 'LaserPwrControl' in the left column.
- 20 3) Select the 'PowerCmd' network variable in the right column.
- 4) Click on the text entry box below the left column (labeled OutputNVValue) and enter 128, which corresponds to 25 milliwatts of laser power. Click on the 'SendToLON'
- 25 button to change the current laser power setting to 25mW.
- 5) Click on the 'CurrMonitor' network variable in the right column.
- 30 6) Check and record the value returned as 'NetVarValue:' below the left column.

7) Divide the returned value recorded in step 6 by the conversion factor of 2.3052. This gives the current in amperes required to deliver the 25 milliwatts of laser power. Record this value.

5 8) If the current required to deliver 25 milliwatts of laser power is greater than 9.5 Amps, the laser is nearing the end of its life. A factory service call should be scheduled to replace the laser.

10 The wafer door interlock should be tested at the same time as the laser lifetime is inspected. To test the interlocks, load a wafer of any kind into the instrument using routine operating procedures. Focus on the surface of the wafer and obtain a laser image on the workstation screen. The user then opens  
15 the wafer door using the LON Diagnostics menu and verifies that the laser image goes away because the wafer door interlock has shuttered the laser. The LON Diagnostics menu can only be accessed by users with service or administrator permissions. It is reached by clicking on 'Maintenance' on the  
20 menu bar at the top of the screen, then clicking on 'LONDiagnostics' on the resultant pull-down menu. In the LON Diagnostics menu is a button labeled 'WaferDoor'; click on this button to open or close the wafer door. Verify that the laser image disappears when the door is open and  
25 re-appears when the door closes.

#### E.3.9 Robot Alignment Check

Testing robot alignment is a qualitative procedure to ensure that the robot is interfacing properly to the cassette,  
30 pre-aligner, and XYZ stage. The only reason that robot

alignment should change is if one of the components in question suffers a mechanical impact of some kind.

#### **E.3.9.1 Robot to Cassette Alignment Check**

5 Bring up the Manual Robot Control window by clicking on  
'Maintenance' on the top menu bar, and then on  
'ManualRobotControl'. With a full cassette of wafers on the  
cassette platform, highlight each of the cassette slots in the  
10 source/destination list and use the 'Robot Get' and 'Robot Put'  
buttons to load and replace wafers from all 25 slots. Use  
extreme care with the Manual Robot Control window because  
the safeguards in the normal operating software are not in  
effect when using this window: for example, the system will  
not prevent the user from attempting to put two wafers in the  
15 same slot from the Manual Robot Control window.  
As the robot's end effector goes in and out of the cassette to  
pick up and replace each wafer, observe that the wafers are  
moved smoothly and that there is no rubbing between the  
wafer and the cassette.

20

#### **E.3.9.2 Robot to Pre-Aligner Alignment Check**

Using the Manual Robot Control window in the same manner  
as described in section 3.9.1, execute a 'Robot Get' operation  
to retrieve a wafer from the cassette, then execute a 'Robot  
25 Put' with the pre-aligner highlighted in the source/destination  
list. Observe the motion of the robot as it enters the  
pre-aligner and be sure that there is no rubbing between the  
robot and any parts of the pre-aligner. Check that the robot  
places the wafer gently on the pre-aligner chuck and picks it  
30 back up in the same manner.

**E.3.9.3 Robot to XYZ Stage Alignment Check**

This operation proceeds in exactly the same manner as section 3.9.2, except that the 'Stage' selection is highlighted in the source/destination list. As the robot enters the optics chamber to drop off or retrieve a wafer, be sure that there is no rubbing between the robot or wafer and the wafer door facia plate or the chuck on the XYZ stage.

**E.3.10 Changing the Halogen Bulb**

Changing the halogen bulb requires taking off the top cover of the machine and manipulating the lamp house on the optics plate. The top cover is interlocked with the laser power supply and should cause the laser to shut off when the cover is removed. It cannot be turned back on until the interlocks are closed and power is cycled. The top cover interlocks are defeatable, but under no circumstances are customer service personnel authorized to defeat the laser interlocks. If at any time, the customer service personnel are aware of laser radiation (intense blue-green light) while the top cover of the instrument is removed, the instrument should be immediately powered off and a service call placed to the Ultrapointe factory.

With the top cover removed, the customer service personnel can have access to components that are at a high voltage. There is no danger of contacting high voltage components unless the operator disassembles components that have nothing to do with the procedure at hand. Use of controls or adjustments, or performance of procedures other than those specified here may result in the customer service personnel being exposed to unnecessary electrical, mechanical, or radiation hazards.

The halogen bulb used for white light illumination gets extremely hot during normal operation. **Take care before handling the bulb or the lamp house that sufficient time has been allowed after powering off the system to allow the lamp components to cool.**

5

The procedure for replacing the halogen bulb is as follows:

- 1) Turn off the machine using the shutdown procedure described in section 3.2.
- 10 2) Release the top cover of the machine by unscrewing the four quarter-turn captive screws, one at each corner of the top cover. Carefully lift off the cover and place it aside in a clean and safe place.
- 3) Wait fifteen minutes for the lamp house to cool down.
- 15 4) Identify the white light lamp house (see figure 2): it is a roughly rectangular housing, off-white in color, with black louvers facing up to allow heat to escape. The lamp house is on the optics plate immediately beneath the top cover. The lamp socket is a cylindrical object that fits into the side of the lamp house and carries the power cord for the lamp.
- 20 5) There are two thumbscrews attached to the lamp socket. One of them protrudes radially from the socket and is used to adjust and lock the angular position of the lamp. The other protrudes axially from the socket and is used to adjust the axial position of the lamp. Loosen the radial thumbscrew and observe that the socket can now be rotated within the lamp house.
- 25 6) Rotate the lamp house clockwise until it contacts a stop, then pull the socket out axially.
- 30

- 7) The bulb is held in its socket by gold contacts. Press down on the long lever arms of the gold contacts to release the bulb; pull the bulb out of its socket while the levers are depressed.
- 5 8) Install a new bulb (generic part number: FCR 12V, Ultrapointe part number 000616) by depressing the socket levers while pressing the bulb into the socket. Take care not to touch the bulb because any contamination on the outside of the bulb will lead to premature failure.
- 10 9) Replace the socket and bulb in the lamp house by reversing the process of step 5.
- 10) Snug down the radial thumbscrew to hold the socket in place and turn on the instrument using the procedure of section 3.1.
- 15 11) Load a patterned wafer into the machine using normal operating procedures. The autofocus operations will not work because the laser interlocks are open.
- 12) Manually focus on the wafer using the joystick and the lowest power objective available.
- 20 13) Loosen the radial thumbscrew on the lamp socket and rotate the lamp socket to obtain the most uniform illumination of the sample, as indicated by the white light image of the wafer. Snug down the radial thumbscrew.
- 25 14) Adjust the axial thumbscrew to obtain the most uniform illumination of the sample. These adjustments are moving the position of the lamp filament with respect to the optics train.
- 30 15) Repeat steps 13 and 14 until the illumination is optimized. Tighten the radial thumbscrew.



- 16) Replace the top cover of the machine and screw in the four quarter turn fasteners that hold it to the machine frame. This operation closes the laser interlocks.
- 5 17) Turn off the machine as described in section 3.2, then restart it as described in section 3.2. Cycling power should bring back the laser back up.

#### E.3.11 Ultrapointe Factory PM

10 Ultrapointe factory PM procedures are to be performed only by factory authorized personnel. The procedures include cleaning and re-alignment of the optics train, which is a sensitive process that is critical to the performance of the instrument. The procedures involved are described separately in internal

15 Ultrapointe documents.

The following form should be used to log regular PM procedures.

	Date	Procedure Description	Initiator	Comment
5				
10				
15				
20				
25				

**E.5. Problem Report Worksheet**

Use the following form to report bugs and problems to the  
Ultrapointe factory. Fill in as much as possible in the sections  
that are not shaded. Mail or fax the worksheet to Ultrapointe  
at:

5

Ultrapointe Corp.  
163 Baypointe Pkwy  
San Jose, CA 95134  
tel: 408 894 7080  
fax: 408 894 7090  
Attn: Bug List Mgr.

10

	<b>Date:</b>		
5	<b>Reported by:</b>		<b>Company:</b>
10	<b>Description of problem:</b>		
15	<b>Operating Conditions:</b>		
20	<b>Software Revision/Date:</b>		
	<b>Screen in use:</b>		
25	<b>How did you recover?</b>		
	<b>Recommendation:</b>		
30	<b>Assigned to:</b>		<b>Date:</b>

What is claimed is:

1. A method for automatic focusing of a confocal microscope, said microscope transmitting a laser beam through a lens and through a pinhole of a spatial  
5 filter to a surface of a target, said method comprising the steps of:

moving a target relative to a lens of said microscope, said target being moved in a predetermined first direction through a  
10 predetermined first distance;

generating an electronic focus signal during movement of said target, the magnitude of said electronic focus signal being a function of the magnitude of light from said laser beam that is  
15 reflected from said surface and passes through said pinhole;

recording a plurality of first values of the magnitude of said electronic focus signal during the movement of said target in said first  
20 direction; and

stopping said movement of said target when said target has travelled through said first predetermined distance, said first predetermined distance being greater than the depth of focus of  
25 said lens, said depth of focus being larger than the distance between two adjacent positions at which said first values are recorded, wherein after stopping, said target is at a position other than a position at which said microscope is  
30 focused.

2. The method of Claim 1, further comprising, after said stopping step the steps of:

calculating a first estimate of a focus  
35 position at which the microscope is focused by using said plurality of first values;

moving said target relative to said lens in a second direction, said second direction being opposite said first direction;

generating said electronic focus signal during movement of said target in said second direction;

recording a plurality of second values of said magnitude of said electronic focus signal during movement of said target in said second direction; and

stopping said movement of said target in said second direction when said target reaches a stop position, wherein said stop position is derived from said first estimate.

3. The method of Claim 2 further comprising: calculating a second estimate of the position at which said microscope is focused by using said plurality of second values; and

moving said target to said second estimate.

4. The method of Claim 2 wherein said calculating step comprises:

summing all of said plurality of first values to generate a first sum;

dividing said full-sum by 2 to generate a half-sum; and

repeatedly summing a sequence of said plurality of first values to generate a plurality of partial sums until one of said plurality of partial sums exceeds said half-sum.

5. The method of Claim 4 further comprising using the elevation at which said partial sum exceeds said half-sum as said first estimate.

6. The method of Claim 4 further comprising determining a starting position and a stopping position, said stopping position and said starting position being on opposite sides of an estimated focus position, wherein said starting position is at a distance of at least a depth of focus of said lens from said estimated focus position and further wherein said stopping position is at said distance.

7. The method of Claim 6 further comprising a start-up move for moving said target to said starting position if said target is away from said starting position.

8. The method of Claim 2, wherein said target is moved through said first distance at a first velocity and through said second distance at a second velocity, said second velocity being less than said first velocity.

9. The method of Claim 1, further comprising calculating a first estimate of the position at which the microscope is focused by using said plurality of first digitized values.

10. The method of Claim 9, further comprising moving said target to said first estimate.

11. The method of Claim 1, wherein said target is moved by a stepper motor or by a piezo electrically driven element.

12. The method of Claim 1, wherein said step of generating further comprises:  
measuring said intensity of light reflected from said surface of said target through said pin

hole in a direction reverse to the incident path  
of said laser beam; and  
transforming said measured intensity into  
said electronic focus signal.

5

13. The method of Claim 12, wherein said step of  
transmitting further comprises moving said laser beam  
to define a line on said surface of said target.

10

14. The method of Claim 1 further comprising the  
step of setting up gain and adjusting the gain if said  
target moves through said first distance without said  
absolute value of said magnitude of said electronic  
focus signal being a predetermined optimal value.

15

15. The method of Claim 1 wherein said recording  
occurs at periodic intervals.

16. The method of Claim 1 wherein the majority of  
said first digitized values are non-zero.

20

17. A method for positioning a target close to a  
focus position in a confocal microscope, said confocal  
microscope generating an electronic focus signal having  
a magnitude proportional to light reflected by said  
target, said electronic focus signal having a peak that  
indicates a focus condition, said method comprising the  
steps of:

25

positioning said target such that said focus  
position is within a first range of motion of said  
target;

30

moving said target a plurality of first steps  
through said first range of motion of said target,  
each of said first steps having a first width;

35

measuring a plurality of first strengths of  
said electronic focus signal at each of said first



steps;

calculating a first estimate of the position  
at which said confocal microscope is focused by  
using said plurality of first strengths;

5 moving said target a plurality of second  
steps through a second range of motion, wherein  
said second range of motion includes second steps  
above and below said first estimate, each of said  
second steps having a second width;

10 measuring a plurality of second strengths of  
said electronic focus signal at each of said  
second steps;

calculating a second estimate of the position  
at which said confocal microscope is focused by  
15 using said plurality of second strengths;

moving said target to said second estimate.

18. The method of Claim 17, wherein said second  
range of motion is shorter than said first range of  
20 motion.

19. The method of Claim 17, wherein said second  
width is shorter than said first width.

25 20. The method of Claim 17, wherein said first  
width is smaller than width of said focus signal.

21. The method of Claim 17, wherein said target  
is moved in a first direction through said plurality of  
30 first steps and in a second direction through said  
plurality of second steps, said second direction being  
identical to said first direction.

22. The method of Claim 17, wherein said target  
35 is moved in a second direction through said plurality  
of second steps and in a second direction through said

plurality of second steps, said second direction being opposite to said first direction.

23. A method for automatically focusing a  
5 microscope that transmits a laser beam through a lens  
and through a pinhole of a spatial filter to a surface  
of a target, said method comprising the steps of:
- moving a target relative to a lens of said  
microscope in a first direction;
  - 10 generating an electronic focus signal for a  
plurality of points in an area of said target at  
least three points of said plurality of points  
defining a plane;
  - recording as a first value the largest value  
15 of said electronic focus signal for said plurality  
of points for each position in a plurality of  
first positions of said target with respect to  
said lens; and
  - stopping said movement of said target when  
20 said target has travelled through a first  
predetermined distance, wherein said first  
predetermined distance is at a position other than  
a position at which said microscope is focused.
24. The method of Claim 23, further comprising,  
25 after said stopping step the steps of:
- calculating a first estimate of the position  
at which the microscope is focused by using a  
plurality of said first values;
  - 30 moving said target relative to said lens in a  
second direction, said second direction being  
opposite said first direction;
  - generating said electronic focus signal for a  
plurality of points in said area of said target,  
35 at least three points of said plurality of points  
defining a plane;

recording as a second value, the largest value of said electronic focus signal for said plurality of points for each position in a plurality of second positions of said target with respect to said lens; and

stopping said movement of said target in said second direction when said target reaches a stop position, wherein said stop position is derived from said first estimate.

25. The method of Claim 23 further comprising: calculating a second estimate of the position at which said microscope is focused by using said plurality of second values; and moving said target to said second estimate.

26. The method of Claim 23 wherein said target is moved through said first distance at a first velocity and through said second distance at a second velocity, said second velocity being less than said first velocity.

27. The method of Claim 23 wherein said second direction is identical to said first direction.

28. The method of Claim 23 wherein said second direction is opposite to said first direction.

29. An apparatus for automatic focusing of a confocal microscope, said microscope having a laser, a lens and a spatial filter, said apparatus comprising: means for moving a target relative to said lens of said microscope, said target being moved in a predetermined first direction through a predetermined first distance; means for generating an electronic focus

signal during movement of said target, the  
magnitude of said electronic focus signal being a  
function of the magnitude of light from said laser  
beam that is reflected from said surface and  
passes through said spatial filter;

means for recording a plurality of first  
values of the magnitude of said electronic focus  
signal during the movement of said target in said  
first direction, said means for recording being  
coupled to said means for generating; and

means for stopping said movement of said  
target when said target has travelled through said  
first predetermined distance, said first  
predetermined distance being greater than the  
depth of focus of said lens, said depth of focus  
being larger than the distance between two  
adjacent positions at which said first values are  
recorded, wherein after stopping, said target is  
at a position other than a position at which said  
microscope is focused.

30. The apparatus of Claim 29 further comprising  
a peak detector coupled to said means for recording,  
wherein said peak detector supplies the maximum value  
of said electronic focus signal that is encountered  
while a laser beam is scanned across an area of said  
target, while said target is held stationary.

31. A method for automatically focusing a  
microscope, comprising the steps of:

continuously moving a target relative to a  
lens of the microscope, the target being moved in  
a first direction through a first distance;

continuously generating an electronic focus  
signal during movement of the target, the  
magnitude of the electronic focus signal being a

function of the magnitude of light reflected from  
a surface of the target through an optical path of  
the microscope;

5 continuously comparing the absolute value of  
the magnitude of the electronic focus signal to a  
threshold value; and

stopping the movement of the target when the  
absolute value of the magnitude of the electronic  
focus signal exceeds the threshold value;

10 wherein the step of continuously  
generating further comprises:

transmitting a laser beam through the lens to  
a surface of the target;

15 measuring the intensity of the light  
reflected from the surface of the target through  
the optical path of the microscope; and

transforming the measured intensity into the  
electronic focus signal;

20 wherein the step of transmitting further  
comprises moving the laser beam to define an  
area on the surface of the target.

## AMENDED CLAIMS

[received by the International Bureau on 26 June 1995 (26.06.95);  
original claims 1,17,23 and 29 amended; remaining claims  
unchanged (8 pages)]

1. A method for automatic focusing of a confocal  
microscope, said microscope transmitting a laser beam  
through a lens and through a pinhole of a spatial  
5 filter to a surface of a target, said method comprising  
the steps of:

moving a target relative to a lens of said  
microscope, said target being moved in a  
predetermined first direction through a  
10 predetermined first distance;

generating an electronic focus signal during  
movement of said target, the magnitude of said  
electronic focus signal being a function of the  
magnitude of light from said laser beam that is  
15 reflected from said surface and passes through  
said pinhole;

recording a plurality of first values of the  
magnitude of said electronic focus signal during  
the movement of said target in said first  
20 direction; and

stopping said movement of said target when  
said target has travelled through said first  
predetermined distance, said first predetermined  
distance being greater than the depth of focus of  
25 said lens, said depth of focus being larger than  
the distance between two adjacent positions at  
which said first values are recorded, wherein  
after stopping, said target is at a position other  
than a position at which said microscope is  
30 focused;

wherein the step of generating  
comprises:

transmitting a laser beam through the lens to  
a surface of the target;

35 measuring the intensity of the light  
reflected from the surface of the target through  
the optical path of the microscope; and

transforming the measured intensity into the electronic focus signal;

wherein the step of transmitting further comprises moving the laser beam to define an area on the surface of the target.

2. The method of Claim 1, further comprising, after said stopping step the steps of:

calculating a first estimate of a focus position at which the microscope is focused by using said plurality of first values;

steps;

calculating a first estimate of the position at which said confocal microscope is focused by using said plurality of first strengths;

5 moving said target a plurality of second steps through a second range of motion, wherein said second range of motion includes second steps above and below said first estimate, each of said second steps having a second width;

10 measuring a plurality of second strengths of said electronic focus signal at each of said second steps;

calculating a second estimate of the position at which said confocal microscope is focused by using said plurality of second strengths; and

15 moving said target to said second estimate; wherein the step of generating comprises:

transmitting a laser beam through the lens to a surface of the target;

20 measuring the intensity of the light reflected from the surface of the target through the optical path of the microscope; and

transforming the measured intensity into the electronic focus signal;

25 wherein the step of transmitting further comprises moving the laser beam to define an area on the surface of the target.

30 18. The method of Claim 17, wherein said second range of motion is shorter than said first range of motion.

35 19. The method of Claim 17, wherein said second width is shorter than said first width.

20. The method of Claim 17, wherein said first



width is smaller than width of said focus signal.

21. The method of Claim 17, wherein said target  
is moved in a first direction through said plurality of  
5 first steps and in a second direction through said  
plurality of second steps, said second direction being  
identical to said first direction.

22. The method of Claim 17, wherein said target  
10 is moved in a second direction through said plurality  
of second steps and in a second direction through said

plurality of second steps, said second direction being opposite to said first direction.

23. A method for automatically focusing a  
5 microscope that transmits a laser beam through a lens and through a pinhole of a spatial filter to a surface of a target, said method comprising the steps of:

moving a target relative to a lens of said  
microscope in a first direction;

10 generating an electronic focus signal for a plurality of points in an area of said target at least three points of said plurality of points defining a plane;

15 recording as a first value the largest value of said electronic focus signal for said plurality of points for each position in a plurality of first positions of said target with respect to said lens; and

20 stopping said movement of said target when said target has travelled through a first predetermined distance, wherein said first predetermined distance is at a position other than a position at which said microscope is focused;

25 wherein the step of generating comprises:

transmitting a laser beam through the lens to a surface of the target;

30 measuring the intensity of the light reflected from the surface of the target through the optical path of the microscope; and

transforming the measured intensity into the electronic focus signal;

35 wherein the step of transmitting further comprises moving the laser beam to define an area on the surface of the target.

24. The method of Claim 23, further comprising,

after said stopping step the steps of:

calculating a first estimate of the position  
at which the microscope is focused by using a  
plurality of said first values;

5 moving said target relative to said lens in a  
second direction, said second direction being  
opposite said first direction;

generating said electronic focus signal for a  
plurality of points in said area of said target,  
10 at least three points of said plurality of points  
defining a plane;

signal during movement of said target, the  
magnitude of said electronic focus signal being a  
function of the magnitude of light from said laser  
beam that is reflected from said surface and  
passes through said spatial filter;

means for recording a plurality of first  
values of the magnitude of said electronic focus  
signal during the movement of said target in said  
first direction, said means for recording being  
coupled to said means for generating; and

means for stopping said movement of said  
target when said target has travelled through said  
first predetermined distance, said first  
predetermined distance being greater than the  
depth of focus of said lens, said depth of focus  
being larger than the distance between two  
adjacent positions at which said first values are  
recorded, wherein after stopping, said target is  
at a position other than a position at which said  
microscope is focused;

wherein the means for generating  
comprises:

means for transmitting a laser beam through  
the lens to a surface of the target;

means for measuring the intensity of the  
light reflected from the surface of the target  
through the optical path of the microscope; and

means for transforming the measured intensity  
into the electronic focus signal;

wherein the means for transmitting  
further comprises means for moving the laser  
beam to define an area on the surface of the  
target.

30. The apparatus of Claim 29 further comprising  
a peak detector coupled to said means for recording,  
wherein said peak detector supplies the maximum value

of said electronic focus signal that is encountered while a laser beam is scanned across an area of said target, while said target is held stationary.

5           31. A method for automatically focusing a microscope, comprising the steps of:

continuously moving a target relative to a lens of the microscope, the target being moved in a first direction through a first distance;

10           continuously generating an electronic focus signal during movement of the target, the magnitude of the electronic focus signal being a




FIG. 1

RECTIFIED SHEET (RULE 91)




FIG. 2A




FIG. 2B




FIG. 2C






FIG. 3A




FIG. 3B




FIG. 3C




FIG 5A-1





FIG 5A-2

409



RECTIFIED SHEET (RULE 91)




FIG 5A-3






FIG. 5B-1




FIG. 5B-2





FIG. 5B-3






FIG. 5B-4




FIG. 5B-5





FIG. 5B-6

KEY TO FIG. 5B

FIG. 5B-1	FIG. 5B-2	FIG. 5B-3
FIG. 5B-4	FIG. 5B-5	FIG. 5B-6

FIG. 5B







FIG. 5C-2











FIG. 5D-2






FIG. 5D-4


FIG. 5E-1

KEY TO FIG. 5E

FIG. 5E-1

FIG. 5E-2

FIG. 5E






FIG. 5E-2





FIG. 6

27 / 76


RECTIFIED SHEET (RULE 91)



**FIG. 7**



28 / 76

RECTIFIED SHEET (RULE 91)



29 / 76 **FIG. 8A**

RECTIFIED SHEET (RULE 91)



30 / 76

203

RECTIFIED SHEET (RULE 91)




FIG. 9

31 / 76

RECTIFIED SHEET (RULE 91)




FIG. 10

RECTIFIED SHEET (RULE 91)






FIG. 12-1

FIG. 12

KEY TO FIG. 12

FIG. 12-1      FIG. 12-2






FIG. 12-2





NOTE: TIME  
SCALES ARE  
DIFFERENT







FIG. 13C







FIG. 14-2






FIG. 14-4

FIG. 15-1




FIG. 15-2





FIG. 15




FIG. 15-3



404

FIG. 15-4






FIG. 16-1





FIG. 17-1






FIG. 17-2

# INTERNATIONAL SEARCH REPORT

International application No.  
PCT/US95/00665

## A. CLASSIFICATION OF SUBJECT MATTER

IPC(6) :G01J 1/20

US CL :250/201.2, 201.3, 204, 216; 359/368, 381

According to International Patent Classification (IPC) or to both national classification and IPC

## B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

U.S. : 250/201.2, 201.3, 204, 216; 359/368, 381

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched  
NONE

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

APS

S (CONFOCAL(2A)MICROSCOP?) AND FOCUS? AND FILTER AND (APERTURE? OR PINHOLE? OR OPENING?)

## C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	US, A, 4,342,905 (FUJII ET AL.) 03 AUGUST 1982, (see entire document)	1-6
Y	US, A, 4,844,617 (KELDERMAN ET AL.) 04 JULY 1989 (see entire document)	1-30
X	US, A, 5,122,648 (COHEN ET AL.) 16 JUNE 1992 (see entire document)	1-10
A,P	US, A, 5,306,902 (GOODMAN) 26 APRIL 1994 (see entire document)	1-32

☐ Further documents are listed in the continuation of Box C.

☐ See patent family annex.

* Special categories of cited documents:	*T	later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
*A* document defining the general state of the art which is not considered to be of particular relevance	*X*	document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
*E* earlier document published on or after the international filing date	*Y*	document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
*L* document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)	*Z*	document member of the same patent family
*O* document referring to an oral disclosure, use, exhibition or other means		
*P* document published prior to the international filing date but later than the priority date claimed		

Date of the actual completion of the international search

27 APRIL 1995

Date of mailing of the international search report

03 MAY 1995

Name and mailing address of the ISA/US  
Commissioner of Patents and Trademarks  
Box PCT  
Washington, D.C. 20231

Facsimile No. (703) 305-3230

Authorized officer

STEPHON B. ALLEN

Telephone No. (703) 305-4828

Form PCT/ISA/210 (second sheet)(July 1992)\*




FIG. 45




FIG. 46




FIG. 47



Defect Map	38102002.map	Operator	bruce	Status <input type="text"/>
Product ID	Production.7b198a	SiteNumber	3	
Lot ID	1242381	Pos,X/Y[um]	48225 -20930	
ProcessID	siemsADD	DieNum, X/Y	-3	
Description:			<input type="text"/>	
Defect code:	<input type="text" value="0"/>	<input type="button" value="List"/>		

Deskew1,XY um	- 56770	-6700	<input type="button" value="DeskewOFF"/>	Deskew <input type="text"/>
Deskew2	49318	-6700	<input type="button" value="CalcMapOffset"/>	
Deskew1 to deskew2	=29.0 dle		<input type="button" value="Deskew"/>	
Dsk3,Registration	-42392	11383		
Defect Map Offset	0	0		

Filename:	38102002.map	WaferSize:	125mm	EditMap <input type="text"/>
Wafer#:	<input type="text" value="48"/>	SlotNum:	6	
LotID:	<input type="text" value="1242381"/>			
Recipe:	<input type="text" value="9012"/>			
FileDescription:	<input type="text"/>			
				<input type="button" value="Accept"/>

DefectList						DefList <input type="text"/>	
2.	40971	-24758	19	73	0		<input type="button" value="GoToSelected"/>
3.	48122	-20819	19	73	0		<input type="button" value="GoToNext"/>
4.	10476	3232	25	73	0		
5.	6951	25166	25	73	0		

Text and Status Windows

**FIG. 44**

7 5 / 7 6

RECTIFIED SHEET (RULE 91)

Library  
Window

**FIG. 43**

The screenshot displays the Ultrapointe file selection interface. At the top, there is a status bar with icons for a window, a disk, and a help button. Below this, a 'FileSelection' window is open, showing a list of files and a 'User Name' field. The files listed are:


- 1211.bmp
- 1241.bmp
- 1291.bmp
- 3561.bmp
- 4101.bmp
- 4101\_1.bmp
- 4661.bmp
- 4661\_1.bmp
- 5161.bmp
- 5161\_1.bmp
- plugbyte.bmp
- plugfloat.bmp

The 'User Name' field contains the text 'bruce'. Below the file list, there are two buttons: 'Load' and 'Quit'. The background shows a desktop with a 'File' menu and a 'Help' button.

74176

**RECTIFIED SHEET (RULE 91)**

FIG. 42






FIG. 41

<div> <input type="checkbox"/> UV         </div>		<div> <input type="checkbox"/> <input type="checkbox"/> </div>		<div> <input type="checkbox"/> </div>	
<div> <div>File</div> <div>Screen</div> <div>Setup</div> <div>User</div> <div>Config</div> <div>Maintenance</div> <div>PrtScreen</div> <div>Help</div> </div>					
<div> <div>2D</div> <div>Library</div> <div>Wafer</div> <div>Cassette</div> </div>		<div> <div>Laser</div> <div>3D</div> <div>2D</div> </div>		<div> <div>3D Surface</div> <div>LJ Control</div> </div>	
<div> <div>2DSlices</div> <div>LJ Control</div> </div>		<div> <div>Microscope</div> <div>Objectives: 20</div> <div>AutoFocus</div> <div>CurrentPos</div> <div>Joystick</div> </div>			
<div> <div>Ultrapointe</div> <div>Last</div> <div>2d</div> <div>Next</div> <div>15:1 Jun 1,93</div> </div>					

Microscope  
Window

FIG. 40

<input type="checkbox"/> UV		File Screen Setup User Config Maintenance PrintScreen Help			
2D	3D	LaserScan Objective: 20			
2D	Library	Current Z			
Wafer	3D	S F			
Cassette		AutoFocus Start LaserScan Stage AF			
Microscope	L	Laserintensity 10 XY Res 25nm LaserLine 488 Slices 90			
2DSlices	L Control	StepSize: x12nm 3			
		Camintensity Camerafilter Field Select Vol Acq Brightfield Dimensions 256x256			

## Laser Scan Window

**FIG. 39**


<input type="checkbox"/> UV <input type="radio"/> <input type="checkbox"/>		File Screen Setup User Config Maintenance <u>P</u> rtScreen		Help	
<input type="checkbox"/> 2D <input type="checkbox"/> 3D <input type="checkbox"/> Cassette <input type="checkbox"/> Laser	<input type="checkbox"/> 3D <input type="checkbox"/> Library <input type="checkbox"/> uScope	cy342m2f.map.Sin,125mm X[um]: <input type="text" value="0"/> Y[um]: <input type="text" value="0"/>	Select GoTo		
2DSlices <input type="checkbox"/> L <input type="checkbox"/> Control		2DSlices <input type="checkbox"/> L <input type="checkbox"/> Control		Ultrapointe	
				Last <input type="text" value="2d"/> Next <input type="text" value="16:48 Apr 26,93"/>	

FIG. 38

**FIG. 37**

68 / 76  
RECTIFIED SHEET (RULE 91)



<input type="checkbox"/> UV <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>		<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
File   Screen   Setup   User   Config   Maintenance   PriScreen   Help			
<div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">2D</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">2D</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Wafer</div> <div style="border: 1px solid black; padding: 2px;">Cassette</div>	<div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">3D</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Library</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">3D</div>	<div style="border: 1px solid black; padding: 5px;"> <div style="display: flex; justify-content: space-between;"> <span>LaserScan</span> <span>L Objective: 20</span> </div> <div style="display: flex; justify-content: space-between; align-items: flex-start;"> <div style="width: 40%;"> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Current Z</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">S   F</div> <div style="border: 1px solid black; padding: 2px;">0</div> </div> <div style="width: 60%;"> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">AutoFocus</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Start</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">LaserScan</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Stage AF</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Laserintensity</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">10</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">XY Res 25nm</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">LaserLine 488</div> <div style="border: 1px solid black; padding: 2px;">Slices</div> </div> </div> <div style="display: flex; justify-content: space-between; align-items: flex-start; margin-top: 10px;"> <div style="width: 40%;"> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Camintensity</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">1</div> </div> <div style="width: 60%;"> <div style="display: flex; justify-content: space-between;"> <span>Camerafilter</span> <span>Field Select</span> <span>Vol Acq</span> </div> <div style="display: flex; justify-content: space-between; align-items: flex-start;"> <div style="width: 40%;"> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">All pass</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">BrightField</div> <div style="border: 1px solid black; padding: 2px;">Dimensions</div> </div> <div style="width: 60%;"> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">StepSize: x12nm</div> <div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">3</div> </div> </div> </div> </div> </div>	
<div style="border: 1px solid black; padding: 5px;"> <div style="display: flex; justify-content: space-between;"> <span>Microscope</span> <span>L</span> </div> <div style="border: 1px solid black; height: 100px; width: 100%;"></div> </div>		<div style="border: 1px solid black; padding: 5px;"> <div style="display: flex; justify-content: space-between;"> <span>2DSlices</span> <span>L Control</span> </div> <div style="border: 1px solid black; height: 100px; width: 100%;"></div> </div>	
ULTRAPOINTE		Last 2D Next 16:44 Apr 26,93	

FIG. 36

<div style="display: flex; align-items: center;"> <input style="width: 20px; height: 20px; margin-right: 5px;" type="checkbox"/> <span>LaserScan Calibration</span> </div>																																		
<p>Scan parameters:</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">Obj1: XY</td> <td style="width: 15%;">1800 3060,</td> <td style="width: 15%;">900 1700,</td> <td style="width: 15%;">450 850,</td> <td style="width: 15%;">225 425</td> </tr> <tr> <td>Obj2:</td> <td>1700 3400,</td> <td>850 1700,</td> <td>425 850,</td> <td>212 425</td> </tr> <tr> <td>Obj3:</td> <td>3400 3400,</td> <td>1700 1700,</td> <td>850 850,</td> <td>425 425</td> </tr> <tr> <td>Obj4:</td> <td>3400 3400,</td> <td>1700 1700,</td> <td>850 850,</td> <td>425 425</td> </tr> <tr> <td>Obj5:</td> <td>3400 3400,</td> <td>1700 1700,</td> <td>850 850,</td> <td>425 425</td> </tr> </table>										Obj1: XY	1800 3060,	900 1700,	450 850,	225 425	Obj2:	1700 3400,	850 1700,	425 850,	212 425	Obj3:	3400 3400,	1700 1700,	850 850,	425 425	Obj4:	3400 3400,	1700 1700,	850 850,	425 425	Obj5:	3400 3400,	1700 1700,	850 850,	425 425
Obj1: XY	1800 3060,	900 1700,	450 850,	225 425																														
Obj2:	1700 3400,	850 1700,	425 850,	212 425																														
Obj3:	3400 3400,	1700 1700,	850 850,	425 425																														
Obj4:	3400 3400,	1700 1700,	850 850,	425 425																														
Obj5:	3400 3400,	1700 1700,	850 850,	425 425																														
<p>Obj:</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;"> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>20</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div> </td> <td style="width: 50%; text-align: center;"> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>100</span> <span>50</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(-45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>25</span> <span>12.5</span> </div> </td> </tr> </table>					<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>20</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div>	<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>100</span> <span>50</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(-45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>25</span> <span>12.5</span> </div>	<p>Resol:</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;"> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>3400</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div> </td> <td style="width: 50%; text-align: center;"> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>3400</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(-45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div> </td> </tr> </table>					<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>3400</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div>	<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>3400</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(-45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div>																					
<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>20</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div>	<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>100</span> <span>50</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(-45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>25</span> <span>12.5</span> </div>																																	
<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>3400</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div>	<div style="display: flex; justify-content: space-between; padding: 5px;"> <span>3400</span> <span>100</span> </div> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(-45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div> <div style="display: flex; justify-content: space-between; padding: 5px;"> <span>0</span> <span>40</span> </div>																																	
<p>X Amplitude:</p> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div>					<p>Y Amplitude:</p> <div style="border: 1px solid black; height: 20px; width: 100%; position: relative;"> <div style="position: absolute; left: 0; top: 0; bottom: 0; width: 100%; background: repeating-linear-gradient(-45deg, transparent, transparent 2px, black 2px, black 4px);"></div> </div>																													
<div style="display: flex; justify-content: space-around; width: 100%;"> <span>Save</span> <span>Quit</span> </div>																																		

**FIG. 35**

Lon Diagnostics			
WhiteLightPwr	TurretPos	FieldSelect	
0	Home	BrightField	
LaserPower(mW)	LaserSource		
0	ALL		
LaserCurrent	ChuckVacuum		
0	None		
LaserOnOff	Laser/WhiteLight		
Off	Whitelight		
LaserRunIdle	LaserShutter		
Idle	Closed		
LaserAtten	WaferDoor		
High	Close		
CameraFilter	PageScanCtrl	FastZPosCmd	OutputNWValue:
All pass	4	0	
PageScan	PMTAmp...	FastZ...	PollNetVar
			Default
			Quit
			SendToLon

Nodes	NetVariables
1 LaserPwrController	1 SW Ver
2 CameraFilterWhite	2 ResetStatus
3 LaserLineFilter	3 ModeCmd
4 LaserAttenFilter	4 LaserStatus
5 Bright/DarkField	5 LaserOn/Off
6 PageScanner	6 Run/Idle
7 FastZ Controller	7 Light/Curr
8 PMT Amp	8 CurrCmd
9 Turret/optics	9 PowerCmd
10 Cassette sw	10 WhiteLpwr
11 Solenoid Ctl	11 SpareCmd
	12 Interlock
	13 CurrCheck
	14 PwrCheck
	15 WLPwrCheck
	16 SpareChk
	17 HeadTempMon
	18 SpareTemp
	19 LaserPwrMon
	20 CurrMonitor

NetVarValue: 0

FIG. 34

**File Transfer**

From: System To: System File Type: Surface

Odptrnch.sur	Odptrnch.sur
124_0147b.sur	121_0147b.sur
123_sur	123.sur
124_0147bb.sur	124_0147bb.sur
127_0147bb.sur	127_0147bb.sur
129_0147b.sur	129_0147b.sur
129_0147g.sur	129_0147g.sur
356_6321b.sur	356_6321b.sur
410_8131b1.sur	410_8131b1.sur
410_8131b2.sur	410_8131b2.sur
466_8131b1.sur	466_8131b1.sur
466_8131b2.sur	466_8131b2.sur
516_9028b1.sur	516_9028b1.sur
9238_10.sur	9238_10.sur
a0833_2.sur	a0833_2.sur
a0833_3.sur	a0833_3.sur
a0833_4.sur	a0833_4.sur
a0833_5.sur	a0833_5.sur
a9238_1.sur	a9238_1.sur
a9238_10.sur	a9238_10.sur
a9238_11.sur	a9238_11.sur
a9238_2.sur	a9238_2.sur
a9238_3.sur	a9238_3.sur
a9238_4.sur	a9238_4.sur

Eject Floppy

Delete

Move

Copy

Quit

FIG. 32

**UserAccounting**

User Names: Copy

bruce Cut Paste

chris

dick

garylum

hans

john

kelly

kenlee

tim

Name: \_\_\_\_\_

Password: \_\_\_\_\_

^ Operator

√ Administrator


Save

Quit

FIG. 33

Manual Robot Control	
Stage	
Flat Finder	<input type="button" value="Initialize"/>
Slot 25	Robot offline
Slot 24	
Slot 23	
Slot 22	
Slot 21	
Slot 20	
Slot 19	<input type="button" value="Robot Get"/>
Slot 18	
Slot 17	<input type="button" value="Robot Put"/>
Slot 16	
Slot 15	
Slot 14	
Slot 13	<input type="button" value="Stop"/>
Slot 12	
Slot 11	<input type="button" value="Quit"/>
Slot 10	
Slot 9	
Slot 8	
Slot 7	
Slot 6	
Slot 5	
Slot 4	
Slot 3	
Slot 2	
Slot 1	

**FIG. 31**



PrinterConfig


Print Orientation:

☒ Portrait  
☐ Landscape

Number of copies:

2

OK

**FIG. 29**


Manual Stage Control

Stage Positions: X Y Z

Axis:	Command:	
X	Enable Joystick	InitStage
Y	Disable Joystick	Stage Offline
Z	Report Pos	0
	Report Enc	
	% Dist Fr ConfigLimit	
	WakeUnit	
	Stop	
	XY To Load/unload	
	XYRandom move	
	Rel Move To	StopAll
	Abs Move To	Exec
	Set Final Speed	
	Set Acc	Quit
	Move to Home	

**FIG. 30**

DefectCodeTable

Defect Code Table	Default
1 Big hole	Copy PasteUp
2 .3 um scratch	Cut PasteDn
3 .5 um scratch	
4 1 um scratch	
5 5 um scratch	Entry:
6 Deep cut	
7 Big pit	
8 Incomplete etch	
9 Large square lump	Open/Save
10 .1 micron stuff	
11 .3 micron stuff	
12 Blurry stuff	
13 Large thin bump	Accept
14 Small thick resist	
15 Dust	Quit
16 Hair	

FIG. 27

Login

User Names: Password:

bruce  
chris  
dick  
garylum  
hans  
john  
kelly  
kenlee  
tim

^

Login

Logout

Quit

FIG. 28




FIG. 25





FIG. 26





**FIG. 24-3**

KEY TO FIG 24

FIG. 24-1	FIG. 24-2	FIG. 24-3
-----------	-----------	-----------





FIG. 24-2







FIG. 23

**FIG. 22**




FIG. 21





FIG. 20-2

FIG. 19






FIG. 18









FIG. 17-4



KEY TO FIG. 17

FIG. 17-1	FIG. 17-2
FIG. 17-3	FIG. 17-4

FIG. 17

FIG. 17-3